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1 Norm ASTM F67 (states min. tensile strength of annealed titanium). **2** Data on file for Straumann cold-worked titanium and Roxolid® Implants, MAT 13336, 20131009. **3** Gottlow J et al. : Evaluation of a new titanium-zirconium dental implant: a biomechanical and histological comparative study in the mini pig. *Journal of Clinical Implant Dentistry and Related Research* 2012; 14: 538-545 **4** Wen B et al. : The osseointegration behavior of titanium-zirconium implants in ovariectomized rabbits. *Clin Oral Implants Res.* 2013 Feb 21. **5** Barter S et al. : A pilot study to evaluate the success and survival rate of titanium-zirconium implants in partially edentulous patients: results after 24 months of follow-up. *Clin Oral Implants Res.* 2012 Jul;23(7):873-81 **6** Buser D et al. : Enhanced bone apposition to a chemically modified SLA titanium surface. *J. Dent. Res.* 2004 Jul;83(7):529-33. **7** Schwarz F et al. : Histological and immunohistochemical analysis of initial and early osseous integration at chemically modified and conventional SLA® titanium implants: Preliminary results of a pilot study in dogs. *Clin. Oral Impl. Res.* 2007;11(4):481-488. **8** Schwarz F et al. : Histological and immunohistochemical analysis of initial and early subepithelial connective tissue attachment at chemically modified and conventional SLA® titanium implants. A pilot study in dogs. *Clin. Oral Impl. Res.* 2007;11(3):245-455. **9** Schwarz F et al. : Effects of surface hydrophilicity and microtopography on early stages of soft and hard tissue integration at non-submerged titanium implants: An immunohistochemical study in dogs. *J. Periodontol.* 2007;78(11):2171-2184. **10** Schwarz F et al. : Bone regeneration in dehiscence-type defects at chemically modified (SLActive) and conventional SLA titanium: A pilot study in dogs. *J. Clin. Periodontol.* 2007;34(1):78-86. **11** Zöllner et al. : Immediate and early non-occlusal loading of Straumann implants with a chemically modified surface (SLActive®) in the posterior mandible and maxilla: interim results from a prospective multicentre randomized-controlled study. *Clinical Oral Implants Research*, 19(5), 442-450, 2008. **12** Nicolau P et al. : Immediate and early loading of chronically modified implants in posterior jaws: 3-year results from a prospective randomized study. *Clin Implant Dent Relat Res.* 2013 Aug;15(4):600-612 **13** Raghavendra S et al. *Int. J. Oral Maxillofac. Implants.* 2005 May-Jun;20(3):425-31. **14** Lang, NP et al. : Early osseointegration to hydrophilic and hydrophobic implant surfaces in humans. *Clin Oral Implants Res* 22.4 (2011): 349-56 **15** Lekholm U et al. : Patient selection and preparation. *Tissue integration prostheses.* Chicago: Quintessence Publishing Co. Inc. 1985; 199-209 **16** Rupp F et al. : Enhancing surface free energy and hydrophilicity through chemical modification of microstructured titanium implant surfaces. *Journal of Biomedical Materials Research A*, 76(2):323-334, 2006. **17** DeWild M : Superhydrophilic SLActive® implants. Straumann document 151.52, 2005 **18** Maniura K : Laboratory for Materials – Biology Interactions Empa, St. Gallen, Switzerland Protein and blood adsorption on Ti and TiZr implants as a model for osseointegration. EAO 22nd Annual Scientific Meeting, October 17 – 19 2013, Dublin **19** Schwarz F et al. : Bone regeneration in dehiscence-type defects and submerged chemically modified (SLActive®) and conventional SLA® titanium implants: an immunohistochemical study in dogs. *J Clin.Periodontol.* 35.1 (2008): 64– 75. **20** Rausch-fan X et al. : Differentiation and cytokine synthesis of human alveolar osteoblasts compared to osteoblast-like cells (MG63) in response to titanium surfaces. *Dental Materials* 2008 Jan;24(1):102-10. Epub 2007 Apr 27. **21** Schwarz F et al. : Histological and immunohistochemical analysis of initial and early osseous integration at chemically modified and conventional SLA® titanium implants: Preliminary results of a pilot study in dogs. *Clinical Oral Implants Research*, 11(4): 481-488, 2007. **22** Lang, NP et al. : Early osseointegration to hydrophilic and hydrophobic implant surfaces in humans. *Clin Oral Implants Res* 22.4 (2011): 349-56. **23** Raghavendra S et al. : *Int. J. Oral Maxillofac. Implants.* 2005 May-Jun;20(3):425-31. **24** Oates TW et al. : Enhanced implant stability with a chemically modified SLA® surface: a randomized pilot study. *Int. J. Oral Maxillofac. Implants.* 2007;22(5):755-760. **25** Schwarz F et al. : Bone regeneration in dehiscence-type defects at chemically modified (SLActive®) and conventional SLA® titanium implants: a pilot study in dogs. *J Clin.Periodontol.* 34.1 (2007): 78-86

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Institut Straumann AG Peter Merian-weg 12