



















# Tracking sheet for Axiom® cutting instruments usage

This table enables to track the number of uses of cutting instruments from Axiom® Multi Level® and Axiom® 2.8 surgical kits.  
Before first use and after each surgery, all components must be decontaminated and sterilised following the manufacturer's recommendations (cf Axiom® Multi Level® and Axiom® 2.8 User guides).  
For high performance and optimal clinical results, it is recommended that all cutting instruments are limited to **20 uses**.

		Article	Article Number	Number of uses																			
				1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
Axiom®		Round bur	INFB20	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		Pointer drill	OPPO15150	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
Axiom® BL / TL - REG / PX		Initial drill Ø 2.0 S	OPFI20S	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		Initial drill Ø 2.0 L	OPFI20L	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		Step drill Ø 2.0/ 2.4 S	OPFE24S	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		Step drill Ø 2.4/ 3.0 S	OPFE30S	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		Step drill Ø 3.0/ 3.6 S	OPFE36S	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		Step drill Ø 3.6/ 4.2 S	OPFE42S	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		Step drill Ø 4.2/ 4.8 S	OPFE48S	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		Step drill Ø 2.0/ 2.4 L	OPFE24L	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		Step drill Ø 2.4/ 3.0 L	OPFE30L	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		Step drill Ø 3.0/ 3.6 L	OPFE36L	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		Step drill Ø 3.6/ 4.2 L	OPFE42L	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		Step drill Ø 4.2/ 4.8 L	OPFE48L	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
Axiom® BL / TL - REG		Tap Ø 3.4 mm L	OPTA34L	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		Tap Ø 4.0 mm L	OPTA40L	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		Tap Ø 4.6 mm L	OPTA46L	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		Tap Ø 5.2 mm L	OPTA52L	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
Axiom® 2.8		Initial drill Ø 2.0 x 25 mm	OPFI20250	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		Step drill Ø 2.6 x 25 mm	OPFE26250	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		Tap 2.8 Ø 2.8 x 25 mm	OPTA28250	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20