



Uncompromised  
hygiene.  
Ready for action.





**100% high-grade stainless steel cassette for uncompromised performance.**

The unique Straumann® ProClean Cassettes are made solely of stainless steel with no aluminium, plastics, silicones or other materials that could degrade due to frequent, intensive machine-based reprocessing cycles.





# Transform your workflow with the Straumann® ProClean Cassette.

For over six decades, Straumann has made a decisive contribution to progress in dental implantology through pioneering innovation, that redefines the boundaries of possibility for dental professionals and patients.

Whether you are a major hospital SSD<sup>1</sup> or a busy practice, the revolutionary ProClean Cassette will transform your workflow by increasing productivity, efficiency and safety. The ProClean Cassette virtually eliminates labor-intensive manual cleaning – simply put the used instruments back into the cassette during operation and let machines<sup>2</sup> do the rest.



## UNCOMPROMISED HYGIENE

Meets the highest hygiene requirements of dental practices and hospitals.



## SIMPLE WORKFLOW

Simply put the used instruments back into the cassette. Validated, automated reprocessing makes instrument mix-ups a thing of the past.



## INCREASED PRODUCTIVITY


Free up more time for value-added clinical activities.

# Don't compromise on hygiene.

## Rest assured with validated cleaning results.


Designed exclusively for automated reprocessing of dental instruments post operation, the cassette is fully compatible with dental washer-disinfectors to simplify and automate your reprocessing workflow delivering hygienic cleaning results you can rely on.

- Outstanding cleaning results and uncompromised hygiene
- Increased staff safety – handling of sharp contaminated instruments is significantly reduced<sup>3</sup>
- A validated user friendly design that ensures clean and dry instruments every time



THE STATE-OF-THE-ART STRAUMANN® ProClean CASSETTE  
MEETS THE HIGHEST REQUIREMENTS OF DENTAL PRACTICES AND  
HOSPITALS THAT DEMAND UNCOMPROMISED HYGIENE.






All parts of the unique Straumann® ProClean Cassette are made solely of high-grade stainless steel with no aluminum, plastics, silicones or other materials that could degrade due to frequent, intensive machine-based reprocessing cycles.


- 1 Compact, cleaning-friendly instrument holders minimize touch points and keep the drills securely in place
- 2 Instruments organized according to drill sequences. Color-coded workflows with instrument pictograms guide you to simplify the surgical procedure
- 3 After surgery the cassette lid safeguards staff<sup>3</sup> from exposure to sharp contaminated instruments during the cleaning process
- 4 Auxiliary instruments are ergonomically arranged and retained in position for optimal cleaning
- 5 Open cassette design for optimal cleaning, drainage and drying

# Uncompromised hygiene and a simple surgical workflow.



Straumann® ProClean Cassette  
Art. No. 040.561

- **Ease-of-use** – contains all the instruments needed to place any Straumann implant
- **Intuitive** – color-coded workflows guide you through the operation
- **Ease of handling** – simply place the used drills back into the cassette during surgery
- **Increased staff safety** – post-surgery, close the cassette lid, which safeguards staff from exposure to sharp contaminated instruments<sup>3</sup>
- **Avoid instrument mix-ups** – no need to sort soiled drills, all instruments are reprocessed in the cassette



A host of intelligent features, such as specific color-coded workflows, implant-line identification and instrument pictograms, provide clear intuitive navigation to the surgeon throughout the operation.





Don't compromise on hygiene anymore and contact your local Straumann representative now or visit [www.straumann.com](http://www.straumann.com).

### **International Headquarters**

Institut Straumann AG

Peter Merian-Weg 12

CH-4002 Basel, Switzerland

Phone +41 (0)61 965 11 11

Fax +41 (0)61 965 11 01

[www.straumann.com](http://www.straumann.com)

© Institut Straumann AG, 2016. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.