

DIGITAL MATERIALS OVERVIEW

MATERIALS ¹

	MILLING				SERVICE / SOLUTION			APPLICATION ¹													
	CARES® C Series	CARES® M series	CARES® D series	Centralized milling	CARES® X-Stream™	CARES® Scan and Shape	Connectivity	Tooth restoration			Others	Direct connection to implants									
								Inlays ² / Onlays ² / Veneers	Copings / Crowns	Bridges / Bars (maximum units)	Splints	CARES® Customized Abutments and Pre-Milled Blanks	etkon® iDent Customized Abutments and Pre-Milled Blanks	CARES® Screw-Retained Bridges and Bars							
Glass-ceramics	Translucency																				
nIce®				●	●	●	●	●	●												
IPS e.max® CAD	● ³	● ³		●	●	●	●	●	●												
VITA® Suprinity				●	●																
VITA® Mark II				●	●																
VITA® TriLuxe / TriLuxe Forte				●	●																
Hybrids	Translucency																				
3M™ ESPE™ Lava™ Ultimate		●					●	●													
VITA® Enamic	●	●																			
Ceramics	Translucency																				
zerion® UTML					●	●	●	●	●	●	3										
zolid® SHT ML		●									3										
zolid® SHT / SHT preshade		●									3										
zerion® ML					●	●	●	●	●	●	16										
3M™ ESPE™ Lava™ Plus Zirconia					●	●	●	●	●	●	16										
zolid® HT / HT preshade		●									16										
zerion® LT					●	●	●	●	●	●	16										
zi		●									16										
zerion® GI					●	●	●	●	●	●	16										
Zirconium Dioxide for CARES®					●	●	●	●	●	●	16			● ⁴							
Metal																					
coron® Cobalt-Chromium					●	●	●	●	●	●	16		●							●	
Cobalt-Chromium for etkon® iDent	●	●																			●
Sintron®		●					●				16										
TAN for CARES®	●	●						●	●					●							
ticon® Titanium Grade 4								●	●	●											●
Titanium for CARES®					●	●	●	●	●	●				● ⁴							
Titanium for etkon® iDent	●	●						●	●	●											●
Polymer																					
polycon® ae					●	●	●	●	●	●	16										
VITA® CAD-Temp	●	●									16										
Temp / Temp multilayer		●									16										
PEEK		●									16										
Splintec		●											●								
polycon® cast								●	●												●
PMMA (burn out polymer)		●																			
Wax		●																			

¹ application and material availability might differ from country to country.

² veneerable inlay/onlay for 3M™ ESPE™ Lava™ Plus Zirconia, zerion®, ticon®, coron®, polycon® ae and polycon® cast materials.

³ material listed and mill-able

⁴ CARES® Abutments available, no Pre-Milled Abutments Blanks available

MO: Medium Opacity LT: Low Translucency HT: High Translucency GI: Gingiva shaded

VITA® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Germany.

3M™, ESPE™ and Lava™ are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada.

JUVORA™ is a trademark of JUVORA Ltd.

IPS e.max® is a registered trademark of Ivoclar Vivadent AG, Liechtenstein.