

产品目录

2018/2019

缩写列表

-04	4支装
ABS	丙烯腈丁二烯苯乙烯
Al	铝
BL	骨水平
BLT	骨水平锥柱状种植体
EDTA	乙二胺四乙酸
EMD	牙釉基质衍生物
HT	高透明度
LDPE	低密度聚乙烯
LT	低透明度
MO 0	色号为MO 0
MO 1	色号为MO 1
NC	窄CrossFit®（十字锁合）=种植体直径3.3毫米，骨水平种植体
Ncm	扭力扳手测量单位
ND	窄直径=种植体肩台直径3.5mm
NN	窄颈=种植体肩台直径3.5毫米，美学种植体
NNC	窄颈窄CrossFit®（十字锁合）=种植体直径3.3mm，种植体肩台直径3.5mm，美学种植体
PA	聚酰胺
PC	聚碳酸酯
PEEK	聚醚醚酮
PMMA	聚甲基丙烯酸甲酯
POM	缩醛树脂
PPS	聚苯硫醚
PS	聚苯乙烯
PSU	聚砜
RC	常规CrossFit®（十字锁合）=种植体直径4.1和4.8毫米，骨水平种植体
RD	常规直径=种植体肩台直径4.8mm
RN	常规颈=种植体肩台直径4.8毫米，标准/美学/锥柱状种植体
S	标准
SCS	螺丝携带系统
SP	美学种植体
TAN	钛铝铌合金
TAV	钛铝钒合金
TE	锥柱状
Ti	钛
Ti-N	氮化钛
TiZr	钛锆合金（Roxolid®）

2018
/ 2019



Simply doing more

超出您的期望



随着大众对种植治疗的期待越来越高，医生也希望种植解决方案可以更加多元化。您接待的每一位患者都是与众不同的，每一次治疗也是独一无二的。士卓曼努力为您提供完美的齿科部件，开发全新的治疗解决方案。确保为您创造更加简化高效的工作流程，帮助您从容应对日益增长的业务需求，从而使更多缺牙患者享受高品质的医疗服务。

“Simply doing more”还体现在我们的全新产品目录当中。“More”是指其结构更合理，色彩更丰富，可以更便捷地引导您为每个适应症选择到合适的产品。

士卓曼原厂部件为您和您的患者提供了“更多”安全性。我们对产品制造精益求精，种植体和修复部件的设计都精确到毫米，最微小的细节也有着严格的质检把控。使用士卓曼原厂部件更加安全可靠，减少您的后顾之忧。

在中国上市了三种形态的**钛锆合金亲水种植体**：骨水平锥柱状（BLT）以及软组织水平平行壁（TL）、骨水平平行壁（BL）。钛锆合金亲水种植体采用了市面仅有的钛锆合金材料Roxolid®制造而成并结合了SLActive®亲水处理表面技术，开创了种植牙微创新时代，并为健康受损患者带来更高的可预期性治疗方案。

士卓曼CARES®作为**综合的数字化解决方案**，旨在为医生和技师提供最简单的方式实现牙科数字化。士卓曼CARES®拥有包括扫描、术前规划、切削、打印和相应服务的全部口腔治疗流程，这些流程都是高度整合，开放和模块化的——用户可以根据自身情况灵活地选择CARES®数字化方案。

致以最诚挚的问候

A stylized, handwritten signature in black ink, appearing to read 'Gadola'.

Marco Gadola
首席执行官

主打产品



骨水平锥柱状种植体，钛锆合金材料与亲水技术表面的完美结合，为您提供更高的初始稳定性。



SLActive®亲水表面处理技术，更快愈合，速度加倍。



士卓曼原厂基台，种植体与基台深度嵌合，为您提供高稳定性的解决方案。



士卓曼CARES®作为综合的数字化解方案涵盖了数字化外科以及修复等方面，给您和您的患者带来全新的诊疗体验。

目 录

	外科产品 4 种植体概览 6 颜色代码说明 10 软组织水平种植体 11 骨水平种植体 23	外科产品
	修复配件 30 软组织水平 NNC 33 软组织水平 RN 43 软组织水平 WN 55 骨水平 NC 69 骨水平 RC 91 Straumann® 经典 115	修复配件
	CADCAM 118 CADCAM 120	CADCAM
	器械和配件 124 外科手术工具 126 修复配件 144	器械和配件
	索引 152 字母索引 155 货号索引 173	索引



查看种植体章节时，首先选择种植体的类型，如软组织水平、骨水平。然后找到您所需的种植体直径。

每页的侧边都用彩色的矩形对此进行了标记。您可以在相应的章节中找到器械和配件。

不仅仅是一个种植体：更是一种信任

每位患者都与众不同，每次治疗也是个性化的。士卓曼将这种个性化体现在其种植体产品之中。不管您选择软组织水平种植体还是骨水平种植体，所有类型的种植体都为牙医提供了多样的选择、简便的操作，以及可靠性和安全性。

简易是指导原则。有了士卓曼®牙种植体系统，您可以在广泛的治疗方案中进行选择。

- 软组织水平和骨水平种植体系统可提供最大化选择
- SLActive® 亲水表面处理技术实现治疗结果的可预期性
- 选择 Roxolid® 瑞锆-钛锆合金种植体，开创微创种植新时代
- 全世界已经植入了超过1500万颗士卓曼种植体


外科产品

种植体概览	6
颜色代码说明	10
软组织水平种植体	11
种植体 Ø 3.3 毫米	12
种植体	12
封闭螺丝和愈合帽	14
种植体 Ø 4.1 毫米	15
种植体	15
封闭螺丝和愈合帽	17
种植体 Ø 4.8 毫米	18
种植体	18
封闭螺丝和愈合帽	21
骨水平种植体	23
种植体 Ø 3.3 毫米	24
种植体	24
封闭螺丝和愈合基台	25
种植体 Ø 4.1 毫米	26
种植体	26
封闭螺丝和愈合基台	27
种植体 Ø 4.8 毫米	28
种植体	28
封闭螺丝和愈合基台	29

Roxolid® SLActive® 种植体 (Loxim™)

外科产品

标准种植体	RN	RN	RN	RN	WN
 2.8 毫米	 Ø 4.8 毫米 Ø 3.3 毫米	 Ø 4.8 毫米 Ø 3.3 毫米	 Ø 4.8 毫米 Ø 4.1 毫米	 Ø 4.8 毫米 Ø 4.8 毫米	 Ø 6.5 毫米 Ø 4.8 毫米
SLActive®	S Ø 3.3 RN Loxim™	S Ø 3.3 RN OTP	S Ø 4.1 RN Loxim™	S Ø 4.8 RN Loxim™	S Ø 4.8 WN Loxim™
6 毫米	—	—	033.530S	033.580S	033.600S
8 毫米	033.501S	033.431S	033.531S	033.581S	033.601S
10 毫米	033.502S	033.432S	033.532S	033.582S	033.602S
12 毫米	033.503S	033.433S	033.533S	033.583S	033.603S
14 毫米	033.504S	033.434S	033.534S	033.584S	—
16 毫米	033.505S	033.435S	033.535S	—	—

美学种植体	NNC	RN	RN	RN	WN
 1.8 毫米	 Ø 3.5 毫米 Ø 3.3 毫米	 Ø 4.8 毫米 Ø 3.3 毫米	 Ø 4.8 毫米 Ø 4.1 毫米	 Ø 4.8 毫米 Ø 4.8 毫米	 Ø 6.5 毫米 Ø 4.8 毫米
SLActive®	SP Ø 3.3 NNC	SP Ø 3.3 RN Loxim™	SP Ø 4.1 RN Loxim™	SP Ø 4.8 RN Loxim™	SP Ø 4.8 WN Loxim™
6 毫米	—	—	033.560S	033.590S	033.610S
8 毫米	033.416S	033.511S	033.561S	033.591S	033.611S
10 毫米	033.417S	033.512S	033.562S	033.592S	033.612S
12 毫米	033.418S	033.513S	033.563S	033.593S	033.613S
14 毫米	033.419S	033.514S	033.564S	033.594S	—

锥柱状种植体		RN	RN		WN
 1.8 毫米		 Ø 4.8 毫米 Ø 3.3 毫米	 Ø 4.8 毫米 Ø 4.1 毫米		 Ø 6.5 毫米 Ø 4.8 毫米
SLActive®		TE Ø 3.3 RN Loxim™	TE Ø 4.1 RN Loxim™		TE Ø 4.8 WN Loxim™
8 毫米		033.521S	033.571S		—
10 毫米		033.522S	033.572S		033.622S
12 毫米		033.523S	033.573S		033.623S
14 毫米		033.524S	033.574S		033.624S

NNC=窄颈窄CrossFit® (十字锁合), RN=常规颈, S=标准, SP=美学种植体, TE=锥柱状, WN=宽颈

Loxim™: 卡扣式携带体, OTP: 螺丝固位式携带体

*注: Roxolid SLActive® 3.3mm RN 标准有两种携带体类型: OTP 和 Loxim™

Roxolid® SLActive® 种植体 (Loxim™)

骨水平种植体	NC	NC	RC	RC	
	 Ø 3.3 毫米 Ø 3.3 毫米	 Ø 3.3 毫米 Ø 3.3 毫米	 Ø 4.1 毫米 Ø 4.1 毫米	 Ø 4.8 毫米 Ø 4.8 毫米	
SLActive®	BL Ø 3.3 NC Loxim™	BL Ø 3.3 NC OTP	BL Ø 4.1 RC Loxim™	BL Ø 4.8 RC Loxim™	
8 毫米	021.2308	021.2108	021.4308	021.6308	
10 毫米	021.2310	021.2110	021.4310	021.6310	
12 毫米	021.2312	021.2112	021.4312	021.6312	
14 毫米	021.2314	021.2114	021.4314	021.6314	

骨水平锥柱状种植体		NC	RC	RC	
		 Ø 3.3 毫米	 Ø 4.1 毫米	 Ø 4.8 毫米	
SLActive®		BLT Ø 3.3 NC Loxim™	BLT Ø 4.1 RC Loxim™	BLT Ø 4.8 RC Loxim™	
8 毫米		021.3308	021.5308	021.7308	
10 毫米		021.3310	021.5310	021.7310	
12 毫米		021.3312	021.5312	021.7312	
14 毫米		021.3314	021.5314	021.7314	
16 毫米		021.3316	021.5316	021.7316	
18 毫米		021.3318	021.5318	021.7318	

Ti SLActive® 种植体 (OTP)

标准种植体		RN	RN	RN	WN
	 2.8 毫米	 Ø 4.8 毫米	 Ø 4.8 毫米	 Ø 4.8 毫米	 Ø 6.5 毫米
		 Ø 3.3 毫米	 Ø 4.1 毫米	 Ø 4.8 毫米	 Ø 4.8 毫米
SLActive®		S Ø 3.3 RN	S Ø 4.1 RN	S Ø 4.8 RN	S Ø 4.8 WN
6 毫米		—	033.030S	033.230S	033.630S
8 毫米		033.131S	033.031S	033.231S	033.631S
10 毫米		033.132S	033.032S	033.232S	033.632S
12 毫米		033.133S	033.033S	033.233S	033.633S
14 毫米		033.134S	033.034S	033.234S	—
16 毫米		033.135S	033.035S	—	—

美学种植体		NN	RN	RN	RN	WN
	 1.8 毫米	 Ø 3.5 毫米	 Ø 4.8 毫米	 Ø 4.8 毫米	 Ø 4.8 毫米	 Ø 6.5 毫米
		 Ø 3.3 毫米	 Ø 3.3 毫米	 Ø 4.1 毫米	 Ø 4.8 毫米	 Ø 4.8 毫米
SLActive®		SP Ø 3.3 NN	SP Ø 3.3 RN	SP Ø 4.1 RN	SP Ø 4.8 RN	SP Ø 4.8 WN
6 毫米		—	—	033.050S	033.250S	033.650S
8 毫米		033.951S	033.151S	033.051S	033.251S	033.651S
10 毫米		033.952S	033.152S	033.052S	033.252S	033.652S
12 毫米		033.953S	033.153S	033.053S	033.253S	033.653S
14 毫米		033.954S	033.154S	033.054S	033.254S	—

锥柱状种植体		RN	RN		WN
	 1.8 毫米	 Ø 4.8 毫米	 Ø 4.8 毫米		 Ø 6.5 毫米
		 Ø 3.3 毫米	 Ø 4.1 毫米		 Ø 4.8 毫米
SLActive®		TE Ø 3.3 RN	TE Ø 4.1 RN		TE Ø 4.8 WN
8 毫米		033.721S	033.761S		—
10 毫米		033.722S	033.762S		033.712S
12 毫米		033.723S	033.763S		033.713S
14 毫米		033.724S	033.764S		033.714S

骨水平种植体		NC	RC	RC	
		 Ø 3.3 毫米	 Ø 4.1 毫米	 Ø 4.8 毫米	
		 Ø 3.3 毫米	 Ø 4.1 毫米	 Ø 4.8 毫米	
SLActive®		BL Ø 3.3 NC	BL Ø 4.1 RC	BL Ø 4.8 RC	
8 毫米		021.2108	021.4108	021.6108	
10 毫米		021.2110	021.4110	021.6110	
12 毫米		021.2112	021.4112	021.6112	
14 毫米		021.2114	021.4114	021.6114	

Ti SLA® 种植体 (OTP)

外科产品

标准种植体		RN	RN	RN	WN
	2.8 毫米	Ø 4.8 毫米	Ø 4.8 毫米	Ø 4.8 毫米	Ø 6.5 毫米
		Ø 3.3 毫米	Ø 4.1 毫米	Ø 4.8 毫米	Ø 4.8 毫米
SLA®		S Ø 3.3 RN	S Ø 4.1 RN	S Ø 4.8 RN	S Ø 4.8 WN
6 毫米		—	043.030S	043.230S	043.630S
8 毫米		043.131S	043.031S	043.231S	043.631S
10 毫米		043.132S	043.032S	043.232S	043.632S
12 毫米		043.133S	043.033S	043.233S	043.633S
14 毫米		043.134S	043.034S	043.234S	—
16 毫米		043.135S	043.035S	—	—

美学种植体	NN	RN	RN	RN	WN
↓ 1.8 毫米 ↑	Ø 3.5 毫米  Ø 3.3 毫米	Ø 4.8 毫米  Ø 3.3 毫米	Ø 4.8 毫米  Ø 4.1 毫米	Ø 4.8 毫米  Ø 4.8 毫米	Ø 6.5 毫米  Ø 4.8 毫米
SLA®	SP Ø 3.3 NN	SP Ø 3.3 RN	SP Ø 4.1 RN	SP Ø 4.8 RN	SP Ø 4.8 WN
6 毫米	—	—	043.050S	043.250S	043.650S
8 毫米	042.930S	043.151S	043.051S	043.251S	043.651S
10 毫米	042.931S	043.152S	043.052S	043.252S	043.652S
12 毫米	042.932S	043.153S	043.053S	043.253S	043.653S
14 毫米	042.933S	043.154S	043.054S	043.254S	—

锥柱状种植体		RN	RN		WN
		 Ø 4.8 毫米 Ø 3.3 毫米	 Ø 4.8 毫米 Ø 4.1 毫米		 Ø 6.5 毫米 Ø 4.8 毫米
SLA®		TE Ø 3.3 RN	TE Ø 4.1 RN		TE Ø 4.8 WN
8 毫米		043.721S	043.761S		—
10 毫米		043.722S	043.762S		043.712S
12 毫米		043.723S	043.763S		043.713S
14 毫米		043.724S	043.764S		043.714S

骨水平种植体		NC	RC	RC	
		Ø 3.3 毫米  Ø 3.3 毫米	Ø 4.1 毫米  Ø 4.1 毫米	Ø 4.8 毫米  Ø 4.8 毫米	
SLA®		BL Ø 3.3 NC	BL Ø 4.1 RC	BL Ø 4.8 RC	
8 毫米		021.2408	021.4408	021.6408	
10 毫米		021.2410	021.4410	021.6410	
12 毫米		021.2412	021.4412	021.6412	
14 毫米		021.2414	021.4414	021.6414	

代码和型号说明

外科产品

颜色代码		
	黄色	骨内种植体直径3.3毫米
	红色	骨内种植体直径4.1毫米
	绿色	骨内种植体直径4.8毫米

种植体类型		
TL	S	标准种植体
	SP	美学种植体
	TE	锥柱状种植体
BL	BL	骨水平种植体
	BLT	骨水平锥柱状种植体

连接类型	
NNC: 窄颈窄 CrossFit® (十字锁合) Ø 3.5 mm	Ø 3.5 毫米 
RN: 常规颈 Ø 4.8 毫米	Ø 4.8 毫米 
WN: 宽颈 Ø 6.5 毫米	Ø 6.5 毫米 
NC: 窄 CrossFit® (十字锁合) Ø 3.3 毫米	Ø 3.3 毫米 
RC: 常规 CrossFit® (十字锁合) Ø 4.1 毫米 和 Ø 4.8 mm	Ø 4.1 毫米  Ø 4.8 毫米 
NN: 窄颈 Ø 3.5 毫米	Ø 3.5 毫米 

软组织水平种植体

Roxolid® SLActive® (Loxim™)

产品	平台	图片	材料	长度	产品编号
Standard Plus NNC NNC 美学种植体 Ø 3.3	NNC – 窄颈 CrossFit® (十字锁合) Ø 3.5 毫米		Roxolid®	8 毫米	033.416S
				10 毫米	033.417S
				12 毫米	033.418S
				14 毫米	033.419S
Standard Plus 美学种植体 Ø 3.3	RN – 常规颈 Ø 4.8 毫米			8 毫米	033.511S
				10 毫米	033.512S
				12 毫米	033.513S
				14 毫米	033.514S
Standard 标准种植体 Ø 3.3 (Loxim™)				8 毫米	033.501S
				10 毫米	033.502S
				12 毫米	033.503S
				14 毫米	033.504S
Standard 标准种植体 Ø 3.3 (OTP)				16 毫米	033.505S
				8 毫米	033.431S
				10 毫米	033.432S
				12 毫米	033.433S
Tapered Effect 锥柱状种植体 Ø 3.3		14 毫米	033.434S		
		16 毫米	033.435S		
		8 毫米	033.521S		
		10 毫米	033.522S		
				12 毫米	033.523S
				14 毫米	033.524S

SLActive® 引导型

Standard Plus 美学种植体 Ø 3.3	RN – 常规颈 Ø 4.8 毫米		Roxolid®	8 毫米	033.451G
Standard 标准种植体 Ø 3.3				10 毫米	033.452G
				12 毫米	033.453G
		8 毫米		033.431G	
		10 毫米		033.432G	
		12 毫米		033.433G	

Ti SLActive® (OTP)

产品	平台	图片	材料	长度 产品编号	
Standard Plus 美学种植体 Ø 3.3	RN – 常规颈 Ø 3.3 毫米		纯钛	8 毫米	033.151S
				10 毫米	033.152S
				12 毫米	033.153S
				14 毫米	033.154S
Standard 标准种植体 Ø 3.3				8 毫米	033.131S
				10 毫米	033.132S
				12 毫米	033.133S
				14 毫米	033.134S
Tapered Effect 锥柱状种植体 Ø 3.3				16 毫米	033.135S
				8 毫米	033.721S
				10 毫米	033.722S
				12 毫米	033.723S
				14 毫米	033.724S

Ti SLA® (OTP)

产品	平台	图片	材料	长度 产品编号	
Standard Plus 美学种植体 Ø 3.3	RN – 常规颈 Ø 3.3 毫米		纯钛	8 毫米	043.151S
				10 毫米	043.152S
				12 毫米	043.153S
				14 毫米	043.154S
Standard 标准种植体 Ø 3.3				8 毫米	043.131S
				10 毫米	043.132S
				12 毫米	043.133S
				14 毫米	043.134S
Tapered Effect 锥柱状种植体 Ø 3.3				16 毫米	043.135S
				8 毫米	043.721S
				10 毫米	043.722S
				12 毫米	043.723S
				14 毫米	043.724S

NNC 封闭螺丝和愈合帽

产品	图片	描述	高度	材料	产品编号
封闭螺丝					
Sterile NNC Closure Cap NNC 封闭螺丝, 非无菌包装		Ø 2.6 毫米, 小型	0 毫米	纯钛	048.324
		Ø 2.6 毫米, 小型, 4 枚装			048.324V4
		Ø 4 毫米, 大型	1.5 毫米		048.325
		Ø 4 毫米, 大型, 4 枚装			048.325V4
愈合帽					
Sterile NNC Healing Cap NNC 愈合帽, 非无菌包装		Ø 4 毫米, 短型	3 毫米	纯钛	048.071
		Ø 4 毫米, 常规高度	4.5 毫米		048.074
		Ø 4 毫米, 唇侧倾斜	2 毫米	纯钛/TAN	048.082

RN 封闭螺丝和愈合帽

产品	图片	描述	高度	材料	产品编号
封闭螺丝					
Sterile RN Closure Cap RN 封闭螺丝, 非无菌包装		Ø 3.5 毫米, 小型	0 毫米	纯钛	048.371
		Ø 3.5 毫米, 小型, 4 枚装			048.371V4
		Ø 6 毫米, 大型	1.5 毫米		048.373
		Ø 6 毫米, 大型, 4 枚装			048.373V4
愈合帽					
Sterile RN Healing Cap RN 愈合帽, 非无菌包装		Ø 5.5 毫米	2 毫米	纯钛	048.033
			3 毫米		048.034
			4.5 毫米		048.037
		Ø 5.5 毫米, 唇侧倾斜, 小型	2 毫米		048.028
		Ø 5.5 毫米, 唇侧倾斜, 大型	3.5 毫米		048.029

Roxolid® SLActive® (Loxim™)

产品	平台	图片	材料	长度	产品编号
Standard Plus 美学种植体 Ø 4.1	RN – 常规颈 Ø 4.8 毫米		Roxolid®	6 毫米	033.560S
				8 毫米	033.561S
				10 毫米	033.562S
				12 毫米	033.563S
				14 毫米	033.564S
Standard 标准种植体 Ø 4.1				6 毫米	033.530S
				8 毫米	033.531S
				10 毫米	033.532S
				12 毫米	033.533S
				14 毫米	033.534S
锥柱状种植体 Ø 4.1				16 毫米	033.535S
				8 毫米	033.571S
				10 毫米	033.572S
				12 毫米	033.573S

SLActive® 引导型					
Standard Plus 美学种植体 Ø 4.1	RN – 常规颈 Ø 4.8 毫米		Roxolid®	8 毫米	033.561G
				10 毫米	033.562G
				12 毫米	033.563G
Standard 标准种植体 Ø 4.1				8 毫米	033.531G
				10 毫米	033.532G
				12 毫米	033.533G
Tapered Effect 锥柱状种植体 Ø 4.1				8 毫米	033.571G
				10 毫米	033.572G
				12 毫米	033.573G

产品	平台	图片	材料	长度	产品编号
Standard Plus 美学种植体 Ø 4.1	RN – 常规颈 Ø 4.8 毫米		纯钛	6 毫米	033.050S
				8 毫米	033.051S
				10 毫米	033.052S
				12 毫米	033.053S
				14 毫米	033.054S
Standard 标准种植体 Ø 4.1				6 毫米	033.030S
				8 毫米	033.031S
				10 毫米	033.032S
				12 毫米	033.033S
				14 毫米	033.034S
Tapered Effect 锥柱状种植体 Ø 4.1				16 毫米	033.035S
				8 毫米	033.761S
				10 毫米	033.762S
				12 毫米	033.763S

Standard Plus 美学种植体 Ø 4.1	RN – 常规颈 Ø 4.8 毫米		纯钛	8 毫米	033.051G		
Standard 标准种植体 Ø 4.1				10 毫米	033.052G		
				12 毫米	033.053G		
				8 毫米	033.031G		
Tapered Effect 锥柱状种植体 Ø 4.1				10 毫米	033.032G		
				12 毫米	033.033G		
				8 毫米	033.761G		
						10 毫米	033.762G
						12 毫米	033.763G

产品	平台	图片	材料	长度	产品编号
Standard Plus 美学种植体 Ø 4.1	RN – 常规颈 Ø 4.8 毫米		纯钛	6 毫米	043.0505
				8 毫米	043.0515
				10 毫米	043.0525
				12 毫米	043.0535
				14 毫米	043.0545
Standard 标准种植体 Ø 4.1				6 毫米	043.0305
				8 毫米	043.0315
				10 毫米	043.0325
				12 毫米	043.0335
				14 毫米	043.0345
Tapered Effect 锥柱状种植体 Ø 4.1				8 毫米	043.0355
				10 毫米	043.7615
				12 毫米	043.7625
				14 毫米	043.7635

RN 封闭螺丝和愈合帽

产品	图片	描述	高度	材料	产品编号
封闭螺丝					
Sterile RN Closure Cap RN 封闭螺丝，非无菌包装		Ø 3.5 毫米，小型	0 毫米	纯钛	048.371
		Ø 3.5 毫米，小型，4 枚装			048.371V4
		Ø 6 毫米，大型	1.5 毫米		048.373
		Ø 6 毫米，大型，4 枚装			048.373V4
愈合帽					
Sterile RN Healing Cap RN 愈合帽，非无菌包装		Ø 5.5 毫米	2 毫米	纯钛	048.033
			3 毫米		048.034
			4.5 毫米		048.037
		Ø 5.5 毫米，舌侧倾斜，小型	2 毫米		048.028
		Ø 5.5 毫米，舌侧倾斜，大型	3.5 毫米		048.029

Roxolid® SLActive® (Loxim™)

产品	平台	图片	材料	长度	产品编号
Standard Plus 美学种植体 Ø 4.8	RN – 常规颈 Ø 4.8 毫米		Roxolid®	6 毫米	033.590S
				8 毫米	033.591S
				10 毫米	033.592S
				12 毫米	033.593S
				14 毫米	033.594S
Standard 标准种植体 Ø 4.8				6 毫米	033.580S
				8 毫米	033.581S
				10 毫米	033.582S
				12 毫米	033.583S
				14 毫米	033.584S
Standard Plus 美学种植体 Ø 4.8	WN – 宽颈 Ø 6.5 毫米			6 毫米	033.610S
				8 毫米	033.611S
				10 毫米	033.612S
				12 毫米	033.613S
Standard 标准种植体 Ø 4.8				6 毫米	033.600S
				8 毫米	033.601S
				10 毫米	033.602S
				12 毫米	033.603S
Tapered Effect 锥柱状种植体 Ø 4.8				10 毫米	033.622S
			12 毫米	033.623S	
			14 毫米	033.624S	
SLActive® 引导型					
Standard Plus 美学种植体 Ø 4.8	RN – 常规颈 Ø 4.8 毫米		Roxolid®	8 毫米	033.591G
				10 毫米	033.592G
				12 毫米	033.593G

Ti SLActive® (OTP)

4.8
外科产品

产品	平台	图片	材料	长度 产品编号	
Standard Plus 美学种植体 Ø 4.8	RN – 常规颈 Ø 4.8 毫米		纯钛	6 毫米	033.250S
				8 毫米	033.251S
				10 毫米	033.252S
				12 毫米	033.253S
				14 毫米	033.254S
Standard 标准种植体 Ø 4.8				6 毫米	033.230S
				8 毫米	033.231S
				10 毫米	033.232S
				12 毫米	033.233S
				14 毫米	033.234S
Standard Plus 美学种植体 Ø 4.8	WN – 宽颈 Ø 6.5 毫米			6 毫米	033.650S
				8 毫米	033.651S
				10 毫米	033.652S
				12 毫米	033.653S
				14 毫米	033.654S
Standard 标准种植体 Ø 4.8				6 毫米	033.630S
				8 毫米	033.631S
				10 毫米	033.632S
				12 毫米	033.633S
				14 毫米	033.634S
Tapered Effect 锥柱状种植体 Ø 4.8				10 毫米	033.712S
				12 毫米	033.713S
				14 毫米	033.714S
SLActive® 引导型					
Standard 标准种植体 Ø 4.8	RN – 常规颈 Ø 4.8 毫米		纯钛	8 毫米	033.251G
				10 毫米	033.252G
				12 毫米	033.253G

Ti SLA® (OTP)

产品	平台	图片	材料	长度 产品编号	
Standard Plus 美学种植体 Ø 4.8	RN – 常规颈 Ø 4.8 毫米		纯钛	6 毫米	043.250S
				8 毫米	043.251S
				10 毫米	043.252S
				12 毫米	043.253S
				14 毫米	043.254S
Standard 标准种植体 Ø 4.8				6 毫米	043.230S
				8 毫米	043.231S
				10 毫米	043.232S
				12 毫米	043.233S
				14 毫米	043.234S
Standard Plus 美学种植体 Ø 4.8	WN – 宽颈 Ø 6.5 毫米			6 毫米	043.650S
				8 毫米	043.651S
				10 毫米	043.652S
				12 毫米	043.653S
Standard 标准种植体 Ø 4.8				6 毫米	043.630S
			8 毫米	043.631S	
			10 毫米	043.632S	
			12 毫米	043.633S	
Tapered Effect 锥柱状种植体 Ø 4.8			10 毫米	043.712S	
			12 毫米	043.713S	
			14 毫米	043.714S	

RN 封闭螺丝和愈合帽

产品	图片	描述	高度	材料	产品编号
封闭螺丝					
Sterile RN Closure Cap RN 封闭螺丝，非无菌包装		Ø 3.5 毫米，小型	0 毫米	纯钛	048.371
		Ø 3.5 毫米，小型，4 枚装			048.371V4
		Ø 6 毫米，大型	1.5 毫米		048.373
		Ø 6 毫米，大型，4 枚装			048.373V4
愈合帽					
Sterile RN Healing Cap RN 愈合帽，非无菌包装		Ø 5.5 毫米	2 毫米	纯钛	048.033
			3 毫米		048.034
			4.5 毫米		048.037
		Ø 5.5 毫米，唇侧倾斜，小型	2 毫米		048.028
		Ø 5.5 毫米，唇侧倾斜，大型	3.5 毫米		048.029

WN 封闭螺丝和愈合帽

产品	图片	描述	高度	材料	产品编号
封闭螺丝					
Sterile WN Closure Cap WN 封闭螺丝，非无菌包装		Ø 4.3 毫米	0 毫米	纯钛	048.375
愈合帽					
Sterile WN Healing Cap WN 愈合帽，非无菌包装		Ø 7.2 毫米	2 毫米	纯钛	048.038
			3 毫米		048.039
			4.5 毫米		048.053
		Ø 7.2 毫米，唇侧倾斜	2 毫米		048.030

骨水平种植体

Roxolid® SLActive® (Loxim™)

产品	平台	图片	材料	长度	产品编号
Bone Level 骨水平种植体 Ø 3.3 (Loxim™)	NC – 窄 CrossFit® (十字锁合)		Roxolid®	8 毫米	021.2308
				10 毫米	021.2310
				12 毫米	021.2312
				14 毫米	021.2314
Bone Level 骨水平种植体 Ø 3.3 (OTP)				8 毫米	021.2108
				10 毫米	021.2110
				12 毫米	021.2112
				14 毫米	021.2114
Bone Level Tapered 骨水平锥柱状种植 体 Ø 3.3 (Loxim™)				8 毫米	021.3308
				10 毫米	021.3310
				12 毫米	021.3312
				14 毫米	021.3314

SLActive® 引导型					
Bone Level 骨水平种植体 Ø3.3	NC – 窄 CrossFit® (十字锁合)		Roxolid®	8 毫米	021.2208G
Bone Level Tapered 骨水平锥柱状种植体 Ø 3.3				10 毫米	021.2210G
				12 毫米	021.2212G
			8 毫米	021.3308G	
			10 毫米	021.3310G	
			12 毫米	021.3312G	

Ti SLActive® (OTP)

产品	平台	图片	材料	长度	产品编号
Bone Level 骨水平种植体 Ø 3.3	NC – 窄 CrossFit® (十字锁合)		纯钛	8 毫米	021.2108
				10 毫米	021.2110
				12 毫米	021.2112
				14 毫米	021.2114

Ti SLA® (OTP)

产品	平台	图片	材料	长度	产品编号
Bone Level 骨水平种植体 Ø 3.3	NC – 窄 CrossFit® (十字锁合)		纯钛	8 毫米	021.2408
				10 毫米	021.2410
				12 毫米	021.2412
				14 毫米	021.2414

NC 封闭螺丝和愈合基台

产品	图片	描述	高度	材料	产品编号
NC 封闭螺丝					
Sterile NC Closure Cap NC 封闭螺丝，非无菌包装		小型	0 毫米	纯钛	024.2100
		4 枚装			024.2100-04
		大型	0.5 毫米		024.2105
		4 枚装			024.2105-04
NC 愈合基台					
Sterile NC Healing Abutment NC 愈合基台，非无菌包装		Ø 3.6 毫米，圆锥形	2 毫米	纯钛	024.2222
			3.5 毫米		024.2224
			5 毫米		024.2226
		Ø 4.8 毫米，圆锥形	2 毫米		024.2242
			3.5 毫米		024.2244
			5 毫米		024.2246
		Ø 3.3 毫米，宝瓶形	3.5 毫米		024.2234
			5 毫米		024.2236
		Ø 5 毫米，个性化制作	7 毫米	PEEK	024.2270

Roxolid® SLActive® (Loxim™)

产品	平台	图片	材料	长度	产品编号
Bone Level 骨水平种植体 Ø 4.1	RC – 常规 CrossFit® (十字锁合)		Roxolid®	8 毫米	021.4308
				10 毫米	021.4310
				12 毫米	021.4312
				14 毫米	021.4314
Bone Level Tapered 骨水平锥柱状种植体 Ø 4.1				8 毫米	021.5308
				10 毫米	021.5310
				12 毫米	021.5312
				14 毫米	021.5314
				16 毫米	021.5316

SLActive® 引导型					
Bone Level 骨水平种植体 Ø 4.1	RC – 常规 CrossFit® (十字锁合)		Roxolid®	8 毫米	021.4308G
				10 毫米	021.4310G
				12 毫米	021.4312G

Ti SLActive® (OTP)

产品	平台	图片	材料	长度	产品编号
Bone Level 骨水平种植体 Ø 4.1	RC – 常规 CrossFit® (十字锁合)		纯钛	8 毫米	021.4108
				10 毫米	021.4110
				12 毫米	021.4112
				14 毫米	021.4114

SLActive® 引导型					
Bone Level 骨水平种植体 Ø 4.1	RC – 常规 CrossFit® (十字锁合)		纯钛	8 毫米	021.4108G
				10 毫米	021.4110G
				12 毫米	021.4112G

Ti SLA® (OTP)

产品	平台	图片	材料	长度	产品编号
Bone Level 骨水平种植体 Ø 4.1	RC – 常规 CrossFit® (十字锁合)		纯钛	8 毫米	021.4408
				10 毫米	021.4410
				12 毫米	021.4412
				14 毫米	021.4414

RC 封闭螺丝和愈合基台

产品	图片	描述	高度	材料	产品编号
RC 封闭螺丝					
Sterile RC Closure Cap RC 封闭螺丝，非无菌包装		小型	0 毫米	纯钛	024.4100
		4 枚装			024.4100-04
		大型	0.5 毫米		024.4105
		4 枚装			024.4105-04
RC 愈合基台					
Sterile RC Healing Abutment RC 愈合基台，非无菌包装		Ø 4.5 毫米，圆锥形	2 毫米	纯钛	024.4222
			4 毫米		024.4224
			6 毫米		024.4226
		Ø 6 毫米，圆锥形	2 毫米		024.4242
			4 毫米		024.4244
			6 毫米		024.4246
		Ø 4.4 毫米，宝瓶形	4 毫米	PEEK	024.4234
		Ø 4.7 毫米，宝瓶形	6 毫米		024.4236
		Ø 7 毫米，个性化制作	7 毫米		024.4270

Roxolid® SLActive® (Loxim™)

产品	平台	图片	材料	长度	产品编号
Bone Level 骨水平种植体 Ø 4.8	RC – 常规 CrossFit® (十字锁合)		Roxolid®	8 毫米	021.6308
				10 毫米	021.6310
				12 毫米	021.6312
				14 毫米	021.6314
Bone Level Tapered 骨水平锥柱状种植体 Ø 4.8				8 毫米	021.7308
				10 毫米	021.7310
				12 毫米	021.7312
				14 毫米	021.7314
				16 毫米	021.7316
SLActive® 引导型					
Bone Level 骨水平种植体 Ø 4.8	RC – 常规 CrossFit® (十字锁合)		Roxolid®	8 毫米	021.6308G
				10 毫米	021.6310G
				12 毫米	021.6312G
Bone Level Tapered 骨水平锥柱状种植体 Ø 4.8				8 毫米	021.7308G
				10 毫米	021.7310G
				12 毫米	021.7312G

Ti SLActive® (OTP)

产品	平台	图片	材料	长度	产品编号
Bone Level 骨水平种植体 Ø 4.8	RC – 常规 CrossFit® (十字锁合)		纯钛	8 毫米	021.6108
				10 毫米	021.6110
				12 毫米	021.6112
				14 毫米	021.6114
SLActive® 引导型					
Bone Level 骨水平种植体 Ø 4.8	RC – 常规 CrossFit® (十字锁合)		纯钛	8 毫米	021.6108G
				10 毫米	021.6110G
				12 毫米	021.6112G

Ti SLA® (OTP)

产品	平台	图片	材料	长度	产品编号
Bone Level 骨水平种植体 Ø 4.8	RC – 常规 CrossFit® (十字锁合)		纯钛	8 毫米	021.6408
				10 毫米	021.6410
				12 毫米	021.6412
				14 毫米	021.6414

RC 封闭螺丝和愈合基台

产品	图片	描述	高度	材料	产品编号
RC 封闭螺丝					
Sterile RC Closure Cap RC 封闭螺丝，非无菌包装		小型	0 毫米	纯钛	024.4100
		4 枚装			024.4100-04
		大型	0.5 毫米		024.4105
		4 枚装			024.4105-04
RC 愈合基台					
Sterile RC Healing Abutment RC 愈合基台，非无菌包装		Ø 4.5 毫米，圆锥形	2 毫米	纯钛	024.4222
			4 毫米		024.4224
			6 毫米		024.4226
		Ø 6 毫米，圆锥形	2 毫米		024.4242
			4 毫米		024.4244
			6 毫米		024.4246
		Ø 4.4 毫米，宝瓶形	4 毫米	PEEK	024.4234
		Ø 4.7 毫米，宝瓶形	6 毫米		024.4236
		Ø 7 毫米，个性化制作	7 毫米		024.4270



查看修复部件时，首先确定软组织水平种植体或骨水平种植体，然后通过修复平台检索。每页的侧边都用彩色的矩形对此进行了标记。从单个单位到多单位修复，所有对应的部件您都可以找到。器械和配件请查询“器械和配件”章节。

不仅仅是最高的精密度，更是完美的解决方案

士卓曼为所有适应症提供了多样化的种植修复解决方案——从单牙修复到全口解决方案。我们提供高品质的解决方案——士卓曼 SynOcta®（八角连接）和 CrossFit®（十字锁合），确保基台在种植体上的安全就位，以实现长期的稳定性。广泛的材料和螺栓固位、粘接固位和混合义齿修复产品的可选性使我们的产品组合更加完整。

原厂基台

只有士卓曼最了解士卓曼种植体的所有尺寸、公差和材料。借助于这些数据，我们对种植体和修复部件的设计进行校准，使之完美契合。制造的高精度和精湛的测试工序确保所有原装士卓曼部件的质量始终如一。

专用于 CAD/CAM 修复体的 Straumann® CARES® 系统

Straumann® CARES® 提供多种修复体解决方案——包括单牙修复体、牙齿支持式桥体和种植体支持式修复体。多种材料可供选择，流程便捷而简单。确保每一个产品均可实现高精度、高有效性和高美观度的修复。

修复配件

软组织水平 NNC	33
窄颈窄 CrossFit®(十字锁合) (NNC)	34
牙冠和牙桥	34
实心基台	38
混合式义齿	39
软组织水平 RN	43
常规颈 (RN) synOcta®	44
牙冠和牙桥	44
实心基台	49
混合式义齿	51
软组织水平 WN	55
宽颈 (WN) synOcta®	56
牙冠和牙桥	56
实心基台	61
混合式义齿	63
修复设计套装 (RN, WN)	66
模拟基台	66
骨水平 NC	69
窄 CrossFit® (十字锁合) (NC)	70
牙冠和牙桥	70
混合式义齿	88
骨水平 RC	91
常规 CrossFit® (十字锁合) (RC)	92
牙冠和牙桥	92
混合式义齿	108
修复设计套装 (NC, RC)	110
模拟基台	111
Straumann® 经典	115
经典套装	117

软组织水平 NNC
窄颈窄CrossFit®（十字锁合）

印模和制作工作模型

产品	图片	描述	材料	产品编号
转移杆和替代体				
NNC Impression Post NNC 转移杆		适用于螺丝固位，开窗印模，高度10.1毫米	TAN	048.128
		适用于螺丝固位，非开窗印模，高度10.9毫米	POM/TAN	048.134
NNC Implant Analog NNC 种植体替代体		长度12.5毫米	纯钛	048.127
数字化印模				
CARES® NNC Mono Scanbody CARES® NNC 扫描杆		包括固定螺丝，Ø 3.5 毫米，高度10毫米	PEEK/TAN	048.173
骀记录转移体				
NNC Bite Registration Aid NNC 骀记录转移体		短型，高度8毫米，4枚装	POM	048.960V4
		长型，高度12毫米，4枚装		048.961V4

临时修复

产品	图片	描述	材料	产品编号
临时基台				
NNC Post NNC 基台		适用于临时修复，冠用，包括中央螺丝048.313，高度 10.3 毫米	TAN	048.501
		适用于临时修复，桥用，包括中央螺丝048.313，高度 10.3 毫米		048.502
替代螺丝				
NNC Basal Screw NNC 中央螺丝		适用于 NNC 钛合金临时基台，金基台以及 Variobase®，长度 7.1 毫米	TAN	048.313

粘接固位修复，粘接基台

产品	图片	描述	材料	产品编号
粘接基台, 0°				
NNC Cementable Abutment NNC 粘接基台		直型，适用于牙冠和牙桥，包括中央螺丝048.314，高度 5.7 毫米	纯钛	048.598
粘接基台, 15°				
NNC Cementable Abutment NNC 粘接基台		成角，角度15°，A 型，包括中央螺丝048.314，高度5.9 毫米	纯钛	048.621
		成角，角度15°，B 型，包括中央螺丝048.314，高度5.9 毫米		048.622
脱蜡基底				
NNC Plastic Coping NNC 塑料基底		适用于粘接基台，直型，冠用，高度 7 毫米	POM	048.256
		适用于粘接基台，直型，桥用，高度 7 毫米		048.257
		适用于粘接基台，成角，角度15°，适用于牙冠和牙桥，高度7.2 毫米		048.258
替代螺丝				
NNC Basal Screw NNC 中央螺丝		适用于NNC 粘接基台，长度 7.1 毫米	TAN	048.314

螺丝固位或粘接固位修复，CAD/CAM

产品	图片	描述	材料	产品编号
Variobase® 基台				
NNC Variobase® for Crown NNC 多能基台, 适用于单冠		包括中央螺丝 048.313, Ø 3.9 毫米, 高度 3.75 毫米	TAN	048.709
		包括中央螺丝 048.313, Ø 3.9 毫米, 高度 5.75 毫米		022.0021
NNC Variobase® C NNC 多能基台		适用于 Dentsply® Sirona®, 包括中央螺丝 048.313, Ø 3.9 毫米, 高度 4.7 毫米		022.0018
CADCAM蜡型				
NNC Wax-up Sleeve NNC 蜡型		长度 12.75 毫米	PEEK	048.137
		4 枚装		048.137V4
替代螺丝				
NNC Basal Screw NNC 中央螺丝		适用于 NNC 钛合金临时基台, 金基台以及 Variobase®, 长度 7.1 毫米	TAN	048.313

- Straumann® CARES® CAD/CAM

 - 定制解决方案，包括种植上部修复、牙齿支持式桥体和种植体支持式修复体
 - 适用范围广泛、可选材料丰富
 - 针对各种适应症的定制解决方案
 - 更多信息请访问：www.straumann.com

螺丝固位或粘接固位修复，金基台

产品	图片	描述	材料	产品编号
金基台				
NNC Gold Abutment NNC 金基台		适用于牙冠，包括中央螺丝 048.313，高度3.6毫米	Ceramicor®/ TAN/POM	048.592
		适用于牙桥，包括中央螺丝 048.313，高度3.6毫米		048.593
替代螺丝				
NNC Basal Screw NNC 中央螺丝		适用于 NNC 钛合金临时基台，金基台以及 Variobase®，长度 7.1 毫米	TAN	048.313

Ceramicor® 是一种非氧化贵金属合金，适用于铸造贵金属合金。

螺丝固位或粘接固位修复，多能基台

产品	图片	描述	材料	产品编号
Variobase® 基台				
NNC Variobase® for Crown NNC 多能基台，适用于单冠		包括中央螺丝 048.313，Ø 3.9 毫米，高度 3.75 毫米	TAN	048.709
		包括中央螺丝 048.313，Ø 3.9 毫米，高度 5.75 毫米		022.0021
脱蜡基底				
NNC Burn-out Coping NNC 脱蜡基底		适用于 Variobase® 牙冠高度3.75 毫米，Ø 3.9 毫米，高度11.1毫米	POM	048.267
		4枚装		048.267V4
		适用于 Variobase® 牙冠高度5.75 毫米，Ø 3.9 毫米，高度11.1毫米		023.0014
		4枚装		023.0014V4
替代螺丝				
NNC Basal Screw NNC 中央螺丝		适用于 NNC 钛合金临时基台，金基台以及 Variobase®，长度 7.1 毫米	TAN	048.313

产品	图片	描述	材料	产品编号
印模和制作工作模型				
NNC Impression Cap NNC 印模帽		高度 8 毫米	POM	048.579
NNC Positioning Cylinder NNC 定位柱		适用于 048.547, 黄色, 高度 10.2 毫米		048.587
		适用于 048.548, 灰色, 高度 10.2 毫米		048.588
		适用于 048.549, 蓝色, 高度 10.2 毫米		048.589
NNC Implant Analog NNC 可复位种植体替代体		适用于实心基台 048.547, 黄色, 高度 14.3 毫米	TAN	048.527
		适用于实心基台 048.548, 灰色, 高度 15.8 毫米		048.528
		适用于实心基台 048.549, 蓝色, 高度 17.3 毫米		048.529
临时修复				
NNC Temporary Coping NNC 临时基底		适用于实心基台, 牙冠用, 高度 8.5 毫米	PMMA	048.658
		适用于实心基台, 牙桥用, 高度 8.5 毫米		048.659
保护帽				
NNC Protective Cap NNC 固定式保护帽		适用于实心基台 048.547, 高度 5.6 毫米	PEEK	048.700
		适用于实心基台 048.548, 高度 7.1 毫米		048.701
		适用于实心基台 048.549, 高度 8.6 毫米		048.702

产品	图片	描述	材料	产品编号
实心基台				
NNC Solid Abutment NNC 实心基台		黄色，高度 4 毫米	纯钛	048.547
		灰色，高度 5.5 毫米		048.548
		蓝色，高度 7 毫米		048.549
脱蜡基底				
NNC Plastic Coping NNC 塑料基底		适用于实心基台，牙冠用，高度 10.2 毫米	POM	048.263
		适用于实心基台，牙桥用，高度 10.2 毫米		048.264
工具和辅助部件				
NNC Solid Abutment Driver NNC 实心基台螺丝刀		短型，长度 13 毫米	不锈钢	048.155
		长型，长度 19 毫米		048.156

LOCATOR®

产品	图片	描述	材料	产品编号
印模和制作工作模型				
LOCATOR® Impression Coping 取模部件		高度4毫米，4枚装	铝外壳，含 LDPE 衬垫	048.197V4
LOCATOR® Female Analog 基台替代体		长度10毫米，Ø 4 毫米，4枚装	铝	025.0120-04
钛基台				
NC LOCATOR® Abutment NC LOCATOR® 基台		高度1毫米	TAV	048.581
		高度2毫米		048.582
		高度3毫米		048.583
		高度4毫米		048.584
		高度5毫米		048.585
		高度6毫米		048.586

 0086 制造商
Zest Anchors, Inc.
Escondido, CA 92029
美国

NNC
修复配件

LOCATOR®

产品	图片	描述	材料	产品编号
技工室套装，可以达到20°的分散补偿				
LOCATOR® Male Processing Package 阳性部件套装		内含：基底帽(∅ 5.5毫米，高度2.5毫米)，含阳性黑色衬垫(高度1.9毫米)，封闭环(048.196)，透明阳性垫片(048.190)，粉色阳性垫片(048.191)，蓝色阳性垫片(048.192)，2枚装	纯钛/尼龙/ LDPE	048.189V2
LOCATOR® Replacement Male 阳性垫片		透明，0°-10°，高度1.7毫米，2.27kg，4枚装	尼龙	048.190V4
		粉色，轻微固位力，0°-10°，高度1.7毫米，1.36kg，4枚装		048.191V4
		蓝色，超低固位力，0°-10°，高度1.7毫米，0.68kg，4枚装		048.192V4
技工室套装，可以达到40°的分散补偿				
LOCATOR® Male Processing Package 阳性部件套装		加大角度补偿范围，内含：基底帽(∅ 5.5毫米，高度2.5毫米)，含阳性黑色衬垫(高度1.9毫米)，封闭环(048.196)，绿色阳性垫片(048.193)，橙色阳性垫片(048.188)，红色阳性垫片(048.194)，2枚装	纯钛/尼龙/ LDPE	048.182V2
LOCATOR® Replacement Male 阳性垫片		绿色，加大角度补偿范围，10°-20°，高度1.7毫米，1.82kg，4枚装	尼龙	048.193V4
		橙色，轻微固位力，加大角度补偿范围，10°-20°，高度1.7毫米，0.91kg，4枚装		048.188V4
		红色，超低固位力，加大角度补偿范围，10°-20°，高度1.7毫米，0.45kg，4枚装		048.194V4
		灰色，加大角度补偿范围，零(0)固位力，10°-20°，高度1.7毫米，4枚装		048.181V4

CE

0086

制造商

Zest Anchors, Inc.

Escondido, CA 92029

美国

LOCATOR® 配件

产品	图片	描述	材料	产品编号
LOCATOR® Parallel Post 平行柱		高度8毫米,4枚装	LDPE	048.199V4
LOCATOR® Processing Male, black 阳性垫片, 黑色		高度1.9毫米,4枚装		048.195V4
LOCATOR® Block-out Spacer 封闭环		厚度0.4毫米,20枚装	硅树脂 橡胶	048.196V20
LOCATOR® Abutment Holder Sleeve 基台固位套		用于将基台紧扣到工具上, 4 枚装	PSU	046.413V4

辅助配件				
LOCATOR® Core Tool 核心工具		三个部件, 长度100毫米	不锈钢	046.415
LOCATOR® Screwdriver 螺丝刀		适用于棘轮扳手, 短型,长度15毫米		046.416
		适用于棘轮扳手, 长型,长度21毫米		046.417
LOCATOR® Angle Measurement Guide 角度尺		长度50毫米, 宽度15毫米		048.200

 0086 制造商
Zest Anchors, Inc.
Escondido, CA 92029
美国

NNC
修复配件

软组织水平 RN
常规颈

印模和制作工作模型

产品	图片	描述	材料	产品编号
印模（种植体水平）				
RN Impression Post for open-tray impression RN 转移杆，适用于开窗式印模		螺丝固位，适用于开窗式印模，高度 11 毫米	铝/纯钛	048.010
		螺丝固位，适用于开窗式印模，高度 21 毫米		048.090
RN Impression Post for closed-tray impression RN 转移杆，适用于非开窗式印模		印模帽，白色，适用于非开窗式印模，高度 8 毫米，4 枚装	POM	048.017V4
		定位柱，红色，适用于非开窗式印模，高度 12 毫米，4 枚装		048.070V4
		螺丝固位，高度 11 毫米，带1枚引导螺丝及2枚印模帽	铝/纯钛/POM	048.862
印模（基台水平）				
RN Impression Post for open-tray impression RN 转移杆，适用于开窗式印模		螺丝固位，高度 12 毫米，适用于 RN synOcta®（八角）1.5 螺丝固位基台	铝/纯钛	048.866
RN Impression Post for closed-tray impression RN 转移杆，适用于非开窗式印模		螺丝固位，高度 11 毫米，适用于 RN synOcta®（八角）1.5 螺丝固位基台	铝/纯钛/POM	048.864
替代体				
RN Analog RN 替代体		灰色带红色条纹，长度 12 毫米	不锈钢	048.124
		适用于 synOcta®（八角）1.5 螺丝固位基台 048.601，直型，成角 0°，Ø 4.8 毫米		048.108
数字化印模				
CARES® RN Mono Scanbody CARES® RN 扫描杆		包括固定螺丝，Ø 4.8毫米，高度 10 毫米	PEEK/TAN	048.168
咬骀记录转移体				
Bite Registration Aid 咬骀记录转移体		短型，高度 8 毫米，Ø 5 毫米，4 枚装	POM	048.940V4
		长型，高度 12 毫米，Ø 5 毫米，4 枚装		048.941V4

临时修复

产品	图片	描述	材料	产品编号
临时基台				
RN synOcta® Temporary Meso abutment RN 临时基台		包含中央螺丝 048.356，高度 10 毫米	PEEK/TAN	048.668
临时基台				
RN synOcta®(八角) Post RN 临时基台		适用于牙桥临时修复，包含整体导向螺丝，高度9毫米	纯钛	048.650
		适用于牙冠临时修复，包含整体导向螺丝，高度9毫米		048.651
替代螺丝				
synOcta®(八角) Basal screw 中央螺丝		适用于 RN/WN Variobase®，金基台和 PEEK 临时基台，长度 6.7 毫米	TAN	048.356
替代螺丝				
SCS Occlusal Screw SCS 骀面螺丝		适用于 RN/WN 钛合金临时基台，长度 7.6 毫米	TAN	049.181

粘接固位修复，粘接基台

产品	图片	描述	材料	产品编号
粘接基台				
RN synOcta®(八角) Cementable Abutment RN 粘接基台		直型，包含中央螺丝，高度 5.5 毫米	纯钛	048.605
脱蜡基底				
RN synOcta®(八角) Plastic Coping RN 塑料基底		适用于常规颈 synOcta®(八角) 粘接基台 048.605，适用于牙冠，高度 7 毫米	POM	048.663
		适用于常规颈 synOcta®(八角) 粘接基台 048.605，适用于牙桥，高度 7 毫米		048.662
辅助配件				
Transfer Aid 转移辅助配件		适用于常规颈 synOcta®(八角) 粘接基台 048.605，高度 6.5 毫米，4枚装	PC	048.059V4

RN

修复配件

粘接固位或螺丝固位修复，角度基台

产品	图片	描述	材料	产品编号
角度基台，15°				
RN synOcta®(八角) Angled Abutment, 15° RN 角度基台，角度15°		短型，A型，高度5.7毫米，内含中央螺丝	纯钛/TAN	048.612
		长型，A型，高度6.7毫米，内含中央螺丝		048.610
		短型，B型，高度5.7毫米，内含中央螺丝		048.613
		长型，B型，高度6.7毫米，内含中央螺丝		048.611
角度基台，20°				
RN synOcta®(八角) Angled Abutment, 20° RN 角度基台，角度20°		短型，A型，高度5.7毫米，内含中央螺丝	纯钛/TAN	048.617
		长型，A型，高度6.7毫米，内含中央螺丝		048.615
		短型，B型，高度5.7毫米，内含中央螺丝		048.618
		长型，B型，高度6.7毫米，内含中央螺丝		048.616
脱蜡基底				
Extension Shell 加长套筒		适用于 synOcta® (八角)短型和长型，角度15° 和 角度20°基台，高度5.5毫米	POM	048.670
Plastic Shoulder 塑料肩台基底		适用于 synOcta® (八角)短型和长型，角度15° 和 角度20°基台，高度4毫米，4枚装		048.676V4
辅助配件和螺丝				
SCS Occlusal Screw SCS 骀面螺丝		用于 RN/WN 固位基底的骀面螺丝，长度4.4毫米	纯钛	048.350
		4枚装		048.350V4
Transfer Aid 转移辅助配件		适用于 (八角)，角度15°和角度20°基台，短型，高度5.5毫米, 4枚装	PC	048.000V4
		适用于 (八角)，角度15°和角度20°基台，长型，高度5.5毫米, 4枚装		048.002V4

粘接固位或螺丝固位修复，金基台

产品	图片	描述	材料	产品编号
金基台				
RN synOcta® (八角) Gold Abutment RN 金基台		适用于牙冠，含中央螺丝048.356，高度4.3毫米	Ceramicor®/ POM	048.642
替代螺丝				
synOcta®(八角) Basal Screw 中央螺丝		适用于 RN/WN Variobase®，金基台和临时基台，长度6.7 毫米	TAN	048.356

Ceramicor® 是一种非氧化贵金属合金，适用于铸造贵金属合金。

粘接固位或螺丝固位修复，多能基台，CAD/CAM

产品	图片	描述	材料	产品编号
Variobase® 基台				
RN Variobase® for Crown RN 多能基台，适用于单冠		含中央螺丝048.356，Ø 5.05毫米，高度4毫米	TAN	048.710
		含中央螺丝048.356，Ø 5.05毫米，高度6毫米		022.0022
RN Variobase® C RN 多能基台		适用于Dentsply® Sirona®，含中央螺丝 022.0045, Ø 5 毫米，高度7.7毫米		022.0019

替代螺丝				
synOcta®(八角) Basal Screw 中央螺丝		适用于 RN/WN Variobase®，金基台和临时基台，长度 6.7 毫米	TAN	048.356
Basal Screw RN/WN Variobase® C RN/WN 中央螺丝，适用于多能基台 C		适用于 Dentsply® Sirona®，RN/WN Variobase® C，长度 7.2 毫米		022.0045

RN

修复配件

粘接固位或螺丝固位修复，多能基台，铸造

产品	图片	描述	材料	产品编号
Variobase® 基台				
RN Variobase® for Crown RN 多能基台，适用于单冠		含中央螺丝048.356，Ø 5.05毫米，高度4毫米	TAN	048.710
		含中央螺丝048.356，Ø 5.05毫米，高度6毫米		022.0022
脱蜡基底				
RN Burn-out Coping RN 脱蜡基底		适用于Variobase®基台，Ø 5.05毫米，高度11.1毫米	POM	048.268
		4枚装		048.268V4
		适用于Variobase®基台，牙冠高度 6 毫米，Ø 5 毫米，高度 11.1 毫米		023.0015
		4枚装		023.0015V4
替代螺丝				
synOcta®(八角) Basal Screw 中央螺丝		适用于 RN/WN Variobase®，金基台和临时基台，长度 6.7 毫米	TAN	048.356

螺丝固位修复，synOcta® (八角) 1.5

产品	图片	描述	材料	产品编号
基台				
RN synOcta® (八角) 1.5 Screw-retained Abutment RN 1.5 螺丝固位基台		内含中央螺丝，高度1.5毫米	TAN	048.601
可铸金基底				
RN synOcta® (八角) Gold Coping RN 可铸金基底		适用于synOcta®(八角)1.5(包括铸模辅件)，牙冠，高度4.25/12毫米	Ceramicor®/ POM	048.633
		适用于synOcta®(八角)1.5(包括铸模辅件)，牙桥，高度4.25/12毫米		048.632
脱蜡基底				
RN synOcta® (八角) Plastic Coping RN 塑料基底		适用于synOcta®(八角)1.5，牙冠，高度10毫米	POM	048.229
		适用于synOcta®(八角)1.5，牙桥和杆卡，高度10毫米		048.227
钛研磨基底				
RN synOcta® (八角) Meso Milling Cylinder RN Meso 式研磨基底		适用于synOcta®(八角)1.5，牙冠，高度10毫米，Ø 8毫米	纯钛	048.571
RN synOcta® (八角) Milling Cylinder RN 研磨基底		适用于synOcta®(八角)1.5，牙冠，高度8毫米，Ø 6毫米		048.240
		适用于synOcta®(八角)1.5，牙桥，高度8毫米，Ø 6毫米		048.241
辅助配件和螺丝				
SCS Occlusal Screw SCS 骀面螺丝		用于 RN/WN 固位基底的骀面螺丝，长度4.4毫米	纯钛	048.350
		4枚装		048.350V4
Polishing Protector 打磨保护杆		适用于 synOcta®(八角)基底，横向骀面螺丝固位，长度15毫米	不锈钢	046.245

Ceramicor® 是一种非氧化贵金属合金，适用于铸造贵金属合金。

印模和制作工作模型，适用于实心基台

产品	图片	描述	材料	产品编号
印模和制作工作模型				
RN Impression Cap RN 印模帽		白色，适用于非开窗式印模，高度 8 毫米， 4 枚装	POM	048.017V4
Positioning Cylinder 定位柱		适用于048.540，黄色，高度10.2毫米，4枚装		048.060V4
		适用于048.541，灰色，高度10.2毫米，4枚装		048.061V4
		适用于048.542，蓝色，高度10.2毫米，4枚装		048.062V4
Implant Analog 替代体		适用于048.540，黄色，长度14.3毫米	铝	048.160
		适用于048.541，灰色，长度15.8毫米		048.161
		适用于048.542，蓝色，长度17.3毫米		048.162
Shoulder Analog 肩台替代体		适用于RN实心基台，带固定杆，长度10毫米，长度17毫米，4枚装	PPS/黄铜	048.117V4
临时修复				
RN Temporary Coping RN 实心基台临时基底		适用于RN实心基台，牙冠，长度8.5毫米	PMMA	048.655
		适用于RN实心基台，牙桥，长度8.5毫米		048.654
保护帽				
RN Protective Cap RN 保护帽		适用于048.540，粘接，高度5.8毫米，4枚装	PEEK	048.047V4
		适用于048.541，粘接，高度7.3毫米，4枚装		048.048V4
		适用于048.542，粘接，高度8.8毫米，4枚装		048.049V4

RN

修复配件

粘接固位修复，适用于实心基台

产品	图片	描述	材料	产品编号
实心基台				
RN Solid Abutment 6° RN 实心基台 6°		黄色，高度4毫米	纯钛	048.540
		灰色，高度5.5毫米		048.541
		蓝色，高度7毫米		048.542
脱蜡基底				
Plastic Coping 塑料基底		适用于RN实心基台，牙冠，高度10毫米	POM	048.245
		适用于RN实心基台，牙桥，高度10毫米		048.246
辅助配件				
Inserting Device 螺丝刀		适用RN实心基台，，短型，长度13毫米	不锈钢	046.067
		适用RN实心基台，，长型，长度19毫米		046.068

RN

修复配件

杆卡结构，金基底

产品	图片	描述	材料	产品编号
基台				
RN synOcta® (八角) 1.5 Screw-retained Abutment RN 1.5 螺丝固位基台		内含中央螺丝，高度1.5毫米	TAN	048.601
杆卡金基底				
RN synOcta® (八角) Gold Coping RN 金基底		适用于synOcta®(八角)1.5，杆卡，高度5.8毫米，合金重量0.789克	Ceramicor®	048.204
辅助配件和螺丝				
SCS Occlusal Screw SCS 骀面螺丝		用于 RN/WN 固位基底的骀面螺丝，长度4.4毫米	纯钛	048.350
		4枚装		048.350V4

杆卡结构，钛基底

产品	图片	描述	材料	产品编号
基台				
RN synOcta® (八角) 1.5 Screw-retained Abutment RN 1.5 螺丝固位基台		内含中央螺丝，高度1.5毫米	TAN	048.601
钛杆卡部件				
RN synOcta® (八角) Titanium Coping RN 钛基底		适用于synOcta®(八角)1.5，杆卡，高度6毫米	纯钛	048.214
辅助配件和螺丝				
SCS Occlusal Screw SCS 骀面螺丝		用于 RN/WN 固位基底的骀面螺丝，长度4.4毫米	纯钛	048.350
		4枚装		048.350V4

杆卡结构，塑料基底

产品	图片	描述	材料	产品编号
基台				
RN synOcta® (八角) 1.5 Screw-retained Abutment RN 1.5 螺丝固位基台		内含中央螺丝，高度1.5毫米	TAN	048.601
脱蜡基底				
RN synOcta®(八角) Plastic Coping RN 塑料基底		适用于synOcta®(八角)1.5，牙桥和杆卡，高度10毫米	POM	048.227
辅助配件和螺丝				
SCS Occlusal Screw SCS 骀面螺丝		用于 RN/WN 固位基底的骀面螺丝，长度4.4毫米	纯钛	048.350
		4枚装		048.350V4

Ceramicor® 是一种非氧化贵金属合金，适用于铸造贵金属合金。

LOCATOR®

产品	图片	描述	材料	产品编号
印模和制作工作模型				
LOCATOR® Impression Coping 印模帽		高度4毫米，4枚装	铝外壳，含 LDPE 衬垫	048.197V4
LOCATOR® Female Analog 基台替代体		长度10毫米，Ø 5毫米，4枚装	铝	048.198V4
钛基台				
RN LOCATOR® Abutment RN LOCATOR® 基台		高度1毫米	TAV	048.175
		高度2毫米		048.176
		高度3毫米		048.177
		高度4毫米		048.178
		高度5毫米		048.179
		高度6毫米		048.180

CE

0086

制造商
Zest Anchors, Inc.
Escondido, CA 92029
美国

LOCATOR®

产品	图片	描述	材料	产品编号
技工室套装，可以达到20°的分散补偿				
LOCATOR® Male Processing Package 阳性部件套装		内含：基底帽(Ø 5.5毫米，高度2.5毫米)，含阳性黑色衬垫(高度1.9毫米)，封闭环(048.196)，透明阳性垫片(048.190)，粉色阳性垫片(048.191)，蓝色阳性垫片(048.192)，2枚装	纯钛/尼龙/ LDPE	048.189V2
LOCATOR® Replacement Male 阳性垫片		透明，0°-10°，高度1.7毫米，2.27kg，4枚装	尼龙	048.190V4
		粉色，轻微固位力，0°-10°，高度1.7毫米，1.36kg，4枚装		048.191V4
		蓝色，超低固位力，0°-10°，高度1.7毫米，0.68kg，4枚装		048.192V4
技工室套装，可以达到40°的分散补偿				
LOCATOR® Male Processing Package 阳性部件套装		加大角度补偿范围，内含：基底帽(Ø 5.5毫米，高度2.5毫米)，含阳性黑色衬垫(高度1.9毫米)，封闭环(048.196)，绿色阳性垫片(048.193)，橙色阳性垫片(048.188)，红色阳性垫片(048.194)，2枚装	纯钛/尼龙/ LDPE	048.182V2
LOCATOR® Replacement Male 阳性垫片		绿色，加大角度补偿范围，10°-20°，高度1.7毫米，1.82kg，4枚装	尼龙	048.193V4
		橙色，轻微固位力，加大角度补偿范围，10°-20°，高度1.7毫米，0.91kg，4枚装		048.188V4
		红色，超低固位力，加大角度补偿范围，10°-20°，高度1.7毫米，0.45kg，4枚装		048.194V4
		灰色，加大角度补偿范围，零(0)固位力，10°-20°，高度1.7毫米，4枚装		048.181V4

CE 0086 制造商
Zest Anchors, Inc.
Escondido, CA 92029
美国

RN

修复配件

RN

修复配件

辅助配件				
LOCATOR® Core Tool 核心工具		三个部件，长度100毫米	不锈钢	046.415
LOCATOR® Screwdriver 螺丝刀		适用于棘轮扳手，短型,长度15毫米		046.416
		适用于棘轮扳手，长型,长度21毫米		046.417
LOCATOR® Angle Measurement Guide 角度尺		长度50毫米，宽度15毫米		048.200



软组织水平 WN
宽颈

WN

修复配件

56

临时修复

产品	图片	描述	材料	产品编号
临时基台				
WN synOcta® (八角) Post WN 临时基台		适用于牙冠临时修复, 包含中央螺丝, 高度9毫米	纯钛	048.233
		适用于牙桥临时修复, 包含中央螺丝, 高度9毫米		048.234
替代螺丝				
synOcta®(八角) Basal Screw 中央螺丝		适用于 RN/WN Variobase®, 金基台和临时基台, 长度6.7 毫米	TAN	048.356
替代螺丝				
SCS Occlusal Screw SCS 骀面螺丝		适用于 RN/WN 钛合金临时基台, 长度7.6毫米	TAN	049.181

WN

修复配件

粘接固位修复, 粘接基台

产品	图片	描述	材料	产品编号
粘接基台				
WN synOcta®(八角) Cementable Abutment WN 粘接基台		包含中央螺丝，高度5.5毫米	纯钛	048.606
脱蜡基底				
WN synOcta®(八角) Plastic Coping WN 塑料基底		适用于 048.606，牙冠，高度7.35毫米	POM	048.243
		适用于 048.606，牙桥，高度7.35毫米		048.244
辅助配件				
Transfer Aid 转移辅助配件		适用于 048.606，高度6.5毫米，4枚装	PC	048.054V4

粘接固位修复，角度基台

产品	图片	描述	材料	产品编号
角度 15° 基台				
WN synOcta®(八角) Angled Abutment, 15° WN 角度基台，角度 15°		A型，包含中央螺丝，高度5.5毫米	纯钛/TAN	048.608
		B型，包含中央螺丝，高度5.5毫米		048.609
Plastic Shoulder 塑料肩台基底		适用于WN,synOcta® (八角)，成角15° ,高度4毫米	POM	048.678
Transfer Aid 转移辅助配件		适用于WN,synOcta® (八角)，成角15° ,高度5毫米	PC	048.032

螺丝固位或粘接固位修复，多能基台，CAD/CAM

产品	图片	描述	材料	产品编号
Variobase® 基台				
WN Variobase® for Crown WN 多能基台，适用于单冠		包含中央螺丝048.356，Ø 7 毫米，高度4.5毫米	TAN	048.711
		包含中央螺丝048.356，Ø 7 毫米，高度6.5毫米		022.0023
WN Variobase® C WN 多能基台		适用于Dentsply®Sirona®，包含中央螺丝022.0045，Ø 6.7 毫米，高度8.2毫米		022.0020
替代螺丝				
synOcta®(八角) Basal Screw 中央螺丝		适用于 RN/WN Variobase®，金基台和临时基台，长度 6.7 毫米	TAN	048.356
Basal Screw RN/WN Variobase® C RN/WN 中央螺丝，适用于多能基台 C		适用于 Dentsply® Sirona®，RN/WN Variobase® C，长度 7.2 毫米		022.0045

螺丝固位或粘接固位修复，金基台

产品	图片	描述	材料	产品编号
金基台				
WN synOcta® (八角) Gold Abutment WN 金基台		适用于牙冠，带螺丝048.356,高度14.1毫米	Ceramicor®/ POM	048.644
替代螺丝				
synOcta® (八角) Basal Screw 中央螺丝		适用于 RN/WN Variobase®, 金基台和临时基台，长度 6.7 毫米	TAN	048.356

WN

修复配件

螺丝固位或粘接固位修复，多能基台，铸造

产品	图片	描述	材料	产品编号
Variobase® 基台				
WN Variobase® for Crown WN 多能基台，适用于单冠		包含中央螺丝 048.356, Ø 7 毫米，高度 4.5 毫米	TAN	048.711
		包含中央螺丝 048.356, Ø 7 毫米，高度 6.5 毫米		022.0023
脱蜡基底				
WN Burn-out Copings WN 脱蜡基底		用于 Variobase® 基台，Ø 7毫米，高度 11.1 毫米	POM	048.269
		适用于 Variobase®基台，适用于牙冠高度 6.5 毫米，Ø 7 毫米，高度 11.1 毫米		023.0016
		4枚装		023.0016V4
替代螺丝				
synOcta®(八角) Basal Screw 中央螺丝		适用于 RN/WN Variobase® ，金基台和临时基台，长度 6.7 毫米	TAN	048.356

Ceramicor® 是一种非氧化贵金属合金，适用于铸造贵金属合金。

螺丝固位修复，synOcta® (八角) 1.5

产品	图片	描述	材料	产品编号
基台				
WN synOcta® 1.5(八角) Screw-retained Abutment WN 1.5 螺丝固位基台		包含中央螺丝，高度1.5毫米	TAN	048.603
金基底				
WN synOcta®(八角) Gold Coping WN 金基底		适用于synOcta® (八角) 1.5 螺丝固位基台(包括铸模辅件)，牙冠，高度4.75/12毫米，合金重量0.898克	Ceramicor®/ POM	048.639
		适用于synOcta® (八角) 1.5 螺丝固位基台(包括铸模辅件)，牙桥，高度4.75/12毫米，合金重量0.879克		048.638
脱蜡基底				
WN synOcta®(八角) Plastic Coping WN 塑料基底		适用于牙冠，高度10毫米	POM	048.667
		适用于牙桥，高度10毫米		048.666
研磨基底				
WN synOcta®(八角) Meso Milling Cylinder WN Meso 式研磨基底		适用于synOcta® (八角) 1.5 螺丝固位基台，牙冠，高度10毫米，Ø 10毫米	纯钛	048.573
螺丝				
SCS Occlusal Screw SCS 骀面螺丝		用于 RN/WN 固位基底的骀面螺丝，长度4.4毫米	纯钛	048.350
		4枚装		048.350V4

Ceramicor® 是一种非氧化贵金属合金，适用于铸造贵金属合金。

适用于印模和制作工作模型，实心基台

产品	图片	描述	材料	产品编号
印模和制作工作模型				
WN Impression Cap WN 印模帽		白色，适用于非开窗式印模，高度8毫米	POM	048.013
Positioning Cylinder 定位柱		适用于实心基台048.545，绿色，高度10毫米		048.065
		适用于实心基台048.546，褐色，高度10毫米		048.066
Implant Analog 替代体		适用于实心基台048.545，长度16毫米	铝	048.165
		适用于实心基台048.546，长度17.5毫米		048.166
Shoulder Analog 肩台替代体		适用于WN实心基台，长度10毫米		
临时修复				
WN Temporary Coping WN 实心临时基台		适用于WN实心基台，适用于牙冠，高度7.3毫米	PMMA	048.657
		适用于WN实心基台，适用于牙桥，高度7.3毫米		048.656
保护帽				
WN Protective Cap WN 保护帽		适用于048.545实心基台，高度6毫米	PEEK	048.051
		适用于048.546实心基台，高度7.5毫米		048.052

WN

修复配件

实心基台粘接固位修复

产品	图片	描述	材料	产品编号
实心基台				
WN Solid Abutment 6° WN 实心基台 6°		绿色，高度4毫米	纯钛	048.545
		褐色，高度5.5毫米		048.546
脱蜡基底				
Plastic Coping 塑料基底		适用于WN实心基台，适用于牙冠，高度7毫米	POM	048.247
		适用于WN实心基台，适用于牙桥，高度7毫米		048.248
辅助器械				
SCS Screwdriver SCS 螺丝刀		适用于棘轮扳手，超短型，长度15毫米	不锈钢	046.400
		适用于棘轮扳手，短型，长度21毫米		046.401
		适用于棘轮扳手，长型，长度27毫米		046.402

WN

修复配件

LOCATOR®

产品	图片	描述	材料	产品编号
印模和制作工作模型				
LOCATOR® Impression Coping 印模帽		高度4毫米，4枚装	铝外壳，含 LDPE 衬垫	048.197V4
LOCATOR® Female Analog 基台替代体		长度10毫米，Ø 5 毫米，4枚装	铝	048.198V4
钛基台				
WN LOCATOR® Abutment WN 基台		高度1毫米	TAV	048.183
		高度2毫米		048.184
		高度3毫米		048.185
		高度4毫米		048.186
		高度5毫米		048.187

 0086 制造商
Zest Anchors, Inc.
Escondido, CA 92029
美国

WN
修复配件

LOCATOR®

产品	图片	描述	材料	产品编号
技工室套装，可以达到20°的分散补偿				
LOCATOR® Male Processing Package 阳性部件套装		内含：基底帽(∅ 5.5毫米，高度2.5毫米)，含阳性黑色衬垫(高度1.9毫米)，封闭环(048.196)，透明阳性垫片(048.190)，粉色阳性垫片(048.191)，蓝色阳性垫片(048.192)， 2枚装	纯钛/尼龙/ LDPE	048.189V2
LOCATOR® Replacement Male 阳性垫片		透明，0°-10°，高度1.7毫米，2.27kg，4枚装	尼龙	048.190V4
		粉色，轻微固位力，0°-10°，高度1.7毫米，1.36kg，4枚装		048.191V4
		蓝色，超低固位力，0°-10°，高度1.7毫米，0.68kg，4枚装		048.192V4
技工室套装，可以达到40°的分散补偿				
LOCATOR® Male Processing Package 阳性部件套装		加大角度补偿范围，内含：基底帽(∅ 5.5毫米，高度2.5毫米)，含阳性黑色衬垫(高度1.9毫米)，封闭环(048.196)，绿色阳性垫片(048.193)，橙色阳性垫片(048.188)，红色阳性垫片(048.194)， 2枚装	纯钛/尼龙/ LDPE	048.182V2
LOCATOR® Replacement Male 阳性垫片		绿色，加大角度补偿范围，10°-20°，高度1.7毫米，1.82kg，4枚装	尼龙	048.193V4
		橙色，轻微固位力，加大角度补偿范围，10°-20°，高度1.7毫米，0.91kg，4枚装		048.188V4
		红色，超低固位力，加大角度补偿范围， 10°-20°，高度1.7毫米，0.45kg，4枚装		048.194V4
		灰色，加大角度补偿范围，零(0)固位力，10°-20°，高度1.7毫米，4枚装		048.181V4

CE

0086

制造商

Zest Anchors, Inc.

Escondido, CA 92029

美国

LOCATOR® 配件

产品	图片	描述	材料	产品编号
LOCATOR® Parallel Post 平行柱		高度8毫米,4枚装	LDPE	048.199V4
LOCATOR® Processing Male, black 阳性垫片，黑色		高度1.9毫米,4枚装		048.195V4
LOCATOR® Block-out Spacer 封闭环		厚度0.4毫米,20枚装	硅树脂 橡胶	048.196V20
LOCATOR® Abutment Holder Sleeve 基台固位套		用于将基台紧扣到工具上， 4 枚装	PSU	046.413V4

辅助配件				
LOCATOR® Core Tool 核心工具		三个部件，长度100毫米	不锈钢	046.415
LOCATOR® Screwdriver 螺丝刀		适用于棘轮扳手，短型,长度15毫米		046.416
		适用于棘轮扳手，长型,长度21毫米		046.417
LOCATOR® Angle Measurement Guide 角度尺		长度50毫米，宽度15毫米		048.200

 0086 制造商
Zest Anchors, Inc.
Escondido, CA 92029
美国

WN
修复配件

模拟基台

产品	图片	描述	材料	产品编号
修复设计套装				
Soft Tissue Level PLAN-Set 软组织水平模拟套装		包括：所有的软组织水平模拟基台，长度278毫米，宽度124毫米，高度33毫米	PC	048.904
Soft Tissue Level PLAN Cassette 软组织水平模拟套装外盒		空盒，长度278毫米，宽度124毫米，高度33毫米		048.905
NNC 模拟基台				
NNC PLAN Cementable Abutment, straight NNC 直型模拟粘接基台		适用于牙冠和牙桥，高度5.7毫米，4枚装	POM	048.943V4
NNC PLAN Cementable Abutment, angled 15° NNC 成角 15° 模拟角度基台		A型，高度5.9毫米，4枚装		048.944V4
		B型，高度5.9毫米，4枚装		048.945V4
NNC PLAN LOCATOR® Abutment* NNC 模拟 LOCATOR®基台		高度 1 毫米，4枚装		048.946V4
		高度 2 毫米，4枚装		048.947V4
		高度 3 毫米，4枚装		048.948V4
		高度 4 毫米，4枚装		048.949V4
		高度 5 毫米，4枚装		048.950V4
NNC PLAN Solid Abutment NNC 模拟实心基台		高度 4 毫米，4枚装		048.952V4
		高度 5.5 毫米，4枚装		048.953V4
		高度 7 毫米，4枚装		048.954V4

* NNC 模拟 Novaloc® 基台(048.951V4) 亦兼容于 LOCATOR® 部件。

模拟基台

产品	图片	描述	材料	产品编号
RN 模拟基台				
RN synOcta®(八角) 1.5 PLAN Abutment RN 1.5 模拟基台		高度1.5毫米,4枚装	POM	048.929V4
RN synOcta®(八角) Cementable PLAN Abutment RN 模拟粘接基台		高度5.5毫米,4枚装		048.930V4
RN synOcta® (八角) Transversal PLAN Abutment RN 模拟横向基台		高度4毫米,4枚装		048.921V4
RN synOcta®(八角) Angled PLAN Abutment, 15° RN 角度 15° 模拟角度基台		A型, 短型, 高度5.7毫米, 4枚装		048.916V4
		B型, 短型, 高度5.7毫米, 4枚装		048.917V4
RN synOcta®(八角) Angled PLAN Abutment, 20° RN 角度 20° 模拟角度基台		A型, 短型, 高度5.7毫米, 4枚装		048.918V4
		B型, 短型, 高度5.7毫米, 4枚装		048.919V4
RN synOcta®(八角) Angled PLAN Abutment, 15° RN 角度 15° 模拟角度基台		A型, 长型, 高度6.7毫米, 4枚装		048.922V4
		B型, 长型, 高度6.7毫米, 4枚装		048.923V4
RN synOcta®(八角) Angled PLAN Abutment, 20° RN 角度 20° 模拟角度基台		A型, 长型, 高度6.7毫米, 4枚装		048.924V4
		B型, 长型, 高度6.7毫米, 4枚装		048.925V4
RN PLAN LOCATOR® Abutment RN 模拟 LOCATOR® 基台		高度1毫米, 4枚装		048.275V4
		高度2毫米, 4枚装		048.276V4
		高度3毫米, 4枚装		048.277V4
		高度4毫米, 4枚装		048.278V4
		高度5毫米, 4枚装		048.279V4
		高度1-6 毫米, 4枚装		048.280V4
RN PLAN Solid Abutment RN 模拟实心基台		高度4毫米, 4枚装		048.926V4
		高度5.5毫米, 4枚装		048.927V4
		高度7毫米, 4枚装		048.928V4

WN

修复配件

模拟基台

产品	图片	描述	材料	产品编号
WN 模拟基台				
WN synOcta®(八角) 1.5 PLAN Abutment WN 1.5 模拟基台		高度1.5毫米，4枚装	POM	048.931V4
WN synOcta®(八角) Cementable PLAN Abutment WN 模拟粘接基台		高度5.5毫米，4枚装		048.932V4
WN synOcta®(八角) Angled PLAN Abutment WN 成角 15° 模拟角 度基台		A型，高度5.5毫米，4枚装		048.933V4
		B型，高度5.5毫米，4枚装		048.934V4
WN PLAN LOCA- TOR® Abutment WN 模拟 LOCATOR® 基台		高度1毫米，4枚装		048.283V4
		高度2毫米，4枚装		048.284V4
		高度3毫米，4枚装		048.285V4
		高度4毫米，4枚装		048.286V4
		高度5毫米，4枚装		048.287V4
WN PLAN Solid Abutment WN 模拟实心基台6°		高度4毫米，4枚装		048.938V4
		高度5.5毫米，4枚装		048.939V4

WN

修复配件

骨水平 NC

窄 CrossFit®(窄十字锁合)

印模和制作工作模型

产品	图片	描述	材料	产品编号
转移杆和替代体				
NC Impression Post NC 转移杆		适用于开窗式印模，短型，包括导向螺丝，长度16.5毫米	TAN	025.2202
		适用于开窗式印模，长型，包括导向螺丝，长度24毫米		025.2205
		适用于非开窗式印模，包括1个导向螺丝和2个印模帽，长度12.3毫米	TAN/POM	025.2201
NC Implant Analog NC 种植体替代体		长度12毫米	纯钛	025.2101
数字化印模				
CARES® NC Mono Scanbody CARES® NC 扫描杆		包括固定螺丝，Ø 3.5 毫米，高度10毫米	PEEK/TAN	025.2915
骀记录转移体				
NC Bite Registration Aid NC 骀记录转移体		短型，高度8毫米，4枚装	POM	025.2208-04
		长型，高度12毫米，4枚装		025.2212-04

临时修复

产品	图片	描述	材料	产品编号
临时基台				
NC Temporary Abutment VITA CAD-Temp® NC 临时基台		高度 12 毫米， Ø 5 毫米	PEEK/TAN	024.2370
临时基台				
NC Temporary Abutment NC 临时基台		适用于牙冠，包括中央螺丝025.2900，高度 11 毫米， Ø 3.5 毫米	TAN	024.2371
		适用于牙桥，包括中央螺丝025.2900，高度 11 毫米， Ø 3.5 毫米		024.2375
替代螺丝				
NC Basal Screw NC 中央螺丝		适用于 NC 解剖基台，NC Variobase® 冠用，NC 金基台，NC 杆卡金基台，长度 7.9 毫米	TAN	025.2900

粘接固位修复，解剖基台

产品	图片	描述	材料	产品编号
解剖式钛基台				
NC Anatomic Abutment NC 解剖基台		直型，包括中央螺丝025.2900，穿龈高度2毫米	纯钛	022.2102
		直型，包括中央螺丝025.2900，穿龈高度3.5毫米		022.2104
		成角，角度15°，包括中央螺丝025.2900，穿龈高度2毫米		022.2152
		成角，角度15°，包括中央螺丝025.2900，穿龈高度3.5毫米		022.2154
替代螺丝				
NC Basal Screw NC 中央螺丝		适用于 NC 解剖基台，NC Variobase® 冠用，NC 金基台，NC 杆卡金基台，长度 7.9 毫米	TAN	025.2900

粘接固位修复，粘接基台

产品	图片	描述	材料	产品编号
印模帽 (基台水平)，适用于单个牙冠修复				
NC Impression Cap NC 印模帽		适用于粘接基台，Ø 3.5毫米，基台高度4毫米，4枚装	POM	025.2203-04
		适用于粘接基台，Ø 5毫米，基台高度4毫米，4枚装		025.2204-04
		适用于粘接基台，Ø 5毫米，基台高度5.5毫米，4枚装		025.2220-04
替代体				
NC Implant Analog NC基台替代体		适用于粘接基台，Ø 3.5毫米，基台高度4毫米	TAN	025.2112
		适用于粘接基台，Ø 3.5毫米，基台高度5.5毫米		025.2116
		适用于粘接基台，Ø 5毫米，基台高度4毫米		025.2122
		适用于粘接基台，Ø 5毫米，基台高度5.5毫米		025.2126
临时修复				
NC Temporary Coping NC 塑料临时基底		适用于粘接基台，Ø 3.5毫米	PMMA	024.2346
		适用于粘接基台，Ø 5毫米		024.2356
粘接式保护帽				
NC Protective Cap NC 粘接式保护帽		适用于粘接基台，Ø 3.5毫米，基台高度4毫米，4枚装	PEEK	024.2312-04
		适用于粘接基台，Ø 3.5毫米，基台高度5.5毫米，4枚装		024.2316-04
		适用于粘接基台，Ø 5毫米，基台高度4毫米，4枚装		024.2322-04
		适用于粘接基台，Ø 5毫米，基台高度5.5毫米，4枚装		024.2326-04
转移辅助配件				
NC Transfer Aid NC 转移辅助配件		适用于粘接基台，Ø 3.5毫米，4枚装	PC	025.2911-04
		适用于粘接基台，Ø 5毫米，4枚装		025.2912-04

粘接固位修复，粘接基台

产品	图片	描述	材料	产品编号
粘接基台 Ø3.5毫米				
NC Cementable Abutment NC 粘接基台		包括中央螺丝025.2908, Ø 3.5毫米, 穿龈高度1毫米, 基台高度4毫米	纯钛	022.2311
		包括中央螺丝025.2908, Ø 3.5毫米, 穿龈高度2毫米, 基台高度4毫米		022.2312
		包括中央螺丝025.2908, Ø 3.5毫米, 穿龈高度3毫米, 基台高度4毫米		022.2313
		包括中央螺丝025.2908, Ø 3.5毫米, 穿龈高度1毫米, 基台高度5.5毫米		022.2315
		包括中央螺丝025.2908, Ø 3.5毫米, 穿龈高度2毫米, 基台高度5.5毫米		022.2316
		包括中央螺丝025.2908, Ø 3.5毫米, 穿龈高度3毫米, 基台高度5.5毫米		022.2317
粘接基台 Ø5毫米				
NC Cementable Abutment NC 粘接基台		包括中央螺丝025.2908, Ø 5毫米, 穿龈高度1毫米, 基台高度4毫米	纯钛	022.2321
		包括中央螺丝025.2908, Ø 5毫米, 穿龈高度2毫米, 基台高度4毫米		022.2322
		包括中央螺丝025.2908, Ø 5毫米, 穿龈高度3毫米, 基台高度4毫米		022.2323
		包括中央螺丝025.2908, Ø 5毫米, 穿龈高度1毫米, 基台高度5.5毫米		022.2325
		包括中央螺丝025.2908, Ø 5毫米, 穿龈高度2毫米, 基台高度5.5毫米		022.2326
		包括中央螺丝025.2908, Ø 5毫米, 穿龈高度3毫米, 基台高度5.5毫米		022.2327
替代螺丝				
NC Basal Screw NC 中央螺丝		适用于 NC 粘接基台, 长度 7.9 毫米	TAN	025.2908
可铸基底, Ø 3.5 毫米				
NC Coping NC 基底		适用于粘接基台, 牙冠用, Ø 3.5毫米	POM	023.2316
		适用于粘接基台, 桥体用, Ø 3.5毫米		023.2310
可铸基底, Ø 5 毫米				
NC Coping NC 基底		适用于粘接基台, 牙冠用, Ø 5毫米	POM	023.2326
		适用于粘接基台, 桥体用, Ø 5毫米		023.2320

NC

修复配件

螺丝固位或粘接固位修复，多能基台，CAD/CAM

产品	图片	描述	材料	产品编号
Variobase® 基台				
NC Variobase® for Crown NC 多能基台，适用于单冠		包括中央螺丝 025.2900，Ø3.8 毫米，高度 3.5 毫米	TAN	025.2921
		包括中央螺丝 025.2900，Ø3.8 毫米，高度 5.5 毫米		022.0027
NC Variobase® C NC 多能基台		适用于Dentsply®Sirona®，包括中央螺丝 025.2900，Ø4.5 毫米，高度 4.7 毫米		022.0043
替代螺丝				
NC Basal Screw NC 中央螺丝		适用于 NC 解剖基台，NC Variobase® 多能基台冠用，Variobase®C多能基台，NC 金基台，NC 杆卡金基台，长度 7.9 毫米	TAN	025.2900

螺丝固位或粘接固位修复，金基台

产品	图片	描述	材料	产品编号
金基台				
NC Gold Abutment NC 金基台		适用于牙冠，包括中央螺丝025.2900，高度3.7毫米，合金重量0.313克	Ceramicor®/ POM	022.2410
替代螺丝				
NC Basal Screw NC 中央螺丝		适用于 NC 解剖基台，NC Variobase® 多能基台冠用，Variobase®C多能基台，NC 金基台，NC 杆卡金基台，长度 7.9 毫米	TAN	025.2900
打磨辅件				
NC Polishing Aid NC 打磨辅件		长度 16 毫米	不锈钢	025.2920
		4枚装		025.2920-04

Ceramicor® 是一种非氧化贵金属合金，适用于铸造贵金属合金。

螺丝固位或粘接固位修复，多能基台，铸造

产品	图片	描述	材料	产品编号
Variobase® 基台				
NC Variobase® for Crown NC 多能基台，适用于单冠		包括中央螺丝 025.2900，Ø3.8 毫米，高度 3.5 毫米	TAN	025.2921
		包括中央螺丝 025.2900，Ø3.8 毫米，高度 5.5 毫米		022.0027
脱蜡基底				
NC Burn-out Coping NC 脱蜡基底		适用于牙冠高度 3.5毫米 Variobase 基台，Ø 3.8毫米，高度12毫米	POM	023.2756
		4枚装		023.2756-04
		适用于牙冠高度 5.5毫米 Variobase 基台，Ø 3.8毫米，高度12毫米		023.0018
		4枚装		023.0018V4
替代螺丝				
NC Basal Screw NC 中央螺丝		适用于 NC 解剖基台，NC Variobase® 冠用多能基台，Variobase®C多能基台，NC 金基台，NC 杆卡金基台，长度 7.9 毫米	TAN	025.2900

NC

修复配件

螺丝固位修复，螺丝固位基台

产品	图片	描述	材料	产品编号
冠用转移杆（基台水平）				
Impression Posts for Open-tray Impression 开窗式转移杆		适用于螺丝固位基台，基台水平，Ø 3.5 毫米	TAN	025.2243
		适用于螺丝固位基台，基台水平，Ø 4.6 毫米		025.2244
Impression Posts for Closed-tray Impression 非开窗式转移杆		适用于螺丝固位基台，基台水平，Ø 3.5 毫米	TAN/POM	025.2245
		适用于螺丝固位基台，基台水平，Ø 4.6 毫米		025.2246
桥用转移杆（基台水平）				
Impression Posts for Open-tray Impression 开窗式转移杆		适用于螺丝固位基台，基台水平，Ø 3.5 毫米	TAN	025.0011
		适用于螺丝固位基台，基台水平，Ø 4.6 毫米		025.0012
Impression Posts for Closed-tray Impression 非开窗式转移杆		适用于螺丝固位基台，基台水平，Ø 3.5 毫米	TAN/POM	025.0013
		适用于螺丝固位基台，基台水平，Ø 4.6 毫米		025.0014
数字化印模				
CARES® Mono Scanbody CARES® 扫描杆		适用于螺丝固位基台，基台水平，包括固定螺丝，Ø 3.5 毫米	PEEK	025.0000
		适用于螺丝固位基台，基台水平，包括固定螺丝，Ø 4.6 毫米		025.0001
替代体				
NC Analog NC 替代体		适用于螺丝固位基台，直型，Ø 3.5 毫米	TAN	023.2754
		适用于螺丝固位基台，直型，无牙颌，Ø 3.5 毫米		025.0049
NC/RC Analog NC/RC 替代体		适用于螺丝固位基台，直型，Ø 4.6 毫米		023.4756
		适用于螺丝固位基台，直型，无牙颌，Ø 4.6 毫米		025.0050
		适用于螺丝固位基台，成角，角度 17°/30°，Ø 4.6 毫米		023.4757
技工辅助配件				
Polishing Aid 抛光保护帽		适用于螺丝固位基台，Ø 3.5 毫米	TAN	025.0004
		4 枚装		025.0004V4
		适用于螺丝固位基台，Ø 4.6 毫米		025.0005
		4 枚装		025.0005V4
Lab Processing Screw 技工辅助螺丝		适用于螺丝固位基台，长度 20 毫米	不锈钢	025.0006
		适用于螺丝固位基台，长度 10 毫米		025.0052

螺丝固位修复，螺丝固位基台

产品	图片	描述	材料	产品编号
保护帽				
NC Protective Cap NC 保护帽		适用于螺丝固位基台，包括殆面螺丝023.4763， Ø 3.5毫米，高度 5 毫米，4枚装	PEEK	024.2323-04
		适用于螺丝固位基台，包括殆面螺丝023.4763， Ø 3.5毫米，高度6.5毫米，4枚装		024.2324-04
		适用于螺丝固位基台，包括殆面螺丝023.4763， Ø 3.5 毫米，高度8毫米，4枚装		024.2325-04
NC/RC Protective Cap NC/RC 保护帽		适用于螺丝固位基台，包括殆面螺丝023.4763， Ø 4.6 毫米，高度 5.1 毫米，4枚装		024.4323-04
		适用于螺丝固位基台，包括殆面螺丝023.4763， Ø 4.6 毫米，高度 6.6 毫米，4枚装		024.4324-04
		适用于螺丝固位基台，包括殆面螺丝023.4763， Ø 4.6 毫米，高度 8.1 毫米，4枚装		024.4325-04
辅助部件				
Straumann® Planning Template 种植体轴向引导器		Straumann® Pro Arch 种植体轴向引导器		026.0016
Transfer And Alignment Pin 转移导向杆		适用于螺丝固位基台	TAN	025.0009

NC

修复配件

螺丝固位修复，螺丝固位基台

产品	图片	描述	材料	产品编号
螺丝固位基台 0°				
NC Screw-retained Abutment NC 螺丝固位基台		直型，包括转移导向杆 025.0009 以及螺丝 023.4749， Ø 3.5毫米，穿龈高度 1 毫米	TAN	022.2745P
		直型，包括转移导向杆 025.0009 以及螺丝 023.4750， Ø 3.5毫米，穿龈高度 2.5 毫米		022.2746P
		直型，包括转移导向杆 025.0009 以及螺丝 023.4760， Ø 3.5毫米，穿龈高度 4 毫米		022.2753P
螺丝固位基台 0°				
NC Screw-retained Abutment NC 螺丝固位基台		直型，包括转移导向杆 025.0009 以及螺丝 023.4749， Ø 4.6 毫米，穿龈高度 1 毫米	TAN	022.2747P
		直型，包括转移导向杆 025.0009 以及螺丝 023.4750， Ø 4.6 毫米，穿龈高度 2.5 毫米		022.2748P
		直型，包括转移导向杆 025.0009 以及螺丝 023.4760， Ø 4.6 毫米，穿龈高度 4 毫米		022.2754P

螺丝固位修复，螺丝固位基台

产品	图片	描述	材料	产品编号
螺丝固位基台 17°				
NC Screw-retained Abutment NC 螺丝固位基台		成角，角度 17°，A 型，包括转移导向杆 025.0009 以及螺丝 025.0002，Ø 4.6 毫米，穿龈高度 2.5 毫米	TAN	022.2749P
		成角，角度 17°，B 型，包括转移导向杆 025.0009 以及螺丝 025.0002，Ø 4.6 毫米，穿龈高度 2.5 毫米		022.2750P
		成角，角度 17°，A 型，包括转移导向杆 025.0009 以及螺丝 025.0002，Ø 4.6 毫米，穿龈高度 4 毫米		022.2755P
		成角，角度 17°，B 型，包括转移导向杆 025.0009 以及螺丝 025.0002，Ø 4.6 毫米，穿龈高度 4 毫米		022.2756P
		成角，角度 17°，A 型，包括转移导向杆 025.0009 以及螺丝 025.0002，Ø 4.6 毫米，穿龈高度 5.5 毫米		022.0010P
		成角，角度 17°，B 型，包括转移导向杆 025.0009 以及螺丝 025.0002，Ø 4.6 毫米，穿龈高度 5.5 毫米		022.0011P
螺丝固位基台 30°				
NC Screw-retained Abutment NC 螺丝固位基台		成角，角度 30°，A 型，包括转移导向杆 025.0009 以及螺丝 025.0002，Ø 4.6 毫米，穿龈高度 2.5 毫米	TAN	022.2751P
		成角，角度 30°，B 型，包括转移导向杆 025.0009 以及螺丝 025.0002，Ø 4.6 毫米，穿龈高度 2.5 毫米		022.2752P
		成角，角度 30°，A 型，包括转移导向杆 025.0009 以及螺丝 025.0002，Ø 4.6 毫米，穿龈高度 4 毫米		022.2757P
		成角，角度 30°，B 型，包括转移导向杆 025.0009 以及螺丝 025.0002，Ø 4.6 毫米，穿龈高度 4 毫米		022.2758P
		成角，角度 30°，A 型，包括转移导向杆 025.0009 以及螺丝 025.0002，Ø 4.6 毫米，穿龈高度 5.5 毫米		022.0012P
		成角，角度 30°，B 型，包括转移导向杆 025.0009 以及螺丝 025.0002，Ø 4.6 毫米，穿龈高度 5.5 毫米		022.0013P
替代螺丝				
NC/RC Screw NC/RC 中央螺丝		适用于 NC/RC 螺丝固位基台 0°，穿龈高度 1 毫米，长度 8.8 毫米	TAN	023.4749
		适用于 NC/RC 螺丝固位基台 0°，穿龈高度 2.5 毫米，长度 10.3 毫米		023.4750
		适用于 NC/RC 螺丝固位基台 0°，穿龈高度 4 毫米，长度 11.8 毫米		023.4760
		适用于 NC/RC 螺丝固位基台 17° 和 30°，长度 7.2 毫米		025.0002

NC

修复配件

螺丝固位修复，螺丝固位基台

产品	图片	描述	材料	产品编号
钛基底				
NC Coping NC 基底		适用于螺丝固位基台，牙冠用，包括殆面螺丝 023.4763， $\varnothing$ 3.5 毫米，高度 11 毫米	纯钛	023.2747
		适用于螺丝固位基台，牙桥用，包括殆面螺丝 023.4763， $\varnothing$ 3.5 毫米，高度 11 毫米		023.2749
NC/RC Coping NC/RC 基底		适用于螺丝固位基台，牙冠用，包括殆面螺丝 023.4763， $\varnothing$ 4.6 毫米，高度 11 毫米		023.4747
		适用于螺丝固位基台，牙桥用，包括殆面螺丝 023.4763， $\varnothing$ 4.6 毫米，高度 11 毫米		023.4751
金基底				
NC/RC Coping NC/RC 基底		适用于螺丝固位基台，牙冠用，包括殆面螺丝 023.4763，合金重量 0.359克， $\varnothing$ 3.5 毫米	Ceramicor®/ POM	023.2751
		适用于螺丝固位基台，牙桥用，包括殆面螺丝 023.4763，合金重量 0.356克， $\varnothing$ 3.5 毫米		023.2752
		适用于螺丝固位基台，牙冠用，包括殆面螺丝 023.4763，合金重量 0.440克， $\varnothing$ 4.6 毫米		023.4753
		适用于螺丝固位基台，牙桥用，包括殆面螺丝 023.4763，合金重量 0.469克， $\varnothing$ 4.6 毫米		023.4754
脱蜡基底				
NC Burn-out Coping NC 脱蜡基底		适用于螺丝固位基台，牙冠用，包括殆面螺丝 023.4763， $\varnothing$ 3.5 毫米，高度 11 毫米	POM	023.2748
		适用于螺丝固位基台，适用于桥和杆卡，包括殆面螺丝 023.4763， $\varnothing$ 3.5 毫米，高度 11 毫米		023.2755
NC/RC Burn-out Coping NC/RC 脱蜡基底		适用于螺丝固位基台，牙冠用，包括殆面螺丝 023.4763， $\varnothing$ 4.6 毫米，高度 11 毫米		023.4748
		适用于螺丝固位基台，适用于桥和杆卡，包括殆面螺丝 023.4763， $\varnothing$ 4.6 毫米，高度 11 毫米		023.4758
替代螺丝				
NC/RC Occlusal Screw NC/RC 殆面螺丝		适用于螺丝固位 NC/RC 钛基台，金基台，脱蜡基底，长度 3.7 毫米	TAN	023.4763

NC

修复配件

Ceramicor® 是一种非氧化贵金属合金，适用于铸造贵金属合金。

应用于螺丝固位的杆卡修复，基底，钛

产品	图片	描述	材料	产品编号
杆卡钛基底 Ø 3.5 mm				
NC Coping NC 基底		适用于螺丝固位基台，包括殆面螺丝 023.4763， Ø 3.5 毫米，高度 5.5 毫米	纯钛	023.2750
杆卡钛基底 Ø 4.6 mm				
NC/RC Coping NC/RC 基底		适用于螺丝固位基台，包括殆面螺丝 023.4763， Ø 4.6 毫米，高度 5.5 毫米	纯钛	023.4752
替代螺丝				
NC/RC Occlusal Screw NC/RC 殆面螺丝		适用于螺丝固位 NC/RC 钛基台，金基台，脱蜡基底， 长度 3.7 毫米	TAN	023.4763

应用于螺丝固位的杆卡修复，基底，金

产品	图片	描述	材料	产品编号
杆卡金基底 Ø 3.5 mm				
NC Coping NC 基底		适用于螺丝固位基台，适用于杆卡，包括殆面螺丝 023.4763，合金重量 0.681克，Ø 3.5 毫米，高度 5.5 毫米	Ceramicor®	023.2753
杆卡金基底 Ø 4.6 mm				
NC/RC Coping NC/RC 基底		适用于螺丝固位基台，适用于杆卡，包括殆面螺丝 023.4763，合金重量 0.744 克，Ø 4.6 毫米，高度 5.5 毫米	Ceramicor®	023.4755
替代螺丝				
NC/RC Occlusal Screw NC/RC 殆面螺丝		适用于螺丝固位 NC/RC 钛基台，金基台，脱蜡基底， 长度 3.7 毫米	TAN	023.4763

应用于螺丝固位的杆卡修复，基底，塑料

产品	图片	描述	材料	产品编号
杆卡脱蜡基底，适用于 Ø 3.5 mm				
NC Burn-out Coping NC 脱蜡基底		适用于螺丝固位基台，适用于桥和杆卡，包括殆面螺丝 023.4763，Ø 3.5 毫米，高度 11 毫米	POM	023.2755
杆卡脱蜡基底，适用于 Ø 4.6 mm				
NC/RC Burn-out Coping NC/RC 脱蜡基底		适用于螺丝固位基台，适用于桥和杆卡，包括殆面螺丝 023.4763，Ø 4.6 毫米，高度 11 毫米	POM	023.4758
替代螺丝				
NC/RC Occlusal Screw NC/RC 殆面螺丝		适用于螺丝固位 NC/RC 钛基台，金基台，脱蜡基底， 长度 3.7 毫米	TAN	023.4763

Ceramicor® 是一种不氧化贵金属合金，适用于铸造贵金属合金。

钛杆卡基台

产品	图片	描述	材料	产品编号
NC Abutment NC 基台		适用于杆卡，包括中央螺丝 025.2900，高度 10 毫米	纯钛	022.2702

替代螺丝				
NC Basal Screw NC 中央螺丝		适用于 NC 解剖基台，NC Variobase® 冠用，NC 金基台，NC 杆卡金基台，NC 杆卡钛基台，长度 7.9 毫米	TAN	025.2900

金杆卡基台

产品	图片	描述	材料	产品编号
NC Abutment, one-part NC 基台		适用于杆卡，金，包括中央螺丝 025.2900，高度 10 毫米，合金重量 1.5 克	Ceramicor®	022.2701

替代螺丝				
NC Basal Screw NC 中央螺丝		适用于 NC 解剖基台，NC Variobase® 冠用，NC 金基台，NC 杆卡金基台，NC 杆卡钛基台，长度 7.9 毫米	TAN	025.2900

Ceramicor® 是一种非氧化贵金属合金，适用于铸造贵金属合金。

NC

修复配件

螺丝固位修复，多基基台

NC

修复配件

产品	图片	描述	材料	产品编号
多基基台				
NC Multi-base abutment straight NC 直型多基基台		Ø 3.5 毫米，穿龈高度 1 毫米	TAN	022.2731
		Ø 3.5 毫米，穿龈高度 2.5 毫米		022.2732
		Ø 3.5 毫米，穿龈高度 4 毫米		022.2734
		Ø 4.5 毫米，穿龈高度 1 毫米		022.2741
		Ø 4.5 毫米，穿龈高度 2.5 毫米		022.2742
		Ø 4.5 毫米，穿龈高度 4 毫米		022.2744
NC Multi-base abutment angled 25° NC 成角 25°多基基台		A 型，Ø 4 毫米，穿龈高度 2.5 毫米		022.2782
		B 型，Ø 4 毫米，穿龈高度 2.5 毫米		022.2786
印模（基台水平）和制作工作模型				
NC Impression cap NC 印模帽		适用于多基基台，直型，Ø 3.5 毫米	POM	025.2231-04
		适用于多基基台，直型，Ø 4.5 毫米		025.2241-04
NC Impression post NC 开窗式转移杆		适用于多基基台，直型，Ø 3.5 毫米	TAN	025.2232
		适用于多基基台，直型，Ø 4.5 毫米		025.2242
		适用于多基基台，成角25°，Ø 4 毫米		025.2206
NC abutment analog NC 基台替代体		适用于多基基台，直型，Ø 3.5 毫米		025.2131
		适用于多基基台，直型，Ø 4.5 毫米		025.2141
		适用于多基基台，成角25°，Ø 4 毫米		025.2106
临时修复				
NC Temporary coping NC 临时基底		适用于多基基台，直型，Ø 3.5 毫米	纯钛	024.2332
		适用于多基基台，直型，Ø 4.5 毫米		024.2342
		适用于多基基台，成角25°，Ø 4 毫米		024.2307
NC Protective cap NC 固定式保护帽		适用于多基基台，直型，Ø 3.5 毫米	PEEK/TAN	024.2331-04
		适用于多基基台，直型，Ø 4.5 毫米		024.2341-04
		适用于多基基台，成角25°，Ø 4 毫米		024.2306-04

螺丝固位修复，多基基台

产品	图片	描述	材料	产品编号
基底，适用于多基基台				
NC Coping NC 塑料合金基底		适用于多基基台上的牙桥，直型，Ø 3.5 毫米	Ceramicor®/ POM	023.2731
		适用于多基基台上的牙桥，直型，Ø 4.5 毫米		023.2741
		适用于多基基台上的牙桥，成角25°，Ø 4 毫米		023.2701
NC Coping NC 全塑料基底		适用于多基基台上的牙桥，直型，Ø 3.5 毫米	POM	023.2734
		适用于多基基台上的牙桥，直型，Ø 4.5 毫米		023.2744
		适用于多基基台上的牙桥，成角25°，Ø 4 毫米		023.2704
辅助配件				
NC/RC Occlusal screw NC/RC 殆面螺丝		仅适用于固位多基基台上的基底，长度 4.5 毫米	TAN	025.0900-04
Analog holder 替代体抓持器		长度 105 毫米	铝/钢	046.239

NC

修复配件

Ceramicor® 是一种非氧化贵金属合金，适用于铸造贵金属合金。

辅助配件

产品	图片	描述	材料	产品编号
NC Polishing Aid NC 打磨辅件		长度 16 毫米	不锈钢	025.2920
		4枚装		025.2920-04
NC Basal Screw NC 中央螺丝		适用于 NC 解剖基台，NC Variobase® 冠用，NC 金基台，NC 杆卡金基台，长度 7.9 毫米	TAN	025.2900
NC/RC Occlusal screw NC/RC 殆面螺丝		仅适用于固位多基基台上的基底，长度 4.5 毫米		025.0900-04
Analog holder 替代体抓持器		长度 105 毫米	铝/钢	046.239

NC

修复配件

杆卡修复，金

产品	图片	描述	材料	产品编号
NC Abutment, one-part NC 杆卡基台		适用于杆卡，金，包括中央螺丝025.2900，高度10毫米，合金重量1.5克	Ceramicor®	022.2701
NC Multi-base abutment straight NC 直型多基基台		Ø 3.5 毫米，穿龈高度 1 毫米	TAN	022.2731
		Ø 3.5 毫米，穿龈高度 2.5 毫米		022.2732
		Ø 3.5 毫米，穿龈高度 4 毫米		022.2734
		Ø 4.5 毫米，穿龈高度 1 毫米		022.2741
		Ø 4.5 毫米，穿龈高度 2.5 毫米		022.2742
		Ø 4.5 毫米，穿龈高度 4 毫米		022.2744
NC Multi-base abutment angled 25° NC 成角 25°多基基台		A 型，Ø 4 毫米，穿龈高度 2.5 毫米		022.2782
		B 型，Ø 4 毫米，穿龈高度 2.5 毫米		022.2786
NC Coping NC 全金基底		Ø 3.5 毫米，适用于多基基台上的杆卡，直型	Ceramicor®	023.2732
		Ø 4.5 毫米，适用于多基基台上的杆卡，直型		023.2742
		Ø 4 毫米，适用于多基基台上的杆卡，成角 25°		023.2702

辅助配件				
Activator Set 加力钥匙		适用于杆卡，三件式，长度50毫米	不锈钢	046.150
Deactivator 减力钥匙		适用于杆卡，迷你型，长度66毫米	Grilon BS/黄铜	046.151
		适用于杆卡，标准型，长度66毫米		046.152
NC Basal Screw NC 中央螺丝		适用于 NC 解剖基台，NC Variobase® 冠用，NC 金基台，NC 杆卡金基台，长度 7.9 毫米	TAN	025.2900
NC/RC Occlusal screw NC/RC 殆面螺丝		仅适用于固位多基基台上的基底，长度 4.5 毫米		025.0900-04
NC/RC Stabilization pin NC/RC 稳定杆		适用于多基基台上的杆卡，长度16毫米	不锈钢	025.0911-04


Ceramicor® 是一种非氧化贵金属合金，适用于铸造贵金属合金。

杆卡修复，钛

产品	图片	描述	材料	产品编号
NC Abutment NC 杆卡基台		适用于杆卡，包括中央螺丝025.2900，高度10毫米	纯钛	022.2702
NC Multi-base abutment straight NC 直型多基基台		Ø 3.5 毫米，穿龈高度 1 毫米	TAN	022.2731
		Ø 3.5 毫米，穿龈高度 2.5 毫米		022.2732
		Ø 3.5 毫米，穿龈高度 4 毫米		022.2734
		Ø 4.5 毫米，穿龈高度 1 毫米		022.2741
		Ø 4.5 毫米，穿龈高度 2.5 毫米		022.2742
		Ø 4.5 毫米，穿龈高度 4 毫米		022.2744
NC Multi-base abutment angled 25° NC 成角 25°多基基台		A 型，Ø 4 毫米，穿龈高度 2.5 毫米		022.2782
		B 型，Ø 4 毫米，穿龈高度 2.5 毫米		022.2786
NC Coping NC 杆卡基底		Ø 3.5 毫米，适用于多基基台上的杆卡，直型	纯钛	023.2736
		Ø 4.5 毫米，适用于多基基台上的杆卡，直型		023.2746
		Ø 4 毫米，适用于多基基台上的杆卡，成角 25°		023.2706

辅助配件				
Activator Set 加力钥匙		适用于杆卡，三件式，长度50毫米	不锈钢	046.150
Deactivator 减力钥匙		适用于杆卡，迷你型，长度66毫米	Grilon BS/黄铜	046.151
		适用于杆卡，标准型，长度66毫米		046.152
NC Basal Screw NC 中央螺丝		适用于 NC 解剖基台，NC Variobase® 冠用，NC 金基台，NC 杆卡金基台，长度 7.9 毫米	TAN	025.2900
NC/RC Occlusal screw NC/RC 骀面螺丝		仅适用于固位多基基台上的基底，长度 4.5 毫米		025.0900-04
NC/RC Stabilization pin NC/RC 稳定杆		适用于多基基台上的杆卡，长度16毫米	不锈钢	025.0911-04

杆卡辅助配件

产品	图片	描述	材料	产品编号
NC Polishing Aid NC 打磨辅件		长度 16 毫米	不锈钢	025.2920
		4枚装		025.2920-04
NC Basal Screw NC 中央螺丝		适用于 NC 解剖基台, NC Variobase® 冠用, NC 金基台, NC 杆卡金基台, 长度 7.9 毫米	TAN	025.2900
NC/RC Occlusal screw NC/RC 骀面螺丝		仅适用于固位多基基台上的基底, 长度 4.5 毫米		025.0900-04
NC/RC Stabilization pin NC/RC 稳定杆		适用于多基基台上的杆卡, 长度16毫米	不锈钢	025.0911-04
Analog holder 替代体抓持器		长度 105 毫米	铝/钢	046.239

NC

修复配件

LOCATOR®

产品	图片	描述	材料	产品编号
印模和制作工作模型				
LOCATOR® Impression Coping 取模部件		高度4毫米，4枚装	铝外壳，含 LDPE 衬垫	048.197V4
LOCATOR® Female Analog 基台替代体		长度10毫米，Ø 4 毫米，4枚装	铝	025.0120-04
钛基台				
NC LOCATOR® Abutment NC LOCATOR® 基台		高度2毫米	TAV	022.2502
		高度3毫米		022.2503
		高度4毫米		022.2504
		高度5毫米		022.2505
		高度6毫米		022.2506

CE

0086

制造商
Zest Anchors, Inc.
Escondido, CA 92029
美国

LOCATOR®

产品	图片	描述	材料	产品编号
技工室套装，可以达到20°的分散补偿				
LOCATOR® Male Processing Package 阳性部件套装		内含：基底帽(Ø 5.5毫米，高度2.5毫米)，含阳性黑色衬垫(高度1.9毫米)，封闭环(048.196)，透明阳性垫片(048.190)，粉色阳性垫片(048.191)，蓝色阳性垫片(048.192)，2枚装	纯钛/尼龙/ LDPE	048.189V2
LOCATOR® Replacement Male 阳性垫片		透明，0°-10°，高度1.7毫米，2.27kg，4枚装	尼龙	048.190V4
		粉色，轻微固位力，0°-10°，高度1.7毫米，1.36kg，4枚装		048.191V4
		蓝色，超低固位力，0°-10°，高度1.7毫米，0.68kg，4枚装		048.192V4
技工室套装，可以达到40°的分散补偿				
LOCATOR® Male Processing Package 阳性部件套装		加大角度补偿范围，内含：基底帽(Ø 5.5毫米，高度2.5毫米)，含阳性黑色衬垫 (高度1.9毫米)，封闭环(048.196)，绿色阳性垫片(048.193)，橙色阳性垫片(048.188)，红色阳性垫片(048.194)，2枚装	纯钛/尼龙/ LDPE	048.182V2
LOCATOR® Replacement Male 阳性垫片		绿色，加大角度补偿范围，10°-20°，高度1.7毫米，1.82kg，4枚装	尼龙	048.193V4
		橙色，轻微固位力，加大角度补偿范围，10°-20°，高度1.7毫米，0.91kg，4枚装		048.188V4
		红色，超低固位力，加大角度补偿范围，10°-20°，高度1.7毫米，0.45kg，4枚装		048.194V4
		灰色，加大角度补偿范围，零(0)固位力，10°-20°，高度1.7毫米，4枚装		048.181V4

CE 0086 制造商
Zest Anchors, Inc.
Escondido, CA 92029
美国

NC

修复配件

LOCATOR® 配件

产品	图片	描述	材料	产品编号
LOCATOR® Parallel Post 平行柱		高度8毫米,4枚装	LDPE	048.199V4
LOCATOR® Processing Male, black 阳性垫片, 黑色		高度1.9毫米,4枚装		048.195V4
LOCATOR® Block-out Spacer 封闭环		厚度0.4毫米,20枚装	硅树脂 橡胶	048.196V20
LOCATOR® Abutment Holder Sleeve 基台固位套		用于将基台紧扣到工具上, 4 枚装	PSU	046.413V4
辅助配件				
LOCATOR® Core Tool 核心工具		三个部件, 长度100毫米	不锈钢	046.415
LOCATOR® Screwdriver 螺丝刀		适用于棘轮扳手, 短型,长度15毫米		046.416
		适用于棘轮扳手, 长型,长度21毫米		046.417
LOCATOR® Angle Measurement Guide 角度尺		长度50毫米, 宽度15毫米		048.200

CE

0086

制造商

Zest Anchors, Inc.

Escondido, CA 92029

美国

骨水平 RC
常规 CrossFit®(十字锁合)

印模和制作工作模型

产品	图片	描述	材料	产品编号
转移杆和替代体				
RC Impression Post RC 转移杆		螺丝固位，适用于开窗式印模，短型，长度 16.5 毫米	TAN	025.4202
		螺丝固位，适用于开窗式印模，长型，长度 30 毫米		025.4205
		螺丝固位，适用于非开窗式印模，包含 2 个印模帽，长度 12.3 毫米	TAN/POM	025.4201
RC Implant Analog RC 种植体替代体		长度 12 毫米	纯钛	025.4101
数字印模				
CARES® RC Mono Scanbody CARES® RC 扫描杆		包含固定螺丝，Ø4.1 毫米，高度 10 毫米	PEEK/TAN	025.4915
骀记录转移体				
RC Bite Registration Aid RC 骀记录转移体		短型，高度 8 毫米，4 枚装	POM	025.4208-04
		长型，高度 12 毫米，4 枚装		025.4212-04

临时修复

产品	图片	描述	材料	产品编号
临时基台				
RC Temporary Abutment RC 临时基台		高度 12 毫米，Ø 7 毫米	PEEK/TAN	024.4370
替换螺丝				
RC Basal Screw RC 中央螺丝		适用于 RC 粘接基台和 VITA CAD-Temp® 临时粘接基台，仅可用于粘接，长度 7.9 毫米，	TAN	025.4908
临时基台				
RC Temporary Abutment RC 临时基台		适用于牙冠，包含中央螺丝 025.4900 ， 高度 11 毫米，Ø4.5 毫米	TAN	024.4371
		适用于牙桥，包含中央螺丝 025.4900 ， 高度 11 毫米，Ø4.5 毫米		024.4375
替换螺丝				
RC Basal Screw RC 中央螺丝		适用于 RC 解剖基台，RC Variobase® 冠用，RC 金基台，RC 杆卡金基台，长度 7.9 毫米	TAN	025.4900

粘接固位修复，解剖基台

产品	图片	描述	材料	产品编号
解剖式钛基台				
RC Anatomic Abutment RC 解剖基台		直型，包含中央螺丝 025.4900，穿龈高度 2 毫米	纯钛	022.4102
		直型，包含中央螺丝 025.4900，穿龈高度 3.5 毫米		022.4104
		成角，角度 15°，包含中央螺丝 025.4900，穿龈高度 2 毫米		022.4152
		成角，角度 15°，包含中央螺丝 025.4900，穿龈高度 3.5 毫米		022.4154
替换螺丝				
RC Basal Screw RC 中央螺丝		适用于 RC 解剖基台，RC Variobase® 冠用，RC 金基台，RC 杆卡金基台，长度 7.9 毫米	TAN	025.4900

粘接固位修复，粘接基台

产品	图片	描述	材料	产品编号
印模(基台水平), 适用于单冠修复				
RC Impression Cap RC 印模帽		适用于粘接基台, Ø5 毫米, 基台高度 4 毫米, 4 枚装	POM	025.4203-04
		适用于粘接基台, Ø5 毫米, 基台高度 5.5 毫米, 4 枚装		025.4220-04
		适用于粘接基台, Ø6.5 毫米, 基台高度 4 毫米, 4 枚装		025.4204-04
		适用于粘接基台, Ø6.5 毫米, 基台高度 5.5 毫米, 4 枚装		025.4230-04
替代体				
RC Implant Analog RC 基台替代体		适用于粘接基台, Ø5 毫米, 基台高度 4 毫米	TAN	025.4132
		适用于粘接基台, Ø5 毫米, 基台高度 5.5 毫米		025.4136
		适用于粘接基台, Ø6.5 毫米, 基台高度 4 毫米		025.4122
		适用于粘接基台, Ø6.5 毫米, 基台高度 5.5 毫米		025.4126
临时修复				
RC Temporary Coping RC 塑料临时基底		适用于粘接基台, Ø5 毫米	PMMA	024.4356
		适用于粘接基台, Ø6.5 毫米		024.4366
保护帽				
RC Protective Cap RC 粘接式保护帽		适用于粘接基台, Ø5 毫米, 基台高度 4 毫米, 4 枚装	PEEK	024.4322-04
		适用于粘接基台, Ø6.5 毫米, 基台高度 4 毫米, 4 枚装		024.4332-04
		适用于粘接基台, Ø5 毫米, 基台高度 5.5 毫米, 4 枚装		024.4326-04
		适用于粘接基台, Ø6.5 毫米, 基台高度 5.5 毫米, 4 枚装		024.4336-04
转移辅助配件				
RC Transfer Aid RC 转移辅助配件		适用于粘接基台, Ø5 毫米, 4 枚装	PC	025.4912-04
		适用于粘接基台, Ø6.5 毫米, 4 枚装		025.4913-04

粘接固位修复，粘接基台

产品	图片	描述	材料	产品编号
粘接基台(基台轮廓Ø 5 毫米)				
RC Cementable Abutment RC 粘接基台		包含中央螺丝 025.4908 , Ø5 毫米, 穿龈高度 1 毫米, 基台高度 4 毫米	纯钛	022.4321
		包含中央螺丝 025.4908 , Ø5 毫米, 穿龈高度 2 毫米, 基台高度 4 毫米		022.4322
		包含中央螺丝 025.4908 , Ø5 毫米, 穿龈高度 3 毫米, 基台高度 4 毫米		022.4323
		包含中央螺丝 025.4908 , Ø5 毫米, 穿龈高度 1 毫米, 基台高度 5.5 毫米		022.4325
		包含中央螺丝 025.4908 , Ø5 毫米, 穿龈高度 2 毫米, 基台高度 5.5 毫米		022.4326
		包含中央螺丝 025.4908 , Ø5 毫米, 穿龈高度 3 毫米, 基台高度 5.5 毫米		022.4327
粘接基台(基台轮廓Ø 6.5 毫米)				
RC Cementable Abutment RC 粘接基台		包含中央螺丝 025.4908 , Ø6.5 毫米, 穿龈高度 1 毫米, 基台高度 4 毫米	纯钛	022.4331
		包含中央螺丝 025.4908 , Ø6.5 毫米, 穿龈高度 2 毫米, 基台高度 4 毫米		022.4332
		包含中央螺丝 025.4908 , Ø6.5 毫米, 穿龈高度 3 毫米, 基台高度 4 毫米		022.4333
		包含中央螺丝 025.4908 , Ø6.5 毫米, 穿龈高度 1 毫米, 基台高度 5.5 毫米		022.4335
		包含中央螺丝 025.4908 , Ø6.5 毫米, 穿龈高度 2 毫米, 基台高度 5.5 毫米		022.4336
		包含中央螺丝 025.4908 , Ø6.5 毫米, 穿龈高度 3 毫米, 基台高度 5.5 毫米		022.4337
替换螺丝				
RC Basal Screw RC 中央螺丝		适用于 RC 粘接基台和 PEEK 临时基台, 长度 7.9 毫米	TAN	025.4908
脱蜡基底(基台轮廓Ø 5 毫米)				
RC Coping RC 基底		适用于粘接基台上的牙冠, Ø5 毫米, 高度 5.6 毫米	POM	023.4326
		适用于粘接基台上的牙桥, Ø5 毫米, 高度 5.6 毫米		023.4320
脱蜡基底(基台轮廓Ø 6.5 毫米)				
RC Coping RC 基底		适用于粘接基台上的牙冠, Ø6.5 毫米, 高度 6.5 毫米	POM	023.4336
		适用于粘接基台上的牙桥, Ø6.5 毫米, 高度 6.5 毫米		023.4330

RC

修复配件

粘接固位或螺丝固位修复，多能基台，CAD/CAM

产品	图片	描述	材料	产品编号
Variobase® 基台				
RC Variobase® for Crown RC 多能基台, 适用于单冠		包含中央螺丝 025.4900 , Ø4.5 毫米, 高度 3.5 毫米	TAN	025.4921
		包含中央螺丝 025.4900 , Ø4.5 毫米, 高度 5.5 毫米		022.0026
RC Variobase® C RC 多能基台		适用于Dentsply®Sirona®, 包含中央螺丝 025.4900 , Ø4.5 毫米, 高度 4.7 毫米		022.0044
替换螺丝				
RC Basal Screw RC 中央螺丝		适用于 RC 解剖基台, RC Variobase® 冠用多能基台, Variobase®C多能基台, RC 金基台, RC 杆卡金基台, 长度 7.9 毫米	TAN	025.4900
数字印模				
CARES® RC Mono Scanbody CARES® RC 扫描杆		包含固定螺丝, Ø4.1 毫米, 高度 10 毫米	PEEK/TAN	025.4915

粘接固位或螺丝固位修复，金基台

产品	图片	描述	材料	产品编号
金基台				
RC Gold Abutment RC 金基台		适用于牙冠，包含中央螺丝 025.4900 ， 高度 3.8 毫米，合金重量 0.562 克	Ceramicor®/ POM	022.4410
		适用于牙桥，包含中央螺丝 025.4902 ， 长度 5.2 毫米，合金重量 0.551 克		022.4430
辅助部件和替换螺丝				
RC Basal Screw RC 中央螺丝		适用于 RC 解剖基台，RC Variobase® 冠用多能基台，Variobase®C多能基台，RC 金基台，RC 杆卡金基台，长度 7.9 毫米	TAN	025.4900
		适用于桥用金基台，长度 8.3 毫米		025.4902

Ceramicor® 是一种非氧化贵金属合金，适用于铸造贵金属合金。

粘接固位或螺丝固位修复，多能基台，铸造

产品	图片	描述	材料	产品编号
Variobase® 基台				
RC Variobase® for Crown RC 多能基台，适用于单冠		包含中央螺丝 025.4900 ， Ø4.5 毫米，高度 3.5 毫米	TAN	025.4921
		包含中央螺丝 025.4900 ， Ø4.5 毫米，高度 5.5 毫米		022.0026
脱蜡基底				
RC Burn-out Coping RC 脱蜡基底		适用于 Ø4.5 毫米，牙冠高度3.5毫米，Variobase®基台，高度 12 毫米	POM	023.4759
		4枚装		023.4759-04
		适用于 Ø4.5 毫米，牙冠高度5.5 毫米，Variobase®基台，高度 12 毫米		023.0017
		4枚装		023.0017V4
替换螺丝				
RC Basal Screw RC 中央螺丝		适用于 RC 解剖基台，RC Variobase® 冠用基台，Variobase®C基台， ， RC 金基台，RC 杆卡金基台，长度 7.9 毫米	TAN	025.4900

RC

修复配件

应用螺丝固位基台进行的螺丝固位修复

产品	图片	描述	材料	产品编号
冠用转移杆（基台水平）				
Impression Posts for Open-tray Impression 开窗式转移杆		适用于螺丝固位基台，基台水平，Ø 4.6 毫米	TAN	025.2244
Impression Posts for Closed-tray Impression 非开窗式转移杆			TAN/POM	025.2246
桥用转移杆（基台水平）				
Impression Posts for Open-tray Impression 开窗式转移杆		适用于螺丝固位基台，基台水平，Ø 4.6 毫米	TAN	025.0012
Impression Posts for Closed-tray Impression 非开窗式转移杆			TAN/POM	025.0014
数字化印模				
CARES® Mono Scanbody CARES® 扫描杆		适用于螺丝固位基台，基台水平，包括固定螺丝，Ø 4.6 毫米	PEEK	025.0001
替代体				
NC/RC Analog NC/RC 替代体		适用于螺丝固位基台，直型，Ø 4.6 毫米	TAN	023.4756
		适用于螺丝固位基台，成角，角度 17°/30°，Ø 4.6 毫米		023.4757
		适用于螺丝固位基台，直型，无牙颌，Ø 4.6 毫米		025.0050
技工辅助配件				
Polishing Aid 抛光保护帽		适用于螺丝固位基台，Ø 4.6 毫米	TAN	025.0005
		4 枚装		025.0005V4
Lab Processing Screw 技工辅助螺丝		适用于螺丝固位基台，长度 20 毫米	不锈钢	025.0006
		适用于螺丝固位基台，长度 10 毫米		025.0052
保护帽				
NC/RC Protective Cap NC/RC 保护帽		适用于螺丝固位基台，包括殆面螺丝023.4763，Ø 4.6 毫米，高度 5.1 毫米，4枚装	PEEK	024.4323-04
		适用于螺丝固位基台，包括殆面螺丝023.4763，Ø 4.6 毫米，高度 6.6 毫米，4枚装		024.4324-04
		适用于螺丝固位基台，包括殆面螺丝023.4763，Ø 4.6 毫米，高度 8.1 毫米，4枚装		024.4325-04
辅助部件				
Straumann® Planning Template 轴向引导器		Straumann® Pro Arch 种植体轴向引导器		026.0016
Transfer And Alignment Pin 转移导向杆		适用于螺丝固位基台	TAN	025.0009

螺丝固位修复，螺丝固位基台

产品	图片	描述	材料	产品编号
螺丝固位基台 0°				
RC Screw-retained Abutment		直型，包括转移杆 025.0009 以及螺丝 023.4749， Ø 4.6毫米，穿龈高度 1 毫米	TAN	022.4745P
RC 螺丝固位基台		直型，包括转移杆 025.0009 以及螺丝 023.4750， Ø 4.6毫米，穿龈高度 2.5 毫米		022.4746P
		直型，包括转移杆 025.0009 以及螺丝 023.4760， Ø 4.6毫米，穿龈高度 4 毫米		022.4751P
螺丝固位基台 17°				
RC Screw-retained Abutment		成角，角度 17°，A 型，包括转移导向杆 025.0009 以及 螺丝 025.0002，Ø 4.6 毫米，穿龈高度 2.5 毫米	TAN	022.4747P
RC 螺丝固位基台		成角，角度 17°，B 型，包括转移导向杆 025.0009 以及 螺丝 025.0002，Ø 4.6 毫米，穿龈高度 2.5 毫米		022.4748P
		成角，角度 17°，A 型，包括转移导向杆 025.0009 以及 螺丝 025.0002，Ø 4.6 毫米，穿龈高度 4 毫米		022.4752P
		成角，角度 17°，B 型，包括转移导向杆 025.0009 以及 螺丝 025.0002，Ø 4.6 毫米，穿龈高度 4 毫米		022.4753P
		成角，角度 17°，A 型，包括转移导向杆 025.0009 以及 螺丝 025.0002，Ø 4.6 毫米，穿龈高度 5.5 毫米		022.0014P
		成角，角度 17°，B 型，包括转移导向杆 025.0009 以及 螺丝 025.0002，Ø 4.6 毫米，穿龈高度 5.5 毫米		022.0015P
螺丝固位基台 30°				
RC Screw-retained Abutment		成角，角度 30°，A 型，包括转移导向杆 025.0009 以及 螺丝 025.0002，Ø 4.6 毫米，穿龈高度 2.5 毫米	TAN	022.4749P
RC 螺丝固位基台		成角，角度 30°，B 型，包括转移导向杆 025.0009 以及 螺丝 025.0002，Ø 4.6 毫米，穿龈高度 2.5 毫米		022.4750P
		成角，角度 30°，A 型，包括转移导向杆 025.0009 以及 螺丝 025.0002，Ø 4.6 毫米，穿龈高度 4 毫米		022.4754P
		成角，角度 30°，B 型，包括转移导向杆 025.0009 以及 螺丝 025.0002，Ø 4.6 毫米，穿龈高度 4 毫米		022.4755P
		成角，角度 30°，A 型，包括转移导向杆 025.0009 以及 螺丝 025.0002，Ø 4.6 毫米，穿龈高度 5.5 毫米		022.0016P
		成角，角度 30°，B 型，包括转移导向杆 025.0009 以及 螺丝 025.0002，Ø 4.6 毫米，穿龈高度 5.5 毫米		022.0017P
替代螺丝				
NC/RC Screw		适用于 NC/RC 螺丝固位基台 0°，穿龈高度 1 毫米，长度 8.8 毫米	TAN	023.4749
NC/RC 中央螺丝		适用于 NC/RC 螺丝固位基台 0°，穿龈高度 2.5 毫米，长 度 10.3 毫米		023.4750
		适用于 NC/RC 螺丝固位基台 0°，穿龈高度 4 毫米，长度 11.8 毫米		023.4760
		适用于 NC/RC 螺丝固位基台 17° 和 30°，长度 7.2 毫米		025.0002

RC

修复配件

螺丝固位修复，基底

产品	图片	描述	材料	产品编号
钛基底				
NC/RC Coping NC/RC 基底		适用于螺丝固位基台，牙冠用，包括殆面螺丝 023.4763，Ø 4.6 毫米，高度 11 毫米	纯钛	023.4747
		适用于螺丝固位基台，牙桥用，包括殆面螺丝023.4763，Ø 4.6 毫米，高度 11 毫米		023.4751
金基底				
NC/RC Coping NC/RC 基底		适用于螺丝固位基台，牙冠用，包括殆面螺丝 023.4763，合金重量 0.440克，Ø 4.6 毫米	Ceramicor®/ POM	023.4753
		适用于螺丝固位基台，牙桥用，包括殆面螺丝 023.4763，合金重量 0.469克，Ø 4.6 毫米		023.4754
脱蜡基底				
NC/RC Burn-out Coping NC/RC 脱蜡基底		适用于螺丝固位基台，牙冠用，包括殆面螺丝 023.4763，Ø 4.6 毫米，高度 11 毫米	POM	023.4748
		适用于螺丝固位基台，适用于桥和杆卡，包括殆面螺丝 023.4763，Ø 4.6 毫米，高度 11 毫米		023.4758
替代螺丝				
NC/RC Occlusal Screw NC/RC 殆面螺丝		适用于螺丝固位 NC/RC 钛基台，金基台，脱蜡基底，长度 3.7 毫米	TAN	023.4763

螺丝固位杆卡修复，金

产品	图片	描述	材料	产品编号
杆卡金基底				
NC/RC Coping NC/RC 基底		适用于螺丝固位基台，适用于杆卡，包括殆面螺丝 023.4763，合金重量 0.744 克，Ø 4.6 毫米，高度 5.5 毫米	Ceramicor®	023.4755
替代螺丝				
NC/RC Occlusal Screw NC/RC 殆面螺丝		适用于螺丝固位 NC/RC 钛基台，金基台，脱蜡基底，长度 3.7 毫米	TAN	023.4763

Ceramicor® 是一种非氧化贵金属合金，适用于铸造贵金属合金。

金杆卡基台

产品	图片	描述	材料	产品编号
RC Abutment RC 杆卡基台		适用于杆卡，金，包括中央螺丝 025.2900，高度 10 毫米，合金重量 1.9 克	Ceramicor®	022.4701
替代螺丝				
RC Basal Screw RC 中央螺丝		适用于 RC 解剖基台，RC Variobase® 冠用，RC Variobase® CEREC®，RC 金基台 冠用，RC 杆卡金基台，RC 杆卡钛基台，长度 7.9 毫米	TAN	025.4900

螺丝固位杆卡修复，钛

产品	图片	描述	材料	产品编号
杆卡钛基底				
NC/RC Coping NC/RC 基底		适用于螺丝固位基台，包括骀面螺丝 023.4763，Ø 4.6 毫米，高度 5.5 毫米	纯钛	023.4752
替代螺丝				
NC/RC Occlusal Screw NC/RC 骀面螺丝		适用于螺丝固位 NC/RC 钛基台，金基台，脱蜡基底，长度 3.7 毫米	TAN	023.4763

钛杆卡基台

产品	图片	描述	材料	产品编号
RC Abutment RC 杆卡基台		适用于杆卡，包括中央螺丝 025.4900，高度 10 毫米	纯钛	022.4702
替代螺丝				
RC Basal Screw RC 中央螺丝		适用于 RC 解剖基台，RC Variobase® 冠用，RC Variobase® CEREC®，RC 金基台 冠用，RC 杆卡金基台，RC 杆卡钛基台，长度 7.9 毫米	TAN	025.4900
脱蜡基底				
NC/RC Burn-out Coping NC/RC 脱蜡基底		适用于螺丝固位基台，适用于桥和杆卡，包括骀面螺丝 023.4763，Ø 4.6 毫米，高度 11 毫米	POM	023.4758
替代螺丝				
NC/RC Occlusal Screw NC/RC 骀面螺丝		适用于螺丝固位 NC/RC 钛基台，金基台，脱蜡基底，长度 3.7 毫米	TAN	023.4763

Ceramicor® 是一种非氧化贵金属合金，适用于铸造贵金属合金。

螺丝固位修复，多基基台

产品	图片	描述	材料	产品编号
多基基台				
RC Multi-base abutment RC 直型多基基台		Ø 4.5 毫米，穿龈高度 1 毫米	TAN	022.4741
		Ø 4.5 毫米，穿龈高度 2.5 毫米		022.4742
		Ø 4.5 毫米，穿龈高度 4 毫米		022.4744
		Ø 6.5 毫米，穿龈高度 1 毫米		022.4761
		Ø 6.5 毫米，穿龈高度 2.5 毫米		022.4762
		Ø 6.5 毫米，穿龈高度 4 毫米		022.4764
RC Multi-base abutment angled 25° RC 成角 25°多基基台		A 型，Ø 4 毫米，穿龈高度 2.5 毫米		022.4782
		B 型，Ø 4 毫米，穿龈高度 2.5 毫米		022.4786
印模（基台水平）和制作工作模型				
RC Impression cap RC 印模帽		适用于多基基台，直型，Ø 4.5 毫米	POM	025.4241-04
		适用于多基基台，直型，Ø 6.5 毫米		025.4261-04
RC Impression post RC 转移杆		适用于多基基台，直型，Ø 4.5 毫米	TAN	025.4242
		适用于多基基台，直型，Ø 6.5 毫米		025.4262
		适用于多基基台，成角25°，Ø 4 毫米		025.4206
RC abutment analog RC 基台替代体		适用于多基基台，直型，Ø 4.5 毫米		025.4141
		适用于多基基台，直型，Ø 6.5 毫米		025.4161
		适用于多基基台，成角25°，Ø 4 毫米		025.4106
临时修复				
RC Temporary coping RC 临时基底		适用于多基基台，直型，Ø 4.5 毫米	纯钛	024.4342
		适用于多基基台，直型，Ø 6.5 毫米		024.4362
		适用于多基基台，成角25°，Ø 4 毫米		024.4307
RC Protective cap RC 固定式保护帽		适用于多基基台，直型，Ø 4.5 毫米	PEEK/TAN	024.4341-04
		适用于多基基台，直型，Ø 6.5 毫米		024.4361-04
		适用于多基基台，成角25°，Ø 4 毫米		024.4306-04


多基基台

产品	图片	描述	材料	产品编号
脱蜡基底(基台轮廓Ø 5 毫米)				
RC Coping RC 基底		适用于多基基台上的牙桥，直型，Ø 4.5 毫米	Ceramicor®/ POM	023.4741
		适用于多基基台上的牙桥，直型，Ø 6.5 毫米		023.4761
		适用于多基基台上的牙桥，成角25°，Ø 4.4 毫米		023.4701
		适用于多基基台上的牙桥，直型，Ø 4.5 毫米	POM	023.4744
		适用于多基基台上的牙桥，直型，Ø 6.5 毫米		023.4764
		适用于多基基台上的牙桥，成角25°，Ø 4.4 毫米		023.4704
辅助配件				
NC/RC Occlusal screw NC/RC 骀面螺丝		仅适用于固位多基基台上的基底，长度 4.5 毫米	TAN	025.0900-04
Analog holder 替代体抓持器		长度 105 毫米	铝/钢	046.239

RC

修复配件

杆卡金基台

产品	图片	描述	材料	产品编号
RC Abutment RC 杆卡基台		适用于杆卡，金，包括中央螺丝025.4900，高度10毫米，合金重量1.5克	Ceramicor®	022.4701
RC Multi-base abutment RC 直型多基基台		Ø 4.5 毫米，穿龈高度 1 毫米	TAN	022.4741
		Ø 4.5 毫米，穿龈高度 2.5 毫米		022.4742
		Ø 4.5 毫米，穿龈高度 4 毫米		022.4744
		Ø 6.5 毫米，穿龈高度 1 毫米		022.4761
		Ø 6.5 毫米，穿龈高度 2.5 毫米		022.4762
		Ø 6.5 毫米，穿龈高度 4 毫米		022.4764
RC Multi-base abutment angled 25° RC 成角 25°多基基台		A 型，Ø 4 毫米，穿龈高度 2.5 毫米		022.4782
		B 型，Ø 4 毫米，穿龈高度 2.5 毫米		022.4786
RC Coping RC 杆卡基底		Ø 4.5 毫米，适用于多基基台上的杆卡，直型	Ceramicor®	023.4742
		Ø 6.5 毫米，适用于多基基台上的杆卡，直型		023.4762
		Ø 4 毫米，适用于多基基台上的杆卡，成角25°		023.4702
辅助配件				
Activator Set 加力钥匙		适用于杆卡，三件式，长度50毫米	不锈钢	046.150
Deactivator 减力钥匙		适用于杆卡，迷你型，长度66毫米	Grilon BS/黄铜	046.151
		适用于杆卡，标准型，长度66毫米		046.152
RC Basal Screw RC 中央螺丝		适用于 RC 解剖基台，RC Variobase® 冠用，RC 金基台，RC 杆卡金基台，长度 7.9 毫米	TAN	025.4900
NC/RC Occlusal screw NC/RC 骀面螺丝		仅适用于固位多基基台上的基底，长度 4.5 毫米	TAN	025.0900-04
NC/RC Stabilization pin NC/RC 稳定杆		适用于多基基台上的杆卡，长度16毫米	不锈钢	025.0911-04

Ceramicor® 是一种非氧化贵金属合金，适用于铸造贵金属合金。

杆卡钛基台


产品	图片	描述	材料	产品编号
RC Abutment RC 杆卡基台		适用于杆卡，包括中央螺丝025.4900，高度10毫米	纯钛	022.4702
RC Multi-base abutment RC 直型多基基台		Ø 4.5 毫米，穿龈高度 1 毫米	TAN	022.4741
		Ø 4.5 毫米，穿龈高度 2.5 毫米		022.4742
		Ø 4.5 毫米，穿龈高度 4 毫米		022.4744
		Ø 6.5 毫米，穿龈高度 1 毫米		022.4761
		Ø 6.5 毫米，穿龈高度 2.5 毫米		022.4762
		Ø 6.5 毫米，穿龈高度 4 毫米		022.4764
RC Multi-base abutment angled 25° RC 成角 25°多基基台		A 型，Ø 4 毫米，穿龈高度 2.5 毫米		022.4782
		B 型，Ø 4 毫米，穿龈高度 2.5 毫米		022.4786
RC Coping RC 杆卡基底		Ø 4.5 毫米，适用于多基基台上的杆卡，直型	纯钛	023.4746
		Ø 6.5 毫米，适用于多基基台上的杆卡，直型		023.4766
		Ø 4 毫米，适用于多基基台上的杆卡，成角 25°		023.4706

辅助配件				
Activator Set 加力钥匙		适用于杆卡，三件式，长度50毫米	不锈钢	046.150
Deactivator 减力钥匙		适用于杆卡，迷你型，长度66毫米	Grilon BS/黄铜	046.151
		适用于杆卡，标准型，长度66毫米		046.152
RC Basal Screw RC 中央螺丝		适用于 RC 解剖基台，RC Variobase® 冠用，RC 金基台，RC 杆卡金基台，长度 7.9 毫米	TAN	025.4900
NC/RC Occlusal screw NC/RC 骀面螺丝		仅适用于固位多基基台上的基底，长度 4.5 毫米	TAN	025.0900-04
NC/RC Stabilization pin NC/RC 稳定杆		适用于多基基台上的杆卡，长度16毫米	不锈钢	025.0911-04

RC

修复配件

杆卡辅助配件

产品	图片	描述	材料	产品编号
RC Polishing Aid RC 打磨辅件		长度 16 毫米	不锈钢	025.4920
		4 枚装		025.4920-04
RC Basal Screw RC 中央螺丝		适用于 RC 解剖基台，RC Variobase® 冠用，RC 金基台，RC 杆卡金基台，长度 7.9 毫米	TAN	025.4900
NC/RC Occlusal screw NC/RC 殆面螺丝		仅适用于固位多基基台上的基底，长度 4.5 毫米		025.0900-04
NC/RC Stabilization pin NC/RC 稳定杆		适用于多基基台上的杆卡，长度16毫米	不锈钢	025.0911-04
Analog holder 替代体抓持器		长度 105 毫米	铝/钢	046.239

RC

修复配件

LOCATOR®

产品	图片	描述	材料	产品编号
印模和制作工作模型				
LOCATOR® Impression Coping 取模部件		高度4毫米，4枚装	铝外壳，含 LDPE 衬垫	048.197V4
LOCATOR® Female Analog 基台替代体		长度 10 毫米，Ø4 毫米，4 枚装	铝	025.0120-04
钛基台				
RC LOCATOR® Abutment RC LOCATOR® 基台		高度1毫米	TAV	022.4501
		高度2毫米		022.4502
		高度3毫米		022.4503
		高度4毫米		022.4504
		高度5毫米		022.4505
		高度6毫米		022.4506

CE0086

制造商
Zest Anchors, Inc.
Escondido, CA 92029
美国

LOCATOR®

产品	图片	描述	材料	产品编号
技工室套装，可以达到20°的分散补偿				
LOCATOR® Male Processing Package 阳性部件套装		内含：基底帽(Ø 5.5毫米，高度2.5毫米)，含阳性黑色衬垫(高度1.9毫米)，封闭环(048.196)，透明阳性垫片(048.190)，粉色阳性垫片(048.191)，蓝色阳性垫片(048.192)，2枚装	纯钛/尼龙/LDPE	048.189V2
LOCATOR® Replacement Male 阳性垫片		透明，0°-10°，高度1.7毫米，2.27kg，4枚装	尼龙	048.190V4
		粉色，轻微固位力，0°-10°，高度1.7毫米，1.36kg，4枚装		048.191V4
		蓝色，超低固位力，0°-10°，高度1.7毫米，0.68kg，4枚装		048.192V4
技工室套装，可以达到40°的分散补偿				
LOCATOR® Male Processing Package 阳性部件套装		加大角度补偿范围，内含：基底帽(Ø 5.5毫米，高度2.5毫米)，含阳性黑色衬垫(高度1.9毫米)，封闭环(048.196)，绿色阳性垫片(048.193)，橙色阳性垫片(048.188)，红色阳性垫片(048.194)，2枚装	纯钛/尼龙/LDPE	048.182V2
LOCATOR® Replacement Male 阳性垫片		绿色，加大角度补偿范围，10°-20°，高度1.7毫米，1.82kg，4枚装	尼龙	048.193V4
		橙色，轻微固位力，加大角度补偿范围，10°-20°，高度1.7毫米，0.91kg，4枚装		048.188V4
		红色，超低固位力，加大角度补偿范围，10°-20°，高度1.7毫米，0.45kg，4枚装		048.194V4
		灰色，加大角度补偿范围，零(0)固位力，10°-20°，高度1.7毫米，4枚装		048.181V4

CE

0086

制造商

Zest Anchors, Inc.

Escondido, CA 92029

美国

LOCATOR® 配件

产品	图片	描述	材料	产品编号
LOCATOR® Parallel Post 平行柱		高度8毫米,4枚装	LDPE	048.199V4
LOCATOR® Processing Male, black 阳性垫片, 黑色		高度1.9毫米,4枚装		048.195V4
LOCATOR® Block-out Spacer 封闭环		厚度0.4毫米,20枚装	硅树脂 橡胶	048.196V20
LOCATOR® Abutment Holder Sleeve 基台固位套		用于将基台紧扣到工具上, 4 枚装	PSU	046.413V4
辅助配件				
LOCATOR® Core Tool 核心工具		三个部件, 长度100毫米	不锈钢	046.415
LOCATOR® Screwdriver 螺丝刀		适用于棘轮扳手, 短型,长度15毫米		046.416
		适用于棘轮扳手, 长型,长度21毫米		046.417
LOCATOR® Angle Measurement Guide 角度尺		长度50毫米, 宽度15毫米		048.200

CE 0086

制造商
Zest Anchors, Inc.
Escondido, CA 92029
美国

RC

修复配件

修复设计辅件

产品	图片	描述	材料	产品编号
模拟套装				
Bone Level PLAN Set 骨水平设计套装		长度 278 毫米，宽度 124 毫米，高度 33 毫米	POM	026.0902
Bone Level PLAN Cassette 骨水平设计套装外盒		空盒，长度 278 毫米，宽度 124 毫米，高度 33 毫米	PC	026.0000
NC 模拟基台				
NC PLAN Anatomic Abutment NC 模拟解剖基台		直型，穿龈高度 2 毫米，4 枚装	POM	025.2601-04
		直型，穿龈高度 3.5 毫米，4 枚装		025.2602-04
		成角，角度 15°，穿龈高度 2 毫米，4 枚装		025.2603-04
		成角，角度 15°，穿龈高度 3.5 毫米，4 枚装		025.2604-04
NC PLAN Gold Abutment NC 模拟金基台		适用于牙冠，高度 3.2 毫米，4 枚装		025.2610-04
NC PLAN Cementable Abutment NC 模拟粘接基台		Ø 3.5 毫米，穿龈高度 1 毫米，基台高度 4 毫米，4 枚装		025.2611-04
		Ø 3.5 毫米，穿龈高度 2 毫米，基台高度 4 毫米，4 枚装		025.2612-04
		Ø 3.5 毫米，穿龈高度 3 毫米，基台高度 4 毫米，4 枚装		025.2613-04
		Ø 3.5 毫米，穿龈高度 1 毫米，基台高度 5.5 毫米，4 枚装		025.2615-04
		Ø 3.5 毫米，穿龈高度 2 毫米，基台高度 5.5 毫米，4 枚装		025.2616-04
		Ø 3.5 毫米，穿龈高度 3 毫米，基台高度 5.5 毫米，4 枚装		025.2617-04
		Ø 5 毫米，穿龈高度 1 毫米，基台高度 4 毫米，4 枚装		025.2621-04
		Ø 5 毫米，穿龈高度 2 毫米，基台高度 4 毫米，4 枚装		025.2622-04
		Ø 5 毫米，穿龈高度 3 毫米，基台高度 4 毫米，4 枚装		025.2623-04
		Ø 5 毫米，穿龈高度 1 毫米，基台高度 5.5 毫米，4 枚装		025.2625-04
		Ø 5 毫米，穿龈高度 2 毫米，基台高度 5.5 毫米，4 枚装		025.2626-04
		Ø 5 毫米，穿龈高度 3 毫米，基台高度 5.5 毫米，4 枚装		025.2627-04

修复设计辅件

产品	图片	描述	材料	产品编号
NC PLAN LOCATOR® Abutment		高度 2 毫米, 4枚装	POM	025.2642-04
NC LOCATOR® 模拟基台		高度 3 毫米, 4枚装		025.2643-04
		高度 4 毫米, 4枚装		025.2644-04
		高度 5 毫米, 4枚装		025.2645-04
		高度 6 毫米, 4枚装		025.2646-04
				025.2648-04
NC PLAN Screw-retained Abutment NC 模拟螺丝固位基台		直型, Ø 3.5 毫米, 穿龈高度 2.5 毫米, 4枚装		025.2650-04
		直型, Ø 4.6 毫米, 穿龈高度 2.5 毫米, 4枚装		025.2655-04
		成角, 角度 17°, A 型, Ø 4.6 毫米, 穿龈高度 2.5 毫米, 4枚装		025.2658-04
		成角, 角度 17°, B 型, Ø 4.6 毫米, 穿龈高度 2.5 毫米, 4枚装		025.2653-04
		成角, 角度 30°, A 型, Ø 4.6 毫米, 穿龈高度 2.5 毫米, 4枚装		025.2660-04
		成角, 角度 30°, B 型, Ø 4.6 毫米, 穿龈高度 2.5 毫米, 4枚装		

RC 模拟基台				
RC PLAN Anatomic Abutment RC 模拟解剖基台		直型, 穿龈高度 2 毫米, 4枚装	POM	025.4601-04
		直型, 穿龈高度 3.5 毫米, 4枚装		025.4602-04
		成角, 角度 15°, 穿龈高度 2 毫米, 4枚装		025.4603-04
		成角, 角度 15°, 穿龈高度 3.5 毫米, 4枚装		025.4604-04
RC PLAN Gold Abutment RC 模拟金基台		适用于牙冠, 高度 3.3 毫米, 4枚装		025.4610-04
RC PLAN Cementable Abutment RC 模拟粘接基台		Ø 5 毫米, 穿龈高度 1 毫米, 基台高度 4 毫米, 4枚装		025.4621-04
		Ø 5 毫米, 穿龈高度 2 毫米, 基台高度 4 毫米, 4枚装		025.4622-04
		Ø 5 毫米, 穿龈高度 3 毫米, 基台高度 4 毫米, 4枚装		025.4623-04
		Ø 5 毫米, 穿龈高度 1 毫米, 基台高度 5.5 毫米, 4枚装		025.4625-04
		Ø 5 毫米, 穿龈高度 2 毫米, 基台高度 5.5 毫米, 4枚装		025.4626-04
		Ø 5 毫米, 穿龈高度 3 毫米, 基台高度 5.5 毫米, 4枚装		025.4627-04
RC PLAN Gold Abutment RC 模拟金基台		适用于牙桥, 高度 5.2 毫米, 4枚装		025.4630-04

RC

修复配件

修复设计辅件

产品	图片	描述	材料	产品编号
RC PLAN Cementable Abutment RC 模拟粘接基台		Ø 6.5 毫米, 穿龈高度 1 毫米, 基台高度 4 毫米, 4枚装	POM	025.4631-04
		Ø 6.5 毫米, 穿龈高度 2 毫米, 基台高度 4 毫米, 4枚装		025.4632-04
		Ø 6.5 毫米, 穿龈高度 3 毫米, 基台高度 4 毫米, 4枚装		025.4633-04
		Ø 6.5 毫米, 穿龈高度 1 毫米, 基台高度 5.5 毫米, 4枚装		025.4635-04
		Ø 6.5 毫米, 穿龈高度 2 毫米, 基台高度 5.5 毫米, 4枚装		025.4636-04
		Ø 6.5 毫米, 穿龈高度 3 毫米, 基台高度 5.5 毫米, 4枚装		025.4637-04
		高度 1 毫米, 4枚装		025.4641-04
		高度 2 毫米, 4枚装		025.4642-04
		高度 3 毫米, 4枚装		025.4643-04
		高度 4 毫米, 4枚装		025.4644-04
		高度 5 毫米, 4枚装		025.4645-04
		高度 1-6 毫米, 4枚装		025.4646-04
RC PLAN Screw-retained Abutment RC 模拟螺丝固位基台		直型, Ø 4.6 毫米, 穿龈高度 2.5 毫米, 4枚装		025.4648-04
		成角, 角度 17°, A 型, Ø 4.6 毫米, 穿龈高度 2.5 毫米, 4枚装		025.4649-04
		成角, 角度 17°, B 型, Ø 4.6 毫米, 穿龈高度 2.5 毫米, 4枚装		025.4650-04
		成角, 角度 30°, A 型, Ø 4.6 毫米, 穿龈高度 2.5 毫米, 4枚装		025.4653-04
		成角, 角度 30°, B 型, Ø 4.6 毫米, 穿龈高度 2.5 毫米, 4枚装		025.4660-04

修复设计辅件

产品	图片	描述	材料	产品编号
模拟基台				
NC PLAN Multi-base abutment NC 模拟多基基台		直型, Ø 3.5 毫米, 穿龈高度 1 毫米, 4枚装	POM	025.2651-04
		直型, Ø 3.5 毫米, 穿龈高度 2.5 毫米, 4枚装		025.2652-04
		直型, Ø 3.5 毫米, 穿龈高度 4 毫米, 4枚装		025.2654-04
		直型, Ø 4.5 毫米, 穿龈高度 1 毫米, 4枚装		025.2656-04
		直型, Ø 4.5 毫米, 穿龈高度 2.5 毫米, 4枚装		025.2657-04
		直型, Ø 4.5 毫米, 穿龈高度 4 毫米, 4枚装		025.2659-04
		成角, 角度 25°, A型, Ø 4 毫米, 穿龈高度 2.5 毫米, 4枚装		025.2662-04
		成角, 角度 25°, B型, Ø 4 毫米, 穿龈高度 2.5 毫米, 4枚装		025.2665-04
RC PLAN Multi-base abutment RC 模拟多基基台		直型, Ø 4.5 毫米, 穿龈高度 1 毫米, 4枚装		025.4651-04
		直型, Ø 4.5 毫米, 穿龈高度 2.5 毫米, 4枚装		025.4652-04
		直型, Ø 4.5 毫米, 穿龈高度 4 毫米, 4枚装		025.4654-04
		直型, Ø 6.5 毫米, 穿龈高度 1 毫米, 4枚装		025.4656-04
		直型, Ø 6.5 毫米, 穿龈高度 2.5 毫米, 4枚装		025.4657-04
		直型, Ø 6.5 毫米, 穿龈高度 4 毫米, 4枚装		025.4659-04
		成角, 角度 25°, A型, Ø 4 毫米, 穿龈高度 2.5 毫米, 4枚装		025.4662-04
		成角, 角度 25°, B型, Ø 4 毫米, 穿龈高度 2.5 毫米, 4枚装		025.4665-04

RC

修复配件

Straumann® 经典

士卓曼® Classic

质保和客户服务增加了产品的价值，增强了信任感。

此外，我们的士卓曼Classic计划确保您可以在旧的种植体上修复上部结构。

您可以重新订购1974年以后生产的所有士卓曼种植体的修复配件。

士卓曼® Classic 部件的产品编号末尾以“C”为标记。

不是所有国家都可以购买到士卓曼® Classic。

请联系士卓曼的当地官方经销商，咨询此项服务的更多信息。

软组织水平 (NN) 窄颈种植体，修复配件

产品	图片	描述	材料	产品编号
印模和制作工作模型				
NN Impression Cap NN 转移杆		螺丝固位，适用于开窗式印模，长度 17.5 毫米	铝/TAN	048.016
NN abutment Analog NN 基台替代体		长度 14 毫米	不锈钢	048.130
临时基台				
NN Temporary Coping NN 临时基台		包含 NN 殆面螺丝，高度 10 毫米	PEEK/纯钛	048.669
金基台				
NN Gold Coping NN 金基底		铸造，不含螺丝，高度 10 毫米，合金重量 0.286 克	Ceramicor®/ POM	048.635
钛基台				
NN Post NN 基台		可磨改，不含螺丝，高度 9 毫米	纯钛	048.505
		弯角，成角 15°，不含螺丝，高度 8.8 毫米		048.550
		弯角，成角 20°，不含螺丝，高度 8.8 毫米		048.551
保护帽				
NN Protective Cap 保护帽		包含中央固定螺丝，高度 3.4 毫米	PEEK/TAN	048.050
螺丝				
NN Occlusal Screw NN 殆面螺丝		黄色，长度 5 毫米，螺纹 Ø1.8 毫米	TAN	049.177



Straumann® CARES® 包含一系列用于高质量、定制化 CAD/CAM 修复体的硬件、软件和服务。可连通手术设计用 CAD/CAM 软件的数据，为定制美观修复体带来新机遇。

数字化工作流程：精简工作流程，提高结果可预测性

数字化正深刻影响着牙科学和牙科技术，其带来的一系列新应用和新工艺正在从根本上改变治疗的设计方式。新型软件和硬件产品可提高牙科技工室的创新能力和工作效率，同时可推动牙科技工室产品的分散式生产。这些趋势需要我们调整商业模式，以满足不断增长的患者期望。

Straumann® CARES® 系统

新型 Straumann® CARES® 扫描仪产品系列完美继承了先前产品的成熟技术，并进行了人体工程学和功能方面的重大改善。我们的扫描仪用途广、速度快、精度高，专门为满足高通量环境下的稳定性和生产需求而设计。整合的处理单元能够确保所有相关程序间的通信无缝、流畅。

Straumann® CARES® Visual Design 软件是一个完整的 CAD-CAM 软件系统，可提供多种服务和解决方案。



Straumann® CARES® 数字化修复

您可以使用 Straumann® CARES® CAD/CAM 系统进行扫描和设计，然后直接从我们的服务订购定制修复体。

多样化的产品组合

我们多样化的产品组合可为各种案例提供合适的解决方案—包括 CARES® 螺丝固位桥、CARES® 活动式杆卡及 CARES® 简单、复杂固定式杆卡等多种个性化种植修复，可选加工材料有钛金属、钴铬金属。

材料领先同行

Straumann® CARES® CAD/CAM 数字化修复服务的材料均从欧洲进口，保证原材料的高品质，从根源上保证工件内部的物理性质的一致性。

精确度高

引进了士卓曼集团在欧洲以及世界其他地方实践多年的规模集中化生产工艺流程。无论是从机器配置上还是从机加工技术（HSC）上都达到了工业级水准。所使用的超微型机床重达数吨，在加工过程中几乎没有震动，具备超高的定位精度和重复定位精度。最大程度上保证上部结构与种植体的良好匹配，避免机械并发症的出现。



数字化解决方案



通过全面整合化繁为简

Straumann 进一步整合了引导式植牙和 CADCAM 系统以及所需的设备、材料和服务。

定制美观修复体只需一步

设计软件 coDiagnostiX® 和 CARES® Visual 的结合为种植体修复治疗开启了一个全新的世界。

Straumann® CARES®

- 创新、灵活、可持续的解决方案
- 可利用整套数字化工作流程
- 提高工作效率、精简流程步骤、缩短处理时间、降低相关成本



在这一章节中您可以查询到用于种植体和修复治疗的全面的器械和配件清单。

不仅仅是配件，更是完美无瑕治疗的基础

在以下各页中，您可以找到钻头、修复辅助部件、螺丝钉、存储系统等所有完善士卓曼手术和修复产品组合的重要器械和演示材料。

器械和配件

外科手术工具	126
手术设计辅助工具	126
种植床预备器械	127
器材和工具	131
引导手术用	133
取出工具	137
骨块固定 1.5 毫米	140
骨挤压器	141
储存系统	143
修复配件	144
修复工具	144
储存系统和工具箱	147
修复设计辅助工具	149
修复用螺丝	150

手术设计辅助工具

产品	图片	描述	材料	产品编号
X 射线参考球				
X-ray Reference Sphere X 射线参考球		∅ 5 毫米， 4 枚装， 一次性使用	不锈钢	049.076V4
X 射线模板				
X-ray Templates X 射线模板		适用于标准、美学和 Straumann® 纯陶瓷种植体	塑料	150.215
		适用于骨水平种植体		150.216
辅助工具和间距尺				
Diagnostic T T 型诊断尺		诊断辅助工具，长度 65 毫米，宽度 33 毫米	纯钛	046.147
Implant Distance Indicator for Standard, Standard Plus and Tapered Effect 标准、美学和锥柱状种植体 间距尺		长度 100 毫米	纯钛/不锈钢	046.148
Bone Level Implant Distance Indicator 骨水平种植体间距尺				026.0901

球钻

产品	图片	描述	材料	产品编号
球钻				
Round Bur 球钻		Ø 1.4 毫米	不锈钢	044.022
		Ø 2.3 毫米		044.003
		Ø 3.1 毫米		044.004

针形定导钻

产品	图片	描述	材料	产品编号
针形定导钻				
Needle Drill 针形定导钻		短型, Ø 1.6 毫米, 长度 33 毫米	不锈钢	026.0054
		长型, Ø 1.6 毫米, 长度 41 毫米		026.0056

先锋钻

产品	图片	描述	材料	产品编号
先锋钻 — 适用于 S SP TE BC NNC				
Pilot Drill 1 先锋钻 1		短型, 长度 33 毫米, Ø 2.2 毫米	不锈钢	044.210
		长型, 长度 41 毫米, Ø 2.2 毫米		044.211
Pilot Drill 2 先锋钻 2		短型, 长度 33 毫米, Ø 2.8 毫米	不锈钢	044.214
		长型, 长度 41 毫米, Ø 2.8 毫米		044.215
BLT Pilot Drill BLT 先锋钻		短型, 长度 33 毫米, Ø 2.2 毫米	不锈钢	026.0001
		长型, 长度 41 毫米, Ø 2.2 毫米		026.0002

扩孔钻

产品	图片	描述	材料	产品编号
扩孔钻 — 适用于 S SP TE BC NNC				
Twist Drill 扩孔钻		短型, 长度 33 毫米, Ø 3.5 毫米	不锈钢	044.250
		长型, 长度 41 毫米, Ø 3.5 毫米		044.251
		短型, 长度 33 毫米, Ø 4.2 毫米		044.254
		长型, 长度 41 毫米, Ø 4.2 毫米		044.255

BLT 扩孔钻

产品	图片	描述	材料	产品编号
骨水平钻				
BLT Drill BLT 扩孔钻		短型, 长度 33 毫米, Ø 2.8 毫米	不锈钢	026.2200
		长型, 长度 41 毫米, Ø 2.8 毫米		026.2201
		短型, 长度 33 毫米, Ø 3.5 毫米		026.4200
		长型, 长度 41 毫米, Ø 3.5 毫米		026.4201
BLT Drill BLT 扩孔钻		短型, 长度 33 毫米, Ø 4.2 毫米	不锈钢	026.6200
		长型, 长度 41 毫米, Ø 4.2 毫米		026.6201

成型钻


产品	图片	描述	材料	产品编号
成型钻，适用于美学种植体				
SP Profile Drill SP 成型钻		短型，长度 25 毫米，Ø3.3 毫米	不锈钢	044.086
		长型，长度 34 毫米，Ø3.3 毫米		044.087
		短型，长度 25 毫米，Ø4.1 毫米		044.088
		长型，长度 34 毫米，Ø4.1 毫米		044.089
		短型，长度 25 毫米，Ø4.8 毫米		044.084
		长型，长度 34 毫米，Ø4.8 毫米		044.085
成型钻，适用于椎柱状种植体				
TE Profile Drill TE 成型钻		短型，长度 28 毫米，Ø3.3 毫米	不锈钢	044.701
		长型，长度 37 毫米，Ø3.3 毫米		044.708
		短型，长度 28 毫米，Ø4.1 毫米		044.705
		长型，长度 37 毫米，Ø4.1 毫米		044.712
		短型，长度 28 毫米，Ø4.8 毫米		044.703
		长型，长度 37 毫米，Ø4.8 毫米		044.710
成型钻，适用于骨水平，NNC 和 Straumann® 纯陶瓷种植体				
BL/NNC Profile Drill BL/NNC 成型钻		短型，长度 26 毫米，Ø3.3 毫米	不锈钢	026.2303
		长型，长度 35 毫米，Ø3.3 毫米		026.2306
BL Profile Drill BL 成型钻		短型，长度 26 毫米，Ø4.1 毫米		026.4303
		长型，长度 35 毫米，Ø4.1 毫米		026.4306
		短型，长度 26 毫米，Ø4.8 毫米		026.6303
		长型，长度 35 毫米，Ø4.8 毫米		026.6306
成型钻，适用于骨水平椎柱状种植体				
BLT Profile Drill BLT 成型钻		短型，长度 25 毫米，Ø 3.3 毫米	不锈钢	026.0003
		长型，长度 33 毫米，Ø 3.3 毫米		026.0004
		短型，长度 25 毫米，Ø 4.1 毫米		026.0005
		长型，长度 33 毫米，Ø 4.1 毫米		026.0006
		短型，长度 25 毫米，Ø 4.8 毫米		026.0007
		长型，长度 33 毫米，Ø 4.8 毫米		026.0008

器械和配件

攻丝钻

产品	图片	描述	材料	产品编号
攻丝钻，适用窄颈窄 CrossFit®（十字锁合）种植体				
BL/TE/NNC Tap BL/TE/NNC 攻丝钻		适用于适配器，长度 23 毫米，Ø3.3 毫米	不锈钢/TAN	026.2310
攻丝钻，适用标准/美学种植体				
S/SP Tap S/SP 攻丝钻		适用于适配器，长度 23 毫米，Ø3.3 毫米	不锈钢/TAN	044.575
		适用于适配器，长度 23 毫米，Ø4.1 毫米		044.577
		适用于适配器，长度 23 毫米，Ø4.8 毫米		044.579
攻丝钻，适用骨水平，锥柱状种植体				
BL/TE/NNC Tap BL/TE/NNC 攻丝钻		适用于适配器，长度 23 毫米，Ø3.3 毫米	不锈钢/TAN	026.2310
BL/TE Tap BL/TE 攻丝钻		适用于适配器，长度 23 毫米，Ø4.1 毫米		026.4310
		适用于适配器，长度 23 毫米，Ø4.8 毫米		026.6310
攻丝钻，适用骨水平锥柱状种植体				
BLT Tap BLT 攻丝钻		长度 25 毫米，Ø3.3 毫米	不锈钢/TAN	026.0009
		长度 25 毫米，Ø4.1 毫米		026.0010
		长度 25 毫米，Ø4.8 毫米		026.0011

BLT 测量杆和测深杆

产品	图片	描述	材料	产品编号
测量杆和测深杆，适用骨水平锥柱状种植体				
BLT Alignment Pin BLT 测量杆		长度 27 毫米，Ø2.2 毫米	纯钛	046.703
BLT Depth Gauge with Distance Indicator BLT 测量杆带间距尺		长度 27 毫米，Ø2.2/2.8 毫米		046.704
BLT Depth Gauge BLT 测深杆		长度 27 毫米，Ø2.8 毫米		046.705
		长度 27 毫米，Ø3.5 毫米		046.706
		长度 27 毫米，Ø4.2 毫米		046.707

测量杆和测深杆 — 适用于 S SP TE BC NNC


产品	图片	描述	材料	产品编号
测量杆和测深杆				
Alignment pin straight 测量杆		长度 28 毫米, Ø 2.2 毫米	纯钛	046.458
Depth gauge with distance indicator 测量杆带间距尺		长度 27 毫米, Ø 2.2/2.8 毫米		046.455
Depth gauge 测深杆		长度 26 毫米, Ø 3.5 毫米		046.450
		长度 26 毫米, Ø 4.2 毫米		046.451

取骨钻

产品	图片	描述	材料	产品编号
取骨钻，适用于软组织水平				
RN Bone Profiler RN 取骨钻		配合 049.082 使用，内径 Ø3.6/4.8 毫米， 外径 Ø5.4 毫米，长度 28 毫米	不锈钢	044.304
WN Bone Profiler WN 取骨钻		配合 049.082 使用，内径 Ø3.6/6.5 毫米， 外径 Ø7.5 毫米，长度 25 毫米		044.305
SCS Guiding Cylinder SCS 取出钻导向套		适用于 Ø3.3 毫米种植体，长度 10.5 毫米， Ø3.6 毫米		049.082
取骨钻，适用于骨水平				
BL Bone Profiler 1 BL 取骨钻 1		长度 23 毫米, Ø 5.2 毫米	不锈钢	026.0022*
BL Bone Profiler 2 BL 取骨钻 2		长度 23 毫米, Ø 6.6 毫米		026.0023*
BL Bone Profiler 3 BL 取骨钻 3		长度 23 毫米, Ø 6 毫米		026.0024*
NC Guiding Cylinder NC 取出钻导向套		适用于 BL 取骨钻，长度 12.5 毫米, Ø 3.25 毫米	TAN	026.0025S*
RC Guiding Cylinder RC 取出钻导向套		适用于 BL 取骨钻，长度 12.5 毫米, Ø 3.7 毫米		026.0026S*

带 * 产品，暂无注册证


手术盒

产品	图片	描述	材料	产品编号
Straumann® Surgical Cassette Straumann® 外科工具盒		长度 290 毫米，宽度 175 毫米，高度 62 毫米	塑料	040.165
Straumann® Basic Surgical Cassette Straumann® 外科工具盒 基础版				040.564

辅助器械

产品	图片	描述	材料	产品编号
镊子				
Instrument Tweezers 器械镊子		长度 140 毫米	不锈钢	046.110
Tweezers 镊子		尖端采用等离子镀钛，长度 105 毫米	不锈钢/纯钛	046.095
		尖端未采用等离子镀钛，长度 105 毫米	不锈钢	046.057
器材				
Ratchet 棘轮扳手		包含拆装工具，长度 84 毫米	不锈钢	046.119
Torque Control Device 扭矩控制器		适用于棘轮扳手，长度 82 毫米		046.049
Holding Key 固定扳手		长度 85 毫米		046.064
Service Instrument 拆装工具		适用于棘轮扳手，长度 50 毫米		046.108
Cleaning Brush 清洁刷		适用于棘轮扳手，长度 100 毫米，Ø4.5 毫米，4 枚装	不锈钢/尼龙	045.111V4
Release Aid N for Loxim™ 窄颈 Loxim™ 携带体解除器		用于 NC 和 NNC平台解除携带体使用，长度 90 毫米	不锈钢	026.2558
Release Aid R/W for Loxim™ 常规颈/宽颈 Loxim™ 携带体解除器		用于 RC, RN 和 WN平台解除携带体使用，长度 90 毫米		026.4558

辅助器械

产品	图片	描述	材料	产品编号
适配器				
Adapter 适配器		适用于棘轮扳手，超短型，长度 11 毫米	不锈钢	046.460
		适用于棘轮扳手，短型，长度 18 毫米		046.461
		适用于棘轮扳手，长型，长度 28 毫米		046.462
		适用于机用件，超短型，长度 19 毫米		046.470
		适用于机用件，短型，长度 26 毫米		046.471
		适用于机用件，长型，长度 34 毫米		046.472
螺丝刀				
SCS Screwdriver SCS 螺丝刀		适用于棘轮扳手，超短型，长度 15 毫米	不锈钢	046.400
		适用于棘轮扳手，短型，长度 21 毫米		046.401
		适用于棘轮扳手，长型，长度 27 毫米		046.402
		适用于机用件，超短型，长度 20 毫米		046.410
		适用于机用件，短型，长度 26 毫米		046.411
		适用于机用件，长型，长度 32 毫米		046.412
Hexagonal Screwdriver 侧向螺丝用螺丝刀		长度 30 毫米		046.421
杂项辅助器械				
Light Curette for Implants 种植体清洁器		5 把/包	塑料	045.921

CE 0086 制造商
Kerr Hawe SA
Zona Strecce
CH-6934 Bioggio
瑞士

引导手术用模板

产品	图片	描述	材料	产品编号
Straumann® 引导手术用套筒				
T-Sleeve T 型套筒		引导手术用, 高度 5 毫米, Ø5 毫米, 4 枚装	不锈钢	034.053V4
		引导手术用, 高度 6 毫米, Ø2.8 毫米, 4 枚装		034.055V4

种植床预备器械

产品	图片	描述	材料	产品编号
引导手术用粘膜冲头				
Mucosa Punch 粘膜冲头		引导手术用, 长度 30 毫米, Ø3.4 毫米	不锈钢	034.010
		引导手术用, 长度 30 毫米, Ø4 毫米		034.011
		引导手术用, 长度 30 毫米, Ø4.7 毫米		034.012
引导手术用切削刀				
Milling Cutter 切削刀		引导手术用, 长度 32.5 毫米, Ø2.8 毫米	不锈钢	034.215
		引导手术用, 长度 32.5 毫米, Ø3.5 毫米		034.415
		引导手术用, 长度 32.5 毫米, Ø4.2 毫米		034.615
引导手术用先锋钻				
Pilot Drill 先锋钻		短型, 引导手术用, 长度 32 毫米, Ø2.2 毫米	不锈钢	034.123
		长型, 引导手术用, 长度 36 毫米, Ø2.2 毫米		034.126
		超长型, 引导手术用, 长度 40 毫米, Ø2.2 毫米		034.129
BLT Pilot Drill BLT 先锋钻		短型, 引导手术用, 长度 33.4 毫米, Ø2.2 毫米	不锈钢	034.257
		中型, 引导手术用, 长度 37.4 毫米, Ø2.2 毫米		034.258
		长型, 引导手术用, 长度 41.4 毫米, Ø2.2 毫米		034.259
引导手术用扩孔钻				
Twist Drill PRO 扩孔钻		短型, 引导手术用, 长度 32 毫米, Ø2.8 毫米	不锈钢	034.223
		长型, 引导手术用, 长度 36 毫米, Ø2.8 毫米		034.226
		超长型, 引导手术用, 长度 40 毫米, Ø2.8 毫米		034.229
		短型, 引导手术用, 长度 32 毫米, Ø3.5 毫米		034.423
		长型, 引导手术用, 长度 36 毫米, Ø3.5 毫米		034.426
		超长型, 引导手术用, 长度 40 毫米, Ø3.5 毫米		034.429
		短型, 引导手术用, 长度 32 毫米, Ø4.2 毫米		034.623
		长型, 引导手术用, 长度 36 毫米, Ø4.2 毫米		034.626
		超长型, 引导手术用, 长度 40 毫米, Ø4.2 毫米		034.629

种植床预备器械

产品	图片	描述	材料	产品编号
引导手术用 BLT 钻				
BLT Drill BLT 钻		短型，引导手术用，长度 33.4 毫米，Ø2.8 毫米	不锈钢	034.260
		中型，引导手术用，长度 37.4 毫米，Ø2.8 毫米		034.261
		长型，引导手术用，长度 41.4 毫米，Ø2.8 毫米		034.262
		短型，引导手术用，长度 33.4 毫米，Ø3.5 毫米		034.263
		中型，引导手术用，长度 37.4 毫米，Ø3.5 毫米		034.264
		长型，引导手术用，长度 41.4 毫米，Ø3.5 毫米		034.265
		短型，引导手术用，长度 33.4 毫米，Ø4.2 毫米		034.266
		中型，引导手术用，长度 37.4 毫米，Ø4.2 毫米		034.267
		长型，引导手术用，长度 41.4 毫米，Ø4.2 毫米		034.268
引导手术用成型钻				
SP Profile Drill SP 成型钻		引导手术用，长度 35 毫米，Ø3.3 毫米	不锈钢	034.235
TE Profile Drill TE 成型钻		引导手术用，长度 37 毫米，Ø3.3 毫米		034.237
BL Profile Drill BL 成型钻				026.2503
SP Profile Drill SP 成型钻		引导手术用，长度 35 毫米，Ø4.1 毫米		034.435
TE Profile Drill TE 成型钻		引导手术用，长度 37 毫米，Ø4.1 毫米		034.437
BL Profile Drill BL 成型钻				026.4503
		引导手术用，长度 37 毫米，Ø4.8 毫米		026.6503
BLT Profile Drill BLT 成型钻		引导手术用，长度 36 毫米，Ø3.3 毫米		034.269
		引导手术用，长度 36 毫米，Ø4.1 毫米		034.270
		引导手术用，长度 36 毫米，Ø4.8 毫米	034.271	
引导手术用攻丝钻				
S/SP Tap S/SP 攻丝钻		适用于机用件，引导手术用，长度 42 毫米，Ø3.3 毫米	不锈钢	034.245
BL/TE Tap BL/TE 攻丝钻				026.2510
S/SP Tap S/SP 攻丝钻		适用于机用件，引导手术用，长度 42 毫米，Ø4.1 毫米		034.445
BL/TE Tap BL/TE 攻丝钻				026.4510
S/SP Tap S/SP 攻丝钻		适用于机用件，引导手术用，长度 42 毫米，Ø4.8 毫米		034.645
BL/TE Tap BL/TE 攻丝钻				026.6510
BLT Tap BLT 攻丝钻		引导手术用，长度 42 毫米，Ø3.3 毫米		034.272
		引导手术用，长度 42 毫米，Ø4.1 毫米		034.273
		引导手术用，长度 42 毫米，Ø4.8 毫米		034.274

引导型适配器

产品	图片	描述	材料	产品编号
引导型适配器 Guided Adapter 引导型适配器		适配器,美学/锥柱状/NNC种植体用, 引导手术用	不锈钢	046.708
		适配器,标准种植体用, 引导手术用		046.709
		适配器,骨水平/骨水平锥柱状种植体用, 引导手术用		026.0083
		棘轮扳手适配器,描述为SP/NNC/TE, 长16mm		046.710
		棘轮扳手适配器,描述为S, 长17mm		046.711
		棘轮扳手适配器,描述为BL/BLT, 长17mm		026.0084

辅助器械

产品	图片	描述	材料	产品编号
牙钻手柄				
Drill Handle 牙钻手柄		止停高度 1/3 毫米， $\varnothing$ 2.2 毫米，长度 104 毫米	不锈钢	034.150
		止停高度 1/3 毫米， $\varnothing$ 2.8 毫米，长度 104 毫米		034.250
		止停高度 1/3 毫米， $\varnothing$ 3.5 毫米，长度 104 毫米		034.450
		止停高度 1/3 毫米， $\varnothing$ 4.2 毫米，长度 104 毫米		034.650
C-Handle H2 C 型牙钻手柄 H2		止停高度 6 毫米，长度 89 毫米	不锈钢	034.750
C-Handle H4 C 型牙钻手柄 H4		止停高度 4 毫米，长度 89 毫米		034.751
C-Handle H6 C 型牙钻手柄 H6		止停高度 2 毫米，长度 89 毫米		034.752
Template Fixation Pin 固位杆		引导手术用， $\varnothing$ 2.8/2.8 毫米	不锈钢	034.256
		引导手术用， $\varnothing$ 5/2.8 毫米		034.255
		引导手术用， $\varnothing$ 5/3.5 毫米		034.455
		引导手术用， $\varnothing$ 5/4.2 毫米		034.655
Connector 适配器		使用于棘轮扳手，长度 14 毫米	不锈钢	034.005
止停器				
Stop Key 止停器			不锈钢	034.006

引导手术工具盒

产品	图片	描述	材料	产品编号
Straumann® Guided Surgery Cassette Straumann® 引导手术外科工具盒		适用于 Straumann® 引导手术，长度 290 毫米，宽度 175 毫米，高度 62 毫米	塑料	034.001
Straumann® Basic Guided Surgery Cassette Straumann® 引导手术外科工具盒基础版				034.281

取出钻

产品	图片	描述	材料	产品编号
取出钻，适用于牙科种植体 Ø 3.3				
Explantation Drill 取出钻		短型，适用于牙科种植体 Ø3.3 毫米和 Ø3.5 毫米，长度 32 毫米，内径 Ø3.6 毫米，外径 Ø4.2 毫米	不锈钢	044.330
		中型，适用于牙科种植体 Ø3.3 毫米和 Ø3.5 毫米，长度 37.5 毫米，内径 Ø3.6 毫米，外径 Ø4.2 毫米		044.331
		长型，适用于牙科种植体 Ø3.3 毫米和 Ø3.5 毫米，长度 40.5 毫米，内径 Ø3.6 毫米，外径 Ø4.2 毫米		044.332
取出钻，适用于牙科种植体 Ø 4.1				
Explantation Drill 取出钻		短型，适用于牙科种植体 Ø4.1 毫米，长度 32 毫米，内径 Ø4.2 毫米，外径 Ø4.8 毫米	不锈钢	044.340
		中型，适用于牙科种植体 Ø4.1 毫米，长度 37.5 毫米，内径 Ø4.2 毫米，外径 Ø4.8 毫米		044.341
		长型，适用于牙科种植体 Ø4.1 毫米，长度 40.5 毫米，内径 Ø4.2 毫米，外径 Ø4.8 毫米		044.342
取出钻，适用于牙科种植体 Ø 4.8				
Explantation Drill 取出钻		短型，适用于牙科种植体 Ø4.8 毫米，长度 32 毫米，内径 Ø4.9 毫米，外径 Ø5.5 毫米	不锈钢	044.344
		中型，适用于牙科种植体 Ø4.8 毫米，长度 37.5 毫米，内径 Ø4.9 毫米，外径 Ø5.5 毫米		044.345
		长型，适用于牙科种植体 Ø4.8 毫米，长度 40.5 毫米，内径 Ø4.9 毫米，外径 Ø5.5 毫米		044.346

取出辅件

产品	图片	描述	材料	产品编号
骨水平种植体(不含基台)取出				
NC Guiding Cylinder NC 取出钻导向套		适用于取出 Ø3.3 毫米骨水平种植体，长度 10.7 毫米，Ø3.6 毫米	不锈钢	026.2801
RC Guiding Cylinder RC 取出钻导向套		适用于取出 Ø4.1 毫米骨水平种植体，长度 10.9 毫米，Ø4.2 毫米		026.4801
		适用于取出 Ø4.8 毫米骨水平种植体，长度 10.9 毫米，Ø4.9 毫米		026.6801
S/SP/TE 种植体(不含基台)取出				
SCS Guiding Cylinder SCS 取出钻导向套		适用于取出 Ø3.3 毫米窄颈种植体，长度 5 毫米，Ø3.6 毫米	纯钛/不锈钢	049.182
		适用于取出 Ø3.3 毫米种植体，长度 10.5 毫米，Ø3.6 毫米	不锈钢	049.082
		适用于取出 Ø4.1 毫米种植体，长度 10.5 毫米，Ø4.2 毫米		049.083
		适用于取出 Ø4.8 毫米种植体，长度 10.5 毫米，Ø4.9 毫米		049.167
S/SP/TE 种植体(含基台)取出				
RN Guiding Sleeve RN 导向套		适用于取出 Ø3.3 毫米牙科种植体，粘结固定和螺丝固定，不含匹配螺丝 048.350，长度 4.5 毫米，Ø3.6 毫米	不锈钢	049.850
		适用于取出 Ø4.1 毫米牙科种植体，粘结固定和螺丝固定，不含匹配螺丝 048.350，长度 5.6 毫米，Ø4.2 毫米		049.071
		适用于取出 Ø4.8 毫米牙科种植体，粘结固定和螺丝固定，不含匹配螺丝 049.181，长度 6.4 毫米，Ø4.9 毫米		049.851
SCS Occlusal Screw SCS 骀面螺丝		适用于 RN/WN 固位基底的骀面螺丝，长度 4.4 毫米	纯钛	048.350
		适用于 RN/WN 钛合金临时基台，长度 7.6 毫米	TAN	049.181

取出器械

产品	图片	描述	材料	产品编号
48h Explantation Device 48 小时取出器械		适用于骨水平 NC 种植体，长度 31.4 毫米，Ø8 毫米	不锈钢	026.2048
		适用于骨水平 RC 种植体，长度 31.2 毫米，Ø8 毫米		026.4048
		适用于软组织水平 RN/WN synOcta® 种植体，短型，长度 20.6 毫米，Ø8 毫米		80103
		适用于软组织水平 RN/WN synOcta® 种植体，长型，长度 26.6 毫米，Ø8 毫米		80087

骨块固定

产品	图片	描述	材料	产品编号
基础套装				
Straumann® Bone Block Fixation Cassette Straumann® 骨块固定工具箱		工具箱，用于保存骨块固定所使用到的钻头、迷你螺丝、螺丝刀，长度 130 毫米，宽度 118毫米，高度 25毫米	热塑塑料	041.032

适用于骨挤压

产品	图片	描述	材料	产品编号
Osteotome Instrument Kit for Bone Condensation 骨挤压器套装适用于骨挤压		包含器材	不锈钢/塑料	040.500
Cassette for Bone Condensation Osteotomes 骨挤压器盒适用于骨挤压		不包含器材		040.550
Osteotome 骨挤压器		直型, 适用于骨挤压, Ø2.2 毫米	不锈钢	046.500
		直型, 适用于骨挤压, Ø2.8 毫米		046.501
		直型, 适用于骨挤压, Ø3.5 毫米		046.502
		直型, 适用于骨挤压, Ø4.2 毫米		046.503
		弯角, 适用于骨挤压, Ø2.2 毫米		046.510
		弯角, 适用于骨挤压, Ø2.8 毫米		046.511
		弯角, 适用于骨挤压, Ø3.5 毫米		046.512
		弯角, 适用于骨挤压, Ø4.2 毫米		046.513
Depth Stop 止停套		适用于骨挤压器, Ø2.2 毫米	纯钛/不锈钢	046.550
		适用于骨挤压器, Ø2.8 毫米		046.551
		适用于骨挤压器, Ø3.5 毫米		046.552
		适用于骨挤压器, Ø4.2 毫米		046.553
SCS Occlusal Screw SCS 骀面螺丝		适用于 RN/WN 固位基底的骀面螺丝, 长度 4.4 毫米	纯钛	048.350
		4 枚装		048.350V4

上颌窦提升

产品	图片	描述	材料	产品编号
Osteotome Instrument Kit 上颌窦提升器套装		适用于部分上颌窦提升	不锈钢/塑料	040.501
Cassette for Sinus Floor Elevation Osteotomes 上颌窦提升器套装盒		不包含器材		040.560
Osteotome 上颌窦提升器		直型，适用于上颌窦提升，Ø2.8 毫米	不锈钢	046.521
		直型，适用于上颌窦提升，Ø3.5 毫米		046.522
		直型，适用于上颌窦提升，Ø4.2 毫米		046.523
		弯角，适用于上颌窦提升，Ø2.8 毫米		046.531
		弯角，适用于上颌窦提升，Ø3.5 毫米		046.532
		弯角，适用于上颌窦提升，Ø4.2 毫米		046.533
Depth Stop 止停套		适用于上颌窦提升器，Ø2.8 毫米	纯钛/不锈钢	046.551
		适用于上颌窦提升器，Ø3.5 毫米		046.552
		适用于上颌窦提升器，Ø4.2 毫米		046.553
SCS Occlusal Screw SCS 骀面螺丝		适用于 RN/WN 固位基底的骀面螺丝，长度4.4毫米	纯钛	048.350
		4 枚装		048.350V4

器械和配件

储存系统

产品	图片	描述	材料	产品编号
储存系统				
Storage System 储存系统		适用于种植体，包括抽屉和隔板，长度 226 毫米，宽度 280 毫米，高度 248 毫米	PS	041.500
替换抽屉颌分隔板				
Drawer for Storage System 储存系统抽屉		不包括分隔板，长度 282 毫米，宽度 220 毫米，深度 47 毫米	塑料	041.520
Divider Panel Set 分隔板		种植系统隔板，长度 184 毫米，高度 44 毫米，5枚装		041.530V5

盒子

产品	图片	描述	材料	产品编号
Ultrasonic Cleaning Cassette 超声清洁盒		超声清洁盒，长度 150 毫米，宽度 65 毫米，高度 25 毫米	塑料/硅树脂	040.175

修复工具

产品	图片	描述	材料	产品编号
工具				
Ratchet 棘轮扳手		包含拆装工具，长度 84 毫米	不锈钢	046.119
Torque Control Device 扭矩控制器		适用于棘轮扳手，长度 82 毫米		046.049
Service Instrument 拆装工具		适用于棘轮扳手，长度 50 毫米		046.108
Cleaning Brush 清洁刷		适用于棘轮扳手，长度 100 毫米，Ø4.5 毫米，4 枚装	不锈钢/尼龙	045.111V4
螺丝刀				
SCS Screwdriver SCS 螺丝刀		适用于棘轮扳手，超短型，长度 15 毫米	不锈钢	046.400
		适用于棘轮扳手，短型，长度 21 毫米		046.401
		适用于棘轮扳手，长型，长度 27 毫米		046.402
		适用于机用件，超短型，长度 20 毫米		046.410
		适用于机用件，短型，长度 26 毫米		046.411
		适用于机用件，长型，长度 32 毫米		046.412
Screwdriver Handle 螺丝刀把手		适用于 046.410/411/412，长度 109 毫米	不锈钢/PEEK	046.085
球帽工具				
Retentive anchor driver 球帽基台螺丝刀		长度 19 毫米	不锈钢	046.069
Screwdriver 螺丝刀		球帽阴型脱位工具，长度 37 毫米		046.154
		球帽阴型就位工具，长度 50 毫米		048.452
实心基台螺丝刀				
NNC Solid Abutment Driver NNC 实心基台螺丝刀		短型，长度 13 毫米	不锈钢	048.155
		长型，长度 19 毫米		048.156
Inserting Device 螺丝刀		适用于 RN 实心基台 6°，短型，长度 13 毫米		046.067
		适用于 RN 实心基台 6°，长型，长度 19 毫米		046.068
LOCATOR® 工具				
LOCATOR® Core Tool 核心工具		三个部件，长度 100 毫米	不锈钢	046.415
LOCATOR® Screwdriver 螺丝刀		适用于棘轮扳手，短型，长度 15 毫米		046.416
		适用于棘轮扳手，长型，长度 21 毫米		046.417

修复工具

产品	图片	描述	材料	产品编号
铰刀				
NNC Reamer NNC 铰刀		适用于 45° 肩台，长度 17 毫米，Ø8 毫米	不锈钢/Ti-N	046.237
RN/WN Reamer RN/WN 铰刀		适用于 45° 肩台，Ø8 毫米		046.243
Guide Pin 导向杆		适用于铰刀 046.243，适用于 RN，synOcta® 048.605，长度 50 毫米	不锈钢	046.246
		适用于铰刀 046.243，适用于 RN/WN synOcta® 1.5 基台，长度 50 毫米		046.241
		适用于铰刀 046.243，适用于 WN synOcta® 048.606，长度 50 毫米		046.247
		适用于铰刀 046.243，适用于 RN，实心基台，长度 50 毫米		046.242
		适用于铰刀 046.243，适用于 WN 实心基台，长度 50 毫米		046.244
		适用于铰刀 046.237，适用于实心基台，长度 53 毫米		046.249
Handle 把手		适用于铰刀和 TS 攻丝钻，长度 65 毫米	不锈钢/黄铜	046.240

 0086 制造商
 Zest Anchors, Inc.
 Escondido, CA 92029
 美国

修复工具

产品	图片	描述	材料	产品编号
攻丝钻				
Taps 攻丝钻		适用于 TS ， M2 × 0.25 毫米	不锈钢	044.570
杆卡工具				
Activator Set 加力钥匙		适用于杆卡，三件式，长度 50 毫米	不锈钢	046.150
Deactivator 减力钥匙		适用于 杆卡，迷你型，长度 66 毫米	Grilon BS/黄铜	046.151
		适用于 杆卡，迷你型，长度 66 毫米		046.152
把手				
Analog Holder 替代体抓持器		长度 105 毫米	铝/钢	046.239

 1250 **制造商**
Diener AG Precision Machining
CH-8424 Embrach
瑞士
Exclusively distributed by Institut Straumann AG, CH-4002 Basel, 瑞士

储存系统

产品	图片	描述	材料	产品编号
储存系统				
Storage System 储存系统		适用于修复手术，包括抽屉和隔板，长度 226 毫米，宽度 280 毫米，高度 248 毫米	PS	041.501

工具箱

产品	图片	描述	材料	产品编号
Prosthetic Cassette 塑料盒		适用于器械和辅助部件，全套(不含隔板)，长度 244 毫米，宽度 115 毫米，高度 20 毫米	阳极氧化铝	040.185
		适用于器械和辅助部件，长度 244 毫米，宽度 115 毫米，高度 20 毫米		041.022
Tray 托盘		适用于盒 041.022，长度 218 毫米，宽度 109 毫米，高度 16 毫米		041.211
Divider Panel 隔板		大型，适用于盘 041.440，长度 47.5 毫米，宽度 9.5 毫米，高度 1 毫米		041.612
		小型，适用于盘 041.440，长度 27.5 毫米，宽度 9.5 毫米，高度 1 毫米		041.615
Dish 盘子		适用于托盘 041.211，长度 56 毫米，宽度 31 毫米，高度 11 毫米		041.440
Calibration Cup 校准环		黄色，Ø80 毫米		041.412
Plexiglass Cassette/ Membranette 树脂玻璃盒		长度 150 毫米，宽度 100 毫米	塑料	140.055

修复设计辅助工具

产品	图片	描述	材料	产品编号
X 射线参考球				
X-ray Reference Sphere X 射线参考球		Ø5 毫米， 4 枚装， 一次性使用	不锈钢	049.076V4
单个钻模加工部件				
Drill Sleeve with Collar 袖口型牙钻套筒		长度 10 毫米， 内径 Ø2.4 毫米， 外径 Ø3.5 毫米， 4 枚装	纯钛	049.810V4
Stepped Pin 阶梯形杆		适用于 049.810V4 ， 长度 16 毫米， Ø2.2/3.5 毫米， 4 枚装		049.818V4
Pin 杆		适用于 049.810V4 ， 长度 16 毫米， Ø2.2 毫米， 4 枚装		049.816V4
		适用于 049.810V4 ， 长度 10 毫米， Ø2.2 毫米， 4 枚装		049.817V4
		适用于 049.810V4 ， 长度 16 毫米， Ø3.5 毫米， 4 枚装		049.819V4

修复用螺丝

产品	图片	描述	材料	产品编号
RN 和 WN				
SCS Occlusal Screw SCS 骀面螺丝		适用于 RN/WN 固位基底的骀面螺丝，长度4.4毫米	纯钛	048.350
		4 枚装		048.350V4
synOcta® (八角) Basal Screw 中央螺丝		适用于 RN/WN Variobase®，金基台和临时基台，长度 6.7 毫米	TAN	048.356
SCS Occlusal Screw SCS 骀面螺丝		适用于 RN/WN 钛合金临时基台，长度 7.6 毫米		049.181
SCS Guide Screw SCS 导向螺丝		适用于螺丝固定 RN/WN 杆卡基底，不能缩短，长度 6 毫米， 4 枚装	纯钛	048.360V4

修复用螺丝

产品	图片	描述	材料	产品编号
NC 和 RC				
NC Basal Screw NC 中央螺丝		适用于 NC 解剖基台，NC Variobase® 冠用，NC 金基台，NC 杆卡金基台，长度 7.9 毫米	TAN	025.2900
NC Basal Screw NC 中央螺丝		适用于 NC 粘接基台，长度 7.9 毫米		025.2908
RC Basal Screw RC 中央螺丝		适用于 RC 解剖基台，RC Variobase® 冠用，RC 金基台，RC 杆卡金基台，长度 7.9 毫米		025.4900
		适用于桥用金基台，长度 8.3 毫米		025.4902
		适用于 RC 粘接基台和临时粘接基台，仅可用于粘接，长度 7.9 毫米，		025.4908



信息索引：

为了使您更便捷地查找到相关产品，您可以使用信息索引功能查找。信息索引功能分为：字母索引（按照货品名称首字母从A-Z排列）以及货号索引。

索引

字母索引	155
------	-----

货号索引	173
------	-----

字母索引

产品名称	货号	页数
48h Explantation device for BL Implant NC, L 31.4mm, stainless steel	026.2048	139
48h Explantation device for BL Implant RC, L 31.2mm, stainless steel	026.4048	139
Activator set, for all bar matrices, three-part, stainless steel	046.150	85,86,104,105,146
Adapter for handpiece, long, L 34mm, stainless steel	046.472	132
Adapter for handpiece, short, L 26mm, stainless steel	046.471	132
Adapter for handpiece, x-short, L 19mm, stainless steel	046.470	132
Adapter for ratchet, long, L 28mm, stainless steel	046.462	132
Adapter for ratchet, short, L 18mm, stainless steel	046.461	132
Adapter for ratchet, x-short, L 11mm, stainless steel	046.460	132
Alignment pin straight, Ø 2.2 mm, length 28.0 mm, titanium	046.458	130
Alignment Pin, Ø 2.2mm, L 27mm, Ti	046.703	129
Analog for individualized RN Solid Abutments, 2 pieces, PPS/brass	048.117V4	49
Analog for RN Bar constructions, w. synOcta® 1.5, L 11.5mm, stainless steel	048.108	44
Analog for RN synOcta®, L 12mm, grey w. red stripe, stainless steel	048.124	44
Analog holder, length 205.0 mm, steel, anodized aluminium	046.239	83,84,87,103,106,146
Basal Screw RN/WN Variobase®, TAN	022.0045	47,58
Bite registration aid, long, Ø 5 mm, H 12 mm, POM	048.941V4	44,56
Bite registration aid, short, Ø 5 mm, H 8 mm, POM	048.940V4	44,56
BL Bone Profiler 1, Ø 5.2mm, L 23mm, stainless steel	026.0022*	130
BL Bone Profiler 2, Ø 6.6mm, L 23mm, stainless steel	026.0023*	130
BL Bone Profiler 3, Ø 6mm, L 23mm, stainless steel	026.0024*	130
BL Profile drill, Ø 4.1mm, guided, L 37mm, stainless steel	026.4503	134
BL/NNC Implant profile drill, Ø 3.3mm, long, stainless steel	026.2306	128
BL/NNC Implant profile drill, Ø 3.3mm, short, stainless steel	026.2303	128
BL/NNC Profile drill, Ø 3.3mm, guided, L 37mm, stainless steel	026.2503	134
BL/TE Tap f. adapter, f. Ø 4.1mm implants, L 23mm, polymer/stainless steel	026.4310	129
BL/TE Tap f. adapter, f. Ø 4.8mm implants, L 23mm, polymer/stainless steel	026.6310	129
BL/TE Tap f. handpiece, Ø 4.1mm, guided, L 42mm, stainless steel	026.4510	134
BL/TE Tap f. handpiece, Ø 4.8mm, guided, L 42mm, stainless steel	026.6510	134
BL/TE/NNC Tap f. adapter, f. Ø 3.3mm implants, L 23mm, PEEK/stainless steel	026.2310	129
BL/TE/NNC Tap f. handpiece, Ø 3.3mm, guided, L 42mm, stainless steel	026.2510	134
BLT Drill, guided, Ø 2.8mm, long, stainless steel	034.262	134
BLT Drill, guided, Ø 2.8mm, medium, stainless steel	034.261	134
BLT Drill, guided, Ø 2.8mm, short, stainless steel	034.260	134
BLT Drill, guided, Ø 3.5mm, long, stainless steel	034.265	134
BLT Drill, guided, Ø 3.5mm, medium, stainless steel	034.264	134
BLT Drill, guided, Ø 3.5mm, short, stainless steel	034.263	134
BLT Drill, guided, Ø 4.2mm, long, stainless steel	034.268	134
BLT Drill, guided, Ø 4.2mm, medium, stainless steel	034.267	134
BLT Drill, guided, Ø 4.2mm, short, stainless steel	034.266	134
BLT Drill, long, Ø 2.8mm, L 41mm, stainless steel	026.2201	127
BLT Drill, long, Ø 3.5mm, L 41mm, stainless steel	026.4201	127
BLT Drill, long, Ø 4.2mm, L 41mm, stainless steel	026.6201	127
BLT Drill, short, Ø 2.8mm, L 33mm, stainless steel	026.2200	127
BLT Drill, short, Ø 3.5mm, L 33mm, stainless steel	026.4200	127
BLT Drill, short, Ø 4.2mm, L 33mm, stainless steel	026.6200	127
BLT Pilot Drill, guided, Ø 2.2mm, long, stainless steel	034.259	133
BLT Pilot Drill, guided, Ø 2.2mm, medium, stainless steel	034.258	133
BLT Pilot Drill, guided, Ø 2.2mm, short, stainless steel	034.257	133
BLT Pilot Drill, long, Ø 2.2mm, L 41mm, stainless steel	026.0002	127
BLT Pilot Drill, short, Ø 2.2mm, L 33mm, stainless steel	026.0001	127
BLT Profile Drill, guided, Ø 3.3mm, L 36mm, stainless steel	034.269	134
BLT Profile Drill, guided, Ø 4.1mm, L 36mm, stainless steel	034.270	134
BLT Profile Drill, guided, Ø 4.8mm, L 36mm, stainless steel	034.271	134
BLT Profile Drill, long, Ø 3.3mm, L 33mm, stainless steel	026.0004	128
BLT Profile Drill, long, Ø 4.1mm, L 33mm, stainless steel	026.0006	128
BLT Profile Drill, long, Ø 4.8mm, L 33mm, stainless steel	026.0008	128
BLT Profile Drill, short, Ø 3.3mm, L 25mm, stainless steel	026.0003	128
BLT Profile Drill, short, Ø 4.1mm, L 25mm, stainless steel	026.0005	128
BLT Profile Drill, short, Ø 4.8mm, L 25mm, stainless steel	026.0007	128
BLT Tap, guided, Ø 3.3mm, L 42mm, stainless steel	034.272	134

产品名称	货号	页数
BLT Tap, guided, Ø 4.1mm, L 42mm, stainless steel	034.273	134
BLT Tap, guided, Ø 4.8mm, L 42mm, stainless steel	034.274	134
BLT Tap, Ø 3.3mm, L 25mm, stainless steel/TAN	026.0009	129
BLT Tap, Ø 4.1mm, L 25mm, stainless steel/TAN	026.0010	129
BLT Tap, Ø 4.8mm, L 25mm, stainless steel/TAN	026.0011	129
Bone Block Fixation cassette for sterilization and storage of drills,mini-screws, screwdriver, 130mmx118mmx25mm	041.032	140
Bone Level implant distance indicator, L 100mm, Ti/stainless steel	026.0901	126
Bone Level Implant Ø 3.3mm NC, SLActive® 10mm, Roxolid® Loxim®	021.2310	24
Bone Level Implant Ø 3.3mm NC, SLActive® 12mm, Roxolid® Loxim®	021.2312	24
Bone Level Implant Ø 3.3mm NC, SLActive® 14mm, Roxolid® Loxim®	021.2314	24
Bone Level Implant Ø 3.3mm NC, SLActive® 8mm, Roxolid® Loxim®	021.2308	24
Bone Level Implant Ø 4.1mm RC, SLActive® 10mm, Roxolid® Loxim®	021.4310	26
Bone Level Implant Ø 4.1mm RC, SLActive® 12mm, Roxolid® Loxim®	021.4312	26
Bone Level Implant Ø 4.1mm RC, SLActive® 14mm, Roxolid® Loxim®	021.4314	26
Bone Level Implant Ø 4.1mm RC, SLActive® 8mm, Roxolid® Loxim®	021.4308	26
Bone Level Implant Ø 4.8mm RC, SLActive® 10mm, Roxolid® Loxim®	021.6310	28
Bone Level Implant Ø 4.8mm RC, SLActive® 12mm, Roxolid® Loxim®	021.6312	28
Bone Level Implant Ø 4.8mm RC, SLActive® 14mm, Roxolid® Loxim®	021.6314	28
Bone Level Implant Ø 4.8mm RC, SLActive® 8mm, Roxolid® Loxim®	021.6308	28
Bone Level Implant profile drill, Ø 4.1mm, long, stainless steel	026.4306	128
Bone Level Implant profile drill, Ø 4.1mm, short, stainless steel	026.4303	128
Bone Level Implant profile drill, Ø 4.8mm, long, stainless steel	026.6306	128
Bone Level Implant profile drill, Ø 4.8mm, short, stainless steel	026.6303	128
Bone Level Implant, Ø 3.3 mm, NC, SLA 10 mm, Ti	021.2410	24
Bone Level Implant, Ø 3.3 mm, NC, SLA 12 mm, Ti	021.2412	24
Bone Level Implant, Ø 3.3 mm, NC, SLA 14 mm, Ti	021.2414	24
Bone Level Implant, Ø 3.3 mm, NC, SLA 8 mm, Ti	021.2408	24
Bone Level Implant, Ø 3.3 mm, NC, SLActive 10 mm, Ti	021.2110	24
Bone Level Implant, Ø 3.3 mm, NC, SLActive 12 mm, Ti	021.2112	24
Bone Level Implant, Ø 3.3 mm, NC, SLActive 14 mm, Ti	021.2114	24
Bone Level Implant, Ø 3.3 mm, NC, SLActive 8 mm, Ti	021.2108	24
Bone Level Implant, Ø 4.1 mm, RC, SLA 10 mm, Ti	021.4410	26
Bone Level Implant, Ø 4.1 mm, RC, SLA 12 mm, Ti	021.4412	26
Bone Level Implant, Ø 4.1 mm, RC, SLA 14 mm, Ti	021.4414	26
Bone Level Implant, Ø 4.1 mm, RC, SLA 8 mm, Ti	021.4408	26
Bone Level Implant, Ø 4.1 mm, RC, SLActive 10 mm, Ti	021.4110	26
Bone Level Implant, Ø 4.1 mm, RC, SLActive 10 mm, Ti, guided	021.4110G	26
Bone Level Implant, Ø 4.1 mm, RC, SLActive 12 mm, Ti	021.4112	26
Bone Level Implant, Ø 4.1 mm, RC, SLActive 12 mm, Ti, guided	021.4112G	26
Bone Level Implant, Ø 4.1 mm, RC, SLActive 14 mm, Ti	021.4114	26
Bone Level Implant, Ø 4.1 mm, RC, SLActive 8 mm, Ti	021.4108	26
Bone Level Implant, Ø 4.1 mm, RC, SLActive 8 mm, Ti, guided	021.4108G	26
Bone Level Implant, Ø 4.8 mm, RC, SLA 10 mm, Ti	021.6410	28
Bone Level Implant, Ø 4.8 mm, RC, SLA 12 mm, Ti	021.6412	28
Bone Level Implant, Ø 4.8 mm, RC, SLA 14 mm, Ti	021.6414	28
Bone Level Implant, Ø 4.8 mm, RC, SLA 8 mm, Ti	021.6408	28
Bone Level Implant, Ø 4.8 mm, RC, SLActive 10 mm, Ti	021.6110	28
Bone Level Implant, Ø 4.8 mm, RC, SLActive 10 mm, Ti, guided	021.6110G	28
Bone Level Implant, Ø 4.8 mm, RC, SLActive 12 mm, Ti	021.6112	28
Bone Level Implant, Ø 4.8 mm, RC, SLActive 12 mm, Ti, guided	021.6112G	28
Bone Level Implant, Ø 4.8 mm, RC, SLActive 14 mm, Ti	021.6114	28
Bone Level Implant, Ø 4.8 mm, RC, SLActive 8 mm, Ti	021.6108	28
Bone Level Implant, Ø 4.8 mm, RC, SLActive 8 mm, Ti, guided	021.6108G	28
Bone Level Implant, Ø 3.3 mm, NC, SLActive 10 mm, Roxolid, guided	021.2210G	24
Bone Level Implant, Ø 3.3 mm, NC, SLActive 12 mm, Roxolid, guided	021.2212G	24
Bone Level Implant, Ø 3.3 mm, NC, SLActive 8 mm, Roxolid, guided	021.2208G	24
Bone Level Implant, Ø 4.1mm RC SLActive® 10mm, Roxolid®, guided	021.4310G	26
Bone Level Implant, Ø 4.1mm RC SLActive® 12mm, Roxolid®, guided	021.4312G	26
Bone Level Implant, Ø 4.1mm RC SLActive® 8mm, Roxolid®, guided	021.4308G	26
Bone Level Implant, Ø 4.8mm RC SLActive® 10mm, Roxolid®, guided	021.6310G	28
Bone Level Implant, Ø 4.8mm RC SLActive® 12mm, Roxolid®, guided	021.6312G	28
Bone Level Implant, Ø 4.8mm RC SLActive® 8mm, Roxolid®, guided	021.6308G	28
Bone Level profile drill, Ø 4.8mm, guided, L 37mm, stainless steel	026.6503	134
Bone Level Tapered Implant, Ø 3.3mm NC, SLActive® 10mm, Roxolid®, guided	021.3310G	24
Bone Level Tapered Implant, Ø 3.3mm NC, SLActive® 10mm, Roxolid®, Loxim®	021.3310	24
Bone Level Tapered Implant, Ø 3.3mm NC, SLActive® 12mm, Roxolid®, guided	021.3312G	24
Bone Level Tapered Implant, Ø 3.3mm NC, SLActive® 12mm, Roxolid®, Loxim®	021.3312	24

产品名称	货号	页数
Bone Level Tapered Implant, Ø 3.3mm NC, SLActive® 14mm, Roxolid®, Loxim®	021.3314	24
Bone Level Tapered Implant, Ø 3.3mm NC, SLActive® 16mm, Roxolid®, Loxim®	021.3316	24
Bone Level Tapered Implant, Ø 3.3mm NC, SLActive® 8mm, Roxolid®, guided	021.3308G	24
Bone Level Tapered Implant, Ø 3.3mm NC, SLActive® 8mm, Roxolid®, Loxim®	021.3308	24
Bone Level Tapered Implant, Ø 4.1mm RC, SLActive® 10mm, Roxolid®, Loxim®	021.5310	26
Bone Level Tapered Implant, Ø 4.1mm RC, SLActive® 12mm, Roxolid®, Loxim®	021.5312	26
Bone Level Tapered Implant, Ø 4.1mm RC, SLActive® 14mm, Roxolid®, Loxim®	021.5314	26
Bone Level Tapered Implant, Ø 4.1mm RC, SLActive® 16mm, Roxolid®, Loxim®	021.5316	26
Bone Level Tapered Implant, Ø 4.1mm RC, SLActive® 8mm, Roxolid®, Loxim®	021.5308	26
Bone Level Tapered Implant, Ø 4.8mm RC, SLActive® 10mm, Roxolid®, guided	021.7310G	28
Bone Level Tapered Implant, Ø 4.8mm RC, SLActive® 10mm, Roxolid®, Loxim®	021.7310	28
Bone Level Tapered Implant, Ø 4.8mm RC, SLActive® 12mm, Roxolid®, guided	021.7312G	28
Bone Level Tapered Implant, Ø 4.8mm RC, SLActive® 12mm, Roxolid®, Loxim®	021.7312	28
Bone Level Tapered Implant, Ø 4.8mm RC, SLActive® 14mm, Roxolid®, Loxim®	021.7314	28
Bone Level Tapered Implant, Ø 4.8mm RC, SLActive® 16mm, Roxolid®, Loxim®	021.7316	28
Bone Level Tapered Implant, Ø 4.8mm RC, SLActive® 8mm, Roxolid®, guided	021.7308G	28
Bone Level Tapered Implant, Ø 4.8mm RC, SLActive® 8mm, Roxolid®, Loxim®	021.7308	28
Bone profiler for RN, L 28mm, stainless steel, needed guiding cylinder: 049.082	044.304	130
Bone profiler for WN, L 25mm, stainless steel, needed guiding cylinder: 049.082	044.305	130
Calibration cup, yellow, Ø 80mm, anodized aluminium	041.412	148
CARES® Mono Scanbody for screw-retained abutment, abutment level, Ø 3.5mm, PEEK/TAN	025.0000	76
CARES® Mono Scanbody for screw-retained abutment, abutment level, Ø 4.6mm, PEEK/TAN	025.0001	76,98
CARES® NC Mono Scanbody, Ø 3.5 mm, H 10 mm, PEEK, TAN	025.2915	70
CARES® NNC Mono Scanbody, Ø 3.5 mm, H 10 mm, PEEK, TAN	048.173	34
CARES® RC Mono Scanbody, Ø 4.1 mm, H 10 mm, PEEK, TAN	025.4915	92,96
CARES® RN Mono Scanbody, Ø 4.8 mm, H 10 mm, PEEK, TAN	048.168	44
CARES® WN Mono Scanbody, Ø 6.5 mm, H 10 mm, PEEK, TAN	048.169	56
Cassette without content, for bone condensation	040.550	141
Cassette without content, for sinus floor elevation	040.560	142
C-handle H2, 6mm stop height, L 89mm, stainless steel	034.750	136
C-handle H4, 4mm stop height, L 89mm, stainless steel	034.751	136
C-handle H6, 2mm stop height, L 89mm, stainless steel	034.752	136
Cleaning brush for ratchet, L 100mm, Ø 4.5mm, nylon/stainless steel	045.111V4	131,144
Connector for ratchet, L 14mm, stainless steel	034.005	136
CrossFit® Plan Cassette, empty, PC	026.0000	110
CrossFit® Plan Set, POM	026.0902	110
Deactivator for Dolder® Bar matrices, mini, L 66mm, Grilon BS/brass	046.151	85,86,104,105,146
Deactivator for Dolder® Bar matrices, regular, L 66mm, Grilon BS/brass	046.152	85,86,104,105,146
Depth gauge Ø 3.5mm, L 26mm, Ti	046.450	130
Depth gauge Ø 4.2mm, L 26mm, Ti for Ø 4.8mm implants	046.451	130
Depth gauge with distance indicator Ø 2.2/2.8mm, L 27mm, Ti	046.455	130
Depth Gauge, Ø 2.8mm, L 27mm, Ti	046.705	129
Depth Gauge, Ø 3.5mm, L 27mm, Ti	046.706	129
Depth Gauge, Ø 4.2 mm, L 27mm, Ti	046.707	129
Depth Gauge, with distance indicator, Ø 2.2/2.8mm, L 27mm, Ti	046.704	129
Depth stop for osteotomes, Ø 2.2mm, Ti/stainless steel	046.550	141
Depth stop for osteotomes, Ø 2.8mm, Ti/stainless steel	046.551	141,142
Depth stop for osteotomes, Ø 3.5mm, Ti/stainless steel	046.552	141,142
Depth stop for osteotomes, Ø 4.2mm, Ti/stainless steel	046.553	141,142
Diagnostic T, L 65mm, W 33mm, Ti	046.147	126
Dish for tray 041.211, 56x31x11 mm, Al	041.440	148
Divider set, flat, for implants, 184 x 44 mm, plastic	041.530V5	143
Divider, large, for tray 041.211, 47.5 x 9.5 x 1mm, anodized aluminium	041.612	148
Divider, small, for dish 041.440, 27.5 x 9.5 x 1mm, anodized aluminium	041.615	148
Drawer for storage system, high, 282 x 220 x 47 mm, plastic	041.520	143
Drill handle, Ø 2.2mm, 1mm/3mm stop, L 104mm, stainless steel	034.150	136
Drill handle, Ø 2.8mm, 1mm/3mm stop, L 104mm, stainless steel	034.250	136
Drill handle, Ø 3.5mm, 1mm/3mm stop, L 104mm, stainless steel	034.450	136
Drill handle, Ø 4.2mm, 1mm/3mm stop, L 104mm, stainless steel	034.650	136
Drill sleeve with collar, inner Ø 2.4mm, outer Ø 3.5mm, L 10mm, Ti	049.810V4	149
Explantation drill, long, for Ø 3.3mm, Ø 3.5mm, L 40.5mm, stainless steel	044.332	137
Explantation drill, long, for Ø 4.1mm, L 40.5mm, stainless steel	044.342	137
Explantation drill, long, for Ø 4.8mm, L 40.5mm, stainless steel	044.346	137
Explantation drill, medium, for Ø 3.3mm, Ø 3.5mm, L 37.5mm, stainless steel	044.331	137
Explantation drill, medium, for Ø 4.1mm, L 37.5mm, stainless steel	044.341	137
Explantation drill, medium, for Ø 4.8mm, L 37.5mm, stainless steel	044.345	137
Explantation drill, short, for Ø 3.3mm, Ø 3.5mm, L 32mm, stainless steel	044.330	137

产品名称	货号	页数
Explantation drill, short, for Ø 4.1mm, L 32mm, stainless steel	044.340	137
Explantation drill, short, for Ø 4.8mm, L 32mm, stainless steel	044.344	137
Explantation instrument for synOcta® implants, long, L 26.6mm, stainless steel	80087	139
Explantation instrument for synOcta® implants, short, L 20.6mm, stainless steel	80103	139
Extension shell for angled RN synOcta® Abutm., H 5.5mm, POM	048.670	46
Guide cylinder for implant Ø 4.8mm, Ø 4.9mm, L 10.5mm, stainless steel	049.167	138
Guide cylinder for implants Ø 3.3mm & Ø 3.5mm, Ø 3.6mm, L 10.5mm, stainless steel	049.082	130,138
Guide cylinder for implants Ø 4.1mm, Ø 4.2mm, L 10.5mm, stainless steel	049.083	138
Guide pin for reamer, for 048.540-542, L 50mm, stainless steel	046.242	145
Guide pin for reamer, for 048.545-546, L 50mm, stainless steel	046.244	145
Guide pin for reamer, for 048.605, L 50mm, stainless steel	046.246	145
Guide pin for reamer, for 048.606, L 50mm, stainless steel	046.247	145
Guide pin for reamer, for synOcta®, L 50mm, stainless steel	046.241	145
Guide sleeve for implants Ø 3.3mm, Ø 3.6mm, L 4.5mm, stainless steel	049.850	138
Guide sleeve for implants Ø 4.1mm & Ø 3.5mm, Ø 4.2mm, L 5.6mm, stainless steel	049.071	138
Guide sleeve for implants Ø 4.8mm, Ø 4.9mm, L 6.4mm, stainless steel	049.851	138
Handle for reamers and Tap for TS, L 65mm, steel/brass	046.240	145
Hexagonal screwdriver, L 30mm, stainless steel	046.421	132
Holding key, L 85mm, stainless steel	046.064	131
Implant distance indicator, L 100mm, Ti/stainless steel	046.148	126
Impression Post f. closed tray f. screw-ret. abut., abut. level, Ø 3.5mm, non-engaging, TAN/POM	025.0013	76
Impression Post f. closed tray f. screw-ret. abut., abut. level, Ø 4.6mm, non-engaging, TAN/POM	025.0014	76,98
Impression Post f. open tray f. screw-ret. abut., abut. level, Ø 3.5mm, non-engaging, TAN	025.0011	76
Impression Post f. open tray f. screw-ret. abut., abut. level, Ø 4.6mm, non-engaging, TAN	025.0012	76,98
Impression post for closed tray, abutment level, for screw-retained abutment, Ø 3.5mm, TAN/POM	025.2245	76
Impression post for closed tray, abutment level, for screw-retained abutment, Ø 4.6mm, TAN/POM	025.2246	76,98
Impression post for open tray, TAN, for screw-ret. abut., abut. level, 0°, Ø 3.5mm	025.2243	76
Impression post for open tray, TAN, for screw-ret. abut., abut. level, 0°, Ø 4.6mm	025.2244	76,98
Instrument tweezers, L 140mm, stainless steel	046.110	131
Lab Processing Screw, for Screw-retained Abutments, L 10mm, stainless steel	025.0052	76,98
Lab Processing Screw, for Screw-retained Abutments, stainless steel	025.0006	76,98
Light curette (pack of 5 pcs)	045.921	132
LOCATOR® abutment holder sleeve	046.413V4	41,54,65,90,109
LOCATOR® angle measurement guide, stainless steal	048.200	41,54,65,90,109
LOCATOR® black processing male	048.195V4	41,54,65,90,109
LOCATOR® block-out spacer, silicon	048.196V20	41,54,65,90,109
LOCATOR® core tool, stainless steel with Au	046.415	41,54,65,90,109,144
LOCATOR® female analog, ø 5.0 mm, aluminium	048.198V4	52,63
LOCATOR® Female Implant Analog, Ø 4.0 mm	025.0120-04	39,88,107
LOCATOR® impression coping, aluminium housing with black male LDPE	048.197V4	39,52,63,88,107
LOCATOR® Male processing package, extended range, titanium/nylon	048.182V2	40,53,64,89,108
LOCATOR® male processing package, titanium/nylon	048.189V2	40,53,64,89,108
LOCATOR® parallel post, LDPE	048.199V4	41,54,65,90,109
LOCATOR® replacement male, clear, 0°–10°, 5 lbs, 2.27 kg,H 1.7mm, nylon	048.190V4	40,53,64,89,108
LOCATOR® replacement male, extended range, green, height 1.7 mm, nylon	048.193V4	40,53,64,89,108
LOCATOR® Replacement male, extended range, orange, height 1.7 mm, nylon	048.188V4	40,53,64,89,108
LOCATOR® replacement male, extra light retention, blue, height 1.7 mm,nylon	048.192V4	40,53,64,89,108
LOCATOR® replacement male, extra light retention, ext. range, red,height 1.7 mm, nylon	048.194V4	40,53,64,89,108
LOCATOR® replacement male, light retention, pink, height 1.7 mm, nylon	048.191V4	40,53,64,89,108
LOCATOR® screwdriver for ratchet, length 21.0 mm	046.417	41,54,65,90,109,144
LOCATOR® screwdriver for ratchet, short, length 15.0 mm, stainless steel	046.416	41,54,65,90,109,144
LOCATOR® zero retention ext range, gray	048.181V4	40,53,64,89,108
Membranette	140.055	148
Milling cutter, Ø 2.8mm, guided, L 32.5mm, stainless steel	034.215	133
Milling cutter, Ø 3.5mm, guided, L 32.5mm, stainless steel	034.415	133
Milling cutter, Ø 4.2mm, guided, L 32.5mm, stainless steel	034.615	133
Mucosa punch Ø 3.4mm, guided, L 30mm, stainless steel	034.010	133
Mucosa punch Ø 4.0mm, guided, L 30mm, stainless steel	034.011	133
Mucosa punch Ø 4.7mm, guided, L 30mm, stainless steel	034.012	133
NC Abutment for bars, gold, H 10mm, Ceramicor®	022.2701	81,85
NC Abutment, for bars, Ti	022.2702	81,86
NC Analog for screw-retained abutment, TAN, straight 0°, Ø 3.5mm	023.2754	76
NC Analog for Screw-retained Abutments, edentulous, straight 0°, Ø 3.5mm, TAN	025.0049	76
NC Anatomic Abutment, angled 15°, GH 2mm, Ti	022.2152	71
NC Anatomic Abutment, angled 15°, GH 3.5mm, Ti	022.2154	71
NC Anatomic Abutment, straight, GH 2mm, Ti	022.2102	71
NC Anatomic Abutment, straight, GH 3.5mm, Ti	022.2104	71

产品名称	货号	页数
NC Basal Screw, for Cementable Abutment only, TAN	025.2908	73,151
NC Basal Screw, TAN	025.2900	70,71,74,75,81,84,85,86,87,151
NC Bite Registration Aid, H 12mm, POM	025.2212-04	70
NC Bite Registration Aid, H 8mm, POM	025.2208-04	70
NC Bone Level Implant Analog, L 12mm, Ti	025.2101	70
NC Burn-out coping f. screw-retained abut., POM, Bridge/Bar, Ø 3.5mm	023.2755	79,80
NC Burn-out coping for screw-retained abutment, POM, Crown, Ø 3.5mm	023.2748	79
NC Burn-out coping, for Variobase® Abutment, POM	023.2756	75
NC Burn-out coping, for Variobase® Abutment, POM	023.2756-04	75
NC Burn-out Coping, for Variobase®, for crown, AH 5.5mm, POM	023.0018	75
NC Burn-out Coping, for Variobase®, for crown, AH 5.5mm, POM	023.0018V4	75
NC Cementable Abutment, Ø 3.5mm, GH 1mm, AH 4mm, Ti	022.2311	73
NC Cementable Abutment, Ø 3.5mm, GH 1mm, AH 5.5mm, Ti	022.2315	73
NC Cementable Abutment, Ø 3.5mm, GH 2mm, AH 4mm, Ti	022.2312	73
NC Cementable Abutment, Ø 3.5mm, GH 2mm, AH 5.5mm, Ti	022.2316	73
NC Cementable Abutment, Ø 3.5mm, GH 3mm, AH 4mm, Ti	022.2313	73
NC Cementable Abutment, Ø 3.5mm, GH 3mm, AH 5.5mm, Ti	022.2317	73
NC Cementable Abutment, Ø 5mm, GH 1mm, AH 4mm, Ti	022.2321	73
NC Cementable Abutment, Ø 5mm, GH 1mm, AH 5.5mm, Ti	022.2325	73
NC Cementable Abutment, Ø 5mm, GH 2mm, AH 4mm, Ti	022.2322	73
NC Cementable Abutment, Ø 5mm, GH 2mm, AH 5.5mm, Ti	022.2326	73
NC Cementable Abutment, Ø 5mm, GH 3mm, AH 4mm, Ti	022.2323	73
NC Cementable Abutment, Ø 5mm, GH 3mm, AH 5.5mm, Ti	022.2327	73
NC Closure Screw, 4 pieces, Height 0 mm, Ti	024.2100-04	25
NC Closure Screw, 4 pieces, Height 0.5 mm, Ti	024.2105-04	25
NC Closure Screw, H 0.5mm, Ti	024.2105	25
NC Closure Screw, H 0mm, Ti	024.2100	25
NC Coping for screw-retained abutment, Bar, Ø 3.5mm, Ceramicor®	023.2753	80
NC Coping for screw-retained abutment, Bridge, Ø 3.5mm, Ceramicor®/POM	023.2752	79
NC Coping for screw-retained abutment, Crown, Ø 3.5mm, Ceramicor®/POM	023.2751	79
NC Coping for screw-retained abutment, Ti, Bar, Ø 3.5mm	023.2750	80
NC Coping for screw-retained abutment, Ti, Bridge, Ø 3.5mm	023.2749	79
NC Coping for screw-retained abutment, Ti, Crown, Ø 3.5mm	023.2747	79
NC Coping, for bar on multi-base Abutment, angled 25°, Ceramicor®	023.2702	85
NC Coping, for bar on Multi-Base Abutment, angled 25°, Ti	023.2706	86
NC Coping, for bar on multi-base Abutment, straight, Ø 3.5mm, Ceramicor®	023.2732	85
NC Coping, for bar on Multi-Base Abutment, straight, Ø 3.5mm, Ti	023.2736	86
NC Coping, for bar on multi-base Abutment, straight, Ø 4.5mm, Ceramicor®	023.2742	85
NC Coping, for bar on Multi-Base Abutment, straight, Ø 4.5mm, Ti	023.2746	86
NC Coping, for bridge on Cementable Abutment, Ø 3.5mm, POM	023.2310	73
NC Coping, for bridge on Cementable Abutment, Ø 5mm, POM	023.2320	73
NC Coping, for bridge on multi-base Abutment, angled 25°, Ceramicor®/POM	023.2701	83
NC Coping, for bridge on Multi-Base Abutment, angled 25°, POM	023.2704	83
NC Coping, for bridge on multi-base Abutment, straight, Ø 3.5mm, Ceramicor®/POM	023.2731	83
NC Coping, for bridge on Multi-Base Abutment, straight, Ø 3.5mm, POM	023.2734	83
NC Coping, for bridge on multi-base Abutment, straight, Ø 4.5mm, Ceramicor®/POM	023.2741	83
NC Coping, for bridge on Multi-Base Abutment, straight, Ø 4.5mm, POM	023.2744	83
NC Coping, for crown on Cementable Abutment, Ø 3.5mm, POM	023.2316	73
NC Coping, for crown on Cementable Abutment, Ø 5mm, POM	023.2326	73
NC Gold Abutment, for crown, H 3.7mm, Ceramicor®/POM	022.2410	74
NC Guiding Cylinder for BL Bone Profiler, Ø 3.25mm, L 12.5mm, TAN	026.00255*	130
NC Guiding cylinder for Ø 3.3mm, BL Implant explantation, stainless steel	026.2801	138
NC Healing Abutment, bottle shape, Ø 3.3mm, H 5mm, Ti	024.2236	25
NC Healing Abutment, bottle shape, Ø 3.3mm, H3.5mm, Ti	024.2234	25
NC Healing Abutment, conical shape, Ø 3.6mm, H 2mm, Ti	024.2222	25
NC Healing Abutment, conical shape, Ø 3.6mm, H 3.5mm, Ti	024.2224	25
NC Healing Abutment, conical shape, Ø 3.6mm, H 5mm, Ti	024.2226	25
NC Healing Abutment, conical shape, Ø 4.8mm, H 2mm, Ti	024.2242	25
NC Healing Abutment, conical shape, Ø 4.8mm, H 3.5mm, Ti	024.2244	25
NC Healing Abutment, conical shape, Ø 4.8mm, H 5mm, Ti	024.2246	25
NC Healing Abutment, customizable, Ø 5mm, PEEK	024.2270	25
NC Implant Analog, for Cementable Abutment, Ø 3.5mm, AH 4mm, TAN	025.2112	72
NC Implant Analog, for Cementable Abutment, Ø 3.5mm, AH 5.5mm, TAN	025.2116	72
NC Implant Analog, for Cementable Abutment, Ø 5mm, AH 4mm, TAN	025.2122	72
NC Implant Analog, for Cementable Abutment, Ø 5mm, AH 5.5mm, TAN	025.2126	72
NC Implant Analog, for Multi-base Abutment, angled 25°Ø 4mm, TAN	025.2106	82
NC Implant Analog, for Multi-base Abutment, straight, Ø 3.5mm, TAN	025.2131	82

产品名称	货号	页数
NC Implant Analog, for Multi-base Abutment, straight, Ø 4.5mm, TAN	025.2141	82
NC Impression Cap, for Cementable Abutment, Ø 3.5mm, AH 4mm, POM	025.2203-04	72
NC Impression Cap, for Cementable Abutment, Ø 5mm, AH 4mm, POM	025.2204-04	72
NC Impression Cap, for Cementable Abutment, Ø 5mm, AH 5.5mm, POM	025.2220-04	72
NC Impression Cap, for Multi-Base Abutment, straight, Ø 3.5mm, POM	025.2231-04	82
NC Impression Cap, for Multi-Base Abutment, straight, Ø 4.5mm, POM	025.2241-04	82
NC Impression Post Closed Tray, with 1 guide screw & 2 caps, TAN/POM	025.2201	70
NC Impression Post Open Tray, short, with guide screw, L 16.5mm, TAN	025.2202	70
NC Impression Post, for Multi-Base Abutment, angled 25°, TAN	025.2206	82
NC Impression Post, for Multi-Base Abutment, straight, Ø 3.5mm, TAN	025.2232	82
NC Impression Post, for Multi-Base Abutment, straight, Ø 4.5mm, TAN	025.2242	82
NC Impression post, for open tray, with guide screw, L 30mm, TAN	025.2205	70
NC LOCATOR® Abutment, height 2 mm, Ti alloy with TiN coating	022.2502	88
NC LOCATOR® Abutment, height 3 mm, Ti alloy with TiN coating	022.2503	88
NC LOCATOR® Abutment, height 4 mm, Ti alloy with TiN coating	022.2504	88
NC LOCATOR® Abutment, height 5 mm, Ti alloy with TiN coating	022.2505	88
NC LOCATOR® Abutment, height 6 mm, Ti alloy with TiN coating	022.2506	88
NC Multi-Base Abutment, straight, Ø 3.5mm, GH 1mm, TAN	022.2731	82,85,86
NC Multi-Base Abutment, straight, Ø 3.5mm, GH 2.5mm, TAN	022.2732	82,85,86
NC Multi-Base Abutment, straight, Ø 3.5mm, GH 4mm, TAN	022.2734	82,85,86
NC Multi-Base Abutment, straight, Ø 4.5mm, GH 1mm, TAN	022.2741	82,85,86
NC Multi-Base Abutment, straight, Ø 4.5mm, GH 2.5mm, TAN	022.2742	82,85,86
NC Multi-Base Abutment, straight, Ø 4.5mm, GH 4mm, TAN	022.2744	82,85,86
NC Multi-Base Abutment, Type A, angled 25°, TAN	022.2782	82,85,86
NC Multi-Base Abutment, Type B, angled 25°, TAN	022.2786	82,85,86
NC Novaloc® Plan Abutment, H 1-6 mm, POM	025.2646-04	111
NC PLAN Anatomic Abutment, angled 15°, GH 2mm, POM	025.2603-04	110
NC PLAN Anatomic Abutment, angled 15°, GH 3.5mm, POM	025.2604-04	110
NC PLAN Anatomic Abutment, straight, GH 2mm, POM	025.2601-04	110
NC PLAN Anatomic Abutment, straight, GH 3.5mm, POM	025.2602-04	110
NC PLAN Cementable Abutment, Ø 3.5mm, GH 1mm, AH 4mm, POM	025.2611-04	110
NC PLAN Cementable Abutment, Ø 3.5mm, GH 1mm, AH 5.5mm, POM	025.2615-04	110
NC PLAN Cementable Abutment, Ø 3.5mm, GH 2mm, AH 4mm, POM	025.2612-04	110
NC PLAN Cementable Abutment, Ø 3.5mm, GH 2mm, AH 5.5mm, POM	025.2616-04	110
NC PLAN Cementable Abutment, Ø 3.5mm, GH 3mm, AH 4mm, POM	025.2613-04	110
NC PLAN Cementable Abutment, Ø 3.5mm, GH 3mm, AH 5.5mm, POM	025.2617-04	110
NC PLAN Cementable Abutment, Ø 5mm, GH 1mm, AH 4mm, POM	025.2621-04	110
NC PLAN Cementable Abutment, Ø 5mm, GH 1mm, AH 5.5mm, POM	025.2625-04	110
NC PLAN Cementable Abutment, Ø 5mm, GH 2mm, AH 4mm, POM	025.2622-04	110
NC PLAN Cementable Abutment, Ø 5mm, GH 2mm, AH 5.5mm, POM	025.2626-04	110
NC PLAN Cementable Abutment, Ø 5mm, GH 3mm, AH 4mm, POM	025.2623-04	110
NC PLAN Cementable Abutment, Ø 5mm, GH 3mm, AH 5.5mm, POM	025.2627-04	110
NC PLAN Gold Abutment, for crown, POM	025.2610-04	110
NC PLAN LOCATOR® Abutment, H 2mm, POM	025.2642-04	111
NC PLAN LOCATOR® Abutment, H 3mm, POM	025.2643-04	111
NC PLAN LOCATOR® Abutment, H 4mm, POM	025.2644-04	111
NC PLAN LOCATOR® Abutment, H 5mm, POM	025.2645-04	111
NC PLAN Multi-Base Abutment, straight, Ø 3.5mm, GH 1mm, POM	025.2651-04	113
NC PLAN Multi-Base Abutment, straight, Ø 3.5mm, GH 2.5mm, POM	025.2652-04	113
NC PLAN Multi-Base Abutment, straight, Ø 3.5mm, GH 4mm, POM	025.2654-04	113
NC PLAN Multi-Base Abutment, straight, Ø 4.5mm, GH 1mm, POM	025.2656-04	113
NC PLAN Multi-Base Abutment, straight, Ø 4.5mm, GH 2.5mm, POM	025.2657-04	113
NC PLAN Multi-Base Abutment, straight, Ø 4.5mm, GH 4mm, POM	025.2659-04	113
NC PLAN Multi-Base Abutment, Type A, angled 25°, POM	025.2662-04	113
NC PLAN Multi-Base Abutment, Type B, angled 25°, POM	025.2665-04	113
NC Plan screw-retained abutment, POM, angled 17°, Ø 4.6mm, GH 2.5mm, Type A	025.2655-04	111
NC Plan screw-retained abutment, POM, angled 17°, Ø 4.6mm, GH 2.5mm, Type B	025.2658-04	111
NC Plan screw-retained abutment, POM, angled 30°, Ø 4.6mm, GH 2.5mm, Type A	025.2653-04	111
NC Plan screw-retained abutment, POM, angled 30°, Ø 4.6mm, GH 2.5mm, Type B	025.2660-04	111
NC Plan screw-retained abutment, POM, straight 0°, Ø 3.5mm, GH 2.5mm	025.2648-04	111
NC Plan screw-retained abutment, POM, straight 0°, Ø 4.6mm, GH 2.5mm	025.2650-04	111
NC Polishing Aid, stainless steel	025.2920	74,84,87
NC Polishing Aid, stainless steel	025.2920-04	74,84,87
NC Protective cap for screw-retained abutment, Ø 3.5mm, H 5mm, PEEK/TAN	024.2323-04	77
NC Protective cap for screw-retained abutment, Ø 3.5mm, H 6.5mm, PEEK/TAN	024.2324-04	77
NC Protective cap for screw-retained abutment, Ø 3.5mm, H 8mm, PEEK/TAN	024.2325-04	77
NC Protective Cap, for Cementable Abutment, Ø 3.5mm, AH 4mm, PEEK	024.2312-04	72

产品名称	货号	页数
NC Protective Cap, for Cementable Abutment, Ø 3.5mm, AH 5.5mm, PEEK	024.2316-04	72
NC Protective Cap, for Cementable Abutment, Ø 5mm, AH 4mm, PEEK	024.2322-04	72
NC Protective Cap, for Cementable Abutment, Ø 5mm, AH 5.5mm, PEEK	024.2326-04	72
NC Protective Cap, for Multi-Base Abutment, angled 25°, PEEK/TAN	024.2306-04	82
NC Protective Cap, for Multi-Base Abutment, straight, Ø 3.5mm, PEEK/TAN	024.2331-04	82
NC Protective Cap, for Multi-Base Abutment, straight, Ø 4.5mm, PEEK/TAN	024.2341-04	82
NC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 2.5mm, Type A, Pin	022.2749P	78
NC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 2.5mm, Type B, Pin	022.2750P	78
NC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 4mm, Type A, Pin	022.2755P	78
NC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 4mm, Type B, Pin	022.2756P	78
NC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 5.5mm, Type A, Pin	022.0010P	78
NC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 5.5mm, Type B, Pin	022.0011P	78
NC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 2.5mm, Type A, Pin	022.2751P	78
NC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 2.5mm, Type B, Pin	022.2752P	78
NC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 4mm, Type A, Pin	022.2757P	78
NC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 4mm, Type B, Pin	022.2758P	78
NC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 5.5mm, Type A, Pin	022.0012P	78
NC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 5.5mm, Type B, Pin	022.0013P	78
NC Screw-retained Abutment, TAN, straight 0°, Ø 3.5mm, GH 1mm, Pin	022.2745P	77
NC Screw-retained Abutment, TAN, straight 0°, Ø 3.5mm, GH 2.5mm, Pin	022.2746P	77
NC Screw-retained Abutment, TAN, straight 0°, Ø 3.5mm, GH 4mm, Pin	022.2753P	77
NC Screw-retained Abutment, TAN, straight 0°, Ø 4.6mm, GH 1mm, Pin	022.2747P	77
NC Screw-retained Abutment, TAN, straight 0°, Ø 4.6mm, GH 2.5mm, Pin	022.2748P	77
NC Screw-retained Abutment, TAN, straight 0°, Ø 4.6mm, GH 4mm, Pin	022.2754P	77
NC Temporary Abutment, Ø 5mm, H 12mm, PEEK/TAN	024.2370	70
NC Temporary Abutment, Bridge, Ø 3.5mm, H 11mm, TAN	024.2375	70
NC Temporary Abutment, Crown, Ø 3.5mm, H 11mm, TAN	024.2371	70
NC Temporary Coping, for Cementable Abutment, Ø 3.5mm, PMMA	024.2346	72
NC Temporary Coping, for Cementable Abutment, Ø 5mm, PMMA	024.2356	72
NC Temporary Coping, for Multi-Base Abutment, angled 25°, Ti	024.2307	82
NC Temporary Coping, for Multi-Base Abutment, straight, Ø 3.5mm, Ti	024.2332	82
NC Temporary Coping, for Multi-Base Abutment, straight, Ø 4.5mm, Ti	024.2342	82
NC Transfer Aid, for Cementable Abutment, Ø 3.5mm, PC	025.2911-04	72
NC Transfer Aid, for Cementable Abutment, Ø 5mm, PC	025.2912-04	72
NC Variobase® Abutment, with screw, Ø 3.8mm, H 3.5mm, TAN	025.2921	74,75
NC Variobase®, for crown, with screw, Ø 3.8mm, H 5.5mm, TAN	022.0027	74,75
NC/RC Analog for screw-retained abutment, TAN, angled 17°/30°, Ø 4.6mm	023.4757	76,98
NC/RC Analog for screw-retained abutment, TAN, straight 0°, Ø 4.6mm	023.4756	76,98
NC/RC Analog for Screw-retained Abutments, edentulous, straight 0°, Ø 4.6mm, TAN	025.0050	76,98
NC/RC Burn-out coping f. screw-retained abut., POM, Bridge/Bar, Ø 4.6mm	023.4758	79,80,100,101
NC/RC Burn-out coping for screw-retained abutment, POM, Crown, Ø 4.6mm	023.4748	79,100
NC/RC Coping for screw-retained abutment, Bar, Ø 4.6mm, Ceramicor®	023.4755	80,100
NC/RC Coping for screw-retained abutment, Bridge, Ø 4.6mm, Ceramicor®/POM	023.4754	79,100
NC/RC Coping for screw-retained abutment, Crown, Ø 4.6mm, Ceramicor®/POM	023.4753	79,100
NC/RC Coping for screw-retained abutment, Ti, Bar, Ø 4.6mm	023.4752	80,101
NC/RC Coping for screw-retained abutment, Ti, Bridge, Ø 4.6mm	023.4751	79,100
NC/RC Coping for screw-retained abutment, Ti, Crown, Ø 4.6mm	023.4747	79,100
NC/RC Occlusal Screw, for coping on, Multi-Base Abutment, TAN	025.0900-04	83,84,85,86,87,103,104,105,106
NC/RC Occlusal screw, TAN, for coping, screw-retained abutment	023.4763	79,80,100,101
NC/RC Protective cap for screw-retained abutment, Ø 4.6mm, H 5.1mm, PEEK/TAN	024.4323-04	77,98
NC/RC Protective cap for screw-retained abutment, Ø 4.6mm, H 6.6mm, PEEK/TAN	024.4324-04	77,98
NC/RC Protective cap for screw-retained abutment, Ø 4.6mm, H 8.1mm, PEEK/TAN	024.4325-04	77,98
NC/RC Screw for screw-retained abutment, TAN, straight 0°, GH 1mm	023.4749	78,99
NC/RC Screw for screw-retained abutment, TAN, straight 0°, GH 2.5mm	023.4750	78,99
NC/RC Screw for screw-retained abutment, TAN, straight 0°, GH 4mm	023.4760	78,99
NC/RC Screw, for screw-retained abutment 17°/30°, TAN	025.0002	78,99
NC/RC Stabilization pin for bar on Multi-base abutment, L 16mm, Stainless steel	025.0911-04	85,86,87,104,105,106
Needle Drill, long, Ø 1.6mm, L 41mm, stainless steel	026.0056	127
Needle Drill, short, Ø 1.6mm, L 33mm, stainless steel	026.0054	127
NN Analog, L 14mm, stainless steel	048.130	117
NN Coping 15° angled, H 8,8 mm, Ti	048.550	117
NN Coping 20° angled, H 8,8 mm, Ti	048.551	117
NN Coping, H 9mm, Ti	048.505	117
NN Gold coping, cast-on, H 10mm, Ceramicor®/POM	048.635	117
NN Impression cap with integral guide screw, L 17,5 mm, Al/TAN	048.016	117
NN Occlusal screw, L 5mm, Ø 1.8mm, TAN	049.177	117
NN Protective cap with integr.occlusal screw, Ø 4mm, H 3.4mm, PEEK/TAN	048.050	117

产品名称	货号	页数
NN Temporary Coping, H 10mm, PEEK/TAN	048.669	117
NNC Analog for Solid abutment, for 048.547, H 14.3mm, TAN	048.527	38
NNC Analog for Solid abutment, for 048.548, H 15.8mm, TAN	048.528	38
NNC Analog for Solid abutment, for 048.549, H 17.3mm, TAN	048.529	38
NNC Basal Screw, for cementable abutments, L 7.1mm, TAN	048.314	35
NNC Basal screw, L 7.1mm, TAN	048.313	34,36,37
NNC Bite registration aid, H 12mm, POM	048.961V4	34
NNC Bite registration aid, H 8mm, POM	048.960V4	34
NNC Burn-out coping, for Variobase® Abutment, POM	048.267	37
NNC Burn-out coping, for Variobase® Abutment, POM	048.267V4	37
NNC Burn-out Coping, for Variobase®, for crown, AH 5.75mm, POM	023.0014	37
NNC Burn-out Coping, for Variobase®, for crown, AH 5.75mm, POM	023.0014V4	37
NNC Cementable abutment, angled 15°, type A, H 5.9mm, Ti	048.621	35
NNC Cementable abutment, angled 15°, type B, H 5.9mm, Ti	048.622	35
NNC Cementable abutment, straight, for crowns/bridges, H 5.7mm, Ti	048.598	35
NNC Closure cap, Ø 2.6mm, H 0mm, Ti	048.324	14
NNC Closure cap, Ø 2.6mm, H 0mm, Ti	048.324V4	14
NNC Closure cap, Ø 4mm, H 1.5mm, Ti	048.325	14
NNC Closure cap, Ø 4mm, H 1.5mm, Ti	048.325V4	14
NNC Gold abutment, bridge, H 3.6mm, Ceramicor® / TAN / POM	048.593	36
NNC Gold abutment, crown, H 3.6mm, Ceramicor® / TAN / POM	048.592	36
NNC Guide pin for reamer 046.237, for Solid abutments, L 53mm, stainl. steel	046.249	145
NNC Healing cap, labial bevel, Ø 4mm, H 2mm, Ti	048.082	14
NNC Healing cap, Ø 4mm, H 3mm, Ti	048.071	14
NNC Healing cap, Ø 4mm, H 4.5mm, Ti	048.074	14
NNC Implant analog, L 12.5mm, Ti	048.127	34
NNC Impression cap, H 8mm, POM	048.579	38
NNC Impression post, closed tray, H 10.9mm, TAN/POM	048.134	34
NNC Impression post, for open tray, H 10.1mm, TAN	048.128	34
NNC LOCATOR® Abutment, H 1mm, Ti alloy with TiN coating	048.581	39
NNC LOCATOR® Abutment, H 2mm, Ti alloy with TiN coating	048.582	39
NNC LOCATOR® Abutment, H 3mm, Ti alloy with TiN coating	048.583	39
NNC LOCATOR® Abutment, H 4mm, Ti alloy with TiN coating	048.584	39
NNC LOCATOR® Abutment, H 5mm, Ti alloy with TiN coating	048.585	39
NNC LOCATOR® Abutment, H 6mm, Ti alloy with TiN coating	048.586	39
NNC PLAN Cementable abutment, angled 15°, type A, H 5.9mm, POM	048.944V4	66
NNC PLAN Cementable abutment, angled 15°, type B, H 5.9mm, POM	048.945V4	66
NNC PLAN Cementable abutment, straight, for crown/bridge, H 5.7mm, POM	048.943V4	66
NNC PLAN LOCATOR® abutment, H 1mm, POM	048.946V4	66
NNC PLAN LOCATOR® abutment, H 2mm, POM	048.947V4	66
NNC PLAN LOCATOR® abutment, H 3mm, POM	048.948V4	66
NNC PLAN LOCATOR® abutment, H 4mm, POM	048.949V4	66
NNC PLAN LOCATOR® abutment, H 5mm, POM	048.950V4	66
NNC PLAN Solid Abutment, H 4mm, POM	048.952V4	66
NNC PLAN Solid Abutment, H 5.5mm, POM	048.953V4	66
NNC PLAN Solid Abutment, H 7mm, POM	048.954V4	66
NNC Plastic coping for Solid abutment, for bridge, H 10.2mm, POM	048.264	38
NNC Plastic coping for Solid abutment, for crown, H 10.2mm, POM	048.263	38
NNC Plastic coping, for 048.598, for bridge, H 7mm, POM	048.257	35
NNC Plastic coping, for 048.598, for crown, H 7mm, POM	048.256	35
NNC Plastic coping, for 048.621/ 048.622, for crowns & bridges, H 7.2mm, POM	048.258	35
NNC Positioning cylinder for 048.547, H 10.2mm, POM	048.587	38
NNC Positioning cylinder for 048.548, H 10.2mm, POM	048.588	38
NNC Positioning cylinder for 048.549, H 10.2mm, POM	048.589	38
NNC Post for temporary restoration, bridge, H 10.3mm, TAN	048.502	34
NNC Post for temporary restoration, crown, H 10.3mm, TAN	048.501	34
NNC Protective cap for Solid abutment, for 048.547, H 5.6mm, PEEK	048.700	38
NNC Protective cap for Solid abutment, for 048.548, H 7.1mm, PEEK	048.701	38
NNC Protective cap for Solid abutment, for 048.549, H 8.6mm, PEEK	048.702	38
NNC Reamer for 45° shoulder, L 17mm, stainless steel / Ti-N	046.237	145
NNC Solid abutment driver, long, L 19mm, stainless steel	048.156	38,144
NNC Solid abutment driver, short, L 13mm, stainless steel	048.155	38,144
NNC Solid abutment, H 4mm, yellow, Ti	048.547	38
NNC Solid abutment, H 5.5mm, gray, Ti	048.548	38
NNC Solid abutment, H 7mm, blue, Ti	048.549	38
NNC Temporary coping for Solid abutment, for bridge, H 8.5mm, PMMA	048.659	38
NNC Temporary coping for Solid abutment, for crown, H 8.5mm, PMMA	048.658	38

产品名称	货号	页数
NNC Variobase® Abutment, with screw, Ø 3.9mm, H 3.75mm, TAN	048.709	36,37
NNC Variobase®, for crown, with screw, Ø 3.9mm, H 5.75mm, TAN	022.0021	36,37
NNC Wax-up sleeve, L 12.75 mm, PEEK	048.137	36
NNC Wax-up sleeve, L 12.75 mm, PEEK	048.137V4	36
Osteotome angled, for bone condensation, Ø 2.2mm, stainless steel	046.510	141
Osteotome angled, for bone condensation, Ø 2.8mm, stainless steel	046.511	141
Osteotome angled, for bone condensation, Ø 3.5mm, stainless steel	046.512	141
Osteotome angled, for bone condensation, Ø 4.2mm, stainless steel	046.513	141
Osteotome angled, for sinus floor elevation, Ø 2.8mm, stainless steel	046.531	142
Osteotome angled, for sinus floor elevation, Ø 3.5mm, stainless steel	046.532	142
Osteotome angled, for sinus floor elevation, Ø 4.2mm, stainless steel	046.533	142
Osteotome instrument kit, for bone condensation, stainless steel	040.500	141
Osteotome instrument kit, for sinus floor elevation, stainless steel	040.501	142
Osteotome straight, for bone condensation, Ø 2.2mm, stainless steel	046.500	141
Osteotome straight, for bone condensation, Ø 2.8mm, stainless steel	046.501	141
Osteotome straight, for bone condensation, Ø 3.5mm, stainless steel	046.502	141
Osteotome straight, for bone condensation, Ø 4.2mm, stainless steel	046.503	141
Osteotome straight, for sinus floor elevation, Ø 2.8mm, stainless steel	046.521	142
Osteotome straight, for sinus floor elevation, Ø 3.5mm, stainless steel	046.522	142
Osteotome straight, for sinus floor elevation, Ø 4.2mm, stainless steel	046.523	142
Pilot drill 1, long, Ø 2.2mm, L 41mm, stainless steel	044.211	127
Pilot drill 1, short, Ø 2.2mm, L 33mm, stainless steel	044.210	127
Pilot drill 2, long, Ø 2.8mm, L 41mm, stainless steel	044.215	127
Pilot drill 2, short, Ø 2.8mm, L 33mm, stainless steel	044.214	127
Pilot drill, Ø 2.2mm, long, 20mm, guided, L 36mm, stainless steel	034.126	133
Pilot drill, Ø 2.2mm, short, 16mm, guided, L 32mm, stainless steel	034.123	133
Pilot drill, Ø 2.2mm, x-long, 24mm, guided, L 40mm, stainless steel	034.129	133
Pin for 049.810, Ø 2.2mm, L 10mm, Ti	049.817V4	149
Pin for 049.810, Ø 2.2mm, L 16mm, Ti	049.816V4	149
Pin for 049.810, Ø 3.5mm, L 16mm, Ti	049.819V4	149
Plastic shoulder for RN synOcta® angled, angled 15° or 20°, H 4mm, POM	048.676V4	46
Plastic shoulder for WN synOcta® angled, angled 15°, H 4mm, POM	048.678	58
Polishing Aid, for Screw-retained Abutments Ø 3.5mm, TAN	025.0004	76
Polishing Aid, for Screw-retained Abutments Ø 3.5mm, TAN	025.0004V4	76
Polishing Aid, for Screw-retained Abutments Ø 4.6mm, TAN	025.0005	76,98
Polishing Aid, for Screw-retained Abutments Ø 4.6mm, TAN	025.0005V4	76,98
Polishing protector for RN synOcta® Coping, L 15 mm, stainless steel	046.245	48
Prosthetic cassette for instruments, 224 x 115 x 20 mm, anodized aluminium	041.022	148
Prosthetic cassette for instruments, complete, 224 x 115 x 20 mm, anodized aluminium	040.185	148
Ratchet, including service instrument, L 84mm, stainless steel	046.119	131,144
RC Abutment for bars, gold, H 10mm, Ceramicor®	022.4701	101,104
RC Abutment, for bars, Ti	022.4702	101,105
RC Anatomic Abutment, angled 15°, GH 2mm, Ti	022.4152	93
RC Anatomic Abutment, angled 15°, GH 3.5mm, Ti	022.4154	93
RC Anatomic Abutment, straight, GH 2mm, Ti	022.4102	93
RC Anatomic Abutment, straight, GH 3.5mm, Ti	022.4104	93
RC Basal Screw, for Cementable Abutment only, TAN	025.4908	92,95,151
RC Basal Screw, for Gold Abutment bridge only, TAN	025.4902	96,151
RC Basal Screw, TAN	025.4900	92,93,96,97,101,104,105,106,151
RC Bite Registration Aid, H 12mm, POM	025.4212-04	92
RC Bite Registration Aid, H 8mm, POM	025.4208-04	92
RC Bone Level Implant Analog, L 12mm, Ti	025.4101	92
RC Burn-out coping, for Variobase® Abutment, POM	023.4759	97
RC Burn-out coping, for Variobase® Abutment, POM	023.4759-04	97
RC Burn-out Coping, for Variobase®, for crown, AH 5.5mm, POM	023.0017	97
RC Burn-out Coping, for Variobase®, for crown, AH 5.5mm, POM	023.0017V4	97
RC Cementable Abutment, Ø 5mm, GH 1mm, AH 4mm, Ti	022.4321	95
RC Cementable Abutment, Ø 5mm, GH 1mm, AH 5.5mm, Ti	022.4325	95
RC Cementable Abutment, Ø 5mm, GH 2mm, AH 4mm, Ti	022.4322	95
RC Cementable Abutment, Ø 5mm, GH 2mm, AH 5.5mm, Ti	022.4326	95
RC Cementable Abutment, Ø 5mm, GH 3mm, AH 4mm, Ti	022.4323	95
RC Cementable Abutment, Ø 5mm, GH 3mm, AH 5.5mm, Ti	022.4327	95
RC Cementable Abutment, Ø 6.5mm, GH 1mm, AH 4mm, Ti	022.4331	95
RC Cementable Abutment, Ø 6.5mm, GH 1mm, AH 5.5mm, Ti	022.4335	95
RC Cementable Abutment, Ø 6.5mm, GH 2mm, AH 4mm, Ti	022.4332	95
RC Cementable Abutment, Ø 6.5mm, GH 2mm, AH 5.5mm, Ti	022.4336	95
RC Cementable Abutment, Ø 6.5mm, GH 3mm, AH 4mm, Ti	022.4333	95

产品名称	货号	页数
RC Cementable Abutment, Ø 6.5mm, GH 3mm, AH 5.5mm, Ti	022.4337	95
RC Closure Screw, 4 pieces, Height 0 mm, Ti	024.4100-04	27,29
RC Closure Screw, 4 pieces, Height 0.5 mm, Ti	024.4105-04	27,29
RC Closure Screw, H 0.5mm, Ti	024.4105	27,29
RC Closure Screw, H 0mm, Ti	024.4100	27,29
RC Coping, for bar on multi-base Abutment, angled 25°, Ceramicor®	023.4702	104
RC Coping, for bar on Multi-Base Abutment, angled 25°, Ti	023.4706	105
RC Coping, for bar on multi-base Abutment, straight, Ø 4.5mm, Ceramicor®	023.4742	104
RC Coping, for bar on Multi-Base Abutment, straight, Ø 4.5mm, Ti	023.4746	105
RC Coping, for bar on multi-base Abutment, straight, Ø 6.5mm, Ceramicor®	023.4762	104
RC Coping, for bar on Multi-Base Abutment, straight, Ø 6.5mm, Ti	023.4766	105
RC Coping, for bridge on Cementable Abutment, Ø 5mm, POM	023.4320	95
RC Coping, for bridge on Cementable Abutment, Ø 6.5mm, POM	023.4330	95
RC Coping, for crown on Cementable Abutment, Ø 5mm, POM	023.4326	95
RC Coping, for crown on Cementable Abutment, Ø 6.5mm, POM	023.4336	95
RC Gold Abutment, for bridge, H 5.2mm, Ceramicor®/POM	022.4430	96
RC Gold Abutment, for crown, H 3.8mm, Ceramicor®/POM	022.4410	96
RC Guiding Cylinder for BL Bone Profiler, Ø 3.7mm, L 12.5mm, TAN	026.0026S*	130
RC Guiding cylinder for Ø 4.1mm, BL Implant explantation, stainless steel	026.4801	138
RC Guiding cylinder for Ø 4.8mm, BL Implant explantation, stainless steel	026.6801	138
RC Healing Abutment, bottle shape, Ø 4.4mm, H 4mm, Ti	024.4234	27,29
RC Healing Abutment, bottle shape, Ø 4.7mm, H 6mm, Ti	024.4236	27,29
RC Healing Abutment, conical shape, Ø 4.5mm, H 2mm, Ti	024.4222	27,29
RC Healing Abutment, conical shape, Ø 4.5mm, H 4mm, Ti	024.4224	27,29
RC Healing Abutment, conical shape, Ø 4.5mm, H 6mm, Ti	024.4226	27,29
RC Healing Abutment, conical shape, Ø 6mm, H 2mm, Ti	024.4242	27,29
RC Healing Abutment, conical shape, Ø 6mm, H 4mm, Ti	024.4244	27,29
RC Healing Abutment, conical shape, Ø 6mm, H 6mm, Ti	024.4246	27,29
RC Healing Abutment, customizable, Ø 7mm, PEEK	024.4270	27,29
RC Implant Analog for Multi-base abutment,angled 25°, Ø 4mm, TAN	025.4106	102
RC Implant Analog for Multi-base abutment,straight, Ø 4.5mm, TAN	025.4141	102
RC Implant Analog for Multi-base abutment,straight, Ø 6.5mm, TAN	025.4161	102
RC Implant Analog, for Cementable Abutment, Ø 5mm, AH 4mm, TAN	025.4132	94
RC Implant Analog, for Cementable Abutment, Ø 5mm, AH 5.5mm, TAN	025.4136	94
RC Implant Analog, for Cementable Abutment, Ø 6.5mm, AH 4mm, TAN	025.4122	94
RC Implant Analog, for Cementable Abutment, Ø 6.5mm, AH 5.5mm, TAN	025.4126	94
RC Impression Cap, for Cementable Abutment, Ø 5mm, AH 4mm, POM	025.4203-04	94
RC Impression Cap, for Cementable Abutment, Ø 5mm, AH 5.5mm, POM	025.4220-04	94
RC Impression Cap, for Cementable Abutment, Ø 6.5mm, AH 4mm, POM	025.4204-04	94
RC Impression Cap, for Cementable Abutment, Ø 6.5mm, AH 5.5mm, POM	025.4230-04	94
RC Impression Cap, for Multi-Base Abutment, straight, Ø 4.5mm, POM	025.4241-04	102
RC Impression Cap, for Multi-Base Abutment, straight, Ø 6.5mm, POM	025.4261-04	102
RC Impression Post Closed Tray, with 1 guide screw & 2 caps, TAN/POM	025.4201	92
RC Impression Post Open Tray, short, with guide screw, L 16.5mm, TAN	025.4202	92
RC Impression Post, for Multi-Base Abutment, angled 25°, TAN	025.4206	102
RC Impression Post, for Multi-Base Abutment, straight, Ø 4.5mm, TAN	025.4242	102
RC Impression Post, for Multi-Base Abutment, straight, Ø 6.5mm, TAN	025.4262	102
RC Impression post, for open tray, with guide screw, L 30mm, TAN	025.4205	92
RC LOCATOR® Abutment, height 1 mm, Ti alloy with TiN coating	022.4501	107
RC LOCATOR® Abutment, height 2 mm, Ti alloy with TiN coating	022.4502	107
RC LOCATOR® Abutment, height 3 mm, Ti alloy with TiN coating	022.4503	107
RC LOCATOR® Abutment, height 4 mm, Ti alloy with TiN coating	022.4504	107
RC LOCATOR® Abutment, height 5 mm, Ti alloy with TiN coating	022.4505	107
RC LOCATOR® Abutment, height 6 mm, Ti alloy with TiN coating	022.4506	107
RC Multi-Base Abutment, straight, Ø 4.5mm, GH 1mm, TAN	022.4741	102,104,105
RC Multi-Base Abutment, straight, Ø 4.5mm, GH 2.5mm, TAN	022.4742	102,104,105
RC Multi-Base Abutment, straight, Ø 4.5mm, GH 4mm, TAN	022.4744	102,104,105
RC Multi-Base Abutment, straight, Ø 6.5mm, GH 1mm, TAN	022.4761	102,104,105
RC Multi-Base Abutment, straight, Ø 6.5mm, GH 2.5mm, TAN	022.4762	102,104,105
RC Multi-Base Abutment, straight, Ø 6.5mm, GH 4mm, TAN	022.4764	102,104,105
RC Multi-Base Abutment, Type A, angled 25°, TAN	022.4782	102,104,105
RC Multi-Base Abutment, Type B, angled 25°, TAN	022.4786	102,104,105
RC PLAN Anatomic Abutment, angled 15°, GH 2mm, POM	025.4603-04	111
RC PLAN Anatomic Abutment, angled 15°, GH 3.5mm, POM	025.4604-04	111
RC PLAN Anatomic Abutment, straight, GH 2mm, POM	025.4601-04	111
RC PLAN Anatomic Abutment, straight, GH 3.5mm, POM	025.4602-04	111
RC PLAN Cementable Abutment, Ø 5mm, GH 1mm, AH 4mm, POM	025.4621-04	111

产品名称	货号	页数
RC PLAN Cementable Abutment, Ø 5mm, GH 1mm, AH 5.5mm, POM	025.4625-04	111
RC PLAN Cementable Abutment, Ø 5mm, GH 2mm, AH 4mm, POM	025.4622-04	111
RC PLAN Cementable Abutment, Ø 5mm, GH 2mm, AH 5.5mm, POM	025.4626-04	111
RC PLAN Cementable Abutment, Ø 5mm, GH 3mm, AH 4mm, POM	025.4623-04	111
RC PLAN Cementable Abutment, Ø 5mm, GH 3mm, AH 5.5mm, POM	025.4627-04	111
RC PLAN Cementable Abutment, Ø 6.5mm, GH 1mm, AH 4mm, POM	025.4631-04	112
RC PLAN Cementable Abutment, Ø 6.5mm, GH 1mm, AH 5.5mm, POM	025.4635-04	112
RC PLAN Cementable Abutment, Ø 6.5mm, GH 2mm, AH 4mm, POM	025.4632-04	112
RC PLAN Cementable Abutment, Ø 6.5mm, GH 2mm, AH 5.5mm, POM	025.4636-04	112
RC PLAN Cementable Abutment, Ø 6.5mm, GH 3mm, AH 4mm, POM	025.4633-04	112
RC PLAN Cementable Abutment, Ø 6.5mm, GH 3mm, AH 5.5mm, POM	025.4637-04	112
RC PLAN Gold Abutment for bridge, POM	025.4630-04	111
RC PLAN Gold Abutment for crown, POM	025.4610-04	111
RC PLAN LOCATOR® Abutment, H 1mm, POM	025.4641-04	112
RC PLAN LOCATOR® Abutment, H 2mm, POM	025.4642-04	112
RC PLAN LOCATOR® Abutment, H 3mm, POM	025.4643-04	112
RC PLAN LOCATOR® Abutment, H 4mm, POM	025.4644-04	112
RC PLAN LOCATOR® Abutment, H 5mm, POM	025.4645-04	112
RC PLAN Multi-Base Abutment, straight, Ø 4.5mm, GH 1mm, POM	025.4651-04	113
RC PLAN Multi-Base Abutment, straight, Ø 4.5mm, GH 2.5mm, POM	025.4652-04	113
RC PLAN Multi-Base Abutment, straight, Ø 4.5mm, GH 4mm, POM	025.4654-04	113
RC PLAN Multi-Base Abutment, straight, Ø 6.5mm, GH 1mm, POM	025.4656-04	113
RC PLAN Multi-Base Abutment, straight, Ø 6.5mm, GH 2.5mm, POM	025.4657-04	113
RC PLAN Multi-Base Abutment, straight, Ø 6.5mm, GH 4mm, POM	025.4659-04	113
RC PLAN Multi-Base Abutment, Type A, angled 25°, POM	025.4662-04	113
RC PLAN Multi-Base Abutment, Type B, angled 25°, POM	025.4665-04	113
RC PLAN NOVALOC® Abutment, H 1-6mm, POM	025.4646-04	112
RC Plan screw-retained abutment, POM, angled 17°, Ø 4.6mm, GH 2.5mm, Type A	025.4649-04	112
RC Plan screw-retained abutment, POM, angled 17°, Ø 4.6mm, GH 2.5mm, Type B	025.4650-04	112
RC Plan screw-retained abutment, POM, angled 30°, Ø 4.6mm, GH 2.5mm, Type A	025.4653-04	112
RC Plan screw-retained abutment, POM, angled 30°, Ø 4.6mm, GH 2.5mm, Type B	025.4660-04	112
RC Plan screw-retained abutment, POM, straight 0°, Ø 4.6mm, GH 2.5mm	025.4648-04	112
RC Polishing Aid, stainless steel	025.4920	106
RC Polishing Aid, stainless steel	025.4920-04	106
RC Protective Cap, for Cementable Abutment, Ø 5mm, AH 4mm, PEEK	024.4322-04	94
RC Protective Cap, for Cementable Abutment, Ø 5mm, AH 5.5mm, PEEK	024.4326-04	94
RC Protective Cap, for Cementable Abutment, Ø 6.5mm, AH 4mm, PEEK	024.4332-04	94
RC Protective Cap, for Cementable Abutment, Ø 6.5mm, AH 5.5mm, PEEK	024.4336-04	94
RC Protective Cap, for Multi-base Abutment, angled 25°, PEEK/TAN	024.4306-04	102
RC Protective Cap, for Multi-base Abutment, straight, Ø 4.5mm, PEEK/TAN	024.4341-04	102
RC Protective Cap, for Multi-base Abutment, straight, Ø 6.5mm, PEEK/TAN	024.4361-04	102
RC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 2.5mm, Type A, Pin	022.4747P	99
RC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 2.5mm, Type B, Pin	022.4748P	99
RC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 4mm, Type A, Pin	022.4752P	99
RC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 4mm, Type B, Pin	022.4753P	99
RC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 5.5mm, Type A, Pin	022.0014P	99
RC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 5.5mm, Type B, Pin	022.0015P	99
RC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 5.5mm, Type B, Pin	022.0017P	99
RC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 2.5mm, Type A, Pin	022.4749P	99
RC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 2.5mm, Type B, Pin	022.4750P	99
RC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 4mm, Type A, Pin	022.4754P	99
RC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 4mm, Type B, Pin	022.4755P	99
RC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 5.5mm, Type A, Pin	022.0016P	99
RC Screw-retained Abutment, TAN, straight 0°, Ø 4.6mm, GH 1mm, Pin	022.4745P	99
RC Screw-retained Abutment, TAN, straight 0°, Ø 4.6mm, GH 2.5mm, Pin	022.4746P	99
RC Screw-retained Abutment, TAN, straight 0°, Ø 4.6mm, GH 4mm, Pin	022.4751P	99
RC Temporary Abutment, Ø 7mm, H 12mm, PEEK/TAN	024.4370	92
RC Temporary Abutment, Bridge, Ø 4.5mm, H 11mm, TAN	024.4375	92
RC Temporary Abutment, Crown, Ø 4.5mm, H 11mm, TAN	024.4371	92
RC Temporary Coping, for Cementable Abutment, Ø 5mm, PMMA	024.4356	94
RC Temporary Coping, for Cementable Abutment, Ø 6.5mm, PMMA	024.4366	94
RC Temporary Coping, for Multi-Base Abutment, angled 25°, Ti	024.4307	102
RC Temporary Coping, for Multi-Base Abutment, straight, Ø 4.5mm, Ti	024.4342	102
RC Temporary Coping, for Multi-Base Abutment, straight, Ø 6.5mm, Ti	024.4362	102
RC Transfer Aid, for Cementable Abutment, Ø 5mm, PC	025.4912-04	94
RC Transfer Aid, for Cementable Abutment, Ø 6.5mm, PC	025.4913-04	94
RC Variobase® Abutment, with screw, Ø 4.5mm, H 3.5mm, TAN	025.4921	96,97

产品名称	货号	页数
RC Variobase®, for crown, with screw, Ø 4.5mm, H 5.5mm, TAN	022.0026	96,97
Reamer for 45° shoulder, Ø 8.0mm, Ti-N/stainless steel	046.243	145
Release aid N for Loxim®	026.2558	131
Release aid R/W for Loxim®	026.4558	131
Retentive anchor driver, L 19mm, stainless steel	046.069	144
RN Analog for Solid Abutment 048.540, L 14.3mm, yellow, Al	048.160	49
RN Analog for Solid Abutment 048.541, L 15.8mm, grey, Al	048.161	49
RN Analog for Solid Abutment 048.542, L 17.3mm, blue, Al	048.162	49
RN Burn-out coping, for Variobase® Abutment, POM	048.268	47
RN Burn-out coping, for Variobase® Abutment, POM	048.268V4	47
RN Burn-out Coping, for Variobase®, for crown, AH 6mm, POM	023.0015	47
RN Burn-out Coping, for Variobase®, for crown, AH 6mm, POM	023.0015V4	47
RN closure screw, large, height 1.5 mm, titanium	048.373	14,17,21
RN closure screw, large, height 1.5 mm, titanium	048.373V4	14,17,21
RN closure screw, small, titanium	048.371	14,17,21
RN closure screw, small, titanium	048.371V4	14,17,21
RN healing cap height 4.5 mm, titanium	048.037	14,17,21
RN healing cap, height 2.0 mm, titanium	048.033	14,17,21
RN healing cap, height 3.0 mm, titanium	048.034	14,17,21
RN Healing cap, labial bevel, large H 3.5mm, Ti	048.029	14,17,21
RN Healing cap, labial bevel, small H 2mm, Ti	048.028	14,17,21
RN Impression Cap, H 8mm, POM	048.017V4	44,49
RN Impression Post Closed Tray, abutment level, for RN synOcta® 1.5 Abutment, Al/Ti/POM	048.864	44
RN Impression Post Closed Tray, with 1 guide screw & 2 caps, Al/Ti/POM	048.862	44
RN Impression Post Open tray, abutment level, for RN synOcta® 1.5 Abutment, Al/Ti	048.866	44
RN LOCATOR® abutment, height 1.0 mm, titanium alloy, TiN coating	048.175	52
RN LOCATOR® abutment, height 2.0 mm, titanium alloy, TiN coating	048.176	52
RN LOCATOR® abutment, height 3.0 mm, titanium alloy, TiN coating	048.177	52
RN LOCATOR® abutment, height 4.0 mm, titanium alloy, TiN coating	048.178	52
RN LOCATOR® abutment, height 5.0 mm, titanium alloy, TiN coating	048.179	52
RN LOCATOR® abutment, height 6.0 mm, titanium alloy, TiN coating	048.180	52
RN Novaloc® Plan Abutment, H 1-6 mm, POM	048.280V4	67
RN PLAN LOCATOR® Abutment, H 1mm, POM	048.275V4	67
RN PLAN LOCATOR® Abutment, H 2mm, POM	048.276V4	67
RN PLAN LOCATOR® Abutment, H 3mm, POM	048.277V4	67
RN PLAN LOCATOR® Abutment, H 4mm, POM	048.278V4	67
RN PLAN LOCATOR® Abutment, H 5mm, POM	048.279V4	67
RN PLAN Solid Abutment, H 4mm, POM	048.926V4	67
RN PLAN Solid Abutment, H 5.5mm, POM	048.927V4	67
RN PLAN Solid Abutment, H 7mm, POM	048.928V4	67
RN Plastic coping for Solid Abutment, bridge, H 10mm, POM	048.246	50
RN Plastic coping for Solid Abutment, crown, H 10mm, POM	048.245	50
RN Positioning Cylinder for Solid Abutment, for 048.540, H 10.2mm, POM	048.060V4	49
RN Positioning Cylinder for Solid Abutment, for 048.541, H 10.2mm, POM	048.061V4	49
RN Positioning Cylinder for Solid Abutment, for 048.542, H 10.2mm, POM	048.062V4	49
RN Protective Cap for Solid Abutment 048.540, H 5.8mm, PEEK	048.047V4	49
RN Protective Cap for Solid Abutment 048.541, H 7.3mm, PEEK	048.048V4	49
RN Protective Cap for Solid Abutment 048.542, H 8.8mm, PEEK	048.049V4	49
RN Solid Abutment 6°, H 4mm, yellow, Ti	048.540	50
RN Solid Abutment 6°, H 5.5mm, grey, Ti	048.541	50
RN Solid Abutment 6°, H 7mm, blue, Ti	048.542	50
RN Straumann® Variobase® for CEREC®, TAN	022.0019	47
RN synOcta® 1.5 screw-retained Abutment, H 1.5mm, TAN	048.601	48,51
RN synOcta® 1.5 screw-retained PLAN Abutment, POM	048.929V4	67
RN synOcta® angled Abutment, long, Type A, angled 15°, H 6.7mm, Ti	048.610	46
RN synOcta® angled Abutment, long, Type A, angled 20°, H 6.7mm, Ti	048.615	46
RN synOcta® angled Abutment, long, Type B, angled 15°, H 6.7mm, Ti	048.611	46
RN synOcta® angled Abutment, long, Type B, angled 20°, H 6.7mm, Ti	048.616	46
RN synOcta® angled Abutment, short, Type A, angled 15°, H 5.7mm, Ti	048.612	46
RN synOcta® angled Abutment, short, Type A, angled 20°, H 5.7mm, Ti	048.617	46
RN synOcta® angled Abutment, short, Type B, angled 15°, H 5.7mm, Ti	048.613	46
RN synOcta® angled Abutment, short, Type B, angled 20°, H 5.7mm, Ti	048.618	46
RN synOcta® angled PLAN Abutment, short, Type A, angled 15°, POM	048.916V4	67
RN synOcta® angled PLAN Abutment, short, Type A, angled 20°, POM	048.918V4	67
RN synOcta® angled PLAN Abutment, short, Type B, angled 15°, POM	048.917V4	67
RN synOcta® angled PLAN Abutment, short, Type B, angled 20°, POM	048.919V4	67
RN synOcta® angled PLAN Abutment, Type A, angled 15°, POM	048.922V4	67

产品名称	货号	页数
RN synOcta® angled PLAN Abutment, Type A, angled 20°, POM	048.924V4	67
RN synOcta® angled PLAN Abutment, Type B, angled 15°, POM	048.923V4	67
RN synOcta® angled PLAN Abutment, Type B, angled 20°, POM	048.925V4	67
RN synOcta® Cementable Abutment, for crowns/bridges, H 5.5mm, Ti	048.605	45
RN synOcta® Cementable PLAN Abutment, POM	048.930V4	67
RN synOcta® gold Abutment, crown, H 4.3mm, Ceramicor®/POM	048.642	46
RN synOcta® gold coping, bar, for synOcta 1.5, H 5.8mm, Ceramicor®	048.204	51
RN synOcta® gold coping, bridge, for synOcta 1.5, H 4.25mm, Ceramicor®/POM	048.632	48
RN synOcta® gold coping, crown, for synOcta 1.5, H 4.25mm, Ceramicor®/POM	048.633	48
RN synOcta® impression cap, screw-retained, with integral guide screw, H 11mm, Al/Ti	048.010	44
RN synOcta® Impression post, with integral guide screw, H 21.0 mm, Al/Ti	048.090	44
RN synOcta® Meso Milling Cylinder, for synOcta® 1.5, crown, H 10mm, Ø 8mm, Ti	048.571	48
RN synOcta® Milling Cylinder, Bridge, for synOcta® 1.5, Ø 6mm, H 8mm, Ti	048.241	48
RN synOcta® Milling Cylinder, Crown, for synOcta® 1.5, Ø 6mm, H 8mm, Ti	048.240	48
RN synOcta® Plastic Coping for 048.605, for bridge, H 7mm, POM	048.662	45
RN synOcta® Plastic Coping for 048.605, for crown, H 7mm, POM	048.663	45
RN synOcta® Plastic Coping for synOcta® 1.5, for bridge/bar, H 10mm, POM	048.227	48,51
RN synOcta® Plastic Coping for synOcta® 1.5, for crown, H 10mm, POM	048.229	48
RN synOcta® Positioning Cylinder, H 12mm, POM	048.070V4	44
RN synOcta® Post for temporary restorations, for bridge, H 9mm, Ti	048.650	45
RN synOcta® Post for temporary restorations, for crown, H 9mm, Ti	048.651	45
RN synOcta® Temporary Meso Abutment, H 10mm, PEEK/TAN	048.668	45
RN synOcta® Ti Coping, Bar for synOcta® 1.5, H 6mm, Ti	048.214	51
RN synOcta® Transversal PLAN Abutment, POM	048.921V4	67
RN Temporary Coping for Solid Abutment, bridge, H 8.5mm, PMMA	048.654	49
RN Temporary Coping for Solid Abutment, crown, H 8.5mm, PMMA	048.655	49
RN Variobase® Abutment, with screw, Ø 5.05mm, H 4mm, TAN	048.710	47
RN Variobase®, for crown, with screw, Ø 5.05mm, H 6mm, TAN	022.0022	47
Round bur, Ø 1.4mm, stainless steel	044.022	127
Round bur, Ø 2.3mm, stainless steel	044.003	127
Round bur, Ø 3.1mm, stainless steel	044.004	127
S/SP Tap f. handpiece, Ø 3.3mm, guided, L 42mm, stainless steel	034.245	134
S/SP Tap f. handpiece, Ø 4.1mm, guided, L 42mm, stainless steel	034.445	134
S/SP Tap f. handpiece, Ø 4.8mm, guided, L 42mm, stainless steel	034.645	134
S/SP Tap for adapter, Ø 3.3mm, L 23mm, PEEK/stainless steel	044.575	129
S/SP Tap for adapter, Ø 4.1mm, L 23mm, PEEK/stainless steel	044.577	129
S/SP Tap for adapter, Ø 4.8mm, L 23mm, PEEK/stainless steel	044.579	129
Screwdriver for titanium matrix, L 50mm, stainless steel	048.452	144
Screwdriver handle for 046.410/411/412, L 109mm, stainless steel	046.085	144
Screwdriver, L 37mm, stainless steel	046.154	144
SCS Guide screw for RN synOcta® Bar coping, cannot be shortened, L 6 mm, Ti	048.360V4	150
SCS Guiding cylinder for NN implant Ø 3.3mm, L 5mm, Ti/stainless steel	049.182	138
SCS Occlusal screw for guiding sleeve, and post for temp. restorations, L 7.6 mm, TAN	049.181	45,57,138,150
SCS Occlusal screw, L 4.4mm, Ti	048.350	46,48,51,60,138,141,142,150
SCS Occlusal screw, L 4.4mm, Ti	048.350V4	46,48,51,60,141,142,150
SCS Screwdriver for handpiece, long, L 32mm, stainless steel	046.412	132,144
SCS Screwdriver for handpiece, short, L 26mm, stainless steel	046.411	132,144
SCS Screwdriver for handpiece, x-short, L 20mm, stainless steel	046.410	132,144
SCS Screwdriver, for ratchet, long, L 27mm, stainless steel	046.402	62,132,144
SCS Screwdriver, for ratchet, short, L 21mm, stainless steel	046.401	62,132,144
SCS Screwdriver, for ratchet, x-short, L 15mm, stainless steel	046.400	62,132,144
Service instrument for ratchet, L 50mm, stainless steel	046.108	131,144
Shoulder Analog for WN Solid Abutment, L 10mm, Al	048.140	61
Solid abut.driver, long, f. RN solid abut.6°, L 19mm, stainless steel	046.068	50,144
Solid abut.driver, short, f. RN solid abut.6°, L 13mm, stainless steel	046.067	50,144
SP Profile drill, Ø 3.3mm, RN, guided, L 35mm, stainless steel	034.235	134
SP Profile drill, Ø 4.1mm, RN, guided, L 35mm, stainless steel	034.435	134
SP Profile drill, RN, long, Ø 3.3mm, L 34mm, stainless steel	044.087	128
SP Profile drill, RN, long, Ø 4.1mm, L 33mm, stainless steel	044.089	128
SP Profile drill, RN, short, Ø 3.3mm, L 25mm, stainless steel	044.086	128
SP Profile drill, RN, short, Ø 4.1mm, L 24mm, stainless steel	044.088	128
SP Profile drill, WN, long, Ø 4.8mm, L 32.8mm, stainless steel	044.085	128
SP Profile drill, WN, short, Ø 4.8mm, L 23.8mm, stainless steel	044.084	128
Standard Implant Ø 4.1mm RN, SLActive® 8mm, Ti guided	033.031G	16
Standard Implant Ø 4.1mm, RN, SLActive® 10mm, Ti	033.032S	16
Standard Implant Ø 4.1mm, RN, SLActive® 10mm, Ti guided	033.032G	16
Standard Implant Ø 4.1mm, RN, SLActive® 12mm, Ti	033.033S	16

产品名称	货号	页数
Standard Implant Ø 4.1mm, RN, SLActive® 12mm, Ti guided	033.033G	16
Standard Implant Ø 4.1mm, RN, SLActive® 14mm, Ti	033.034S	16
Standard Implant Ø 4.1mm, RN, SLActive® 16mm, Ti	033.035S	16
Standard Implant Ø 4.1mm, RN, SLActive® 6mm, Ti	033.030S	16
Standard Implant Ø 4.1mm, RN, SLActive® 8mm, Ti	033.031S	16
Standard Implant Ø 4.8mm RN, SLActive® 10mm, Ti	033.232S	19
Standard Implant Ø 4.8mm RN, SLActive® 12mm, Ti	033.233S	19
Standard Implant Ø 4.8mm RN, SLActive® 6mm, Ti	033.230S	19
Standard Implant Ø 4.8mm RN, SLActive® 8mm, Ti	033.231S	19
Standard Implant Ø 4.8mm, RN, SLActive® 14mm, Ti	033.234S	19
Standard Implant Ø 4.8mm, WN, SLActive® 10mm, Ti	033.632S	19
Standard Implant Ø 4.8mm, WN, SLActive® 12mm, Ti	033.633S	19
Standard Implant Ø 4.8mm, WN, SLActive® 6mm, Ti	033.630S	19
Standard Implant Ø 4.8mm, WN, SLActive® 8mm, Ti	033.631S	19
Standard Implant Ø 3.3mm RN, SLActive® 10mm, Roxolid® guided	033.432G	12
Standard Implant Ø 3.3mm RN, SLActive® 12mm, Roxolid® guided	033.433G	12
Standard Implant Ø 3.3mm RN, SLActive® 8mm, Roxolid® guided	033.431G	12
Standard Implant, Ø 3.3 mm, RN, SLActive 10 mm, Roxolid	033.432S	12
Standard Implant, Ø 3.3 mm, RN, SLActive 12 mm, Roxolid	033.433S	12
Standard Implant, Ø 3.3 mm, RN, SLActive 14 mm, Roxolid	033.434S	12
Standard Implant, Ø 3.3 mm, RN, SLActive 16 mm, Roxolid	033.435S	12
Standard Implant, Ø 3.3 mm, RN, SLActive 8 mm, Roxolid	033.431S	12
Standard Implant, Ø 3.3mm RN SLActive® 10mm, Ti	033.132S	13
Standard Implant, Ø 3.3mm RN SLActive® 12mm, Ti	033.133S	13
Standard Implant, Ø 3.3mm RN SLActive® 14mm, Ti	033.134S	13
Standard Implant, Ø 3.3mm RN SLActive® 16mm, Ti	033.135S	13
Standard Implant, Ø 3.3mm RN SLActive® 8mm, Ti	033.131S	13
Standard Implant, Ø 3.3mm, RN, SLA 10mm, Ti	043.132S	13
Standard Implant, Ø 3.3mm, RN, SLA 12mm, Ti	043.133S	13
Standard Implant, Ø 3.3mm, RN, SLA 14mm, Ti	043.134S	13
Standard Implant, Ø 3.3mm, RN, SLA 16mm, Ti	043.135S	13
Standard Implant, Ø 3.3mm, RN, SLA 8mm, Ti	043.131S	13
Standard Implant, Ø 4.1mm, RN, SLA 8mm, Ti	043.031S	16
Standard Implant, Ø 4.1mm, RN, SLA 10mm, Ti	043.032S	16
Standard Implant, Ø 4.1mm, RN, SLA 14mm, Ti	043.034S	16
Standard Implant, Ø 4.1mm, RN, SLA 16mm, Ti	043.035S	16
Standard Implant, Ø 4.1mm, RN, SLA 6mm, Ti	043.030S	16
Standard Implant, Ø 4.8mm, RN, SLA 10mm, Ti	043.232S	20
Standard Implant, Ø 4.8mm, RN, SLA 12mm, Ti	043.233S	20
Standard Implant, Ø 4.8mm, RN, SLA 14mm, Ti	043.234S	20
Standard Implant, Ø 4.8mm, RN, SLA 6mm, Ti	043.230S	20
Standard Implant, Ø 4.8mm, RN, SLA 8mm, Ti	043.231S	20
Standard Implant, Ø 4.8mm, WN, SLA 10mm, Ti	043.632S	20
Standard Implant, Ø 4.8mm, WN, SLA 12mm, Ti	043.633S	20
Standard Implant, Ø 4.8mm, WN, SLA 6mm, Ti	043.630S	20
Standard Implant, Ø 4.8mm, WN, SLA 8mm, Ti	043.631S	20
Standard Implant, Ø 3.3mm RN, SLActive® 10mm, Roxolid®, Loxim®	033.502S	12
Standard Implant, Ø 3.3mm RN, SLActive® 12mm, Roxolid®, Loxim®	033.503S	12
Standard Implant, Ø 3.3mm RN, SLActive® 14mm, Roxolid®, Loxim®	033.504S	12
Standard Implant, Ø 3.3mm RN, SLActive® 16mm, Roxolid®, Loxim®	033.505S	12
Standard Implant, Ø 3.3mm RN, SLActive® 8mm, Roxolid®, Loxim®	033.501S	12
Standard Implant, Ø 4.1 RN, SLActive® 10mm, Roxolid®, guided	033.532G	15
Standard Implant, Ø 4.1 RN, SLActive® 12mm, Roxolid®, guided	033.533G	15
Standard Implant, Ø 4.1 RN, SLActive® 8mm, Roxolid®, guided	033.531G	15
Standard Implant, Ø 4.1mm RN, SLActive® 10mm, Roxolid®, Loxim®	033.532S	15
Standard Implant, Ø 4.1mm RN, SLActive® 12mm, Roxolid®, Loxim®	033.533S	15
Standard Implant, Ø 4.1mm RN, SLActive® 14mm, Roxolid®, Loxim®	033.534S	15
Standard Implant, Ø 4.1mm RN, SLActive® 16mm, Roxolid®, Loxim®	033.535S	15
Standard Implant, Ø 4.1mm RN, SLActive® 6mm, Roxolid®, Loxim®	033.530S	15
Standard Implant, Ø 4.1mm RN, SLActive® 8mm, Roxolid®, Loxim®	033.531S	15
Standard Implant, Ø 4.8mm RN, SLActive® 10mm, Roxolid®, Loxim®	033.582S	18
Standard Implant, Ø 4.8mm RN, SLActive® 12mm, Roxolid®, Loxim®	033.583S	18
Standard Implant, Ø 4.8mm RN, SLActive® 14mm, Roxolid®, Loxim®	033.584S	18
Standard Implant, Ø 4.8mm RN, SLActive® 6mm, Roxolid®, Loxim®	033.580S	18
Standard Implant, Ø 4.8mm RN, SLActive® 8mm, Roxolid®, Loxim®	033.581S	18
Standard Implant, Ø 4.8mm WN, SLActive® 10mm, Roxolid®, Loxim®	033.602S	18

产品名称	货号	页数
Standard Implant, Ø 4.8mm WN, SLActive® 12mm, Roxolid®, Loxim®	033.603S	18
Standard Implant, Ø 4.8mm WN, SLActive® 6mm, Roxolid®, Loxim®	033.600S	18
Standard Implant, Ø 4.8mm WN, SLActive® 8mm, Roxolid®, Loxim®	033.601S	18
Standard Implant; Ø 4.1mm, RN, SLA 12mm, Ti	043.033S	16
Standard Plus Implant Ø 4.1mm RN SLActive® 10mm, Ti	033.052S	16
Standard Plus Implant Ø 4.1mm RN, SLActive® 10mm, Ti guided	033.052G	16
Standard Plus Implant Ø 4.1mm RN, SLActive® 12mm, Ti guided	033.053G	16
Standard Plus Implant Ø 4.1mm RN, SLActive® 8mm, Ti guided	033.051G	16
Standard Plus Implant Ø 4.1mm, RN, SLActive® 6mm, Ti	033.050S	16
Standard Plus Implant Ø 4.1mm, RN, SLActive® 8mm, Ti	033.051S	16
Standard Plus Implant Ø 4.8mm RN, SLActive® 10mm, Ti guided	033.252G	19
Standard Plus Implant Ø 4.8mm RN, SLActive® 12mm, Ti guided	033.253G	19
Standard Plus Implant Ø 4.8mm RN, SLActive® 8mm, Ti guided	033.251G	19
Standard Plus Implant Ø 4.8mm, RN, SLActive® 10mm, Ti	033.252S	19
Standard Plus Implant Ø 4.8mm, RN, SLActive® 12mm, Ti	033.253S	19
Standard Plus Implant Ø 4.8mm, RN, SLActive® 14mm, Ti	033.254S	19
Standard Plus Implant Ø 4.8mm, RN, SLActive® 6mm, Ti	033.250S	19
Standard Plus Implant Ø 4.8mm, RN, SLActive® 8mm, Ti	033.251S	19
Standard Plus Implant Ø 4.8mm, WN, SLActive® 10mm, Ti	033.652S	19
Standard Plus Implant Ø 4.8mm, WN, SLActive® 12mm, Ti	033.653S	19
Standard Plus Implant Ø 4.8mm, WN, SLActive® 6mm, Ti	033.650S	19
Standard Plus Implant Ø 4.8mm, WN, SLActive® 8mm, Ti	033.651S	19
Standard Plus Implant Ø 3.3mm RN, SLActive® 10mm, Roxolid® guided	033.452G	12
Standard Plus Implant Ø 3.3mm RN, SLActive® 12mm, Roxolid® guided	033.453G	12
Standard Plus Implant Ø 3.3mm RN, SLActive® 8mm, Roxolid® guided	033.451G	12
Standard Plus Implant, Ø 3.3mm RN SLActive® 10mm, Ti	033.152S	13
Standard Plus Implant, Ø 3.3mm RN SLActive® 12mm, Ti	033.153S	13
Standard Plus Implant, Ø 3.3mm RN SLActive® 14mm, Ti	033.154S	13
Standard Plus Implant, Ø 3.3mm RN SLActive® 8mm, Ti	033.151S	13
Standard Plus Implant, Ø 3.3mm, RN, SLA 10mm, Ti	043.152S	13
Standard Plus Implant, Ø 3.3mm, RN, SLA 12mm, Ti	043.153S	13
Standard Plus Implant, Ø 3.3mm, RN, SLA 8mm, Ti	043.151S	13
Standard Plus Implant, Ø 4.1mm RN SLActive® 12mm, Ti	033.053S	16
Standard Plus Implant, Ø 4.1mm RN SLActive® 14mm, Ti	033.054S	16
Standard Plus Implant, Ø 4.1mm, RN, SLA 10mm, Ti	043.052S	16
Standard Plus Implant, Ø 4.1mm, RN, SLA 12mm, Ti	043.053S	16
Standard Plus Implant, Ø 4.1mm, RN, SLA 14mm, Ti	043.054S	16
Standard Plus Implant, Ø 4.1mm, RN, SLA 6mm, Ti	043.050S	16
Standard Plus Implant, Ø 4.1mm, RN, SLA 8mm, Ti	043.051S	16
Standard Plus Implant, Ø 4.8mm, RN, SLA 10mm, Ti	043.252S	20
Standard Plus Implant, Ø 4.8mm, RN, SLA 12mm, Ti	043.253S	20
Standard Plus Implant, Ø 4.8mm, RN, SLA 6mm, Ti	043.250S	20
Standard Plus Implant, Ø 4.8mm, RN, SLA 8mm, Ti	043.251S	20
Standard Plus Implant, Ø 4.8mm, WN, SLA 10mm, Ti	043.652S	20
Standard Plus Implant, Ø 4.8mm, WN, SLA 12mm, Ti	043.653S	20
Standard Plus Implant, Ø 4.8mm, WN, SLA 6mm, Ti	043.650S	20
Standard Plus Implant, Ø 4.8mm, WN, SLA 8mm, Ti	043.651S	20
Standard Plus Implant, Ø 3.3mm NNC, SLActive® 10mm, Roxolid®	033.417S	12
Standard Plus Implant, Ø 3.3mm NNC, SLActive® 12mm, Roxolid®	033.418S	12
Standard Plus Implant, Ø 3.3mm NNC, SLActive® 14mm, Roxolid®	033.419S	12
Standard Plus Implant, Ø 3.3mm NNC, SLActive® 8mm, Roxolid®	033.416S	12
Standard Plus Implant, Ø 3.3mm RN, SLActive® 10mm, Roxolid®, Loxim®	033.512S	12
Standard Plus Implant, Ø 3.3mm RN, SLActive® 12mm, Roxolid®, Loxim®	033.513S	12
Standard Plus Implant, Ø 3.3mm RN, SLActive® 14mm, Roxolid®, Loxim®	033.514S	12
Standard Plus Implant, Ø 3.3mm RN, SLActive® 8mm, Roxolid®, Loxim®	033.511S	12
Standard Plus Implant, Ø 4.1 RN, SLActive® 10mm, Roxolid®, guided	033.562G	15
Standard Plus Implant, Ø 4.1 RN, SLActive® 12mm, Roxolid®, guided	033.563G	15
Standard Plus Implant, Ø 4.1 RN, SLActive® 8mm, Roxolid®, guided	033.561G	15
Standard Plus Implant, Ø 4.1mm RN, SLActive® 10mm, Roxolid®, Loxim®	033.562S	15
Standard Plus Implant, Ø 4.1mm RN, SLActive® 12mm, Roxolid®, Loxim®	033.563S	15
Standard Plus Implant, Ø 4.1mm RN, SLActive® 14mm, Roxolid®, Loxim®	033.564S	15
Standard Plus Implant, Ø 4.1mm RN, SLActive® 6mm, Roxolid®, Loxim®	033.560S	15
Standard Plus Implant, Ø 4.1mm RN, SLActive® 8mm, Roxolid®, Loxim®	033.561S	15
Standard Plus Implant, Ø 4.8 RN, SLActive® 10mm, Roxolid®, guided	033.592G	18
Standard Plus Implant, Ø 4.8 RN, SLActive® 12mm, Roxolid®, guided	033.593G	18
Standard Plus Implant, Ø 4.8 RN, SLActive® 8mm, Roxolid®, guided	033.591G	18
Standard Plus Implant, Ø 4.8mm RN, SLActive® 10mm, Roxolid®, Loxim®	033.592S	18

产品名称	货号	页数
Standard Plus Implant, Ø 4.8mm RN, SLActive® 12mm, Roxolid®, Loxim®	033.593S	18
Standard Plus Implant, Ø 4.8mm RN, SLActive® 14mm, Roxolid®, Loxim®	033.594S	18
Standard Plus Implant, Ø 4.8mm RN, SLActive® 6mm, Roxolid®, Loxim®	033.590S	18
Standard Plus Implant, Ø 4.8mm RN, SLActive® 8mm, Roxolid®, Loxim®	033.591S	18
Standard Plus Implant, Ø 4.8mm WN, SLActive® 10mm, Roxolid®, Loxim®	033.612S	18
Standard Plus Implant, Ø 4.8mm WN, SLActive® 12mm, Roxolid®, Loxim®	033.613S	18
Standard Plus Implant, Ø 4.8mm WN, SLActive® 6mm, Roxolid®, Loxim®	033.610S	18
Standard Plus Implant, Ø 4.8mm WN, SLActive® 8mm, Roxolid®, Loxim®	033.611S	18
Standard Plus Implant,Ø 3.3mm, RN, SLA 14mm, Ti	043.154S	13
Standard Plus Implant,Ø 4.8mm, RN, SLA 14mm, Ti	043.254S	20
Stepped pin for 049.810, Ø 2.2mm/3.5mm, L 16mm, Ti	049.818V4	149
Stop key for guided implants, S/SP/TE/BL, stainless steel	034.006	136
Storage System for implants (SLA & SLActive), 280 x 226 x 248 mm, plastic	041.500	143
Storage System for Prosthetics, 280 x 226 x 248mm, plastic	041.501	147
Straumann® Basic Surgical Cassette	040.564	131
Straumann® Basic Surgical Cassette, guided	034.281	136
Straumann® Planning Guide	026.0016	77,98
Straumann® surgical cassette	040.165	131
Surgical cassette, for guided surgery	034.001	136
synOcta® Basal screw, L 6,7 mm, TAN	048.356	45,46,47,57,58,59,150
Tap for transversal screw retention, M2 x 0.25mm, stainless steel	044.570	146
Tapered Effect Implant Ø 3.3mm RN, SLActive® 10mm, Ti	033.722S	13
Tapered Effect Implant Ø 3.3mm RN, SLActive® 12mm, Ti	033.723S	13
Tapered Effect Implant Ø 3.3mm RN, SLActive® 14mm, Ti	033.724S	13
Tapered Effect Implant Ø 3.3mm RN, SLActive® 8mm, Ti	033.721S	13
Tapered Effect Implant Ø 4.1mm RN, SLActive® 10mm, Ti	033.762S	16
Tapered Effect Implant Ø 4.1mm RN, SLActive® 10mm, Ti guided	033.762G	16
Tapered Effect Implant Ø 4.1mm RN, SLActive® 12mm, Ti	033.763S	16
Tapered Effect Implant Ø 4.1mm RN, SLActive® 12mm, Ti guided	033.763G	16
Tapered Effect Implant Ø 4.1mm RN, SLActive® 14mm, Ti	033.764S	16
Tapered Effect Implant Ø 4.1mm RN, SLActive® 8mm, Ti	033.761S	16
Tapered Effect Implant Ø 4.1mm RN, SLActive® 8mm, Ti guided	033.761G	16
Tapered Effect Implant Ø 4.8mm WN, SLActive® 10mm, Ti	033.712S	19
Tapered Effect Implant Ø 4.8mm WN, SLActive® 12mm, Ti	033.713S	19
Tapered Effect Implant Ø 4.8mm WN, SLActive® 14mm, Ti	033.714S	19
Tapered Effect Implant, Ø 3.3mm, RN, SLA 10mm, Ti	043.722S	13
Tapered Effect Implant, Ø 3.3mm, RN, SLA 12mm, Ti	043.723S	13
Tapered Effect Implant, Ø 3.3mm, RN, SLA 14mm, Ti	043.724S	13
Tapered Effect Implant, Ø 3.3mm, RN, SLA 8mm, Ti	043.721S	13
Tapered Effect Implant, Ø 4.1mm, RN, SLA 10mm, Ti	043.762S	16
Tapered Effect Implant, Ø 4.1mm, RN, SLA 12mm, Ti	043.763S	16
Tapered Effect Implant, Ø 4.1mm, RN, SLA 14mm, Ti	043.764S	16
Tapered Effect Implant, Ø 4.1mm, RN, SLA 8mm, Ti	043.761S	16
Tapered Effect Implant, Ø 4.8mm, WN, SLA 10mm, Ti	043.712S	20
Tapered Effect Implant, Ø 4.8mm, WN, SLA 12mm, Ti	043.713S	20
Tapered Effect Implant, Ø 4.8mm, WN, SLA 14mm, Ti	043.714S	20
Tapered Effect Implant, Ø 3.3mm RN, SLActive® 10mm, Roxolid®, Loxim®	033.522S	12
Tapered Effect Implant, Ø 3.3mm RN, SLActive® 12mm, Roxolid®, Loxim®	033.523S	12
Tapered Effect Implant, Ø 3.3mm RN, SLActive® 14mm, Roxolid®, Loxim®	033.524S	12
Tapered Effect Implant, Ø 3.3mm RN, SLActive® 8mm, Roxolid®, Loxim®	033.521S	12
Tapered Effect Implant, Ø 4.1 RN, SLActive® 10mm, Roxolid®, guided	033.572G	15
Tapered Effect Implant, Ø 4.1 RN, SLActive® 12mm, Roxolid®, guided	033.573G	15
Tapered Effect Implant, Ø 4.1 RN, SLActive® 8mm, Roxolid®, guided	033.571G	15
Tapered Effect Implant, Ø 4.1mm RN, SLActive® 10mm, Roxolid®, Loxim®	033.572S	15
Tapered Effect Implant, Ø 4.1mm RN, SLActive® 12mm, Roxolid®, Loxim®	033.573S	15
Tapered Effect Implant, Ø 4.1mm RN, SLActive® 14mm, Roxolid®, Loxim®	033.574S	15
Tapered Effect Implant, Ø 4.1mm RN, SLActive® 8mm, Roxolid®, Loxim®	033.571S	15
Tapered Effect Implant, Ø 4.8mm WN, SLActive® 10mm, Roxolid®, Loxim®	033.622S	18
Tapered Effect Implant, Ø 4.8mm WN, SLActive® 12mm, Roxolid®, Loxim®	033.623S	18
Tapered Effect Implant, Ø 4.8mm WN, SLActive® 14mm, Roxolid®, Loxim®	033.624S	18
TE Profile drill, Ø 3.3mm, RN, guided, L 37mm, stainless steel	034.237	134
TE Profile drill, Ø 4.1mm, RN, guided, L 37mm, stainless steel	034.437	134
TE Profile drill, RN, long, Ø 3.3mm, L 37mm, stainless steel	044.708	128
TE Profile drill, RN, long, Ø 4.1mm, L 37mm, stainless steel	044.712	128
TE Profile drill, RN, short, Ø 3.3mm, L 28mm, stainless steel	044.701	128
TE Profile drill, RN, short, Ø 4.1mm, L 28mm, stainless steel	044.705	128
TE Profile drill, WN, long, Ø 4.8mm, L 37mm, stainless steel	044.710	128

产品名称	货号	页数
TE Profile drill, WN, short, Ø 4.8mm, L 28mm, stainless steel	044.703	128
Template fixation pin, Ø 2.8mm/2.8mm, guided, L 25mm, stainless steel	034.256	136
Template fixation pin, Ø 5.0mm/2.8mm, guided, L 24mm, stainless steel	034.255	136
Template fixation pin, Ø 5.0mm/3.5mm, guided, L 24mm, stainless steel	034.455	136
Template fixation pin, Ø 5mm/4.2mm, guided, L 24mm, stainless steel	034.655	136
Tissue Level PLAN cassette, empty, 278x124x33mm, polymer	048.905	66
Tissue Level PLAN Set, 278x124x33 mm, PC	048.904	66
Torque control device for ratchet, stainless steel	046.049	131,144
Transfer aid for 048.605, H 6,5 mm, PC	048.059V4	45
Transfer aid for 048.606, H 6,5 mm, PC	048.054V4	57
Transfer aid for RN synOcta®, angled, long, H 4 mm, PC	048.002V4	46
Transfer aid for RN synOcta®, angled, short, H 5.5 mm, PC	048.000V4	46
Transfer aid for WN synOcta®, angled, H 5 mm, PC	048.032	58
Transfer and Alignment Pin, for Screw-retained Abutments, TAN	025.0009	77,98
Tray for cassette 041.022, 218x109x16 mm, Al	041.211	148
T-sleeve, Ø 2.8mm, H 6mm, guided, stainless steel	034.055V4	133
T-sleeve, Ø 5mm, H 5mm, guided, stainless steel	034.053V4	133
Tweezers with titanium plasma-coated tips, L 105mm, Ti/stainless steel	046.095	131
Tweezers without titanium plasma-coated tips, L 105mm, stainless steel	046.057	131
Twist drill PRO, long, Ø 3.5mm, L 41mm, stainless steel	044.251	127
Twist drill PRO, long, Ø 4.2mm, L 41mm, stainless steel	044.255	127
Twist drill PRO, Ø 2.8mm, long, 20mm, guided, L 36mm, stainless steel	034.226	133
Twist drill PRO, Ø 2.8mm, short, 16mm, guided, L 32mm, stainless steel	034.223	133
Twist drill PRO, Ø 2.8mm, x-long, 24mm, guided, L 40mm, stainless steel	034.229	133
Twist drill PRO, Ø 3.5mm, long, 20mm, guided, L 36mm, stainless steel	034.426	133
Twist drill PRO, Ø 3.5mm, short, 16mm, guided, L 32mm, stainless steel	034.423	133
Twist drill PRO, Ø 3.5mm,x-long, 24mm, guided, L 40mm, stainless steel	034.429	133
Twist drill PRO, Ø 4.2mm, long, 20mm, guided, L 36mm, stainless steel	034.626	133
Twist drill PRO, Ø 4.2mm, short, 16mm, guided, L 32mm, stainless steel	034.623	133
Twist drill PRO, Ø 4.2mm,x-long, 24mm, guided, L 40mm, stainless steel	034.629	133
Twist drill PRO, short, Ø 3.5mm, L 33mm, stainless steel	044.250	127
Twist drill PRO, short, Ø 4.2mm, L 33mm, stainless steel	044.254	127
Ultrasonic Cleaning Cassette, 150 x 65 x 25 mm, plastic/silicone	040.175	143
WN Analog for Solid Abutment 048.545, L 16mm, green, Al	048.165	61
WN Analog for Solid Abutment 048.546, L 17.5mm, brown, Al	048.166	61
WN Burn-out coping, for Variobase® Abutment, POM	048.269	59
WN Burn-out Coping, for Variobase®, for crown, AH 6.5mm, POM	023.0016	59
WN Burn-out Coping, for Variobase®, for crown, AH 6.5mm, POM	023.0016V4	59
WN closure screw, small, titanium	048.375	21
WN Healing cap, H 2mm, Ti	048.038	21
WN Healing cap, H 3mm, Ti	048.039	21
WN healing cap, height 4.5 mm, titanium	048.053	21
WN healing cap, with labial bevel, height 2.0 mm, titanium	048.030	21
WN Implant Analog for synOcta® Abutment, abutment level, straight 0°, Ø6.5mm, stainless steel	048.872	56
WN Impression Cap, H 8mm, POM	048.013	56,61
WN Impression post closed tray, abutment level, for WN synOcta® 1.5 Abutment, Al/Ti/POM	048.871	56
WN Impression Post Closed Tray, with 1 guide screw & 2 caps, Al/Ti/POM	048.870	56
WN Impression Post Open Tray, abutment level, for WN synOcta® 1.5 Abutment, Al/Ti	048.869	56
WN LOCATOR®, abutment, tissue cuff height 4.0 mm, titanium alloy, TiNcoating	048.186	63
WN LOCATOR®, abutment, tissue cuff height 5.0 mm, titanium alloy, TiNcoating	048.187	63
WN LOCATOR®, solid abutment, tissue cuff height 1.0 mm, titanium alloy, TiN coating	048.183	63
WN LOCATOR®, solid abutment, tissue cuff height 3.0 mm, titanium alloy, TiN coating	048.185	63
WN LOCATOR®, solid abutment, tissue cuff, height 2.0mm, titanium alloy, TiN coating	048.184	63
WN PLAN LOCATOR® Abutment, H 1mm, POM	048.283V4	68
WN PLAN LOCATOR® Abutment, H 2mm, POM	048.284V4	68
WN PLAN LOCATOR® Abutment, H 3mm, POM	048.285V4	68
WN PLAN LOCATOR® Abutment, H 4mm, POM	048.286V4	68
WN PLAN LOCATOR® Abutment, H 5mm, POM	048.287V4	68
WN PLAN Solid Abutment 6°, H 4mm, green, POM	048.938V4	68
WN PLAN Solid Abutment 6°, H 5.5mm, brown, POM	048.939V4	68
WN Plastic coping for Solid Abutment, bridge, H 7mm, POM	048.248	62
WN Plastic coping for Solid Abutment, crown, H 7mm, POM	048.247	62
WN Positioning Cylinder for Solid Abutment, for 048.545, H 10mm, POM	048.065	61
WN Positioning Cylinder for Solid Abutment, for 048.546, H 10mm, POM	048.066	61
WN Protective Cap for Solid Abutment 048.545, H 6mm, PEEK	048.051	61
WN Protective Cap for Solid Abutment 048.546, H 7.5mm, PEEK	048.052	61
WN Solid Abutment 6°, H 4mm, green, Ti	048.545	62

产品名称	货号	页数
WN Solid Abutment 6°, H 5.5mm, brown, Ti	048.546	62
WN Straumann® Variobase® for CEREC®, TAN	022.0020	58
WN synOcta® 1.5 screw-retained Abutment, H 1.5mm, TAN	048.603	60
WN synOcta® 1.5 screw-retained PLAN Abutment, POM	048.931V4	68
WN synOcta® Analog, L 12mm, stainless steel	048.171	56
WN synOcta® angled Abutment, Type A, angled 15°, H 5.5mm, Ti	048.608	58
WN synOcta® angled Abutment, Type B, angled 15°, H 5.5mm, Ti	048.609	58
WN synOcta® angled PLAN Abutment, Type A, angled 15°, POM	048.933V4	68
WN synOcta® angled PLAN Abutment, Type B, angled 15°, POM	048.934V4	68
WN synOcta® Cementable Abutment, H 5.5mm, Ti	048.606	57
WN synOcta® Cementable PLAN Abutment, POM	048.932V4	68
WN synOcta® gold Abutment, crown, H 14.1mm, Ceramicor®/POM	048.644	59
WN synOcta® gold coping, bridge, for synOcta 1.5, H 4.75mm, Ceramicor®/POM	048.638	60
WN synOcta® gold coping, crown, for synOcta 1.5, H 4.75mm, Ceramicor®/POM	048.639	60
WN synOcta® impression cap, screw-retained, with integral guide screw, H 11mm, Al	048.091	56
WN synOcta® Meso Milling Cylinder, for synOcta® 1.5, crown, H 10mm, Ø 10mm, Ti	048.573	60
WN synOcta® Plastic Coping for 048.606, for bridge, H 7.35mm, POM	048.244	57
WN synOcta® Plastic Coping for 048.606, for crown, H 7.35mm, POM	048.243	57
WN synOcta® Plastic Coping, Bridge, H 10mm, POM	048.666	60
WN synOcta® Plastic Coping, Crown, H 10mm, POM	048.667	60
WN synOcta® Positioning Cylinder, H 12mm, POM	048.095	56
WN synOcta® Post for temp. restoration, bridge, H 9mm, Ti	048.234	57
WN synOcta® Post, for temporary restoration, crown, H 9mm, Ti	048.233	57
WN Temporary Coping for Solid Abutment, bridge, H 7.3mm, PMMA	048.656	61
WN Temporary Coping for Solid Abutment, crown, H 7.3mm, PMMA	048.657	61
WN Variobase® Abutment, with screw, Ø 7mm, H 4.5mm, TAN	048.711	58,59
WN Variobase®, for crown, with screw, Ø 7mm, H 6.5mm, TAN	022.0023	58,59
X-ray reference sphere, Ø 5.0mm, stainless steel	049.076V4	126,149
x-ray template for Bone Level Implants	150.216	126
x-ray template for Standard and Standard Plus Implants	150.215	126

货号索引

货号	产品名称	页数
021.2108	Bone Level Implant, Ø 3.3 mm, NC, SLActive 8 mm, Ti	24
021.2110	Bone Level Implant, Ø 3.3 mm, NC, SLActive 10 mm, Ti	24
021.2112	Bone Level Implant, Ø 3.3 mm, NC, SLActive 12 mm, Ti	24
021.2114	Bone Level Implant, Ø 3.3 mm, NC, SLActive 14 mm, Ti	24
021.2208G	Bone Level Implant, Ø 3.3 mm, NC, SLActive 8 mm, Roxolid, guided	24
021.2210G	Bone Level Implant, Ø 3.3 mm, NC, SLActive 10 mm, Roxolid, guided	24
021.2212G	Bone Level Implant, Ø 3.3 mm, NC, SLActive 12 mm, Roxolid, guided	24
021.2308	Bone Level Implant Ø 3.3mm NC, SLActive® 8mm, Roxolid® Loxim®	24
021.2310	Bone Level Implant Ø 3.3mm NC, SLActive® 10mm, Roxolid® Loxim®	24
021.2312	Bone Level Implant Ø 3.3mm NC, SLActive® 12mm, Roxolid® Loxim®	24
021.2314	Bone Level Implant Ø 3.3mm NC, SLActive® 14mm, Roxolid® Loxim®	24
021.2408	Bone Level Implant, Ø 3.3 mm, NC, SLA 8 mm, Ti	24
021.2410	Bone Level Implant, Ø 3.3 mm, NC, SLA 10 mm, Ti	24
021.2412	Bone Level Implant, Ø 3.3 mm, NC, SLA 12 mm, Ti	24
021.2414	Bone Level Implant, Ø 3.3 mm, NC, SLA 14 mm, Ti	24
021.3308	Bone Level Tapered Implant, Ø 3.3mm NC, SLActive® 8mm, Roxolid®, Loxim®	24
021.3308G	Bone Level Tapered Implant, Ø 3.3mm NC, SLActive® 8mm, Roxolid®, guided	24
021.3310	Bone Level Tapered Implant, Ø 3.3mm NC, SLActive® 10mm, Roxolid®, Loxim®	24
021.3310G	Bone Level Tapered Implant, Ø 3.3mm NC, SLActive® 10mm, Roxolid®, guided	24
021.3312	Bone Level Tapered Implant, Ø 3.3mm NC, SLActive® 12mm, Roxolid®, Loxim®	24
021.3312G	Bone Level Tapered Implant, Ø 3.3mm NC, SLActive® 12mm, Roxolid®, guided	24
021.3314	Bone Level Tapered Implant, Ø 3.3mm NC, SLActive® 14mm, Roxolid®, Loxim®	24
021.3316	Bone Level Tapered Implant, Ø 3.3mm NC, SLActive® 16mm, Roxolid®, Loxim®	24
021.4108	Bone Level Implant, Ø 4.1 mm, RC, SLActive 8 mm, Ti	26
021.4108G	Bone Level Implant, Ø 4.1 mm, RC, SLActive 8 mm, Ti, guided	26
021.4110	Bone Level Implant, Ø 4.1 mm, RC, SLActive 10 mm, Ti	26
021.4110G	Bone Level Implant, Ø 4.1 mm, RC, SLActive 10 mm, Ti, guided	26
021.4112	Bone Level Implant, Ø 4.1 mm, RC, SLActive 12 mm, Ti	26
021.4112G	Bone Level Implant, Ø 4.1 mm, RC, SLActive 12 mm, Ti, guided	26
021.4114	Bone Level Implant, Ø 4.1 mm, RC, SLActive 14 mm, Ti	26
021.4308	Bone Level Implant Ø 4.1mm RC, SLActive® 8mm, Roxolid® Loxim®	26
021.4308G	Bone Level Implant, Ø 4.1mm RC SLActive® 8mm, Roxolid®, guided	26
021.4310	Bone Level Implant Ø 4.1mm RC, SLActive® 10mm, Roxolid® Loxim®	26
021.4310G	Bone Level Implant, Ø 4.1mm RC SLActive® 10mm, Roxolid®, guided	26
021.4312	Bone Level Implant Ø 4.1mm RC, SLActive® 12mm, Roxolid® Loxim®	26
021.4312G	Bone Level Implant, Ø 4.1mm RC SLActive® 12mm, Roxolid®, guided	26
021.4314	Bone Level Implant Ø 4.1mm RC, SLActive® 14mm, Roxolid® Loxim®	26
021.4408	Bone Level Implant, Ø 4.1 mm, RC, SLA 8 mm, Ti	26
021.4410	Bone Level Implant, Ø 4.1 mm, RC, SLA 10 mm, Ti	26
021.4412	Bone Level Implant, Ø 4.1 mm, RC, SLA 12 mm, Ti	26
021.4414	Bone Level Implant, Ø 4.1 mm, RC, SLA 14 mm, Ti	26
021.5308	Bone Level Tapered Implant, Ø 4.1mm RC, SLActive® 8mm, Roxolid®, Loxim®	26
021.5310	Bone Level Tapered Implant, Ø 4.1mm RC, SLActive® 10mm, Roxolid®, Loxim®	26
021.5312	Bone Level Tapered Implant, Ø 4.1mm RC, SLActive® 12mm, Roxolid®, Loxim®	26
021.5314	Bone Level Tapered Implant, Ø 4.1mm RC, SLActive® 14mm, Roxolid®, Loxim®	26
021.5316	Bone Level Tapered Implant, Ø 4.1mm RC, SLActive® 16mm, Roxolid®, Loxim®	26
021.6108	Bone Level Implant, Ø 4.8 mm, RC, SLActive 8 mm, Ti	28
021.6108G	Bone Level Implant, Ø 4.8 mm, RC, SLActive 8 mm, Ti, guided	28
021.6110	Bone Level Implant, Ø 4.8 mm, RC, SLActive 10 mm, Ti	28
021.6110G	Bone Level Implant, Ø 4.8 mm, RC, SLActive 10 mm, Ti, guided	28
021.6112	Bone Level Implant, Ø 4.8 mm, RC, SLActive 12 mm, Ti	28
021.6112G	Bone Level Implant, Ø 4.8 mm, RC, SLActive 12 mm, Ti, guided	28
021.6114	Bone Level Implant, Ø 4.8 mm, RC, SLActive 14 mm, Ti	28
021.6308	Bone Level Implant Ø 4.8mm RC, SLActive® 8mm, Roxolid® Loxim®	28
021.6308G	Bone Level Implant, Ø 4.8mm RC SLActive® 8mm, Roxolid®, guided	28
021.6310	Bone Level Implant Ø 4.8mm RC, SLActive® 10mm, Roxolid® Loxim®	28
021.6310G	Bone Level Implant, Ø 4.8mm RC SLActive® 10mm, Roxolid®, guided	28
021.6312	Bone Level Implant Ø 4.8mm RC, SLActive® 12mm, Roxolid® Loxim®	28
021.6312G	Bone Level Implant, Ø 4.8mm RC SLActive® 12mm, Roxolid®, guided	28
021.6314	Bone Level Implant Ø 4.8mm RC, SLActive® 14mm, Roxolid® Loxim®	28
021.6408	Bone Level Implant, Ø 4.8 mm, RC, SLA 8 mm, Ti	28

货号	产品名称	页数
021.6410	Bone Level Implant, Ø 4.8 mm, RC, SLA 10 mm, Ti	28
021.6412	Bone Level Implant, Ø 4.8 mm, RC, SLA 12 mm, Ti	28
021.6414	Bone Level Implant, Ø 4.8 mm, RC, SLA 14 mm, Ti	28
021.7308	Bone Level Tapered Implant, Ø 4.8mm RC, SLActive® 8mm, Roxolid®, Loxim®	28
021.7308G	Bone Level Tapered Implant, Ø 4.8mm RC, SLActive® 8mm, Roxolid®, guided	28
021.7310	Bone Level Tapered Implant, Ø 4.8mm RC, SLActive® 10mm, Roxolid®, Loxim®	28
021.7310G	Bone Level Tapered Implant, Ø 4.8mm RC, SLActive® 10mm, Roxolid®, guided	28
021.7312	Bone Level Tapered Implant, Ø 4.8mm RC, SLActive® 12mm, Roxolid®, Loxim®	28
021.7312G	Bone Level Tapered Implant, Ø 4.8mm RC, SLActive® 12mm, Roxolid®, guided	28
021.7314	Bone Level Tapered Implant, Ø 4.8mm RC, SLActive® 14mm, Roxolid®, Loxim®	28
021.7316	Bone Level Tapered Implant, Ø 4.8mm RC, SLActive® 16mm, Roxolid®, Loxim®	28
022.0010P	NC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 5.5mm, Type A, Pin	78
022.0011P	NC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 5.5mm, Type B, Pin	78
022.0012P	NC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 5.5mm, Type A, Pin	78
022.0013P	NC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 5.5mm, Type B, Pin	78
022.0014P	RC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 5.5mm, Type A, Pin	99
022.0015P	RC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 5.5mm, Type B, Pin	99
022.0016P	RC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 5.5mm, Type A, Pin	99
022.0017P	RC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 5.5mm, Type B, Pin	99
022.0019	RN Straumann® Variobase® for CEREC®, TAN	47
022.0020	WN Straumann® Variobase® for CEREC®, TAN	58
022.0021	NNC Variobase®, for crown, with screw, Ø 3.9mm, H 5.75mm, TAN	36,37
022.0022	RN Variobase®, for crown, with screw, Ø 5.05mm, H 6mm, TAN	47
022.0023	WN Variobase®, for crown, with screw, Ø 7mm, H 6.5mm, TAN	58,59
022.0026	RC Variobase®, for crown, with screw, Ø 4.5mm, H 5.5mm, TAN	96,97
022.0027	NC Variobase®, for crown, with screw, Ø 3.8mm, H 5.5mm, TAN	74,75
022.0045	Basal Screw RN/WN Variobase®, TAN	47,58
022.2102	NC Anatomic Abutment, straight, GH 2mm, Ti	71
022.2104	NC Anatomic Abutment, straight, GH 3.5mm, Ti	71
022.2152	NC Anatomic Abutment, angled 15°, GH 2mm, Ti	71
022.2154	NC Anatomic Abutment, angled 15°, GH 3.5mm, Ti	71
022.2311	NC Cementable Abutment, Ø 3.5mm, GH 1mm, AH 4mm, Ti	73
022.2312	NC Cementable Abutment, Ø 3.5mm, GH 2mm, AH 4mm, Ti	73
022.2313	NC Cementable Abutment, Ø 3.5mm, GH 3mm, AH 4mm, Ti	73
022.2315	NC Cementable Abutment, Ø 3.5mm, GH 1mm, AH 5.5mm, Ti	73
022.2316	NC Cementable Abutment, Ø 3.5mm, GH 2mm, AH 5.5mm, Ti	73
022.2317	NC Cementable Abutment, Ø 3.5mm, GH 3mm, AH 5.5mm, Ti	73
022.2321	NC Cementable Abutment, Ø 5mm, GH 1mm, AH 4mm, Ti	73
022.2322	NC Cementable Abutment, Ø 5mm, GH 2mm, AH 4mm, Ti	73
022.2323	NC Cementable Abutment, Ø 5mm, GH 3mm, AH 4mm, Ti	73
022.2325	NC Cementable Abutment, Ø 5mm, GH 1mm, AH 5.5mm, Ti	73
022.2326	NC Cementable Abutment, Ø 5mm, GH 2mm, AH 5.5mm, Ti	73
022.2327	NC Cementable Abutment, Ø 5mm, GH 3mm, AH 5.5mm, Ti	73
022.2410	NC Gold Abutment, for crown, H 3.7mm, Ceramicor®/POM	74
022.2502	NC LOCATOR® Abutment, height 2 mm, Ti alloy with TiN coating	88
022.2503	NC LOCATOR® Abutment, height 3 mm, Ti alloy with TiN coating	88
022.2504	NC LOCATOR® Abutment, height 4 mm, Ti alloy with TiN coating	88
022.2505	NC LOCATOR® Abutment, height 5 mm, Ti alloy with TiN coating	88
022.2506	NC LOCATOR® Abutment, height 6 mm, Ti alloy with TiN coating	88
022.2701	NC Abutment for bars, gold, H 10mm, Ceramicor®	81,85
022.2702	NC Abutment, for bars, Ti	81,86
022.2731	NC Multi-Base Abutment, straight, Ø 3.5mm, GH 1mm, TAN	82,85,86
022.2732	NC Multi-Base Abutment, straight, Ø 3.5mm, GH 2.5mm, TAN	82,85,86
022.2734	NC Multi-Base Abutment, straight, Ø 3.5mm, GH 4mm, TAN	82,85,86
022.2741	NC Multi-Base Abutment, straight, Ø 4.5mm, GH 1mm, TAN	82,85,86
022.2742	NC Multi-Base Abutment, straight, Ø 4.5mm, GH 2.5mm, TAN	82,85,86
022.2744	NC Multi-Base Abutment, straight, Ø 4.5mm, GH 4mm, TAN	82,85,86
022.2745P	NC Screw-retained Abutment, TAN, straight 0°, Ø 3.5mm, GH 1mm, Pin	77
022.2746P	NC Screw-retained Abutment, TAN, straight 0°, Ø 3.5mm, GH 2.5mm, Pin	77
022.2747P	NC Screw-retained Abutment, TAN, straight 0°, Ø 4.6mm, GH 1mm, Pin	77
022.2748P	NC Screw-retained Abutment, TAN, straight 0°, Ø 4.6mm, GH 2.5mm, Pin	77
022.2749P	NC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 2.5mm, Type A, Pin	78
022.2750P	NC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 2.5mm, Type B, Pin	78
022.2751P	NC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 2.5mm, Type A, Pin	78
022.2752P	NC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 2.5mm, Type B, Pin	78
022.2753P	NC Screw-retained Abutment, TAN, straight 0°, Ø 3.5mm, GH 4mm, Pin	77
022.2754P	NC Screw-retained Abutment, TAN, straight 0°, Ø 4.6mm, GH 4mm, Pin	77

货号	产品名称	页数
022.2755P	NC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 4mm, Type A, Pin	78
022.2756P	NC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 4mm, Type B, Pin	78
022.2757P	NC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 4mm, Type A, Pin	78
022.2758P	NC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 4mm, Type B, Pin	78
022.2782	NC Multi-Base Abutment, Type A, angled 25°, TAN	82,85,86
022.2786	NC Multi-Base Abutment, Type B, angled 25°, TAN	82,85,86
022.4102	RC Anatomic Abutment, straight, GH 2mm, Ti	93
022.4104	RC Anatomic Abutment, straight, GH 3.5mm, Ti	93
022.4152	RC Anatomic Abutment, angled 15°, GH 2mm, Ti	93
022.4154	RC Anatomic Abutment, angled 15°, GH 3.5mm, Ti	93
022.4321	RC Cementable Abutment, Ø 5mm, GH 1mm, AH 4mm, Ti	95
022.4322	RC Cementable Abutment, Ø 5mm, GH 2mm, AH 4mm, Ti	95
022.4323	RC Cementable Abutment, Ø 5mm, GH 3mm, AH 4mm, Ti	95
022.4325	RC Cementable Abutment, Ø 5mm, GH 1mm, AH 5.5mm, Ti	95
022.4326	RC Cementable Abutment, Ø 5mm, GH 2mm, AH 5.5mm, Ti	95
022.4327	RC Cementable Abutment, Ø 5mm, GH 3mm, AH 5.5mm, Ti	95
022.4331	RC Cementable Abutment, Ø 6.5mm, GH 1mm, AH 4mm, Ti	95
022.4332	RC Cementable Abutment, Ø 6.5mm, GH 2mm, AH 4mm, Ti	95
022.4333	RC Cementable Abutment, Ø 6.5mm, GH 3mm, AH 4mm, Ti	95
022.4335	RC Cementable Abutment, Ø 6.5mm, GH 1mm, AH 5.5mm, Ti	95
022.4336	RC Cementable Abutment, Ø 6.5mm, GH 2mm, AH 5.5mm, Ti	95
022.4337	RC Cementable Abutment, Ø 6.5mm, GH 3mm, AH 5.5mm, Ti	95
022.4410	RC Gold Abutment, for crown, H 3.8mm, Ceramicor®/POM	96
022.4430	RC Gold Abutment, for bridge, H 5.2mm, Ceramicor®/POM	96
022.4501	RC LOCATOR® Abutment, height 1 mm, Ti alloy with TiN coating	107
022.4502	RC LOCATOR® Abutment, height 2 mm, Ti alloy with TiN coating	107
022.4503	RC LOCATOR® Abutment, height 3 mm, Ti alloy with TiN coating	107
022.4504	RC LOCATOR® Abutment, height 4 mm, Ti alloy with TiN coating	107
022.4505	RC LOCATOR® Abutment, height 5 mm, Ti alloy with TiN coating	107
022.4506	RC LOCATOR® Abutment, height 6 mm, Ti alloy with TiN coating	107
022.4701	RC Abutment for bars, gold, H 10mm, Ceramicor®	101,104
022.4702	RC Abutment, for bars, Ti	101,105
022.4741	RC Multi-Base Abutment, straight, Ø 4.5mm, GH 1mm, TAN	102,104,105
022.4742	RC Multi-Base Abutment, straight, Ø 4.5mm, GH 2.5mm, TAN	102,104,105
022.4744	RC Multi-Base Abutment, straight, Ø 4.5mm, GH 4mm, TAN	102,104,105
022.4745P	RC Screw-retained Abutment, TAN, straight 0°, Ø 4.6mm, GH 1mm, Pin	99
022.4746P	RC Screw-retained Abutment, TAN, straight 0°, Ø 4.6mm, GH 2.5mm, Pin	99
022.4747P	RC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 2.5mm, Type A, Pin	99
022.4748P	RC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 2.5mm, Type B, Pin	99
022.4749P	RC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 2.5mm, Type A, Pin	99
022.4750P	RC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 2.5mm, Type B, Pin	99
022.4751P	RC Screw-retained Abutment, TAN, straight 0°, Ø 4.6mm, GH 4mm, Pin	99
022.4752P	RC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 4mm, Type A, Pin	99
022.4753P	RC Screw-retained Abutment, TAN, angled 17°, Ø 4.6mm, GH 4mm, Type B, Pin	99
022.4754P	RC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 4mm, Type A, Pin	99
022.4755P	RC Screw-retained Abutment, TAN, angled 30°, Ø 4.6mm, GH 4mm, Type B, Pin	99
022.4761	RC Multi-Base Abutment, straight, Ø 6.5mm, GH 1mm, TAN	102,104,105
022.4762	RC Multi-Base Abutment, straight, Ø 6.5mm, GH 2.5mm, TAN	102,104,105
022.4764	RC Multi-Base Abutment, straight, Ø 6.5mm, GH 4mm, TAN	102,104,105
022.4782	RC Multi-Base Abutment, Type A, angled 25°, TAN	102,104,105
022.4786	RC Multi-Base Abutment, Type B, angled 25°, TAN	102,104,105
023.0014	NNC Burn-out Coping, for Variobase®, for crown, AH 5.75mm, POM	37
023.0014V4	NNC Burn-out Coping, for Variobase®, for crown, AH 5.75mm, POM	37
023.0015	RN Burn-out Coping, for Variobase®, for crown, AH 6mm, POM	47
023.0015V4	RN Burn-out Coping, for Variobase®, for crown, AH 6mm, POM	47
023.0016	WN Burn-out Coping, for Variobase®, for crown, AH 6.5mm, POM	59
023.0016V4	WN Burn-out Coping, for Variobase®, for crown, AH 6.5mm, POM	59
023.0017	RC Burn-out Coping, for Variobase®, for crown, AH 5.5mm, POM	97
023.0017V4	RC Burn-out Coping, for Variobase®, for crown, AH 5.5mm, POM	97
023.0018	NC Burn-out Coping, for Variobase®, for crown, AH 5.5mm, POM	75
023.0018V4	NC Burn-out Coping, for Variobase®, for crown, AH 5.5mm, POM	75
023.2310	NC Coping, for bridge on Cementable Abutment, Ø 3.5mm, POM	73
023.2316	NC Coping, for crown on Cementable Abutment, Ø 3.5mm, POM	73
023.2320	NC Coping, for bridge on Cementable Abutment, Ø 5mm, POM	73
023.2326	NC Coping, for crown on Cementable Abutment, Ø 5mm, POM	73
023.2701	NC Coping, for bridge on multi-base Abutment, angled 25°, Ceramicor®/POM	83
023.2702	NC Coping, for bar on multi-base Abutment, angled 25°, Ceramicor®	85

货号	产品名称	页数
023.2704	NC Coping, for bridge on Multi-Base Abutment, angled 25°, POM	83
023.2706	NC Coping, for bar on Multi-Base Abutment, angled 25°, Ti	86
023.2731	NC Coping, for bridge on multi-base Abutment, straight, Ø 3.5mm, Ceramicor®/POM	83
023.2732	NC Coping, for bar on multi-base Abutment, straight, Ø 3.5mm, Ceramicor®	85
023.2734	NC Coping, for bridge on Multi-Base Abutment, straight, Ø 3.5mm, POM	83
023.2736	NC Coping, for bar on Multi-Base Abutment, straight, Ø 3.5mm, Ti	86
023.2741	NC Coping, for bridge on multi-base Abutment, straight, Ø 4.5mm, Ceramicor®/POM	83
023.2742	NC Coping, for bar on multi-base Abutment, straight, Ø 4.5mm, Ceramicor®	85
023.2744	NC Coping, for bridge on Multi-Base Abutment, straight, Ø 4.5mm, POM	83
023.2746	NC Coping, for bar on Multi-Base Abutment, straight, Ø 4.5mm, Ti	86
023.2747	NC Coping for screw-retained abutment, Ti, Crown, Ø 3.5mm	79
023.2748	NC Burn-out coping for screw-retained abutment, POM, Crown, Ø 3.5mm	79
023.2749	NC Coping for screw-retained abutment, Ti, Bridge, Ø 3.5mm	79
023.2750	NC Coping for screw-retained abutment, Ti, Bar, Ø 3.5mm	80
023.2751	NC Coping for screw-retained abutment, Crown, Ø 3.5mm, Ceramicor®/POM	79
023.2752	NC Coping for screw-retained abutment, Bridge, Ø 3.5mm, Ceramicor®/POM	79
023.2753	NC Coping for screw-retained abutment, Bar, Ø 3.5mm, Ceramicor®	80
023.2754	NC Analog for screw-retained abutment, TAN, straight 0°, Ø 3.5mm	76
023.2755	NC Burn-out coping f. screw-retained abut., POM, Bridge/Bar, Ø 3.5mm	79,80
023.2756	NC Burn-out coping, for Variobase® Abutment, POM	75
023.2756-04	NC Burn-out coping, for Variobase® Abutment, POM	75
023.4320	RC Coping, for bridge on Cementable Abutment, Ø 5mm, POM	95
023.4326	RC Coping, for crown on Cementable Abutment, Ø 5mm, POM	95
023.4330	RC Coping, for bridge on Cementable Abutment, Ø 6.5mm, POM	95
023.4336	RC Coping, for crown on Cementable Abutment, Ø 6.5mm, POM	95
023.4702	RC Coping, for bar on multi-base Abutment, angled 25°, Ceramicor®	104
023.4706	RC Coping, for bar on Multi-Base Abutment, angled 25°, Ti	105
023.4742	RC Coping, for bar on multi-base Abutment, straight, Ø 4.5mm, Ceramicor®	104
023.4746	RC Coping, for bar on Multi-Base Abutment, straight, Ø 4.5mm, Ti	105
023.4747	NC/RC Coping for screw-retained abutment, Ti, Crown, Ø 4.6mm	79,100
023.4748	NC/RC Burn-out coping for screw-retained abutment, POM, Crown, Ø 4.6mm	79,100
023.4749	NC/RC Screw for screw-retained abutment, TAN, straight 0°, GH 1mm	78,99
023.4750	NC/RC Screw for screw-retained abutment, TAN, straight 0°, GH 2.5mm	78,99
023.4751	NC/RC Coping for screw-retained abutment, Ti, Bridge, Ø 4.6mm	79,100
023.4752	NC/RC Coping for screw-retained abutment, Ti, Bar, Ø 4.6mm	80,101
023.4753	NC/RC Coping for screw-retained abutment, Crown, Ø 4.6mm, Ceramicor®/POM	79,100
023.4754	NC/RC Coping for screw-retained abutment, Bridge, Ø 4.6mm, Ceramicor®/POM	79,100
023.4755	NC/RC Coping for screw-retained abutment, Bar, Ø 4.6mm, Ceramicor®	80,100
023.4756	NC/RC Analog for screw-retained abutment, TAN, straight 0°, Ø 4.6mm	76,98
023.4757	NC/RC Analog for screw-retained abutment, TAN, angled 17°/30°, Ø 4.6mm	76,98
023.4758	NC/RC Burn-out coping f. screw-retained abut., POM, Bridge/Bar, Ø 4.6mm	79,80,100,101
023.4759	RC Burn-out coping, for Variobase® Abutment, POM	97
023.4759-04	RC Burn-out coping, for Variobase® Abutment, POM	97
023.4760	NC/RC Screw for screw-retained abutment, TAN, straight 0°, GH 4mm	78,99
023.4762	RC Coping, for bar on multi-base Abutment, straight, Ø 6.5mm, Ceramicor®	104
023.4763	NC/RC Occlusal screw, TAN, for coping, screw-retained abutment	79,80,100,101
023.4766	RC Coping, for bar on Multi-Base Abutment, straight, Ø 6.5mm, Ti	105
024.2100	NC Closure Screw, H 0mm, Ti	25
024.2100-04	NC Closure Screw, 4 pieces, Height 0 mm, Ti	25
024.2105	NC Closure Screw, H 0.5mm, Ti	25
024.2105-04	NC Closure Screw, 4 pieces, Height 0.5 mm, Ti	25
024.2222	NC Healing Abutment, conical shape, Ø 3.6mm, H 2mm, Ti	25
024.2224	NC Healing Abutment, conical shape, Ø 3.6mm, H 3.5mm, Ti	25
024.2226	NC Healing Abutment, conical shape, Ø 3.6mm, H 5mm, Ti	25
024.2234	NC Healing Abutment, bottle shape, Ø 3.3mm, H3.5mm, Ti	25
024.2236	NC Healing Abutment, bottle shape, Ø 3.3mm, H 5mm, Ti	25
024.2242	NC Healing Abutment, conical shape, Ø 4.8mm, H 2mm, Ti	25
024.2244	NC Healing Abutment, conical shape, Ø 4.8mm, H 3.5mm, Ti	25
024.2246	NC Healing Abutment, conical shape, Ø 4.8mm, H 5mm, Ti	25
024.2270	NC Healing Abutment, customizable, Ø 5mm, PEEK	25
024.2306-04	NC Protective Cap, for Multi-Base Abutment, angled 25°, PEEK/TAN	82
024.2307	NC Temporary Coping, for Multi-Base Abutment, angled 25°, Ti	82
024.2312-04	NC Protective Cap, for Cementable Abutment, Ø 3.5mm, AH 4mm, PEEK	72
024.2316-04	NC Protective Cap, for Cementable Abutment, Ø 3.5mm, AH 5.5mm, PEEK	72
024.2322-04	NC Protective Cap, for Cementable Abutment, Ø 5mm, AH 4mm, PEEK	72
024.2323-04	NC Protective cap for screw-retained abutment, Ø 3.5mm, H 5mm, PEEK/TAN	77
024.2324-04	NC Protective cap for screw-retained abutment, Ø 3.5mm, H 6.5mm, PEEK/TAN	77

货号	产品名称	页数
024.2325-04	NC Protective cap for screw-retained abutment, Ø 3.5mm, H 8mm, PEEK/TAN	77
024.2326-04	NC Protective Cap, for Cementable Abutment, Ø 5mm, AH 5.5mm, PEEK	72
024.2331-04	NC Protective Cap, for Multi-Base Abutment, straight, Ø 3.5mm, PEEK/TAN	82
024.2332	NC Temporary Coping, for Multi-Base Abutment, straight, Ø 3.5mm, Ti	82
024.2341-04	NC Protective Cap, for Multi-Base Abutment, straight, Ø 4.5mm, PEEK/TAN	82
024.2342	NC Temporary Coping, for Multi-Base Abutment, straight, Ø 4.5mm, Ti	82
024.2346	NC Temporary Coping, for Cementable Abutment, Ø 3.5mm, PMMA	72
024.2356	NC Temporary Coping, for Cementable Abutment, Ø 5mm, PMMA	72
024.2370	NC Temporary Abutment, Ø 5mm, H 12mm, PEEK/TAN	70
024.2371	NC Temporary Abutment, Crown, Ø 3.5mm, H 11mm, TAN	70
024.2375	NC Temporary Abutment, Bridge, Ø 3.5mm, H 11mm, TAN	70
024.4100	RC Closure Screw, H 0mm, Ti	27,29
024.4100-04	RC Closure Screw, 4 pieces, Height 0 mm, Ti	27,29
024.4105	RC Closure Screw, H 0.5mm, Ti	27,29
024.4105-04	RC Closure Screw, 4 pieces, Height 0.5 mm, Ti	27,29
024.4222	RC Healing Abutment, conical shape, Ø 4.5mm, H 2mm, Ti	27,29
024.4224	RC Healing Abutment, conical shape, Ø 4.5mm, H 4mm, Ti	27,29
024.4226	RC Healing Abutment, conical shape, Ø 4.5mm, H 6mm, Ti	27,29
024.4234	RC Healing Abutment, bottle shape, Ø 4.4mm, H 4mm, Ti	27,29
024.4236	RC Healing Abutment, bottle shape, Ø 4.7mm, H 6mm, Ti	27,29
024.4242	RC Healing Abutment, conical shape, Ø 6mm, H 2mm, Ti	27,29
024.4244	RC Healing Abutment, conical shape, Ø 6mm, H 4mm, Ti	27,29
024.4246	RC Healing Abutment, conical shape, Ø 6mm, H 6mm, Ti	27,29
024.4270	RC Healing Abutment, customizable, Ø 7mm, PEEK	27,29
024.4306-04	RC Protective Cap, for Multi-base Abutment, angled 25°, PEEK/TAN	102
024.4307	RC Temporary Coping, for Multi-Base Abutment, angled 25°, Ti	102
024.4322-04	RC Protective Cap, for Cementable Abutment, Ø 5mm, AH 4mm, PEEK	94
024.4323-04	NC/RC Protective cap for screw-retained abutment, Ø 4.6mm, H 5.1mm, PEEK/TAN	77,98
024.4324-04	NC/RC Protective cap for screw-retained abutment, Ø 4.6mm, H 6.6mm, PEEK/TAN	77,98
024.4325-04	NC/RC Protective cap for screw-retained abutment, Ø 4.6mm, H 8.1mm, PEEK/TAN	77,98
024.4326-04	RC Protective Cap, for Cementable Abutment, Ø 5mm, AH 5.5mm, PEEK	94
024.4332-04	RC Protective Cap, for Cementable Abutment, Ø 6.5mm, AH 4mm, PEEK	94
024.4336-04	RC Protective Cap, for Cementable Abutment, Ø 6.5mm, AH 5.5mm, PEEK	94
024.4341-04	RC Protective Cap, for Multi-base Abutment, straight, Ø 4.5mm, PEEK/TAN	102
024.4342	RC Temporary Coping, for Multi-Base Abutment, straight, Ø 4.5mm, Ti	102
024.4356	RC Temporary Coping, for Cementable Abutment, Ø 5mm, PMMA	94
024.4361-04	RC Protective Cap, for Multi-base Abutment, straight, Ø 6.5mm, PEEK/TAN	102
024.4362	RC Temporary Coping, for Multi-Base Abutment, straight, Ø 6.5mm, Ti	102
024.4366	RC Temporary Coping, for Cementable Abutment, Ø 6.5mm, PMMA	94
024.4370	RC Temporary Abutment, Ø 7mm, H 12mm, PEEK/TAN	92
024.4371	RC Temporary Abutment, Crown, Ø 4.5mm, H 11mm, TAN	92
024.4375	RC Temporary Abutment, Bridge, Ø 4.5mm, H 11mm, TAN	92
025.0000	CARES® Mono Scanbody for screw-retained abutment, abutment level, Ø 3.5mm, PEEK/TAN	76
025.0001	CARES® Mono Scanbody for screw-retained abutment, abutment level, Ø 4.6mm, PEEK/TAN	76,98
025.0002	NC/RC Screw, for screw-retained abutment 17°/30°, TAN	78,99
025.0004	Polishing Aid, for Screw-retained Abutments Ø 3.5mm, TAN	76
025.0004V4	Polishing Aid, for Screw-retained Abutments Ø 3.5mm, TAN	76
025.0005	Polishing Aid, for Screw-retained Abutments Ø 4.6mm, TAN	76,98
025.0005V4	Polishing Aid, for Screw-retained Abutments Ø 4.6mm, TAN	76,98
025.0006	Lab Processing Screw, for Screw-retained Abutments, stainless steel	76,98
025.0009	Transfer and Alignment Pin, for Screw-retained Abutments, TAN	77,98
025.0011	Impression Post f. open tray f. screw-ret. abut., abut. level, Ø 3.5mm, non-engaging, TAN	76
025.0012	Impression Post f. open tray f. screw-ret. abut., abut. level, Ø 4.6mm, non-engaging, TAN	76,98
025.0013	Impression Post f. closed tray f. screw-ret. abut., abut. level, Ø 3.5mm, non-engaging, TAN/POM	76
025.0014	Impression Post f. closed tray f. screw-ret. abut., abut. level, Ø 4.6mm, non-engaging, TAN/POM	76,98
025.0049	NC Analog for Screw-retained Abutments, edentulous, straight 0°, Ø 3.5mm, TAN	76
025.0050	NC/RC Analog for Screw-retained Abutments, edentulous, straight 0°, Ø 4.6mm, TAN	76,98
025.0052	Lab Processing Screw, for Screw-retained Abutments, L 10mm, stainless steel	76,98
025.0120-04	LOCATOR® Female Implant Analog, Ø 4.0 mm	39,88,107
025.0900-04	NC/RC Occlusal Screw, for coping on, Multi-Base Abutment, TAN	83,84,85,86,87,103,104,105,106
025.0911-04	NC/RC Stabilization pin for bar on Multi-base abutment, L 16mm, Stainless steel	85,86,87,104,105,106
025.2101	NC Bone Level Implant Analog, L 12mm, Ti	70
025.2106	NC Implant Analog, for Multi-base Abutment, angled 25°Ø 4mm, TAN	82
025.2112	NC Implant Analog, for Cementable Abutment, Ø 3.5mm, AH 4mm, TAN	72
025.2116	NC Implant Analog, for Cementable Abutment, Ø 3.5mm, AH 5.5mm, TAN	72
025.2122	NC Implant Analog, for Cementable Abutment, Ø 5mm, AH 4mm, TAN	72
025.2126	NC Implant Analog, for Cementable Abutment, Ø 5mm, AH 5.5mm, TAN	72

索引 货号索引		
货号	产品名称	页数
025.2131	NC Implant Analog, for Multi-base Abutment, straight, Ø 3.5mm, TAN	82
025.2141	NC Implant Analog, for Multi-base Abutment, straight, Ø 4.5mm, TAN	82
025.2201	NC Impression Post Closed Tray, with 1 guide screw & 2 caps, TAN/POM	70
025.2202	NC Impression Post Open Tray, short, with guide screw, L 16.5mm, TAN	70
025.2203-04	NC Impression Cap, for Cementable Abutment, Ø 3.5mm, AH 4mm, POM	72
025.2204-04	NC Impression Cap, for Cementable Abutment, Ø 5mm, AH 4mm, POM	72
025.2205	NC Impression post, for open tray, with guide screw, L 30mm, TAN	70
025.2206	NC Impression Post, for Multi-Base Abutment, angled 25°, TAN	82
025.2208-04	NC Bite Registration Aid, H 8mm, POM	70
025.2212-04	NC Bite Registration Aid, H 12mm, POM	70
025.2220-04	NC Impression Cap, for Cementable Abutment, Ø 5mm, AH 5.5mm, POM	72
025.2231-04	NC Impression Cap, for Multi-Base Abutment, straight, Ø 3.5mm, POM	82
025.2232	NC Impression Post, for Multi-Base Abutment, straight, Ø 3.5mm, TAN	82
025.2241-04	NC Impression Cap, for Multi-Base Abutment, straight, Ø 4.5mm, POM	82
025.2242	NC Impression Post, for Multi-Base Abutment, straight, Ø 4.5mm, TAN	82
025.2243	Impression post for open tray, TAN, for screw-ret. abut., abut. level, 0°, Ø 3.5mm	76
025.2244	Impression post for open tray, TAN, for screw-ret. abut., abut. level, 0°, Ø 4.6mm	76,98
025.2245	Impression post for closed tray, abutment level, for screw-retained abutment, Ø 3.5mm, TAN/POM	76
025.2246	Impression post for closed tray, abutment level, for screw-retained abutment, Ø 4.6mm, TAN/POM	76,98
025.2601-04	NC PLAN Anatomic Abutment, straight, GH 2mm, POM	110
025.2602-04	NC PLAN Anatomic Abutment, straight, GH 3.5mm, POM	110
025.2603-04	NC PLAN Anatomic Abutment, angled 15°, GH 2mm, POM	110
025.2604-04	NC PLAN Anatomic Abutment, angled 15°, GH 3.5mm, POM	110
025.2610-04	NC PLAN Gold Abutment, for crown, POM	110
025.2611-04	NC PLAN Cementable Abutment, Ø 3.5mm, GH 1mm, AH 4mm, POM	110
025.2612-04	NC PLAN Cementable Abutment, Ø 3.5mm, GH 2mm, AH 4mm, POM	110
025.2613-04	NC PLAN Cementable Abutment, Ø 3.5mm, GH 3mm, AH 4mm, POM	110
025.2615-04	NC PLAN Cementable Abutment, Ø 3.5mm, GH 1mm, AH 5.5mm, POM	110
025.2616-04	NC PLAN Cementable Abutment, Ø 3.5mm, GH 2mm, AH 5.5mm, POM	110
025.2617-04	NC PLAN Cementable Abutment, Ø 3.5mm, GH 3mm, AH 5.5mm, POM	110
025.2621-04	NC PLAN Cementable Abutment, Ø 5mm, GH 1mm, AH 4mm, POM	110
025.2622-04	NC PLAN Cementable Abutment, Ø 5mm, GH 2mm, AH 4mm, POM	110
025.2623-04	NC PLAN Cementable Abutment, Ø 5mm, GH 3mm, AH 4mm, POM	110
025.2625-04	NC PLAN Cementable Abutment, Ø 5mm, GH 1mm, AH 5.5mm, POM	110
025.2626-04	NC PLAN Cementable Abutment, Ø 5mm, GH 2mm, AH 5.5mm, POM	110
025.2627-04	NC PLAN Cementable Abutment, Ø 5mm, GH 3mm, AH 5.5mm, POM	110
025.2642-04	NC PLAN LOCATOR® Abutment, H 2mm, POM	111
025.2643-04	NC PLAN LOCATOR® Abutment, H 3mm, POM	111
025.2644-04	NC PLAN LOCATOR® Abutment, H 4mm, POM	111
025.2645-04	NC PLAN LOCATOR® Abutment, H 5mm, POM	111
025.2646-04	NC Novaloc® Plan Abutment, H 1-6 mm, POM	111
025.2648-04	NC Plan screw-retained abutment, POM, straight 0°, Ø 3.5mm, GH 2.5mm	111
025.2650-04	NC Plan screw-retained abutment, POM, straight 0°, Ø 4.6mm, GH 2.5mm	111
025.2651-04	NC PLAN Multi-Base Abutment, straight, Ø 3.5mm, GH 1mm, POM	113
025.2652-04	NC PLAN Multi-Base Abutment, straight, Ø 3.5mm, GH 2.5mm, POM	113
025.2653-04	NC Plan screw-retained abutment, POM, angled 30°, Ø 4.6mm, GH 2.5mm, Type A	111
025.2654-04	NC PLAN Multi-Base Abutment, straight, Ø 3.5mm, GH 4mm, POM	113
025.2655-04	NC Plan screw-retained abutment, POM, angled 17°, Ø 4.6mm, GH 2.5mm, Type A	111
025.2656-04	NC PLAN Multi-Base Abutment, straight, Ø 4.5mm, GH 1mm, POM	113
025.2657-04	NC PLAN Multi-Base Abutment, straight, Ø 4.5mm, GH 2.5mm, POM	113
025.2658-04	NC Plan screw-retained abutment, POM, angled 17°, Ø 4.6mm, GH 2.5mm, Type B	111
025.2659-04	NC PLAN Multi-Base Abutment, straight, Ø 4.5mm, GH 4mm, POM	113
025.2660-04	NC Plan screw-retained abutment, POM, angled 30°, Ø 4.6mm, GH 2.5mm, Type B	111
025.2662-04	NC PLAN Multi-Base Abutment, Type A, angled 25°, POM	113
025.2665-04	NC PLAN Multi-Base Abutment, Type B, angled 25°, POM	113
025.2900	NC Basal Screw, TAN	70,71,74,75,81,84,85,86,87,151
025.2908	NC Basal Screw, for Cementable Abutment only, TAN	73,151
025.2911-04	NC Transfer Aid, for Cementable Abutment, Ø 3.5mm, PC	72
025.2912-04	NC Transfer Aid, for Cementable Abutment, Ø 5mm, PC	72
025.2915	CARES® NC Mono Scanbody, Ø 3.5 mm, H 10 mm, PEEK, TAN	70
025.2920	NC Polishing Aid, stainless steel	74,84,87
025.2920-04	NC Polishing Aid, stainless steel	74,84,87
025.2921	NC Variobase® Abutment, with screw, Ø 3.8mm, H 3.5mm, TAN	74,75
025.4101	RC Bone Level Implant Analog, L 12mm, Ti	92
025.4106	RC Implant Analog for Multi-base abutment,angled 25°, Ø 4mm, TAN	102
025.4122	RC Implant Analog, for Cementable Abutment, Ø 6.5mm, AH 4mm, TAN	94
025.4126	RC Implant Analog, for Cementable Abutment, Ø 6.5mm, AH 5.5mm, TAN	94

货号	产品名称	页数
025.4132	RC Implant Analog, for Cementable Abutment, Ø 5mm, AH 4mm, TAN	94
025.4136	RC Implant Analog, for Cementable Abutment, Ø 5mm, AH 5.5mm, TAN	94
025.4141	RC Implant Analog for Multi-base abutment, straight, Ø 4.5mm, TAN	102
025.4161	RC Implant Analog for Multi-base abutment, straight, Ø 6.5mm, TAN	102
025.4201	RC Impression Post Closed Tray, with 1 guide screw & 2 caps, TAN/POM	92
025.4202	RC Impression Post Open Tray, short, with guide screw, L 16.5mm, TAN	92
025.4203-04	RC Impression Cap, for Cementable Abutment, Ø 5mm, AH 4mm, POM	94
025.4204-04	RC Impression Cap, for Cementable Abutment, Ø 6.5mm, AH 4mm, POM	94
025.4205	RC Impression post, for open tray, with guide screw, L 30mm, TAN	92
025.4206	RC Impression Post, for Multi-Base Abutment, angled 25°, TAN	102
025.4208-04	RC Bite Registration Aid, H 8mm, POM	92
025.4212-04	RC Bite Registration Aid, H 12mm, POM	92
025.4220-04	RC Impression Cap, for Cementable Abutment, Ø 5mm, AH 5.5mm, POM	94
025.4230-04	RC Impression Cap, for Cementable Abutment, Ø 6.5mm, AH 5.5mm, POM	94
025.4241-04	RC Impression Cap, for Multi-Base Abutment, straight, Ø 4.5mm, POM	102
025.4242	RC Impression Post, for Multi-Base Abutment, straight, Ø 4.5mm, TAN	102
025.4261-04	RC Impression Cap, for Multi-Base Abutment, straight, Ø 6.5mm, POM	102
025.4262	RC Impression Post, for Multi-Base Abutment, straight, Ø 6.5mm, TAN	102
025.4601-04	RC PLAN Anatomic Abutment, straight, GH 2mm, POM	111
025.4602-04	RC PLAN Anatomic Abutment, straight, GH 3.5mm, POM	111
025.4603-04	RC PLAN Anatomic Abutment, angled 15°, GH 2mm, POM	111
025.4604-04	RC PLAN Anatomic Abutment, angled 15°, GH 3.5mm, POM	111
025.4610-04	RC PLAN Gold Abutment for crown, POM	111
025.4621-04	RC PLAN Cementable Abutment, Ø 5mm, GH 1mm, AH 4mm, POM	111
025.4622-04	RC PLAN Cementable Abutment, Ø 5mm, GH 2mm, AH 4mm, POM	111
025.4623-04	RC PLAN Cementable Abutment, Ø 5mm, GH 3mm, AH 4mm, POM	111
025.4625-04	RC PLAN Cementable Abutment, Ø 5mm, GH 1mm, AH 5.5mm, POM	111
025.4626-04	RC PLAN Cementable Abutment, Ø 5mm, GH 2mm, AH 5.5mm, POM	111
025.4627-04	RC PLAN Cementable Abutment, Ø 5mm, GH 3mm, AH 5.5mm, POM	111
025.4630-04	RC PLAN Gold Abutment for bridge, POM	111
025.4631-04	RC PLAN Cementable Abutment, Ø 6.5mm, GH 1mm, AH 4mm, POM	112
025.4632-04	RC PLAN Cementable Abutment, Ø 6.5mm, GH 2mm, AH 4mm, POM	112
025.4633-04	RC PLAN Cementable Abutment, Ø 6.5mm, GH 3mm, AH 4mm, POM	112
025.4635-04	RC PLAN Cementable Abutment, Ø 6.5mm, GH 1mm, AH 5.5mm, POM	112
025.4636-04	RC PLAN Cementable Abutment, Ø 6.5mm, GH 2mm, AH 5.5mm, POM	112
025.4637-04	RC PLAN Cementable Abutment, Ø 6.5mm, GH 3mm, AH 5.5mm, POM	112
025.4641-04	RC PLAN LOCATOR® Abutment, H 1mm, POM	112
025.4642-04	RC PLAN LOCATOR® Abutment, H 2mm, POM	112
025.4643-04	RC PLAN LOCATOR® Abutment, H 3mm, POM	112
025.4644-04	RC PLAN LOCATOR® Abutment, H 4mm, POM	112
025.4645-04	RC PLAN LOCATOR® Abutment, H 5mm, POM	112
025.4646-04	RC PLAN NOVALOC® Abutment, H 1-6mm, POM	112
025.4648-04	RC Plan screw-retained abutment, POM, straight 0°, Ø 4.6mm, GH 2.5mm	112
025.4649-04	RC Plan screw-retained abutment, POM, angled 17°, Ø 4.6mm, GH 2.5mm, Type A	112
025.4650-04	RC Plan screw-retained abutment, POM, angled 17°, Ø 4.6mm, GH 2.5mm, Type B	112
025.4651-04	RC PLAN Multi-Base Abutment, straight, Ø 4.5mm, GH 1mm, POM	113
025.4652-04	RC PLAN Multi-Base Abutment, straight, Ø 4.5mm, GH 2.5mm, POM	113
025.4653-04	RC Plan screw-retained abutment, POM, angled 30°, Ø 4.6mm, GH 2.5mm, Type A	112
025.4654-04	RC PLAN Multi-Base Abutment, straight, Ø 4.5mm, GH 4mm, POM	113
025.4656-04	RC PLAN Multi-Base Abutment, straight, Ø 6.5mm, GH 1mm, POM	113
025.4657-04	RC PLAN Multi-Base Abutment, straight, Ø 6.5mm, GH 2.5mm, POM	113
025.4659-04	RC PLAN Multi-Base Abutment, straight, Ø 6.5mm, GH 4mm, POM	113
025.4660-04	RC Plan screw-retained abutment, POM, angled 30°, Ø 4.6mm, GH 2.5mm, Type B	112
025.4662-04	RC PLAN Multi-Base Abutment, Type A, angled 25°, POM	113
025.4665-04	RC PLAN Multi-Base Abutment, Type B, angled 25°, POM	113
025.4900	RC Basal Screw, TAN	92,93,96,97,101,104,105,106,151
025.4902	RC Basal Screw, for Gold Abutment bridge only, TAN	96,151
025.4908	RC Basal Screw, for Cementable Abutment only, TAN	92,95,151
025.4912-04	RC Transfer Aid, for Cementable Abutment, Ø 5mm, PC	94
025.4913-04	RC Transfer Aid, for Cementable Abutment, Ø 6.5mm, PC	94
025.4915	CARES® RC Mono Scanbody, Ø 4.1 mm, H 10 mm, PEEK, TAN	92,96
025.4920	RC Polishing Aid, stainless steel	106
025.4920-04	RC Polishing Aid, stainless steel	106
025.4921	RC Variobase® Abutment, with screw, Ø 4.5mm, H 3.5mm, TAN	96,97
026.0000	CrossFit® Plan Cassette, empty, PC	110
026.0001	BLT Pilot Drill, short, Ø 2.2mm, L 33mm, stainless steel	127
026.0002	BLT Pilot Drill, long, Ø 2.2mm, L 41mm, stainless steel	127

货号	产品名称	页数
026.0003	BLT Profile Drill, short, Ø 3.3mm, L 25mm, stainless steel	128
026.0004	BLT Profile Drill, long, Ø 3.3mm, L 33mm, stainless steel	128
026.0005	BLT Profile Drill, short, Ø 4.1mm, L 25mm, stainless steel	128
026.0006	BLT Profile Drill, long, Ø 4.1mm, L 33mm, stainless steel	128
026.0007	BLT Profile Drill, short, Ø 4.8mm, L 25mm, stainless steel	128
026.0008	BLT Profile Drill, long, Ø 4.8mm, L 33mm, stainless steel	128
026.0009	BLT Tap, Ø 3.3mm, L 25mm, stainless steel/TAN	129
026.0010	BLT Tap, Ø 4.1mm, L 25mm, stainless steel/TAN	129
026.0011	BLT Tap, Ø 4.8mm, L 25mm, stainless steel/TAN	129
026.0016	Straumann® Planning Guide	77,98
026.0022*	BL Bone Profiler 1, Ø 5.2mm, L 23mm, stainless steel	130
026.0023*	BL Bone Profiler 2, Ø 6.6mm, L 23mm, stainless steel	130
026.0024*	BL Bone Profiler 3, Ø 6mm, L 23mm, stainless steel	130
026.00255*	NC Guiding Cylinder for BL Bone Profiler, Ø 3.25mm, L 12.5mm, TAN	130
026.00265*	RC Guiding Cylinder for BL Bone Profiler, Ø 3.7mm, L 12.5mm, TAN	130
026.0054	Needle Drill, short, Ø 1.6mm, L 33mm, stainless steel	127
026.0056	Needle Drill, long, Ø 1.6mm, L 41mm, stainless steel	127
026.0901	Bone Level implant distance indicator, L 100mm, Ti/stainless steel	126
026.0902	CrossFit® Plan Set, POM	110
026.2048	48h Explantation device for BL Implant NC, L 31.4mm, stainless steel	139
026.2200	BLT Drill, short, Ø 2.8mm, L 33mm, stainless steel	127
026.2201	BLT Drill, long, Ø 2.8mm, L 41mm, stainless steel	127
026.2303	BL/NNC Implant profile drill, Ø 3.3mm, short, stainless steel	128
026.2306	BL/NNC Implant profile drill, Ø 3.3mm, long, stainless steel	128
026.2310	BL/TE/NNC Tap f. adapter, f. Ø 3.3mm implants, L 23mm, PEEK/stainless steel	129
026.2503	BL/NNC Profile drill, Ø 3.3mm, guided, L 37mm, stainless steel	134
026.2510	BL/TE/NNC Tap f. handpiece, Ø 3.3mm, guided, L 42mm, stainless steel	134
026.2558	Release aid N for Loxim®	131
026.2801	NC Guiding cylinder for Ø 3.3mm, BL Implant explantation, stainless steel	138
026.4048	48h Explantation device for BL Implant RC, L 31.2mm, stainless steel	139
026.4200	BLT Drill, short, Ø 3.5mm, L 33mm, stainless steel	127
026.4201	BLT Drill, long, Ø 3.5mm, L 41mm, stainless steel	127
026.4303	Bone Level Implant profile drill, Ø 4.1mm, short, stainless steel	128
026.4306	Bone Level Implant profile drill, Ø 4.1mm, long, stainless steel	128
026.4310	BL/TE Tap f. adapter, f. Ø 4.1mm implants, L 23mm, polymer/stainless steel	129
026.4503	BL Profile drill, Ø 4.1mm, guided, L 37mm, stainless steel	134
026.4510	BL/TE Tap f. handpiece, Ø 4.1mm, guided, L 42mm, stainless steel	134
026.4558	Release aid R/W for Loxim®	131
026.4801	RC Guiding cylinder for Ø 4.1mm, BL Implant explantation, stainless steel	138
026.6200	BLT Drill, short, Ø 4.2mm, L 33mm, stainless steel	127
026.6201	BLT Drill, long, Ø 4.2mm, L 41mm, stainless steel	127
026.6303	Bone Level Implant profile drill, Ø 4.8mm, short, stainless steel	128
026.6306	Bone Level Implant profile drill, Ø 4.8mm, long, stainless steel	128
026.6310	BL/TE Tap f. adapter, f. Ø 4.8mm implants, L 23mm, polymer/stainless steel	129
026.6503	Bone Level profile drill, Ø 4.8mm, guided, L 37mm, stainless steel	134
026.6510	BL/TE Tap f. handpiece, Ø 4.8mm, guided, L 42mm, stainless steel	134
026.6801	RC Guiding cylinder for Ø 4.8mm, BL Implant explantation, stainless steel	138
033.030S	Standard Implant Ø 4.1mm, RN, SLActive® 6mm, Ti	16
033.031G	Standard Implant Ø 4.1mm RN, SLActive® 8mm, Ti guided	16
033.031S	Standard Implant Ø 4.1mm, RN, SLActive® 8mm, Ti	16
033.032G	Standard Implant Ø 4.1mm, RN, SLActive® 10mm, Ti guided	16
033.032S	Standard Implant Ø 4.1mm, RN, SLActive® 10mm, Ti	16
033.033G	Standard Implant Ø 4.1mm, RN, SLActive® 12mm, Ti guided	16
033.033S	Standard Implant Ø 4.1mm, RN, SLActive® 12mm, Ti	16
033.034S	Standard Implant Ø 4.1mm, RN, SLActive® 14mm, Ti	16
033.035S	Standard Implant Ø 4.1mm, RN, SLActive® 16mm, Ti	16
033.050S	Standard Plus Implant Ø 4.1mm, RN, SLActive® 6mm, Ti	16
033.051G	Standard Plus Implant Ø 4.1mm RN, SLActive® 8mm, Ti guided	16
033.051S	Standard Plus Implant Ø 4.1mm, RN, SLActive® 8mm, Ti	16
033.052G	Standard Plus Implant Ø 4.1mm RN, SLActive® 10mm, Ti guided	16
033.052S	Standard Plus Implant Ø 4.1mm RN SLActive® 10mm, Ti	16
033.053G	Standard Plus Implant Ø 4.1mm RN, SLActive® 12mm, Ti guided	16
033.053S	Standard Plus Implant, Ø 4.1mm RN SLActive® 12mm, Ti	16
033.054S	Standard Plus Implant, Ø 4.1mm RN SLActive® 14mm, Ti	16
033.131S	Standard Implant, Ø 3.3mm RN SLActive® 8mm, Ti	13
033.132S	Standard Implant, Ø 3.3mm RN SLActive® 10mm, Ti	13
033.133S	Standard Implant, Ø 3.3mm RN SLActive® 12mm, Ti	13

货号	产品名称	页数
033.134S	Standard Implant, Ø 3.3mm RN SLActive® 14mm, Ti	13
033.135S	Standard Implant, Ø 3.3mm RN SLActive® 16mm, Ti	13
033.151S	Standard Plus Implant, Ø 3.3mm RN SLActive® 8mm, Ti	13
033.152S	Standard Plus Implant, Ø 3.3mm RN SLActive® 10mm, Ti	13
033.153S	Standard Plus Implant, Ø 3.3mm RN SLActive® 12mm, Ti	13
033.154S	Standard Plus Implant, Ø 3.3mm RN SLActive® 14mm, Ti	13
033.230S	Standard Implant Ø 4.8mm RN, SLActive® 6mm, Ti	19
033.231S	Standard Implant Ø 4.8mm RN, SLActive® 8mm, Ti	19
033.232S	Standard Implant Ø 4.8mm RN, SLActive® 10mm, Ti	19
033.233S	Standard Implant Ø 4.8mm RN, SLActive® 12mm, Ti	19
033.234S	Standard Implant Ø 4.8mm, RN, SLActive® 14mm, Ti	19
033.250S	Standard Plus Implant Ø 4.8mm, RN, SLActive® 6mm, Ti	19
033.251G	Standard Plus Implant Ø 4.8mm RN, SLActive® 8mm, Ti guided	19
033.251S	Standard Plus Implant Ø 4.8mm, RN, SLActive® 8mm, Ti	19
033.252G	Standard Plus Implant Ø 4.8mm RN, SLActive® 10mm, Ti guided	19
033.252S	Standard Plus Implant Ø 4.8mm, RN, SLActive® 10mm, Ti	19
033.253G	Standard Plus Implant Ø 4.8mm RN, SLActive® 12mm, Ti guided	19
033.253S	Standard Plus Implant Ø 4.8mm, RN, SLActive® 12mm, Ti	19
033.254S	Standard Plus Implant Ø 4.8mm, RN, SLActive® 14mm, Ti	19
033.416S	Standard Plus Implant, Ø 3.3mm NNC, SLActive® 8mm, Roxolid®	12
033.417S	Standard Plus Implant, Ø 3.3mm NNC, SLActive® 10mm, Roxolid®	12
033.418S	Standard Plus Implant, Ø 3.3mm NNC, SLActive® 12mm, Roxolid®	12
033.419S	Standard Plus Implant, Ø 3.3mm NNC, SLActive® 14mm, Roxolid®	12
033.431G	Standard Implant Ø 3.3mm RN, SLActive® 8mm, Roxolid® guided	12
033.431S	Standard Implant, Ø 3.3 mm, RN, SLActive 8 mm, Roxolid	12
033.432G	Standard Implant Ø 3.3mm RN, SLActive® 10mm, Roxolid® guided	12
033.432S	Standard Implant, Ø 3.3 mm, RN, SLActive 10 mm, Roxolid	12
033.433G	Standard Implant Ø 3.3mm RN, SLActive® 12mm, Roxolid® guided	12
033.433S	Standard Implant, Ø 3.3 mm, RN, SLActive 12 mm, Roxolid	12
033.434S	Standard Implant, Ø 3.3 mm, RN, SLActive 14 mm, Roxolid	12
033.435S	Standard Implant, Ø 3.3 mm, RN, SLActive 16 mm, Roxolid	12
033.451G	Standard Plus Implant Ø 3.3mm RN, SLActive® 8mm, Roxolid® guided	12
033.452G	Standard Plus Implant Ø 3.3mm RN, SLActive® 10mm, Roxolid® guided	12
033.453G	Standard Plus Implant Ø 3.3mm RN, SLActive® 12mm, Roxolid® guided	12
033.501S	Standard Implant, Ø 3.3mm RN, SLActive® 8mm, Roxolid®, Loxim®	12
033.502S	Standard Implant, Ø 3.3mm RN, SLActive® 10mm, Roxolid®, Loxim®	12
033.503S	Standard Implant, Ø 3.3mm RN, SLActive® 12mm, Roxolid®, Loxim®	12
033.504S	Standard Implant, Ø 3.3mm RN, SLActive® 14mm, Roxolid®, Loxim®	12
033.505S	Standard Implant, Ø 3.3mm RN, SLActive® 16mm, Roxolid®, Loxim®	12
033.511S	Standard Plus Implant, Ø 3.3mm RN, SLActive® 8mm, Roxolid®, Loxim®	12
033.512S	Standard Plus Implant, Ø 3.3mm RN, SLActive® 10mm, Roxolid®, Loxim®	12
033.513S	Standard Plus Implant, Ø 3.3mm RN, SLActive® 12mm, Roxolid®, Loxim®	12
033.514S	Standard Plus Implant, Ø 3.3mm RN, SLActive® 14mm, Roxolid®, Loxim®	12
033.521S	Tapered Effect Implant, Ø 3.3mm RN, SLActive® 8mm, Roxolid®, Loxim®	12
033.522S	Tapered Effect Implant, Ø 3.3mm RN, SLActive® 10mm, Roxolid®, Loxim®	12
033.523S	Tapered Effect Implant, Ø 3.3mm RN, SLActive® 12mm, Roxolid®, Loxim®	12
033.524S	Tapered Effect Implant, Ø 3.3mm RN, SLActive® 14mm, Roxolid®, Loxim®	12
033.530S	Standard Implant, Ø 4.1mm RN, SLActive® 6mm, Roxolid®, Loxim®	15
033.531G	Standard Implant, Ø 4.1 RN, SLActive® 8mm, Roxolid®, guided	15
033.531S	Standard Implant, Ø 4.1mm RN, SLActive® 8mm, Roxolid®, Loxim®	15
033.532G	Standard Implant, Ø 4.1 RN, SLActive® 10mm, Roxolid®, guided	15
033.532S	Standard Implant, Ø 4.1mm RN, SLActive® 10mm, Roxolid®, Loxim®	15
033.533G	Standard Implant, Ø 4.1 RN, SLActive® 12mm, Roxolid®, guided	15
033.533S	Standard Implant, Ø 4.1mm RN, SLActive® 12mm, Roxolid®, Loxim®	15
033.534S	Standard Implant, Ø 4.1mm RN, SLActive® 14mm, Roxolid®, Loxim®	15
033.535S	Standard Implant, Ø 4.1mm RN, SLActive® 16mm, Roxolid®, Loxim®	15
033.560S	Standard Plus Implant, Ø 4.1mm RN, SLActive® 6mm, Roxolid®, Loxim®	15
033.561G	Standard Plus Implant, Ø 4.1 RN, SLActive® 8mm, Roxolid®, guided	15
033.561S	Standard Plus Implant, Ø 4.1mm RN, SLActive® 8mm, Roxolid®, Loxim®	15
033.562G	Standard Plus Implant, Ø 4.1 RN, SLActive® 10mm, Roxolid®, guided	15
033.562S	Standard Plus Implant, Ø 4.1mm RN, SLActive® 10mm, Roxolid®, Loxim®	15
033.563G	Standard Plus Implant, Ø 4.1 RN, SLActive® 12mm, Roxolid®, guided	15
033.563S	Standard Plus Implant, Ø 4.1mm RN, SLActive® 12mm, Roxolid®, Loxim®	15
033.564S	Standard Plus Implant, Ø 4.1mm RN, SLActive® 14mm, Roxolid®, Loxim®	15
033.571G	Tapered Effect Implant, Ø 4.1 RN, SLActive® 8mm, Roxolid®, guided	15
033.571S	Tapered Effect Implant, Ø 4.1mm RN, SLActive® 8mm, Roxolid®, Loxim®	15
033.572G	Tapered Effect Implant, Ø 4.1 RN, SLActive® 10mm, Roxolid®, guided	15

货号	产品名称	页数
033.572S	Tapered Effect Implant, Ø 4.1mm RN, SLActive® 10mm, Roxolid®, Loxim®	15
033.573G	Tapered Effect Implant, Ø 4.1 RN, SLActive® 12mm, Roxolid®, guided	15
033.573S	Tapered Effect Implant, Ø 4.1mm RN, SLActive® 12mm, Roxolid®, Loxim®	15
033.574S	Tapered Effect Implant, Ø 4.1mm RN, SLActive® 14mm, Roxolid®, Loxim®	15
033.580S	Standard Implant, Ø 4.8mm RN, SLActive® 6mm, Roxolid®, Loxim®	18
033.581S	Standard Implant, Ø 4.8mm RN, SLActive® 8mm, Roxolid®, Loxim®	18
033.582S	Standard Implant, Ø 4.8mm RN, SLActive® 10mm, Roxolid®, Loxim®	18
033.583S	Standard Implant, Ø 4.8mm RN, SLActive® 12mm, Roxolid®, Loxim®	18
033.584S	Standard Implant, Ø 4.8mm RN, SLActive® 14mm, Roxolid®, Loxim®	18
033.590S	Standard Plus Implant, Ø 4.8mm RN, SLActive® 6mm, Roxolid®, Loxim®	18
033.591G	Standard Plus Implant, Ø 4.8 RN, SLActive® 8mm, Roxolid®, guided	18
033.591S	Standard Plus Implant, Ø 4.8mm RN, SLActive® 8mm, Roxolid®, Loxim®	18
033.592G	Standard Plus Implant, Ø 4.8 RN, SLActive® 10mm, Roxolid®, guided	18
033.592S	Standard Plus Implant, Ø 4.8mm RN, SLActive® 10mm, Roxolid®, Loxim®	18
033.593G	Standard Plus Implant, Ø 4.8 RN, SLActive® 12mm, Roxolid®, guided	18
033.593S	Standard Plus Implant, Ø 4.8mm RN, SLActive® 12mm, Roxolid®, Loxim®	18
033.594S	Standard Plus Implant, Ø 4.8mm RN, SLActive® 14mm, Roxolid®, Loxim®	18
033.600S	Standard Implant, Ø 4.8mm WN, SLActive® 6mm, Roxolid®, Loxim®	18
033.601S	Standard Implant, Ø 4.8mm WN, SLActive® 8mm, Roxolid®, Loxim®	18
033.602S	Standard Implant, Ø 4.8mm WN, SLActive® 10mm, Roxolid®, Loxim®	18
033.603S	Standard Implant, Ø 4.8mm WN, SLActive® 12mm, Roxolid®, Loxim®	18
033.610S	Standard Plus Implant, Ø 4.8mm WN, SLActive® 6mm, Roxolid®, Loxim®	18
033.611S	Standard Plus Implant, Ø 4.8mm WN, SLActive® 8mm, Roxolid®, Loxim®	18
033.612S	Standard Plus Implant, Ø 4.8mm WN, SLActive® 10mm, Roxolid®, Loxim®	18
033.613S	Standard Plus Implant, Ø 4.8mm WN, SLActive® 12mm, Roxolid®, Loxim®	18
033.622S	Tapered Effect Implant, Ø 4.8mm WN, SLActive® 10mm, Roxolid®, Loxim®	18
033.623S	Tapered Effect Implant, Ø 4.8mm WN, SLActive® 12mm, Roxolid®, Loxim®	18
033.624S	Tapered Effect Implant, Ø 4.8mm WN, SLActive® 14mm, Roxolid®, Loxim®	18
033.630S	Standard Implant Ø 4.8mm, WN, SLActive® 6mm, Ti	19
033.631S	Standard Implant Ø 4.8mm, WN, SLActive® 8mm, Ti	19
033.632S	Standard Implant Ø 4.8mm, WN, SLActive® 10mm, Ti	19
033.633S	Standard Implant Ø 4.8mm, WN, SLActive® 12mm, Ti	19
033.650S	Standard Plus Implant Ø 4.8mm, WN, SLActive® 6mm, Ti	19
033.651S	Standard Plus Implant Ø 4.8mm, WN, SLActive® 8mm, Ti	19
033.652S	Standard Plus Implant Ø 4.8mm, WN, SLActive® 10mm, Ti	19
033.653S	Standard Plus Implant Ø 4.8mm, WN, SLActive® 12mm, Ti	19
033.712S	Tapered Effect Implant Ø 4.8mm WN, SLActive® 10mm, Ti	19
033.713S	Tapered Effect Implant Ø 4.8mm WN, SLActive® 12mm, Ti	19
033.714S	Tapered Effect Implant Ø 4.8mm WN, SLActive® 14mm, Ti	19
033.721S	Tapered Effect Implant Ø 3.3mm RN, SLActive® 8mm, Ti	13
033.722S	Tapered Effect Implant Ø 3.3mm RN, SLActive® 10mm, Ti	13
033.723S	Tapered Effect Implant Ø 3.3mm RN, SLActive® 12mm, Ti	13
033.724S	Tapered Effect Implant Ø 3.3mm RN, SLActive® 14mm, Ti	13
033.761G	Tapered Effect Implant Ø 4.1mm RN, SLActive® 8mm, Ti guided	16
033.761S	Tapered Effect Implant Ø 4.1mm RN, SLActive® 8mm, Ti	16
033.762G	Tapered Effect Implant Ø 4.1mm RN, SLActive® 10mm, Ti guided	16
033.762S	Tapered Effect Implant Ø 4.1mm RN, SLActive® 10mm, Ti	16
033.763G	Tapered Effect Implant Ø 4.1mm RN, SLActive® 12mm, Ti guided	16
033.763S	Tapered Effect Implant Ø 4.1mm RN, SLActive® 12mm, Ti	16
033.764S	Tapered Effect Implant Ø 4.1mm RN, SLActive® 14mm, Ti	16
034.001	Surgical cassette, for guided surgery	136
034.005	Connector for ratchet, L 14mm, stainless steel	136
034.006	Stop key for guided implants, S/SP/TE/BL, stainless steel	136
034.010	Mucosa punch Ø 3.4mm, guided, L 30mm, stainless steel	133
034.011	Mucosa punch Ø 4.0mm, guided, L 30mm, stainless steel	133
034.012	Mucosa punch Ø 4.7mm, guided, L 30mm, stainless steel	133
034.053V4	T-sleeve, Ø 5mm, H 5mm, guided, stainless steel	133
034.055V4	T-sleeve, Ø 2.8mm, H 6mm, guided, stainless steel	133
034.123	Pilot drill, Ø 2.2mm, short, 16mm, guided, L 32mm, stainless steel	133
034.126	Pilot drill, Ø 2.2mm, long, 20mm, guided, L 36mm, stainless steel	133
034.129	Pilot drill, Ø 2.2mm, x-long, 24mm, guided, L 40mm, stainless steel	133
034.150	Drill handle, Ø 2.2mm, 1mm/3mm stop, L 104mm, stainless steel	136
034.215	Milling cutter, Ø 2.8mm, guided, L 32.5mm, stainless steel	133
034.223	Twist drill PRO, Ø 2.8mm, short, 16mm, guided, L 32mm, stainless steel	133
034.226	Twist drill PRO, Ø 2.8mm, long, 20mm, guided, L 36mm, stainless steel	133
034.229	Twist drill PRO, Ø 2.8mm, x-long, 24mm, guided, L 40mm, stainless steel	133
034.235	SP Profile drill, Ø 3.3mm, RN, guided, L 35mm, stainless steel	134

货号	产品名称	页数
034.237	TE Profile drill, Ø 3.3mm, RN, guided, L 37mm, stainless steel	134
034.245	S/SP Tap f. handpiece, Ø 3.3mm, guided, L 42mm, stainless steel	134
034.250	Drill handle, Ø 2.8mm, 1mm/3mm stop, L 104mm, stainless steel	136
034.255	Template fixation pin, Ø 5.0mm/2.8mm, guided, L 24mm, stainless steel	136
034.256	Template fixation pin, Ø 2.8mm/2.8mm, guided, L 25mm, stainless steel	136
034.257	BLT Pilot Drill, guided, Ø 2.2mm, short, stainless steel	133
034.258	BLT Pilot Drill, guided, Ø 2.2mm, medium, stainless steel	133
034.259	BLT Pilot Drill, guided, Ø 2.2mm, long, stainless steel	133
034.260	BLT Drill, guided, Ø 2.8mm, short, stainless steel	134
034.261	BLT Drill, guided, Ø 2.8mm, medium, stainless steel	134
034.262	BLT Drill, guided, Ø 2.8mm, long, stainless steel	134
034.263	BLT Drill, guided, Ø 3.5mm, short, stainless steel	134
034.264	BLT Drill, guided, Ø 3.5mm, medium, stainless steel	134
034.265	BLT Drill, guided, Ø 3.5mm, long, stainless steel	134
034.266	BLT Drill, guided, Ø 4.2mm, short, stainless steel	134
034.267	BLT Drill, guided, Ø 4.2mm, medium, stainless steel	134
034.268	BLT Drill, guided, Ø 4.2mm, long, stainless steel	134
034.269	BLT Profile Drill, guided, Ø 3.3mm, L 36mm, stainless steel	134
034.270	BLT Profile Drill, guided, Ø 4.1mm, L 36mm, stainless steel	134
034.271	BLT Profile Drill, guided, Ø 4.8mm, L 36mm, stainless steel	134
034.272	BLT Tap, guided, Ø 3.3mm, L 42mm, stainless steel	134
034.273	BLT Tap, guided, Ø 4.1mm, L 42mm, stainless steel	134
034.274	BLT Tap, guided, Ø 4.8mm, L 42mm, stainless steel	134
034.281	Straumann® Basic Surgical Cassette, guided	136
034.415	Milling cutter, Ø 3.5mm, guided, L 32.5mm, stainless steel	133
034.423	Twist drill PRO, Ø 3.5mm, short, 16mm, guided, L 32mm, stainless steel	133
034.426	Twist drill PRO, Ø 3.5mm, long, 20mm, guided, L 36mm, stainless steel	133
034.429	Twist drill PRO, Ø 3.5mm,x-long, 24mm, guided, L 40mm, stainless steel	133
034.435	SP Profile drill, Ø 4.1mm, RN, guided, L 35mm, stainless steel	134
034.437	TE Profile drill, Ø 4.1mm, RN, guided, L 37mm, stainless steel	134
034.445	S/SP Tap f. handpiece, Ø 4.1mm, guided, L 42mm, stainless steel	134
034.450	Drill handle, Ø 3.5mm, 1mm/3mm stop, L 104mm, stainless steel	136
034.455	Template fixation pin, Ø 5.0mm/3.5mm, guided, L 24mm, stainless steel	136
034.615	Milling cutter, Ø 4.2mm, guided, L 32.5mm, stainless steel	133
034.623	Twist drill PRO, Ø 4.2mm, short, 16mm, guided, L 32mm, stainless steel	133
034.626	Twist drill PRO, Ø 4.2mm, long, 20mm, guided, L 36mm, stainless steel	133
034.629	Twist drill PRO, Ø 4.2mm,x-long, 24mm, guided, L 40mm, stainless steel	133
034.645	S/SP Tap f. handpiece, Ø 4.8mm, guided, L 42mm, stainless steel	134
034.650	Drill handle, Ø 4.2mm, 1mm/3mm stop, L 104mm, stainless steel	136
034.655	Template fixation pin, Ø 5mm/4.2mm, guided, L 24mm, stainless steel	136
034.750	C-handle H2, 6mm stop height, L 89mm, stainless steel	136
034.751	C-handle H4, 4mm stop height, L 89mm, stainless steel	136
034.752	C-handle H6, 2mm stop height, L 89mm, stainless steel	136
040.165	Straumann® surgical cassette	131
040.175	Ultrasonic Cleaning Cassette, 150 x 65 x 25 mm, plastic/silicone	143
040.185	Prosthetic cassette for instruments, complete, 224 x 115 x 20 mm, anodized aluminium	148
040.500	Osteotome instrument kit, for bone condensation, stainless steel	141
040.501	Osteotome instrument kit, for sinus floor elevation, stainless steel	142
040.550	Cassette without content, for bone condensation	141
040.560	Cassette without content, for sinus floor elevation	142
040.564	Straumann® Basic Surgical Cassette	131
041.022	Prosthetic cassette for instruments, 224 x 115 x 20 mm, anodized aluminium	148
041.032	Bone Block Fixation cassette for sterilization and storage of drills,mini-screws, screwdriver, 130mmx118mmx25mm	140
041.211	Tray for cassette 041.022, 218x109x16 mm, Al	148
041.412	Calibration cup, yellow, Ø 80mm, anodized aluminium	148
041.440	Dish for tray 041.211, 56x31x11 mm, Al	148
041.500	Storage System for implants (SLA & SLActive), 280 x 226 x 248 mm, plastic	143
041.501	Storage System for Prosthetics, 280 x 226 x 248mm, plastic	147
041.520	Drawer for storage system, high, 282 x 220 x 47 mm, plastic	143
041.530V5	Divider set, flat, for implants, 184 x 44 mm, plastic	143
041.612	Divider, large, for tray 041.211, 47.5 x 9.5 x 1mm, anodized aluminium	148
041.615	Divider, small, for dish 041.440, 27.5 x 9.5 x 1mm, anodized aluminium	148
043.030S	Standard Implant, Ø 4.1mm, RN, SLA 6mm, Ti	16
043.031S	Standard Implant, Ø 4.1mm, RN, SLA 8mm, Ti	16
043.032S	Standard Implant, Ø 4.1mm, RN, SLA 10mm, Ti	16
043.033S	Standard Implant, Ø 4.1mm, RN, SLA 12mm, Ti	16
043.034S	Standard Implant, Ø 4.1mm, RN, SLA 14mm, Ti	16

货号	产品名称	页数
043.035S	Standard Implant, Ø 4.1mm, RN, SLA 16mm, Ti	16
043.050S	Standard Plus Implant, Ø 4.1mm, RN, SLA 6mm, Ti	16
043.051S	Standard Plus Implant, Ø 4.1mm, RN, SLA 8mm, Ti	16
043.052S	Standard Plus Implant, Ø 4.1mm, RN, SLA 10mm, Ti	16
043.053S	Standard Plus Implant, Ø 4.1mm, RN, SLA 12mm, Ti	16
043.054S	Standard Plus Implant, Ø 4.1mm, RN, SLA 14mm, Ti	16
043.131S	Standard Implant, Ø 3.3mm, RN, SLA 8mm, Ti	13
043.132S	Standard Implant, Ø 3.3mm, RN, SLA 10mm, Ti	13
043.133S	Standard Implant, Ø 3.3mm, RN, SLA 12mm, Ti	13
043.134S	Standard Implant, Ø 3.3mm, RN, SLA 14mm, Ti	13
043.135S	Standard Implant, Ø 3.3mm, RN, SLA 16mm, Ti	13
043.151S	Standard Plus Implant, Ø 3.3mm, RN, SLA 8mm, Ti	13
043.152S	Standard Plus Implant, Ø 3.3mm, RN, SLA 10mm, Ti	13
043.153S	Standard Plus Implant, Ø 3.3mm, RN, SLA 12mm, Ti	13
043.154S	Standard Plus Implant, Ø 3.3mm, RN, SLA 14mm, Ti	13
043.230S	Standard Implant, Ø 4.8mm, RN, SLA 6mm, Ti	20
043.231S	Standard Implant, Ø 4.8mm, RN, SLA 8mm, Ti	20
043.232S	Standard Implant, Ø 4.8mm, RN, SLA 10mm, Ti	20
043.233S	Standard Implant, Ø 4.8mm, RN, SLA 12mm, Ti	20
043.234S	Standard Implant, Ø 4.8mm, RN, SLA 14mm, Ti	20
043.250S	Standard Plus Implant, Ø 4.8mm, RN, SLA 6mm, Ti	20
043.251S	Standard Plus Implant, Ø 4.8mm, RN, SLA 8mm, Ti	20
043.252S	Standard Plus Implant, Ø 4.8mm, RN, SLA 10mm, Ti	20
043.253S	Standard Plus Implant, Ø 4.8mm, RN, SLA 12mm, Ti	20
043.254S	Standard Plus Implant, Ø 4.8mm, RN, SLA 14mm, Ti	20
043.630S	Standard Implant, Ø 4.8mm, WN, SLA 6mm, Ti	20
043.631S	Standard Implant, Ø 4.8mm, WN, SLA 8mm, Ti	20
043.632S	Standard Implant, Ø 4.8mm, WN, SLA 10mm, Ti	20
043.633S	Standard Implant, Ø 4.8mm, WN, SLA 12mm, Ti	20
043.650S	Standard Plus Implant, Ø 4.8mm, WN, SLA 6mm, Ti	20
043.651S	Standard Plus Implant, Ø 4.8mm, WN, SLA 8mm, Ti	20
043.652S	Standard Plus Implant, Ø 4.8mm, WN, SLA 10mm, Ti	20
043.653S	Standard Plus Implant, Ø 4.8mm, WN, SLA 12mm, Ti	20
043.712S	Tapered Effect Implant, Ø 4.8mm, WN, SLA 10mm, Ti	20
043.713S	Tapered Effect Implant, Ø 4.8mm, WN, SLA 12mm, Ti	20
043.714S	Tapered Effect Implant, Ø 4.8mm, WN, SLA 14mm, Ti	20
043.721S	Tapered Effect Implant, Ø 3.3mm, RN, SLA 8mm, Ti	13
043.722S	Tapered Effect Implant, Ø 3.3mm, RN, SLA 10mm, Ti	13
043.723S	Tapered Effect Implant, Ø 3.3mm, RN, SLA 12mm, Ti	13
043.724S	Tapered Effect Implant, Ø 3.3mm, RN, SLA 14mm, Ti	13
043.761S	Tapered Effect Implant, Ø 4.1mm, RN, SLA 8mm, Ti	16
043.762S	Tapered Effect Implant, Ø 4.1mm, RN, SLA 10mm, Ti	16
043.763S	Tapered Effect Implant, Ø 4.1mm, RN, SLA 12mm, Ti	16
043.764S	Tapered Effect Implant, Ø 4.1mm, RN, SLA 14mm, Ti	16
044.003	Round bur, Ø 2.3mm, stainless steel	127
044.004	Round bur, Ø 3.1mm, stainless steel	127
044.022	Round bur, Ø 1.4mm, stainless steel	127
044.084	SP Profile drill, WN, short, Ø 4.8mm, L 23.8mm, stainless steel	128
044.085	SP Profile drill, WN, long, Ø 4.8mm, L 32.8mm, stainless steel	128
044.086	SP Profile drill, RN, short, Ø 3.3mm, L 25mm, stainless steel	128
044.087	SP Profile drill, RN, long, Ø 3.3mm, L 34mm, stainless steel	128
044.088	SP Profile drill, RN, short, Ø 4.1mm, L 24mm, stainless steel	128
044.089	SP Profile drill, RN, long, Ø 4.1mm, L 33mm, stainless steel	128
044.210	Pilot drill 1, short, Ø 2.2mm, L 33mm, stainless steel	127
044.211	Pilot drill 1, long, Ø 2.2mm, L 41mm, stainless steel	127
044.214	Pilot drill 2, short, Ø 2.8mm, L 33mm, stainless steel	127
044.215	Pilot drill 2, long, Ø 2.8mm, L 41mm, stainless steel	127
044.250	Twist drill PRO, short, Ø 3.5mm, L 33mm, stainless steel	127
044.251	Twist drill PRO, long, Ø 3.5mm, L 41mm, stainless steel	127
044.254	Twist drill PRO, short, Ø 4.2mm, L 33mm, stainless steel	127
044.255	Twist drill PRO, long, Ø 4.2mm, L 41mm, stainless steel	127
044.304	Bone profiler for RN, L 28mm, stainless steel, needed guiding cylinder: 049.082	130
044.305	Bone profiler for WN, L 25mm, stainless steel, needed guiding cylinder: 049.082	130
044.330	Explantation drill, short, for Ø 3.3mm, Ø 3.5mm, L 32mm, stainless steel	137
044.331	Explantation drill, medium, for Ø 3.3mm, Ø 3.5mm, L 37.5mm, stainless steel	137
044.332	Explantation drill, long, for Ø 3.3mm, Ø 3.5mm, L 40.5mm, stainless steel	137
044.340	Explantation drill, short, for Ø 4.1mm, L 32mm, stainless steel	137

货号	产品名称	页数
044.341	Explantation drill, medium, for Ø 4.1mm, L 37.5mm, stainless steel	137
044.342	Explantation drill, long, for Ø 4.1mm, L 40.5mm, stainless steel	137
044.344	Explantation drill, short, for Ø 4.8mm, L 32mm, stainless steel	137
044.345	Explantation drill, medium, for Ø 4.8mm, L 37.5mm, stainless steel	137
044.346	Explantation drill, long, for Ø 4.8mm, L 40.5mm, stainless steel	137
044.570	Tap for transversal screw retention, M2 x 0.25mm, stainless steel	146
044.575	S/SP Tap for adapter, Ø 3.3mm, L 23mm, PEEK/stainless steel	129
044.577	S/SP Tap for adapter, Ø 4.1mm, L 23mm, PEEK/stainless steel	129
044.579	S/SP Tap for adapter, Ø 4.8mm, L 23mm, PEEK/stainless steel	129
044.701	TE Profile drill, RN, short, Ø 3.3mm, L 28mm, stainless steel	128
044.703	TE Profile drill, WN, short, Ø 4.8mm, L 28mm, stainless steel	128
044.705	TE Profile drill, RN, short, Ø 4.1mm, L 28mm, stainless steel	128
044.708	TE Profile drill, RN, long, Ø 3.3mm, L 37mm, stainless steel	128
044.710	TE Profile drill, WN, long, Ø 4.8mm, L 37mm, stainless steel	128
044.712	TE Profile drill, RN, long, Ø 4.1mm, L 37mm, stainless steel	128
045.111V4	Cleaning brush for ratchet, L 100mm, Ø 4.5mm, nylon/stainless steel	131,144
045.921	Light curette (pack of 5 pcs)	132
046.049	Torque control device for ratchet, stainless steel	131,144
046.057	Tweezers without titanium plasma-coated tips, L 105mm, stainless steel	131
046.064	Holding key, L 85mm, stainless steel	131
046.067	Solid abut.driver, short, f. RN solid abut.6°, L 13mm, stainless steel	50,144
046.068	Solid abut.driver, long, f. RN solid abut.6°, L 19mm, stainless steel	50,144
046.069	Retentive anchor driver, L 19mm, stainless steel	144
046.085	Screwdriver handle for 046.410/411/412, L 109mm, stainless steel	144
046.095	Tweezers with titanium plasma-coated tips, L 105mm, Ti/stainless steel	131
046.108	Service instrument for ratchet, L 50mm, stainless steel	131,144
046.110	Instrument tweezers, L 140mm, stainless steel	131
046.119	Ratchet, including service instrument, L 84mm, stainless steel	131,144
046.147	Diagnostic T, L 65mm, W 33mm, Ti	126
046.148	Implant distance indicator, L 100mm, Ti/stainless steel	126
046.150	Activator set, for all bar matrices, three-part, stainless steel	85,86,104,105,146
046.151	Deactivator for Dolder® Bar matrices, mini, L 66mm, Grilon BS/brass	85,86,104,105,146
046.152	Deactivator for Dolder® Bar matrices, regular, L 66mm, Grilon BS/brass	85,86,104,105,146
046.154	Screwdriver, L 37mm, stainless steel	144
046.237	NNC Reamer for 45° shoulder, L 17mm, stainless steel / Ti-N	145
046.239	Analog holder, length 205.0 mm, steel, anodized aluminium	83,84,87,103,106,146
046.240	Handle for reamers and Tap for TS, L 65mm, steel/brass	145
046.241	Guide pin for reamer, for synOcta®, L 50mm, stainless steel	145
046.242	Guide pin for reamer, for 048.540-542, L 50mm, stainless steel	145
046.243	Reamer for 45° shoulder, Ø 8.0mm, Ti-N/stainless steel	145
046.244	Guide pin for reamer, for 048.545-546, L 50mm, stainless steel	145
046.245	Polishing protector for RN synOcta® Coping, L 15 mm, stainless steel	48
046.246	Guide pin for reamer, for 048.605, L 50mm, stainless steel	145
046.247	Guide pin for reamer, for 048.606, L 50mm, stainless steel	145
046.249	NNC Guide pin for reamer 046.237, for Solid abutments, L 53mm, stainl. steel	145
046.400	SCS Screwdriver, for ratchet, x-short, L 15mm, stainless steel	62,132,144
046.401	SCS Screwdriver, for ratchet, short, L 21mm, stainless steel	62,132,144
046.402	SCS Screwdriver, for ratchet, long, L 27mm, stainless steel	62,132,144
046.410	SCS Screwdriver for handpiece, x-short, L 20mm, stainless steel	132,144
046.411	SCS Screwdriver for handpiece, short, L 26mm, stainless steel	132,144
046.412	SCS Screwdriver for handpiece, long, L 32mm, stainless steel	132,144
046.413V4	LOCATOR® abutment holder sleeve	41,54,65,90,109
046.415	LOCATOR® core tool, stainless steel with Au	41,54,65,90,109,144
046.416	LOCATOR® screwdriver for ratchet, short, length 15.0 mm, stainless steel	41,54,65,90,109,144
046.417	LOCATOR® screwdriver for ratchet, length 21.0 mm	41,54,65,90,109,144
046.421	Hexagonal screwdriver, L 30mm, stainless steel	132
046.450	Depth gauge Ø 3.5mm, L 26mm, Ti	130
046.451	Depth gauge Ø 4.2mm, L 26mm, Ti for Ø 4.8mm implants	130
046.455	Depth gauge with distance indicator Ø 2.2/2.8mm, L 27mm, Ti	130
046.458	Alignment pin straight, Ø 2.2 mm, length 28.0 mm, titanium	130
046.460	Adapter for ratchet, x-short, L 11mm, stainless steel	132
046.461	Adapter for ratchet, short, L 18mm, stainless steel	132
046.462	Adapter for ratchet, long, L 28mm, stainless steel	132
046.470	Adapter for handpiece, x-short, L 19mm, stainless steel	132
046.471	Adapter for handpiece, short, L 26mm, stainless steel	132
046.472	Adapter for handpiece, long, L 34mm, stainless steel	132
046.500	Osteotome straight, for bone condensation, Ø 2.2mm, stainless steel	141

货号	产品名称	页数
046.501	Osteotome straight, for bone condensation, Ø 2.8mm, stainless steel	141
046.502	Osteotome straight, for bone condensation, Ø 3.5mm, stainless steel	141
046.503	Osteotome straight, for bone condensation, Ø 4.2mm, stainless steel	141
046.510	Osteotome angled, for bone condensation, Ø 2.2mm, stainless steel	141
046.511	Osteotome angled, for bone condensation, Ø 2.8mm, stainless steel	141
046.512	Osteotome angled, for bone condensation, Ø 3.5mm, stainless steel	141
046.513	Osteotome angled, for bone condensation, Ø 4.2mm, stainless steel	141
046.521	Osteotome straight, for sinus floor elevation, Ø 2.8mm, stainless steel	142
046.522	Osteotome straight, for sinus floor elevation, Ø 3.5mm, stainless steel	142
046.523	Osteotome straight, for sinus floor elevation, Ø 4.2mm, stainless steel	142
046.531	Osteotome angled, for sinus floor elevation, Ø 2.8mm, stainless steel	142
046.532	Osteotome angled, for sinus floor elevation, Ø 3.5mm, stainless steel	142
046.533	Osteotome angled, for sinus floor elevation, Ø 4.2mm, stainless steel	142
046.550	Depth stop for osteotomes, Ø 2.2mm, Ti/stainless steel	141
046.551	Depth stop for osteotomes, Ø 2.8mm, Ti/stainless steel	141,142
046.552	Depth stop for osteotomes, Ø 3.5mm, Ti/stainless steel	141,142
046.553	Depth stop for osteotomes, Ø 4.2mm, Ti/stainless steel	141,142
046.703	Alignment Pin, Ø 2.2mm, L 27mm, Ti	129
046.704	Depth Gauge, with distance indicator, Ø 2.2/2.8mm, L 27mm, Ti	129
046.705	Depth Gauge, Ø 2.8mm, L 27mm, Ti	129
046.706	Depth Gauge, Ø 3.5mm, L 27mm, Ti	129
046.707	Depth Gauge, Ø 4.2 mm, L 27mm, Ti	129
048.000V4	Transfer aid for RN synOcta®, angled, short, H 5.5 mm, PC	46
048.002V4	Transfer aid for RN synOcta®, angled, long, H 4 mm, PC	46
048.010	RN synOcta® impression cap, screw-retained, with integral guide screw, H 11mm, Al/Ti	44
048.013	WN Impression Cap, H 8mm, POM	56,61
048.016	NN Impression cap with integral guide screw, L 17.5 mm, Al/TAN	117
048.017V4	RN Impression Cap, H 8mm, POM	44,49
048.028	RN Healing cap, labial bevel, small H 2mm, Ti	14,17,21
048.029	RN Healing cap, labial bevel, large H 3.5mm, Ti	14,17,21
048.030	WN healing cap, with labial bevel, height 2.0 mm, titanium	21
048.032	Transfer aid for WN synOcta®, angled, H 5 mm, PC	58
048.033	RN healing cap, height 2.0 mm, titanium	14,17,21
048.034	RN healing cap, height 3.0 mm, titanium	14,17,21
048.037	RN healing cap height 4.5 mm, titanium	14,17,21
048.038	WN Healing cap, H 2mm, Ti	21
048.039	WN Healing cap, H 3mm, Ti	21
048.047V4	RN Protective Cap for Solid Abutment 048.540, H 5.8mm, PEEK	49
048.048V4	RN Protective Cap for Solid Abutment 048.541, H 7.3mm, PEEK	49
048.049V4	RN Protective Cap for Solid Abutment 048.542, H 8.8mm, PEEK	49
048.050	NN Protective cap with integr.occlusal screw, Ø 4mm, H 3.4mm, PEEK/TAN	117
048.051	WN Protective Cap for Solid Abutment 048.545, H 6mm, PEEK	61
048.052	WN Protective Cap for Solid Abutment 048.546, H 7.5mm, PEEK	61
048.053	WN healing cap, height 4.5 mm, titanium	21
048.054V4	Transfer aid for 048.606, H 6,5 mm, PC	57
048.059V4	Transfer aid for 048.605, H 6,5 mm, PC	45
048.060V4	RN Positioning Cylinder for Solid Abutment, for 048.540, H 10.2mm, POM	49
048.061V4	RN Positioning Cylinder for Solid Abutment, for 048.541, H 10.2mm, POM	49
048.062V4	RN Positioning Cylinder for Solid Abutment, for 048.542, H 10.2mm, POM	49
048.065	WN Positioning Cylinder for Solid Abutment, for 048.545, H 10mm, POM	61
048.066	WN Positioning Cylinder for Solid Abutment, for 048.546, H 10mm, POM	61
048.070V4	RN synOcta® Positioning Cylinder, H 12mm, POM	44
048.071	NNC Healing cap, Ø 4mm, H 3mm, Ti	14
048.074	NNC Healing cap, Ø 4mm, H 4.5mm, Ti	14
048.082	NNC Healing cap, labial bevel, Ø 4mm, H 2mm, Ti	14
048.090	RN synOcta® Impression post, with integral guide screw, H 21.0 mm, Al/Ti	44
048.091	WN synOcta® impression cap, screw-retained, with integral guide screw, H 11mm, Al	56
048.095	WN synOcta® Positioning Cylinder, H 12mm, POM	56
048.108	Analog for RN Bar constructions, w. synOcta® 1.5, L 11.5mm, stainless steel	44
048.117V4	Analog for individualized RN Solid Abutments, 2 pieces, PPS/brass	49
048.124	Analog for RN synOcta®, L 12mm, grey w. red stripe, stainless steel	44
048.127	NNC Implant analog, L 12.5mm, Ti	34
048.128	NNC Impression post, for open tray, H 10.1mm, TAN	34
048.130	NN Analog, L 14mm, stainless steel	117
048.134	NNC Impression post, closed tray, H 10.9mm, TAN/POM	34
048.137	NNC Wax-up sleeve, L 12.75 mm, PEEK	36
048.137V4	NNC Wax-up sleeve, L 12.75 mm, PEEK	36

货号	产品名称	页数
048.140	Shoulder Analog for WN Solid Abutment, L 10mm, Al	61
048.155	NNC Solid abutment driver, short, L 13mm, stainless steel	38,144
048.156	NNC Solid abutment driver, long, L 19mm, stainless steel	38,144
048.160	RN Analog for Solid Abutment 048.540, L 14.3mm, yellow, Al	49
048.161	RN Analog for Solid Abutment 048.541, L 15.8mm, grey, Al	49
048.162	RN Analog for Solid Abutment 048.542, L 17.3mm, blue, Al	49
048.165	WN Analog for Solid Abutment 048.545, L 16mm, green, Al	61
048.166	WN Analog for Solid Abutment 048.546, L 17.5mm, brown, Al	61
048.168	CARES® RN Mono Scanbody, Ø 4.8 mm, H 10 mm, PEEK, TAN	44
048.169	CARES® WN Mono Scanbody, Ø 6.5 mm, H 10 mm, PEEK, TAN	56
048.171	WN synOcta® Analog, L 12mm, stainless steel	56
048.173	CARES® NNC Mono Scanbody, Ø 3.5 mm, H 10 mm, PEEK, TAN	34
048.175	RN LOCATOR® abutment, height 1.0 mm, titanium alloy, TiN coating	52
048.176	RN LOCATOR® abutment, height 2.0 mm, titanium alloy, TiN coating	52
048.177	RN LOCATOR® abutment, height 3.0 mm, titanium alloy, TiN coating	52
048.178	RN LOCATOR® abutment, height 4.0 mm, titanium alloy, TiN coating	52
048.179	RN LOCATOR® abutment, height 5.0 mm, titanium alloy, TiN coating	52
048.180	RN LOCATOR® abutment, height 6.0 mm, titanium alloy, TiN coating	52
048.181V4	LOCATOR® zero retention ext range, gray	40,53,64,89,108
048.182V2	LOCATOR® Male processing package, extended range, titanium/nylon	40,53,64,89,108
048.183	WN LOCATOR®, solid abutment, tissue cuff height 1.0 mm, titanium alloy, TiN coating	63
048.184	WN LOCATOR®, solid abutment, tissue cuff, height 2.0mm, titanium alloy, TiN coating	63
048.185	WN LOCATOR®, solid abutment, tissue cuff height 3.0 mm, titanium alloy, TiN coating	63
048.186	WN LOCATOR®, abutment, tissue cuff height 4.0 mm, titanium alloy, TiNcoating	63
048.187	WN LOCATOR®, abutment, tissue cuff height 5.0 mm, titanium alloy, TiNcoating	63
048.188V4	LOCATOR® Replacement male, extended range, orange, height 1.7 mm, nylon	40,53,64,89,108
048.189V2	LOCATOR® male processing package, titanium/nylon	40,53,64,89,108
048.190V4	LOCATOR® replacement male, clear, 0°–10°, 5 lbs, 2.27 kg, H 1.7mm, nylon	40,53,64,89,108
048.191V4	LOCATOR® replacement male, light retention, pink, height 1.7 mm, nylon	40,53,64,89,108
048.192V4	LOCATOR® replacement male, extra light retention, blue, height 1.7 mm, nylon	40,53,64,89,108
048.193V4	LOCATOR® replacement male, extended range, green, height 1.7 mm, nylon	40,53,64,89,108
048.194V4	LOCATOR® replacement male, extra light retention, ext. range, red, height 1.7 mm, nylon	40,53,64,89,108
048.195V4	LOCATOR® black processing male	41,54,65,90,109
048.196V20	LOCATOR® block-out spacer, silicon	41,54,65,90,109
048.197V4	LOCATOR® impression coping, aluminium housing with black male LDPE	39,52,63,88,107
048.198V4	LOCATOR® female analog, ø 5.0 mm, aluminium	52,63
048.199V4	LOCATOR® parallel post, LDPE	41,54,65,90,109
048.200	LOCATOR® angle measurement guide, stainless steal	41,54,65,90,109
048.204	RN synOcta® gold coping, bar, for synOcta 1.5, H 5.8mm, Ceramicor®	51
048.214	RN synOcta® Ti Coping, Bar for synOcta® 1.5, H 6mm, Ti	51
048.227	RN synOcta® Plastic Coping for synOcta® 1.5, for bridge/bar, H 10mm, POM	48,51
048.229	RN synOcta® Plastic Coping for synOcta® 1.5, for crown, H 10mm, POM	48
048.233	WN synOcta® Post, for temporary restoration, crown, H 9mm, Ti	57
048.234	WN synOcta® Post for temp. restoration, bridge, H 9mm, Ti	57
048.240	RN synOcta® Milling Cylinder, Crown, for synOcta® 1.5, Ø 6mm, H 8mm, Ti	48
048.241	RN synOcta® Milling Cylinder, Bridge, for synOcta® 1.5, Ø 6mm, H 8mm, Ti	48
048.243	WN synOcta® Plastic Coping for 048.606, for crown, H 7.35mm, POM	57
048.244	WN synOcta® Plastic Coping for 048.606, for bridge, H 7.35mm, POM	57
048.245	RN Plastic coping for Solid Abutment, crown, H 10mm, POM	50
048.246	RN Plastic coping for Solid Abutment, bridge, H 10mm, POM	50
048.247	WN Plastic coping for Solid Abutment, crown, H 7mm, POM	62
048.248	WN Plastic coping for Solid Abutment, bridge, H 7mm, POM	62
048.256	NNC Plastic coping, for 048.598, for crown, H 7mm, POM	35
048.257	NNC Plastic coping, for 048.598, for bridge, H 7mm, POM	35
048.258	NNC Plastic coping, for 048.621/ 048.622, for crowns & bridges, H 7.2mm, POM	35
048.263	NNC Plastic coping for Solid abutment, for crown, H 10.2mm, POM	38
048.264	NNC Plastic coping for Solid abutment, for bridge, H 10.2mm, POM	38
048.267	NNC Burn-out coping, for Variobase® Abutment, POM	37
048.267V4	NNC Burn-out coping, for Variobase® Abutment, POM	37
048.268	RN Burn-out coping, for Variobase® Abutment, POM	47
048.268V4	RN Burn-out coping, for Variobase® Abutment, POM	47
048.269	WN Burn-out coping, for Variobase® Abutment, POM	59
048.275V4	RN PLAN LOCATOR® Abutment, H 1mm, POM	67
048.276V4	RN PLAN LOCATOR® Abutment, H 2mm, POM	67
048.277V4	RN PLAN LOCATOR® Abutment, H 3mm, POM	67
048.278V4	RN PLAN LOCATOR® Abutment, H 4mm, POM	67
048.279V4	RN PLAN LOCATOR® Abutment, H 5mm, POM	67

货号	产品名称	页数
048.280V4	RN Novaloc® Plan Abutment, H 1-6 mm, POM	67
048.283V4	WN PLAN LOCATOR® Abutment, H 1mm, POM	68
048.284V4	WN PLAN LOCATOR® Abutment, H 2mm, POM	68
048.285V4	WN PLAN LOCATOR® Abutment, H 3mm, POM	68
048.286V4	WN PLAN LOCATOR® Abutment, H 4mm, POM	68
048.287V4	WN PLAN LOCATOR® Abutment, H 5mm, POM	68
048.313	NNC Basal screw, L 7.1mm, TAN	34,36,37
048.314	NNC Basal Screw, for cementable abutments, L 7.1mm, TAN	35
048.324	NNC Closure cap, Ø 2.6mm, H 0mm, Ti	14
048.324V4	NNC Closure cap, Ø 2.6mm, H 0mm, Ti	14
048.325	NNC Closure cap, Ø 4mm, H 1.5mm, Ti	14
048.325V4	NNC Closure cap, Ø 4mm, H 1.5mm, Ti	14
048.350	SCS Occlusal screw, L 4.4mm, Ti	46,48,51,60,138,141,142,150
048.350V4	SCS Occlusal screw, L 4.4mm, Ti	46,48,51,60,141,142,150
048.356	synOcta® Basal screw, L 6,7 mm, TAN	45,46,47,57,58,59,150
048.360V4	SCS Guide screw for RN synOcta® Bar coping, cannot be shortened, L 6 mm, Ti	150
048.371	RN closure screw, small, titanium	14,17,21
048.371V4	RN closure screw, small, titanium	14,17,21
048.373	RN closure screw, large, height 1.5 mm, titanium	14,17,21
048.373V4	RN closure screw, large, height 1.5 mm, titanium	14,17,21
048.375	WN closure screw, small, titanium	21
048.452	Screwdriver for titanium matrix, L 50mm, stainless steel	144
048.501	NNC Post for temporary restoration, crown, H 10.3mm, TAN	34
048.502	NNC Post for temporary restoration, bridge, H 10.3mm, TAN	34
048.505	NN Coping, H 9mm, Ti	117
048.527	NNC Analog for Solid abutment, for 048.547, H 14.3mm, TAN	38
048.528	NNC Analog for Solid abutment, for 048.548, H 15.8mm, TAN	38
048.529	NNC Analog for Solid abutment, for 048.549, H 17.3mm, TAN	38
048.540	RN Solid Abutment 6°, H 4mm, yellow, Ti	50
048.541	RN Solid Abutment 6°, H 5.5mm, grey, Ti	50
048.542	RN Solid Abutment 6°, H 7mm, blue, Ti	50
048.545	WN Solid Abutment 6°, H 4mm, green, Ti	62
048.546	WN Solid Abutment 6°, H 5.5mm, brown, Ti	62
048.547	NNC Solid abutment, H 4mm, yellow, Ti	38
048.548	NNC Solid abutment, H 5.5mm, gray, Ti	38
048.549	NNC Solid abutment, H 7mm, blue, Ti	38
048.550	NN Coping 15° angled, H 8,8 mm, Ti	117
048.551	NN Coping 20° angled, H 8,8 mm, Ti	117
048.571	RN synOcta® Meso Milling Cylinder, for synOcta® 1.5, crown, H 10mm, Ø 8mm, Ti	48
048.573	WN synOcta® Meso Milling Cylinder, for synOcta® 1.5, crown, H 10mm, Ø 10mm, Ti	60
048.579	NNC Impression cap, H 8mm, POM	38
048.581	NNC LOCATOR® Abutment, H 1mm, Ti alloy with TiN coating	39
048.582	NNC LOCATOR® Abutment, H 2mm, Ti alloy with TiN coating	39
048.583	NNC LOCATOR® Abutment, H 3mm, Ti alloy with TiN coating	39
048.584	NNC LOCATOR® Abutment, H 4mm, Ti alloy with TiN coating	39
048.585	NNC LOCATOR® Abutment, H 5mm, Ti alloy with TiN coating	39
048.586	NNC LOCATOR® Abutment, H 6mm, Ti alloy with TiN coating	39
048.587	NNC Positioning cylinder for 048.547, H 10.2mm, POM	38
048.588	NNC Positioning cylinder for 048.548, H 10.2mm, POM	38
048.589	NNC Positioning cylinder for 048.549, H 10.2mm, POM	38
048.592	NNC Gold abutment, crown, H 3.6mm, Ceramicor® / TAN / POM	36
048.593	NNC Gold abutment, bridge, H 3.6mm, Ceramicor®/ TAN / POM	36
048.598	NNC Cementable abutment, straight, for crowns/bridges, H 5.7mm, Ti	35
048.601	RN synOcta® 1.5 screw-retained Abutment, H 1.5mm, TAN	48,51
048.603	WN synOcta® 1.5 screw-retained Abutment, H 1.5mm, TAN	60
048.605	RN synOcta® Cementable Abutment, for crowns/bridges, H 5.5mm, Ti	45
048.606	WN synOcta® Cementable Abutment, H 5.5mm, Ti	57
048.608	WN synOcta® angled Abutment, Type A, angled 15°, H 5.5mm, Ti	58
048.609	WN synOcta® angled Abutment, Type B, angled 15°, H 5.5mm, Ti	58
048.610	RN synOcta® angled Abutment, long, Type A, angled 15°, H 6.7mm, Ti	46
048.611	RN synOcta® angled Abutment, long, Type B, angled 15°, H 6.7mm, Ti	46
048.612	RN synOcta® angled Abutment, short, Type A, angled 15°, H 5.7mm, Ti	46
048.613	RN synOcta® angled Abutment, short, Type B, angled 15°, H 5.7mm, Ti	46
048.615	RN synOcta® angled Abutment, long, Type A, angled 20°, H 6.7mm, Ti	46
048.616	RN synOcta® angled Abutment, long, Type B, angled 20°, H 6.7mm, Ti	46
048.617	RN synOcta® angled Abutment, short, Type A, angled 20°, H 5.7mm, Ti	46
048.618	RN synOcta® angled Abutment, short, Type B, angled 20°, H 5.7mm, Ti	46

货号	产品名称	页数
048.621	NNC Cementable abutment, angled 15°, type A, H 5.9mm, Ti	35
048.622	NNC Cementable abutment, angled 15°, type B, H 5.9mm, Ti	35
048.632	RN synOcta® gold coping, bridge, for synOcta 1.5, H 4.25mm, Ceramicor®/POM	48
048.633	RN synOcta® gold coping, crown, for synOcta 1.5, H 4.25mm, Ceramicor®/POM	48
048.635	NN Gold coping, cast-on, H 10mm, Ceramicor®/POM	117
048.638	WN synOcta® gold coping, bridge, for synOcta 1.5, H 4.75mm, Ceramicor®/POM	60
048.639	WN synOcta® gold coping, crown, for synOcta 1.5, H 4.75mm, Ceramicor®/POM	60
048.642	RN synOcta® gold Abutment, crown, H 4.3mm, Ceramicor®/POM	46
048.644	WN synOcta® gold Abutment, crown, H 14.1mm, Ceramicor®/POM	59
048.650	RN synOcta® Post for temporary restorations, for bridge, H 9mm, Ti	45
048.651	RN synOcta® Post for temporary restorations, for crown, H 9mm, Ti	45
048.654	RN Temporary Coping for Solid Abutment, bridge, H 8.5mm, PMMA	49
048.655	RN Temporary Coping for Solid Abutment, crown, H 8.5mm, PMMA	49
048.656	WN Temporary Coping for Solid Abutment, bridge, H 7.3mm, PMMA	61
048.657	WN Temporary Coping for Solid Abutment, crown, H 7.3mm, PMMA	61
048.658	NNC Temporary coping for Solid abutment, for crown, H 8.5mm, PMMA	38
048.659	NNC Temporary coping for Solid abutment, for bridge, H 8.5mm, PMMA	38
048.662	RN synOcta® Plastic Coping for 048.605, for bridge, H 7mm, POM	45
048.663	RN synOcta® Plastic Coping for 048.605, for crown, H 7mm, POM	45
048.666	WN synOcta® Plastic Coping, Bridge, H 10mm, POM	60
048.667	WN synOcta® Plastic Coping, Crown, H 10mm, POM	60
048.668	RN synOcta® Temporary Meso Abutment, H 10mm, PEEK/TAN	45
048.669	NN Temporary Coping, H 10mm, PEEK/TAN	117
048.670	Extension shell for angled RN synOcta® Abutment, H 5.5mm, POM	46
048.676V4	Plastic shoulder for RN synOcta® angled, angled 15° or 20°, H 4mm, POM	46
048.678	Plastic shoulder for WN synOcta® angled, angled 15°, H 4mm, POM	58
048.700	NNC Protective cap for Solid abutment, for 048.547, H 5.6mm, PEEK	38
048.701	NNC Protective cap for Solid abutment, for 048.548, H 7.1mm, PEEK	38
048.702	NNC Protective cap for Solid abutment, for 048.549, H 8.6mm, PEEK	38
048.709	NNC Variobase® Abutment, with screw, Ø 3.9mm, H 3.75mm, TAN	36,37
048.710	RN Variobase® Abutment, with screw, Ø 5.05mm, H 4mm, TAN	47
048.711	WN Variobase® Abutment, with screw, Ø 7mm, H 4.5mm, TAN	58,59
048.862	RN Impression Post Closed Tray, with 1 guide screw & 2 caps, Al/Ti/POM	44
048.864	RN Impression Post Closed Tray, abutment level, for RN synOcta® 1.5 Abutment, Al/Ti/POM	44
048.866	RN Impression Post Open tray, abutment level, for RN synOcta® 1.5 Abutment, Al/Ti	44
048.869	WN Impression Post Open Tray, abutment level, for WN synOcta® 1.5 Abutment, Al/Ti	56
048.870	WN Impression Post Closed Tray, with 1 guide screw & 2 caps, Al/Ti/POM	56
048.871	WN Impression post closed tray, abutment level, for WN synOcta® 1.5 Abutment, Al/Ti/POM	56
048.872	WN Implant Analog for synOcta® Abutment, abutment level, straight 0°, Ø6.5mm, stainless steel	56
048.904	Tissue Level PLAN Set, 278x124x33 mm, PC	66
048.905	Tissue Level PLAN cassette, empty, 278x124x33mm, polymer	66
048.916V4	RN synOcta® angled PLAN Abutment, short, Type A, angled 15°, POM	67
048.917V4	RN synOcta® angled PLAN Abutment, short, Type B, angled 15°, POM	67
048.918V4	RN synOcta® angled PLAN Abutment, short, Type A, angled 20°, POM	67
048.919V4	RN synOcta® angled PLAN Abutment, short, Type B, angled 20°, POM	67
048.921V4	RN synOcta® Transversal PLAN Abutment, POM	67
048.922V4	RN synOcta® angled PLAN Abutment, Type A, angled 15°, POM	67
048.923V4	RN synOcta® angled PLAN Abutment, Type B, angled 15°, POM	67
048.924V4	RN synOcta® angled PLAN Abutment, Type A, angled 20°, POM	67
048.925V4	RN synOcta® angled PLAN Abutment, Type B, angled 20°, POM	67
048.926V4	RN PLAN Solid Abutment, H 4mm, POM	67
048.927V4	RN PLAN Solid Abutment, H 5.5mm, POM	67
048.928V4	RN PLAN Solid Abutment, H 7mm, POM	67
048.929V4	RN synOcta® 1.5 screw-retained PLAN Abutment, POM	67
048.930V4	RN synOcta® Cementable PLAN Abutment, POM	67
048.931V4	WN synOcta® 1.5 screw-retained PLAN Abutment, POM	68
048.932V4	WN synOcta® Cementable PLAN Abutment, POM	68
048.933V4	WN synOcta® angled PLAN Abutment, Type A, angled 15°, POM	68
048.934V4	WN synOcta® angled PLAN Abutment, Type B, angled 15°, POM	68
048.938V4	WN PLAN Solid Abutment 6°, H 4mm, green, POM	68
048.939V4	WN PLAN Solid Abutment 6°, H 5.5mm, brown, POM	68
048.940V4	Bite registration aid, short, Ø 5 mm, H 8 mm, POM	44,56
048.941V4	Bite registration aid, long, Ø 5 mm, H 12 mm, POM	44,56
048.943V4	NNC PLAN Cementable abutment, straight, for crown/bridge, H 5.7mm, POM	66
048.944V4	NNC PLAN Cementable abutment, angled 15°, type A, H 5.9mm, POM	66
048.945V4	NNC PLAN Cementable abutment, angled 15°, type B, H 5.9mm, POM	66
048.946V4	NNC PLAN LOCATOR® abutment, H 1mm, POM	66

货号	产品名称	页数
048.947V4	NNC PLAN LOCATOR® abutment, H 2mm, POM	66
048.948V4	NNC PLAN LOCATOR® abutment, H 3mm, POM	66
048.949V4	NNC PLAN LOCATOR® abutment, H 4mm, POM	66
048.950V4	NNC PLAN LOCATOR® abutment, H 5mm, POM	66
048.952V4	NNC PLAN Solid Abutment, H 4mm, POM	66
048.953V4	NNC PLAN Solid Abutment, H 5.5mm, POM	66
048.954V4	NNC PLAN Solid Abutment, H 7mm, POM	66
048.960V4	NNC Bite registration aid, H 8mm, POM	34
048.961V4	NNC Bite registration aid, H 12mm, POM	34
049.071	Guide sleeve for implants Ø 4.1mm & Ø 3.5mm, Ø 4.2mm, L 5.6mm, stainless steel	138
049.076V4	X-ray reference sphere, Ø 5.0mm, stainless steel	126,149
049.082	Guide cylinder for implants Ø 3.3mm & Ø 3.5mm, Ø 3.6mm, L 10.5mm, stainless steel	130,138
049.083	Guide cylinder for implants Ø 4.1mm, Ø 4.2mm, L 10.5mm, stainless steel	138
049.167	Guide cylinder for implant Ø 4.8mm, Ø 4.9mm, L 10.5mm, stainless steel	138
049.177	NN Occlusal screw, L 5mm, Ø 1.8mm, TAN	117
049.181	SCS Occlusal screw for guiding sleeve, and post for temp. restorations, L 7.6 mm, TAN	45,57,138,150
049.182	SCS Guiding cylinder for NN implant Ø 3.3mm, L 5mm, Ti/stainless steel	138
049.810V4	Drill sleeve with collar, inner Ø 2.4mm, outer Ø 3.5mm, L 10mm, Ti	149
049.816V4	Pin for 049.810, Ø 2.2mm, L 16mm, Ti	149
049.817V4	Pin for 049.810, Ø 2.2mm, L 10mm, Ti	149
049.818V4	Stepped pin for 049.810, Ø 2.2mm/3.5mm, L 16mm, Ti	149
049.819V4	Pin for 049.810, Ø 3.5mm, L 16mm, Ti	149
049.850	Guide sleeve for implants Ø 3.3mm, Ø 3.6mm, L 4.5mm, stainless steel	138
049.851	Guide sleeve for implants Ø 4.8mm, Ø 4.9mm, L 6.4mm, stainless steel	138
140.055	Membranette	148
150.215	x-ray template for Standard and Standard Plus Implants	126
150.216	x-ray template for Bone Level Implants	126
80087	Explantation instrument for synOcta® implants, long, L 26.6mm, stainless steel	139
80103	Explantation instrument for synOcta® implants, short, L 20.6mm, stainless steel	139

Valid from July 01, 2018, until publication of the new catalog.

Trademarks

bonebuilder®, cerabone®, collacone®, collprotect®, Jason®, mucoder®m, maxgraft®, maxresorb® are registered trademarks of botiss Medical AG, Germany.
Ceramicor®, Dalbo®-Rotex® Matrix, Elitor®, Estetico®r and Pekkton® are registered trademarks of Cendres & Métaux SA, Switzerland.

CEREC® is a registered trademark of Sirona Dental Systems GmbH, Germany.

Cronidur® 30 is a registered trademark of Energietechnik Essen GmbH, Germany.

Dental Wings® is a registered trademark of Dental Wings Inc., Canada.

Dolder® is a registered trademark of Prof. Eugen Dolder, exclusively licensed to Cendres + Metaux Holding SA, Switzerland.

IPS e.max® is a registered trademark of Ivoclar Vivadent AG, Liechtenstein.

LOCATOR® is a registered trademark of Zest IP Holdings LLC, USA.

Novaloc® is a registered trademark of Valoc AG, Switzerland.

PlanMill® is a registered trademark of Planmeca Oy, Finland.

VITA CAD-Temp® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Germany.



请扫描上方二维码，关注士卓曼微信公众号
海量信息，尽在掌握
www.straumann.cn

士卓曼（北京）医疗器械贸易有限公司

地址：北京市朝阳区东三环北路27号嘉铭中心B座3层303号
100020
电话：+86 10 5775 6555
传真：+86 10 5775 6556

STCN452.200 - 201806

